

Instructions for Use - **Safe-T-Flush™ Pre-Filled Saline Syringe**

Indications for Use

- The Safe-T-Flush™ Pre-Filled Saline Syringe is intended to be used for the flushing of the central and peripheral vascular access devices.

Product Description

- The Safe-T-Flush™ Pre-Filled Saline Syringe is a three-piece, sterile, single use syringe with a 6% conical taper luer connector, prefilled with 0.9% Sodium Chloride (0.9% NaCl) USP, and sealed with a nozzle cap.
- The Safe-T-Flush™ Pre-Filled Saline Syringe is provided with a sterile fluid path, sterilized via steam or dry heat.

Product Components

- Components:- barrel, plunger, piston, nozzle cap and 0.9% Sodium Chloride (0.9% NaCl), conforming with USP standards.

Product Specification

- 10ml.

Sterilization Method

- Steam or dry heat sterilization.

Shelf Life

- 3 years.



















Storage

- Device should be stored between the temperature of 59-86°F / 15-30°C.

Instructions for Use

- Step 1: Remove the Safe-T-Flush™ Pre-Filled Saline Syringe from the packaging.
- Step 2: Push the plunger upward to release the resistance between the piston and the barrel. Note: during this step do not unscrew the nozzle cap.
- Step 3: If the product is connected to an intravenous infusion system with a needleless infusion connector, please connect it to the needleless infusion connector directly. If the product is connected to a traditional infusion connector, that requires a steel needle connection, please connect it to a standard or safety needle according to your establishment's protocol.
- Step 4: Connect the product to the connector valve, or needleless system, and flush relevant to the catheter design.
- Step 5: Used Safe-T-Flush™ Pre-Filled Saline Syringe should be disposed into an approved sharps container.

Symbols used for Medical Device Labels

	RX Only		Sterilized using Steam or Dry Heat
	Single Use		Do Not Resterilize
	Latex Free		Temperature Limit
	Do Not Use if the Package is Damaged		Consult Instructions For Use
	Lot Number		Reference
	Keep Dry		Non-Pyrogenic
	Manufacturer		Date of Manufacture
	Use By Date		Stacking Limit by Mass
	This Way Up		Keep Away from Sunlight