Bard Access Systems

BARD

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Attention, See Instructions For Use

Attention, voir le mode d'emploi
Achtung! Gebrauchsanweisung beachten
Attenzione, leggere attentamente le istruzioni per l'uso
Atención, consulte las instrucciones de uso
Raadpleeg de gebruiksaanwijzing
Atenção, Consulte as Instruções de Utilização
Προσοχή, συμβουλευθείτε τις Οδηγίες Χρήσης
OBS! Se brugsanvisningen
Obs! Se bruksanvisningen
Huomio! Lue käyttöohjeet.
Forsiktig, se bruksanvisningen
Figyelem, lásd a használati utasítást
Uwaga, przeczytać instrukcję przez użyciem
Pozor, čtěte návod k použití.
Dikkat, Kullanım Yönergelerini İnceleyin
Внимание! См. инструкции по применению
주의, 사용 지침 참조
注意，請參閱使用說明書

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TYPE

ΤÍP

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Non-Coring Needle & Non-Coring Needle Winged Infusion Sets

Instructions For Use

Description

Non-Coring Needle:

A stainless steel needle designed to be introduced into a silicone septum of a subcutaneously placed implanted port for infusion of medications, I.V. fluids, parenteral nutrition solutions, blood products and for the withdrawal of blood samples. The Non-Coring needle is designed not to core the silicone septum.

Winged Infusion Set:

Contains a Non-Coring needle with a pre-attached extension set with clamp and with or without an injection Y site.

Indications for Use

The Bard Access Systems Non-Coring Needle & Non-Coring Needle Winged Infusion Sets are indicated for the accessing of subcutaneously placed implanted ports with silicone septums. See **Bard Access Systems Implanted Port Instructions for Use** on specific indications, contraindications, precautions, cautions, warnings and procedures for implanted ports.

Warnings

- Intended for Single Use. **DO NOT REUSE.** Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.
- Resterilization of incompletely cleaned devices may not be effective. Any device that has been contaminated by blood should not be reused or resterilized.
- After use, this product may be a potential biohazard. Handle and discard in accordance with accepted medical practice and applicable local, state and federal laws and regulations.
- Fully tighten all connections before use. Failure to attach caps can result in an embolism or bleeding.

Cautions

- Rx Only
- Carefully read and follow all instructions prior to use.
 - Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
 - Only qualified healthcare practitioners should insert, manipulate and remove these devices.
 - It is necessary to know the maximum psi (kPa) ratings for the implanted port which you are accessing. High pressures or use with power injectors may cause leakage or damage.

Precautions

- Follow all contraindications, warnings, cautions, precautions and instructions for all infusates as specified by its manufacturer.
- Examine package carefully before opening to confirm its integrity and that the expiration date has not passed. The device is supplied in a sterile package and is non-pyrogenic. Do not use if package is damaged, opened or the expiration date has passed. Sterilized by ethylene oxide. **DO NOT RESTERILIZE.**

Site Preparation and Port Access

Always inspect and aseptically prepare the injection site prior to accessing the port.

Equipment:

- Alcohol wipe
- Antiseptic swabs
- Sterile gloves
- Non-Coring Needle or Non-Coring Needle Winged Infusion Set
- 10 ml syringe filled with normal sterile saline
- Extension set with clamp (if not using Winged Infusion Set)

Procedure:

1. Explain procedure to patient. Warn of needle prick sensation.
2. Wash hands thoroughly.
3. Put on sterile gloves.
4. Attach non-coring needle to extension set (if necessary) and 10 ml syringe filled with sterile normal saline. Expel all air and clamp extension.
5. Paint area with alcohol wipe starting at the port and working outward in a spiral motion over an area 10 - 13 cm in diameter.
6. Repeat step 5 with antiseptic swabs three times.
7. Locate port septum by palpation.
 - a. Locate base of port with non-dominant hand.
 - b. Triangulate port between thumb and first two fingers of non-dominant hand.
8. Insert needle perpendicular to port septum. Advance needle through the skin and septum until reaching bottom of reservoir.
9. Verify correct needle placement by blood aspiration.

Note: Follow established hospital or institutional policy for changing I.V. tubing and accessing cannula. Center for Disease Control (CDC) guidelines recommend that the I.V.

tubing be changed every 48 hours. Oncology Nursing Society (ONS) guidelines recommend implanted port access needles be changed every 7 days.

Deaccessing

To reduce the potential for blood back-flow into the catheter tip and possible catheter clotting, always remove a non-coring needle slowly, while injecting the last 0.5 ml of solution. Stabilize the port with two fingers during needle withdrawal.

Bolus Injection

Equipment:

- Non-Coring Needle or Non-Coring Needle Winged Infusion Set
- 10 ml syringe filled with sterile normal saline
- Extension set with clamp (if not using Winged Infusion Set)

Procedure:

1. Prepare injection site and aseptically locate and access port as defined above.
2. Flush port with 10 ml sterile normal saline. Clamp the extension set and remove syringe.
3. Connect syringe containing the drug to extension set. Release clamp and administer injection.
4. Examine the injection site for signs of extravasation; if noted, immediately discontinue the injection and initiate appropriate intervention.
5. When the injection is completed, clamp the extension set.
6. Flush port after each injection with 10 ml of sterile normal saline to help prevent interaction between incompatible drugs.
7. Perform heparin lock procedure on ports with **Hickman*** (open-ended) catheters and saline lock procedure on ports with **Groshong*** (valved) catheters.

Note: The needle hub should not be left open to air while it is in the port. Do not manipulate the needle once it is in the septum.

Continuous Infusion

Equipment:

- Non-Coring Needle or Non-Coring Needle Winged Infusion Set
- Prescribed I.V. solution
- Extension set with clamp (if not using Winged Infusion Set)
- 10 ml syringe filled with sterile normal saline
- I.V. pole
- I.V. pump (if ordered)
- Transparent dressing
- Antibacterial ointment
- 2 in x 2 in (5 cm x 5 cm) gauze pads

Procedure:

1. Prepare injection site and aseptically locate and access port as defined above.
2. Flush port with 10 ml normal sterile saline. Clamp the extension set and remove syringe.
3. Apply antibacterial ointment to injection site and place a rolled gauze pad under needle hub. Secure needle with transparent dressing to help prevent inadvertent dislodgement.
4. Open clamp and flush port with sterile normal saline. Clamp extension set and remove syringe.
5. Connect fluid delivery system.

Note: If a pump is used for infusion, turn it on before opening clamp. To provide additional security during pump infusion, tape all tubing connections. Pumps must incorporate a functional automatic pressure limiting switch which will shut pump off before pressure exceeds 25 psi (172 kPa).

-
6. Release clamp and initiate infusion. Examine the infusion site for signs of extravasation; if noted, or if patient experiences pain, immediately discontinue infusion and initiate appropriate intervention.
 7. When infusion is completed, clamp extension set and then remove the fluid delivery system.
 8. Flush port after each injection with 10 ml of sterile normal saline to help prevent interaction between incompatible drugs.
 9. Perform heparin lock procedure on ports with **Hickman*** (open-ended) catheters and saline lock procedure on ports with **Groshong*** (valved) catheters.

For other infusion uses, please refer to the implanted port instructions for use.

An issued or revision date for these instructions is included for the user's information. In the event two years elapsed between this date and product use, the user should contact **Bard Access Systems** to see if additional product information is available. Revised date: August 2013

***Bard**, **Groshong**, and **Hickman** are trademarks and/or registered trademarks of C. R. Bard, Inc.

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All rights reserved.
© 2013 C. R. Bard, Inc.
Bard Access Systems, Inc.