### Catheter Stabilization Device

Read carefully before use

Safety and Efficacy Considerations:

- ② Single use only.
- ⊘ Do not sterilize.

Rx Only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

**STERILE** **EO** Sterile by ethylene oxide. Do not use if package is damaged. Fluid path (if present) is nonpyrogenic.

**REF** Product code, **LOT** lot number and expiration date information – see packaging.

Do not alter the StatLock® stabilization device or components. Procedure must be performed by trained personnel with knowledge of anatomical landmarks, safe technique and potential complications.

LATEX-FREE

#### Indications for Use:

The StatLock® stabilization device is for compatible medical tubes and catheters. Multiple pad and retainer designs are available.

Choose the StatLock® stabilization device indicated for the tube or catheter to be secured.

#### Contraindications:

Known tape or adhesive allergies.

#### Warnings and Precautions:

1. Do not use the StatLock® stabilization device where loss of adherence could occur, such as with a confused patient, diaphoretic or nonadherent skin, or when the access device is not monitored daily.
2. Observe universal blood and body fluid precautions and infection control procedures, during application and removal of the StatLock® stabilization device.
3. **SUTURE THE STATLOCK® STABILIZATION DEVICE PAD TO SKIN IF DESIRED OR DEEMED NECESSARY.**
4. Avoid contact with alcohol or acetone, both can weaken bonding of components and pad adherence.
5. Minimize catheter/tube manipulation during application and removal of the StatLock® stabilization device.
6. A StatLock® stabilization device Luer-lock connector must be used to secure central venous and arterial catheters.
7. Always apply adhesive strip to central venous and arterial catheters at or near insertion site when using StatLock® stabilization device.
8. Remove oil and moisturizer from targeted skin area.
9. StatLock® stabilization device pad adherence and catheter/tube position should be routinely inspected.
10. Orient StatLock® stabilization device so arrow points toward catheter tip.
11. StatLock® stabilization device should be replaced at least every 7 days.
12. This is a single use device. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.
13. Do not re-sterilize. The sterility of the single use device is not guaranteed following re-sterilization because of an indeterminate degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Re-sterilization may compromise the structural integrity, essential material and/or design characteristics and may lead to an unpredictable loss of functionality and/or device failure.

Use sterile technique for StatLock® stabilization device application and removal. See illustrations for specific details.

Extension Set: If present, Prime/Flush/Connect per accompanying instructions.

#### StatLock® stabilization Device Application:

1. **PREP:** Prepare targeted area using alcohol, to remove oil and moisturizer, and skin prep for skin protection and enhanced adhesion. **ALLOW TO DRY COMPLETELY AS CONTACT WITH ALCOHOL CAN WEAKEN BONDING OF COMPONENTS AND PAD ADHERENCE. ALWAYS CONNECT THE STATLOCK® STABILIZATION DEVICE TO CATHETER BEFORE ADHERING PAD TO SKIN.**
2. **PRESS:** Before applying the StatLock® stabilization device pad to skin, secure tube or catheter to StatLock® retainer.
3. **PEEL:** Peel off the StatLock® stabilization device paper backing, and
4. **PLACE:** Place the StatLock® stabilization device on targeted skin one side at a time.

#### Note:

- The StatLock® stabilization device pad can be adhered after X-ray confirmation of proper catheter/tube position.
- Apply adhesive strip to central venous and arterial catheters for added security during StatLock® stabilization device dressing changes.

#### The StatLock® stabilization device general removal:


Carefully remove overlying dressing using stretch technique and sutures (if applicable).

1. **DISENGAGE** catheter/tube from StatLock® stabilization device retainer.
2. **DISSOLVE** StatLock® stabilization device adhesive with alcohol swab while gently lifting StatLock® stabilization device pad.
3. **DOCUMENT** StatLock® stabilization device dressing change in the patient chart.
4. **DISPOSE** of all equipment in appropriate containers.

#### Extension Set

#### Indications for Use:

For use as an extension to the IV fluid delivery system. StatLock® stabilization device extension set male Luer adapter is specifically designed to secure the extension set when used as part of the StatLock® stabilization device system.

#### Directions for Use:

Use aseptic technique.

1. Tighten all connections.
2. Remove Protector Caps (A) as required.
3. Connect Female Port (B) to fluid delivery line or device and prime/flush.
4. Connect Male Port (C) to catheter hub.
5. Open Slide Clamp (D) for unimpeded flow.
6. Swab septum of injection site and SmartSite™ device (if present) with preferred antiseptic.

#### General Precautions:

1. Avoid air trapped in set. Puncturing tube can cause air embolism or leakage.
2. Avoid contact with alcohol, as contact with alcohol can weaken bonding of components and pad adherence.

#### Options / Precautions:

Needleless Valve (E) and Y-Site w/ Needleless Valve (F)

1. Connect fluid line or device to center of valve using Luer Lock connections only.
2. Do not use with a blunt tip cannula or conventional needle. If a needle must be used, attach a PRN adapter and use as a standard injection site.

T-Connector (G)

Stopcock (H)

1. Position handle over port you wish to close.
2. Open port for fluid delivery.
3. Close and cap available port when not in use.

Standard Y-Site (I)

1. Use with a conventional needle only. Do not use with a blunt tip cannula.

#### NOTES

1. Tubing is non-DEHP PVC
2. Dispose used set in appropriate container.
3. Replace extension set and ports per CDC guidelines.

Micro	No Valve	W/Valve	Fill Volume (mL)
12" (304.8mm) ext set	0.40	NA	
4" (101.6mm) ext set	0.18	NA	
6" (152.4mm) ext set with SmartSite®	0.25	0.34	
6" (152.4mm) ext set with T-Connector and SmartSite®	0.18	0.27	
8" (203.2mm) ext set with 3-way stopcock	0.36	0.55	
8" (203.2mm) ext set with bonded 3-way stopcock	NA	0.55	
9" (228.6mm) bifurcated ext set with SmartSite®	0.57	0.75	

#### Macro

6" (152.4mm) ext set with Smart Y® and SmartSite®	0.90	1.06
6" (152.4mm) ext set with SmartSite®	0.90	0.99

Bard and StatLock are trademarks and/or registered trademarks of C. R. Bard, Inc.

All other trademarks are the property of their respective owners.  
© 2014 C. R. Bard, Inc. All Rights Reserved.

Manufacturer:  
Bard Access Systems, Inc.  
605 North 5600 West  
Salt Lake City, UT 84116 USA  
800-443-5505, 801-522-5000  
www.bardaccess.com  
Assembled in Mexico

**EC REP** Authorized Representative in the European Community  
Bard Limited  
Forest House, Brighton Road  
Crawley, West Sussex  
RH11 9BP UK

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PAD SHAPES <sup>1</sup>	StatLock <sup>®</sup>	AN APPLICATION TECHNIQUE <sup>2</sup> (USE SKIN PREPARATION PAD)					THE CATHETER SECURED <sup>3</sup>					A REMOVAL TECHNIQUE <sup>4</sup> (USE ALCOHOL)					TRANSLATIONS		NOTIFICATION TO CLINICIANS <sup>5</sup>	
	IV Plus																TRADUCCIÓN KÄÄNKES TRADUCTION TRADUČIO ÖVERSÄTTNING ÜBERSETZUNG METAPPAZEIZ VERTALING TRADUKZIONE ÖVERSÄTTELSE TLUMACZENIA FORDÍTÁSOK PŘEKLADY ÇEVİRILIR	S. OR O ELLER OU OU OU TAI ODER H OF OPPURE ELLER LUB VAGY NEBO VEVA	<p><b>*FOR USA ONLY. SÓLO PARA LOS EE.UU. JNK SIESTA USA ANIDA. POUR LES ÉTATS-UNIS UNIQUEMENT. APENAS PARA EUA. FOR USA ENDA. NUR FÜR DIE USA. MOVO YO HIA. ALLEEN VOOR DE V.S. PER USA SOLTANTO. KUN USA BARE. DOTYCZY WYŁĄCZNIE WYROBÓW UŻYTKOWANYCH W USA. KIZAROLAG AZ USA ESETEBEN. POUIZE PRO USA. SADECE A.B.D. İÇİN.</b></p> <p><b>NOTIFICATION TO CLINICIANS ON SHARPS INJURY PREVENTION</b></p> <p>Healthcare workers who use or may be exposed to needles are at increased risk of needlestick injury. Such injuries can lead to serious or fatal infections with bloodborne pathogens such as hepatitis B virus, hepatitis C virus, or human immunodeficiency virus (HIV).</p> <p>In order to prevent bloodborne pathogens exposure, clinicians are advised to utilize both Primary and Secondary sharps injury prevention techniques.</p> <p>Primary prevention practices and technologies eliminate the need to introduce sharps into the workplace, thereby reducing the total number of sharps used.</p> <p>Secondary prevention practices and technologies render safer those sharps that must be introduced into the workplace.</p> <p>The combination of Primary and Secondary Prevention maximizes healthcare worker safety with respect to blood borne pathogens exposure from sharps injury.</p> <p>[Supporting excerpts from NIOSH and OSHA documents]:</p> <p>"Employers of health care workers should implement the use of improved engineering controls to reduce needlestick injuries: Eliminate the use of needles where safe and effective alternatives are available..."</p> <p>"Health care workers should take the following steps to protect themselves and their fellow workers from needlestick injuries: Avoid the use of needles where safe and effective alternative are available..."</p> <p>– NIOSH Alert: "Preventing Needlestick Injuries in Health Care Settings" NOVEMBER 1999</p> <p>"Where engineering controls will reduce employee exposure either by removing, eliminating or isolating the hazard, they must be used."</p> <p>– OSHA Compliance Directive on Occupational Exposure to Bloodborne Pathogens NOVEMBER 1999</p>	
	IV																VERSCHIEDENE POLSTERFORMEN ZKRMÁVATA STROBAMATIN VORMEN VAN CESTOPRUGSSENS FORME DEI CUSONETTI ADESIVI PUDEMODELLER KISZALTALY FODRUKADEK PÁRNA FORMÁK TVÁRY TAMPONŮ PED ŞEKİLLERİ	1. PAD SHAPES FORMAS DEL PARCHÉ KUODFORMER FORMES DE TAMPONS FORMAS DA COMPRESSA TYNYNY MUODOT		
	IV T-Connector																2. AN APPLICATION TECHNIQUE (USE SKIN PREPARATION PAD) UNA TÉCNICA DE APLICACIÓN (USAR EL PARCHÉ DE PREPARACIÓN DE LA PIEL) EN APLICACIÓNSMETOD (ANVÁND HUDPREPARERINGSKUDDE) TECHNIQUE D'APPLICATION (UTILISER LE TAMPON DE PRÉPARATION DE LA PEAU) UMA TÉCNICA DE APLICACIÃO (USAR A COMPRESSA DE PREPARAÇÃO DA PELE) KINN TYSMENETELMÄ (KÄYTÄ IHON VALMISTELUTYNYÄ) EINE BEFESTIGUNGSTECHNIK (HAUTVORBEREITUNGSWATTE VERWENDEN) MAA TEKNIKI EN APMETODE (PHEKIMONHOHETE TAMPONIN AEPMATIKHE PPOETOMAZIAE) METHODE VAN ANBRENGEN (DOESKIE MET HUIDPREPARATIEMIDDEL GEBRUIKEN) UNA TÉCNICA DI APPLICAZIONE (CON L'USO DEL TAMPONE PER LA PREPARAZIONE DELLA PELLE) PÁSETNINGSTECHNIK (BRUG SERVIET TIL KLÆDNING AF HUD) TECHNIKA ZAKŁADANIA (WYKORZYSTANIE WACIKA DO ODKAZANIA SKÓRY) EGY ALKALMAZOTT TECHNIKA (HASZNÁLJON BŐRÉLŐKESZÍTŐ PÁRNÁT) TECHNIKA APLIKACE (POUŽÍTE TAMPON PRO PŘÍPRAVU KŮŽE) BIR UYGULAMA TEKNİĞİ (CILT HAZIRLAMA PEDI KULLANINIZ)	2. AN APPLICATION TECHNIQUE (USE SKIN PREPARATION PAD) UNA TÉCNICA DE APLICACIÓN (USAR EL PARCHÉ DE PREPARACIÓN DE LA PIEL) EN APLICACIÓNSMETOD (ANVÁND HUDPREPARERINGSKUDDE) TECHNIQUE D'APPLICATION (UTILISER LE TAMPON DE PRÉPARATION DE LA PEAU) UMA TÉCNICA DE APLICACIÃO (USAR A COMPRESSA DE PREPARAÇÃO DA PELE) KINN TYSMENETELMÄ (KÄYTÄ IHON VALMISTELUTYNYÄ) EINE BEFESTIGUNGSTECHNIK (HAUTVORBEREITUNGSWATTE VERWENDEN) MAA TEKNIKI EN APMETODE (PHEKIMONHOHETE TAMPONIN AEPMATIKHE PPOETOMAZIAE) METHODE VAN ANBRENGEN (DOESKIE MET HUIDPREPARATIEMIDDEL GEBRUIKEN) UNA TÉCNICA DI APPLICAZIONE (CON L'USO DEL TAMPONE PER LA PREPARAZIONE DELLA PELLE) PÁSETNINGSTECHNIK (BRUG SERVIET TIL KLÆDNING AF HUD) TECHNIKA ZAKŁADANIA (WYKORZYSTANIE WACIKA DO ODKAZANIA SKÓRY) EGY ALKALMAZOTT TECHNIKA (HASZNÁLJON BŐRÉLŐKESZÍTŐ PÁRNÁT) TECHNIKA APLIKACE (POUŽÍTE TAMPON PRO PŘÍPRAVU KŮŽE) BIR UYGULAMA TEKNİĞİ (CILT HAZIRLAMA PEDI KULLANINIZ)		
	CV Plus/PICC Plus																3. THE CATHETER SECURED EL CATERER FIJADO KATETERN FASTSATT LE CATHÈTER EST FERMEMENT FIXÉ FIXAÇÃO DO CATERER KINN TETTY KATETRI DER STABILISIERTE KATHETER ETEPEDMENEDE KAGETHPAE VASTGEZETTE KATHETER IL CATERERE STABILIZZATO DET SIKREDE KATETER ZABEZPIECZONY CEWNIK A KATERER RÖZG TETT ZAJISTENÝ KATETR KATETER SABITLENMİŞ	3. THE CATHETER SECURED EL CATERER FIJADO KATETERN FASTSATT LE CATHÈTER EST FERMEMENT FIXÉ FIXAÇÃO DO CATERER KINN TETTY KATETRI DER STABILISIERTE KATHETER ETEPEDMENEDE KAGETHPAE VASTGEZETTE KATHETER IL CATERERE STABILIZZATO DET SIKREDE KATETER ZABEZPIECZONY CEWNIK A KATERER RÖZG TETT ZAJISTENÝ KATETR KATETER SABITLENMİŞ		
	CV/PICC																4. A REMOVAL TECHNIQUE (USE ALCOHOL) UNA TÉCNICA DE RETIRADA (USAR ALCOHOL) EN BORTTÄGNINGSMETOD (ANVÁND ALKOHOL) TECHNIQUE DE RETRAIT (UTILISER DE L'ALCOOL) UMA TÉCNICA DE REMOÇÃO (UTILIZE ALCOOL) POSTMENE TELMÄ (KÄYTÄ ALKOHOLIA) EINE ENTFERNUNGSTECHNIK (ALKOHOL VERWENDEN) MAA TEKNIKI ARAPIEZHKE (PHEKIMONHOHETE AKKOHOLIN) METHODE VAN VERWIJDEREN (ALCOHOL GEBRUIKEN) UNA TÉCNICA DI RIMOZIONE (CON L'USO DI ALCOOL) FJERNELSESTECHNIK (BRUG SPRIT) TECHNIKA ZDEJMOWANIA (WYKORZYSTANIE ALKOHOLI) ELTÁVOL TÓ TECHNIKA (HASZNÁLJON ALKOHOLT) TECHNIKA ODSTRANĚNÍ (POUŽÍTE ALKOHOL) BIR ÇIKARTMA TEKNİĞİ (ALKOL KULLANINIZ)	4. A REMOVAL TECHNIQUE (USE ALCOHOL) UNA TÉCNICA DE RETIRADA (USAR ALCOHOL) EN BORTTÄGNINGSMETOD (ANVÁND ALKOHOL) TECHNIQUE DE RETRAIT (UTILISER DE L'ALCOOL) UMA TÉCNICA DE REMOÇÃO (UTILIZE ALCOOL) POSTMENE TELMÄ (KÄYTÄ ALKOHOLIA) EINE ENTFERNUNGSTECHNIK (ALKOHOL VERWENDEN) MAA TEKNIKI ARAPIEZHKE (PHEKIMONHOHETE AKKOHOLIN) METHODE VAN VERWIJDEREN (ALCOHOL GEBRUIKEN) UNA TÉCNICA DI RIMOZIONE (CON L'USO DI ALCOOL) FJERNELSESTECHNIK (BRUG SPRIT) TECHNIKA ZDEJMOWANIA (WYKORZYSTANIE ALKOHOLI) ELTÁVOL TÓ TECHNIKA (HASZNÁLJON ALKOHOLT) TECHNIKA ODSTRANĚNÍ (POUŽÍTE ALKOHOL) BIR ÇIKARTMA TEKNİĞİ (ALKOL KULLANINIZ)		
	CV Ultra																RO	RO		
	PICC II																RO	RO		
	first StatLock																RO	RO		
	B. Braun																RO	RO		
	Winged Catheters																RO	RO		
	Dialysis																RO	RO		
	Dialysis II																RO	RO		
	Arterial/Arterial Plus																RO	RO		
	Arrow Arterial																RO	RO		
	Epidural																RO	RO		
	Arrow Epidural																RO	RO		
	Universal Plus																RO	RO		
	Universal																RO	RO		
	Multi Purpose/General																RO	RO		
	J-Loop																RO	RO		

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Bard ve StatLock C. R. Bard, Inc.'ın veya bağlı kuruluşunun ticari markaları ve/veya tescilli ticari markalarıdır.
© 2014 C. R. Bard, Inc. Tüm Hakları Saklıdır.