INDICATIONS FOR USE: The **Huber Plus™** Non-coring Needle Safety Infusion Set is a safety IV administration set with a non-coring right angle Huber needle used to access surgically implanted vascular ports. The **Huber Plus™** Set is used to administer fluids or to withdraw blood. The **Huber Plus™** Set facilitates safe removal of the needle by encapsulating the needle within the wings to help prevent needlestick injuries when using the device for vascular port access.

CAUTIONS AND WARNINGS:

- Intended for Single Use. DO NOT REUSE. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.Ⓢ
- When handling the **Huber Plus™** Set, always follow universal (infection control) precautions as specified by the CDC (U.S.A) or local equivalent.
- Leakage may occur with Luer fitting connections over time or when disconnecting components—use universal precautions to prevent exposure to bloodborne pathogens.
- Verify that clamp is open and that there are no kinks in the IV line before priming or initiating flow.
- If the folding wing safety feature fails to completely cover and lock the contaminated needle during removal, use care in handling the infusion set and immediately dispose of it into a sharps container.
- Do not clean, re-sterilize or reuse the **Huber Plus™** Set.
- Rx Only** CAUTION: Federal law restricts this device to sale by or on the order of a physician.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

INSTRUCTIONS FOR USE (refer to figures 1—4):

Always inspect for integrity of patient's skin over the port septum prior to accessing site. Avoid accessing port in non-intact skin.

INSERTION:

- Prepare the port site for sterile accessing following institution protocol.
- Prime infusion set using aseptic technique following institution protocol.
- Insertion method:
 - Grasp needle guard with non-dominant hand.
 - Grasp vertical fin between thumb and middle finger of dominant hand while placing tip of the index finger onto the front finger pad located directly over the front of the needle (figure 1).
 - Remove needle guard with non-dominant hand.
 - Locate port septum with non-dominant hand and stabilize.
 - Insert needle perpendicular into vascular access port by pressing down with index finger while aiming for the center of the port septum. Advance needle through the skin and septum until reaching the bottom of reservoir.
- Orient the wings in single or dual ports as shown in (figure 2).
- Verify correct needle placement by aspirating for blood.
- Flush per institution protocol and close clamp on the extension tubing.
- Secure needle to site following institution protocol. Do not manipulate the needle once it is in the septum.

REMOVAL:

- Follow institution protocol for irrigation and prepping prior to preparing for needle removal.
- With non-dominant hand, stabilize the port.
- Do not grasp vertical fin during removal. This could result in inactivation of the safety feature and possible injury. Grasp the outer edges of the wing with thumb and side knuckle of index finger (figure 4A) or thumb and index finger tip (figure 4B) and lift and squeeze fingers together.
- Wings will fold together and lock around the needle as it is removed from the patient (figure 5). A snapping sound should be heard indicating that the wings are fully closed and locked around the needle.
- After removal, handle the **Huber Plus™** Set by holding the extension tubing behind the wing assembly. Carefully discard the **Huber Plus™** Set into a sharps container. Change set according to Centers for Disease Control (CDC) guidelines (U.S.A) and/ or local or institution guidelines.

SmartSite™ NEEDLE-FREE INJECTION SITE:

“Y” versions of the **Huber Plus™** Sets are supplied with **SmartSite™** Y needle-free injection sites which should only be used with standard male Luer fittings on syringes or administration sets. Aseptically prepare the **SmartSite™** Injection Site entry port per institutional protocol before attaching a male Luer fitting. Do not insert hypodermic or blunt needles into the **SmartSite™** Injection Site.

Bard and HuberPlus are trademarks and/or registered trademarks of C. R. Bard, Inc. All other trademarks are the property of their respective owners.

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Bard et HuberPlus sont des marques commerciales et/ou des marques commerciales déposées de C. R. Bard, Inc. Toutes les autres marques appartiennent à leurs propriétaires respectifs. © 2017 C. R. Bard, Inc. Tous droits réservés.
© 2017 C. R. Bard, Inc. Všechna práva vyhrazena.