SafeStep*

HUBER NEEDLE SET

English

Instructions for Use

6

Français

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8

Deutsch

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10

Italiano

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12

Español

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14

Nederlands

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16

Português

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18

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
HUBER NEEDLE SET

Indications For Use:

The SafeStep* Huber Needle Set is a device intended for insertion into the septum of a subcutaneously implanted port and for the infusion of fluids into the port. The SafeStep* Huber Needle safety feature is manually activated during needle removal, and is designed to aid in the prevention of accidental needlesticks.

Warnings:

- Verify needle length is correct based on port reservoir depth, tissue thickness and the thickness of any dressing beneath the bend of the needle; if too long, needle and/or port may be damaged at insertion; if too short, needle may not completely pierce port septum, and medication may be delivered into surrounding tissue and/or needle may be blocked.
- Fully tighten all connections, Y-site end caps, or needleless connectors before use. Failure to attach an end cap after removing a male Luer locking end cap can result in an embolism or bleeding.
- Needleless Y-injection site valve may leak slightly at certain pressures resulting in a small bead of fluid on the valve.
- The textured handle and base are not connected. The base will move along the length of the needle.
- Failure to use the device correctly, when removing needle from port site could result in the needle tip re-emerging from the base, resulting in an accidental needle stick with a contaminated needle. A needle stick with a contaminated needle may cause infectious disease.
- Handle and discard in accordance with accepted medical practice and all applicable regulations. After use, this product is a potential biohazard.
- Intended for single use. DO NOT REUSE. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.
- Do not alter the device.

Cautions:

- Care must be taken to avoid accidental needle sticks.
- Only qualified health care practitioners should insert, manipulate and remove these devices.
- Carefully read all instructions prior to use.
- Follow all instructions, contraindications, warnings, cautions and precautions for all infusates, ports, IV pumps, IV sets and needleless systems, as specified by the CDC (USA) or local equivalent.
- Do not use this product if it appears damaged or if the package has been previously opened or damaged.
- Federal law restricts the device to sale by or on the order of a physician.
- High pressure or use with power injectors may cause leakage or damage.
- Confirm correct needle placement in the port reservoir by aspiration of blood before infusion of any substance. If there is doubt regarding proper needle placement, perform a radiographic dye procedure to confirm placement per institutional protocol.
- Do not remove and reinsert the needle into the port.
- Avoid excessive manipulation once the needle is in the port.
- Do not access the needleless Y-injection site valve with a needle. Valve puncture may result in an embolism.

Instructions For Use:

1. Prepare port site for sterile needle insertion. (1)
2. Attach a 10 mL syringe containing normal saline to the proximal Luer lock connection of the safety Huber needle infusion set. Grasp the textured handle and remove the needle cover. Prime and flush the infusion set. (2)
3. Insert needle perpendicular into the port septum. (3) Note: The textured handle and base are not connected. You may position the base on the skin as needed.
4. Verify patency. Dress and secure site per institutional protocol. (4) Begin injection, infusion, or blood aspiration per institutional procedure.
5. **PORT DEACCESS:** After therapy completion, flush port per institutional protocol. Close the clamp while injecting the last 0.5 mL of solution. Stabilize the port by **securely holding** the base down. (5A)
6. **Firmly pull** the textured handle up until you feel a firm stop (5B) **and the needle is locked in the safety position.** Dispose of set in a sharps container.

MR Information:

Non-clinical testing has demonstrated the device is MR conditional. It can be scanned safely under:

- Static magnetic field of 3 tesla or less
- Maximum specific absorption rate (SAR) of 4 W/kg for 15 minutes of scanning
- Spatial gradient field of 360 Gauss/cm or less

In non-clinical testing, the SafeStep* Huber Needle Set produced a temperature rise of less than 1.6 C at a maximum whole body average specific absorption rate (SAR) of 4 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 3 Tesla Siemens Trio using software version VB15.


For minimal Image Artifact

- MR image quality may be compromised if the area of interest is within the distance, specified in Table 1, to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

Huber Needle Set	Image Artifact Distance
SafeStep* Huber Needle Set	7.3 cm

Table 1

* Bard and SafeStep are trademarks and/or registered trademarks of C. R. Bard, Inc.
© 2011 C. R. Bard, Inc. All rights reserved.

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© 2011 C.R. Bard, Inc. All Rights Reserved.
Bard Limited
Bard Access Systems, Inc.