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
EN Contains Phthalates: Bis(2-ethylhexyl)phthalate (DEHP)

Attention: If applicable, the presence and type of Phthalate will be identified on the individual product label.

Directions For Use (Port Identification Key)

- A. Vacuum Port.
- B. Patient Port.
- C. Tandem Port (Orthopedic).
- D. Accessory Port.

Setup and Use

1. **Guardian™ Canisters** (blue lid with rigid canister) and **Guardian™ PPV Canisters** (blue lid with yellow basket underneath and rigid canister): Attach the lid to the canister, making sure to engage all locking tabs (Figure 1). If applicable, place canister/lid assembly into the reusable mounting bracket (Figure 2). **Flex Advantage™ Liners** (blue lid with white filter and flexible liner): Remove from packaging; then twist and fully extend the liner (Figure 3). Insert the liner into the reusable canister. Press down firmly on the lid to ensure proper seal (Figure 4). **CRD™ Liners:** Snap the lid firmly to the liner, making sure to engage all locking tabs (Figure 5). Insert the liner into the reusable hard canister (Figure 6). Press down firmly on the lid to ensure proper seal.

2. Connect the vacuum tubing to Port A.
3. Connect the patient tubing to Port B.
4. Turn On the vacuum supply.
If a tandem setup is required, refer to the instructions supplied with the tandem tubing. It is recommended that the last canister in the tandem setup **not** be filled completely to allow for fluid overflow from the other canisters.

Disposal

Note: Canister should be kept under vacuum until steps 1-3 are complete.

1. Disconnect the patient tubing from Port B and cap the port. If using canisters in tandem, refer to the instructions supplied with the tandem tubing.

2. Firmly seal all ports other than vacuum port (Port A).

3. Disconnect the vacuum tube from Port A and cap the port (Figure 7).

4. Turn Off the vacuum.
5. Ensure all CRD™ ring tabs are engaged, securing the lid into position.

6. Remove canister by lifting upward on the rim of the lid (Figure 8) or by the handle (121 size only).

Note: Canisters are disposable and should not be reused. Dispose of in compliance with hospital protocol.

Caution: Product contents are considered potentially hazardous.

In the Event of Vacuum Loss: Check to ensure the lid on the canister is sealed properly. Check that all connections are firmly secured and that the vacuum is turned On. Note that the liner may collapse if air leaks into the canister (Flex Advantage™ and CRD™ liners only). Check for system damage and replace affected unit.

Warnings

- 1. It is not recommended to subject canister to vacuum when not in use.

2. Premature shut off of the vacuum may occur under the following conditions for the mechanical shut off system:

- The canister is bumped or moved suddenly and designed to stop the vacuum under these conditions. Allow for enough patient tubing during use to avoid tilting the canister when the tubing is pulled and moved during a procedure.
- Suctioning gas-releasing agents, such as hydrogen peroxide or povidone-iodine, can produce high amounts of foam that can activate the shut off valve or cause the fluid to forcefully eject out of the canister. If flow is stopped prematurely, disconnect the vacuum tube from Port A and allow the shut off valve to move back into the open position.

• If a Guardian™, CRD™ or Flex Advantage™ canister is being used by a new canister immediately. Reconnect the vacuum tube to Port A.

3. It is recommended that any solidification agent be introduced into the canister prior to switching off the vacuum source. Failure to do so may cause liquid to flow out of the accessory port when the cap is removed.

4. For suction canister tandem setups, it is recommended that all used and unused disposable canisters be discarded and subsequently replaced following each procedure.

5. Not intended for thoracic or pleural drainage.

6. Flex Advantage™ and CRD™ liners must only be used with the equivalently sized reusable hard canister.

7. Medi-Vac™ canisters must only be used with Medi-Vac™ hardware.

Risks Associated With Reuse: This device has not been validated for reuse. The design of the device may not facilitate cleaning and if inadequately decontaminated may result in risk of infection. The device may not perform as intended by the manufacturer if it is reused and/or decontaminated and may result in material degradation resulting in failure of the device.

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