



# avoset™

## Quick Guide

### Clinican PCA Delivery




# Table of Content

---


## PCA Infusion – Basic Functions

Review Preprogrammed Infusion Parameters .....	3
Change Infusion Parameters .....	4
Review Infusion Parameters During Infusion .....	5
Deliver Patient Bolus .....	6
Deliver Clinician Dose .....	7
Quit Clinician Dose .....	8
Start Loading Dose .....	9
Quit Loading Dose .....	10

 This quick guide provides guidelines for the basic functionalities of the Avoset infusion pump, please refer to the user manual and the training materials terms of use on our website: [www.eitanmedical.com/term\\_training\\_materials/](http://www.eitanmedical.com/term_training_materials/)




# Review Preprogrammed Infusion Parameters



Press the *REVIEW* multi-function button to see the infusion parameters



1st parameter (optional):  
The drug concentration




2nd parameter:  
The total volume to be infused during the infusion



3rd parameter:  
The infusion rate between boluses



4th parameter:  
The volume to be infused in a single patient bolus



5th parameter:  
The time between boluses in which the user cannot initiate a bolus



Back to main display  
(*REVIEW* option is available again)

**i** Based on the infusion programming, additional parameters may be displayed. Refer to the User Manual for further information.



# Change Infusion Parameters


Editing a PCA infusion is possible before the infusion started or when the continuous rate is paused. The changes apply to the current infusion and future repeated infusions.



1. If the infusion is running, press **START/PAUSE** to pause the infusion



2. Press the **EDIT** multi-function button




3. The pump prompts you to enter password. Use the **+** / **-** keys to set the three digits and press the **OK** multi-function button



4. To edit the displayed parameter, use the **+** / **-** keys, and press **OK**. To keep the parameter unedited, press **SKIP**



5. Press **NO** to finish editing, or press **YES** to edit additional parameters



6. Review the edited parameters, to confirm press **CONFIRM**, to change press **EDIT**



7. Press **START/PAUSE** to resume the infusion



8. Back to running continuous rate

**i** The parameters that can be edited on the pump are: Continuous Dose Rate, Bolus Amount and Lockout Duration, to edit other parameters send a new program from the PC programming tool.




# Review Infusion Parameters During Infusion



Press the **INFO** multi-function button to see the infusion parameters




**1st parameter:**  
Time remaining until a bolus can be delivered.  
\*When a bolus is available, this parameter is not presented



**2nd parameter:**  
The total volume infused so far



**3rd parameter:**  
The remaining volume to be infused




**4th parameter:**  
The time elapsed since bolus history was last cleared

**5th parameter:**  
The number of requested boluses and given boluses



**6th parameter:**  
Continuous Rate



**7th parameter:**  
The volume to be infused in a single bolus



**8th parameter (optional):**  
The drug concentration  
Press **DONE** to finish reviewing the parameters




Back to main display (**INFO** option is available again)

**i** Based on the infusion programming, additional parameters may be displayed. Refer to the User Manual for further information.  
**Tip:** long press on the **NEXT** multi-function button will return to home screen



# Deliver Patient Bolus



1. Press  to initiate a bolus



2. Bolus delivery indication

After delivering a bolus, a lockout period between boluses prevents the user from requesting additional bolus until the lockout period elapses.

Pressing the  key will prompt this notification




# Deliver Clinician Dose

Clinician dose is an optional way for the clinician to add dose for the patient during a PCA treatment.



1. Press the **CLINICIAN DOSE** multi-function button to start setting the clinician dose



2. The pump prompts you to enter password. Use the **+** / **-** keys to set the three digits and press the **OK** multi-function button



3. Use the **+** / **-** keys to change the value. Press the **OK** multi-function button to approve and proceed



4. Press the **YES** multi-function button to approve and start delivering clinician dose



5. Clinician dose delivery indication



# Quit Clinician Dose



1. Press  to pause the Clinician dose



2. Press the **QUIT** multi-function button to quit the Clinician dose



3. Press the **YES** multi-function button to quit the Clinician dose



4. Press the **YES** multi-function button to resume the continuous rate



5. Back to running continuous rate





# Start Loading Dose

A Loading Dose is an optional clinician dose given when starting the first infusion right after programming. Loading Dose is part of the program and is available only when the first infusion is delivered. When repeating the same infusion, the infusion will start with the continuous rate.



1. Press  to start the Loading Dose



2. Loading Dose delivery indication



# Quit Loading Dose



1. Press **START PAUSE** to pause the Loading Dose



2. Press the **QUIT** multi-function button to quit Loading Dose



3. Press the **YES** multi-function button



4. Back to continuous rate



This page is left intentionally blank





Avoset English Clinician Quick Guide PCA | 15132-048-0005-UM, Rev.02



MedNet EC-REP IIb GmbH  
Borkstrasse 10  
48163 Münster, Germany



Eitan Medical Ltd.  
29 Yad Haruzim St.  
P.O. Box 8639  
Netanya 4250529, Israel



CE  
0123

