## EN Instructions for Use - Spike Sets

### Warnings:

- This set is not intended for delivery of cellular blood products.
- This set is not intended for gravity infusions.
- To prevent potential air delivery to patient, remove all air from the bag and tubing prior to initiating infusion when not hanging on an IV pole. Air delivery could result in an embolism, which can cause serious injury or death.
- When using a vial or rigid medication reservoir, a vented spike adapter or an administration set with a vented spike must be used.
- When using an administration set with a vented spike, always hang the medication reservoir from an IV pole throughout the duration of the infusion to prevent air delivery to the patient.
- To prevent route misconnections, trace all lines from the medication reservoir to the patient's access site. Verify that all attachments are correct before making any connections or reconnections, administering medications, solutions, or other products.
- Set use should not exceed 96 hours or 9 liters (whichever comes first). Pump performance may be affected beyond these limits.
- Inspect before use. Do not use if package seal is broken, if the caps are removed from the tubing, if the white slide clamp(s) (C) is closed, or if the break-away tab is missing from the yellow "Flow-stop" (B).
- Recapping of the bag spike is not recommended. If recapping is necessary, avoid touching the spike or scraping the inside of the cap with the spike.
- Avoid leaving tubing and cords where infants, children, or those deemed at high risk for medical line entanglement, can become entangled. Be aware that if these items become wrapped around a patient's neck, there is an increased risk of strangulation or death. Caregivers of patients who are at risk for entanglement should discuss with their health care provider how to properly manage their lines and properly monitor patients based on their needs while in use.
- Do not prime while the tubing is connected to the patient's access device. Doing so could result in overdose or air embolism and may cause injury.

**Warning:** This product can expose you to chemicals including vinyl chloride, which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

### Before you begin (Precautions):

- Read all instructions prior to use. Refer to the CURLIN 8000 Pump User Manual (LS-61765), CURLIN 6000 Pump User Manual (LS-360-9009) or the CURLIN PainSmart Pump User Manual (LS-360-9043) for further instructions, cautions, and warnings pertaining to the system.
- Before use, patients and caregivers must be

educated on the use of the administration set by a qualified health care provider. A qualified health care provider should monitor ongoing use.

- Sterile technique is required for all fluid pathway connections.
- When the tubing is in the pump, always close the white slide clamp (C) prior to opening the pump door.
- For sets with an NRFit connector, the neuroaxial hubs are narrower than standard luer connectors to prevent mis-connection. Align hubs, press, and twist connectors together.
- The administration set contains small parts such as caps, which could be a choking hazard. Keep all small parts away from children.
- When possible, bring the medication/solution to room temperature prior to starting the infusion to minimize "air-in-line" alarms.
- Refer to facility, drug manufacturer, and other applicable standards such as INS Standards of Practice and CDC Guidelines for additional recommendations on topics such as set change frequency, use of add-on devices, needleless connectors, etc.
- Filter information: 0.2 micron filters have been tested for bacterial retention up to 96 hours. 1.2 micron filters have been tested for fungus retention up to 24 hours.
- Delivery accuracy may be affected by in-line components such as filters. See pump user manual for more information on pump delivery accuracy.
- Follow facility protocols for the safe handling and disposal of the administration set in accordance with local, state, and federal regulations.
- For disposal in home settings, contact your care provider for instruction as medications may impact method of disposal.

### Priming and Air Removal Instructions for End Users (Clinician/Patient/Caregiver):

#### Refer to figures on front of pouch and below for additional instructions.

The air must be removed from the administration set (tubing) prior to connecting to the patient. This is called "priming the set." Most sets can be primed before or after inserting into the pump. Follow your healthcare provider's instructions on which method of priming you should use.

Sets with anti-siphon valves must be primed using the pump. If priming has not been completed, follow the instructions below:

#### Priming BEFORE inserting the administration set into the pump:

To maintain sterility during setup, blue cap may be left on luer while priming.

Note: If an add-on "R-lock check valve" is provided, connect it to the luer connector before priming.

- With the slide clamp closed, **twist (do**

**not snap)** and remove the breakaway tab from the yellow "Flow-Stop" (B).

- Remove protector cover from outlet port on medication reservoir. Hold port with your non-dominant hand and insert the spike into port until it will go no further. Do not touch the spike while inserting into the port. Do not use if spike punctures medication reservoir or port.
- If the bag **will not** be hung from an IV pole throughout the infusion, all air in the bag must be removed. To do so, turn the bag upside down during priming so that air is at the top and can exit through the tubing port.
- Open the white slide clamp(s) (C).
- "Pinch to Prime": To remove air from the bag, gently squeeze the inverted bag while simultaneously pinching the top of the yellow "Flow-Stop" (B). Continue priming until all air has been removed from the administration set (and bag if applicable). Fluid should be visible throughout the length of the tubing from the medication reservoir to the luer connector.

*Note: Administration Sets with Filters: The set you are using may contain an air-eliminating filter. When using sets with filters, always follow the priming instructions included with the set you are using. While priming, filters should be inverted so arrow points upward. This position will allow the filter to fill with fluid on both sides. Clamp immediately after priming is complete to prevent emptying of line.*

- Do not place tape over the filter. This could block the filter vent and prevent it from working properly.
- Avoid raising the filter above the height of the patient's catheter entry point throughout the infusion (refer to User Manual for more details).
- Release the yellow "Flow-Stop" (B) and manually close the white slide clamp (C).
- Follow the instructions on the back of pouch label for installing the administration set into the pump.

### Priming AFTER inserting the administration set into the Pump:

To remove air from the administration set using the pump's prime feature, refer to pump's "User Manual" for instructions.

### Additional instructions for Sets with a Y-site:

- The Y-site should remain capped at all times. Do not leave the Y-site open to the air.
- Consider replacing the cap on the Y-site with a needleless connector.
- Prime the Y-site using a syringe and preferred priming solution prior to priming the full administration set.
- Close the white slide clamp (C) on the Y-site after priming and when not infusing through the Y-site.
- Close the white slide clamp (C) between the pump and Y-site to prevent retrograde flow only when infusing through the Y-site.
- When infusing through both lines, the slower flow rate should be infused through the Y-site.

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