

LTV™ and PTV ventilators

Vyaire recommends the following steps to avoid infection and contamination when using our ventilators.

- Please follow the [WHO guidelines](#) to prevent epidemic and pandemic infections from spreading
- Please follow the instructions for use for Vyaire (LTV and PTV) ventilators at www.vyaire.com/Covid-19
- Please follow the cleaning and disinfection instructions for www.vyaire.com/Covid-19
- Always use an inspiratory bacteria filter (www.vyaire.com/Covid-19) to protect the inspiratory side from contamination
- If bacterial filters are used in conjunction with the LTV and PTV ventilators, comply with all procedures as specified by the filter manufacturer
- Use a heat moisture exchanger with HEPA Filter (HMEF) at the patient wye to protect against contamination. Please consider that resistance and apparatus dead space will increase and follow the recommendations of the manufacturer of the HMEF for changing intervals
- When changing the HMEF please follow your PPE policy and procedure
- Single-use consumables, such as breathing circuits, nebulizers, airway adapters and filters, can help to mitigate to cross contamination when the ventilator is reprocessed for use of the next patient. A list of consumables suitable for Vyaire ventilators is available on www.vyaire.com/Covid-19
- When reusable materials are in use follow the hospital, manufacturers or the [WHO guidelines](#) for reprocessing these items
- For suctioning patients please use a closed suction system like the [Verso](#) to prevent cross contamination
- Disinfect the outer surfaces of the ventilators regularly during ventilation or after treatment of a patient with a suitable registered disinfectant. Consult the hygiene specialist in your facility regarding the appropriate disinfectant or follow the instructions for use of the ventilator



PTV ventilator



LTV™ 1200

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA

vyaire.com

REONLY For U.S. distribution.

© 2020 Vyaire. Vyaire, the Vyaire logo and all other trademarks or registered trademarks are property of Vyaire Medical, Inc., or one of its affiliates. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-GBL-2000111

vyaire[™]
M E D I C A L