



# Reference Guide for MicroClave<sup>®</sup> Clear

Neutral Displacement Connector

## Directions for Use



### Disinfect to Protect

- › When placing a new MicroClave on a catheter, disinfect the catheter hub and prime the MicroClave.
- › Before accessing a MicroClave, always disinfect the injection site with the approved antiseptic per facility protocol.
- › Scrub the injection site in accordance with facility protocol for appropriate scrubbing and dry times.



### Administer or Aspirate

- › Attach IV tubing, syringe or blood tube holder to MicroClave by inserting the luer and twisting 1/4 turn, or until a friction fit is achieved.
- › Do not over-tighten a luer beyond the friction fit as this may damage both the luer and the MicroClave.
- › To disconnect, grasp MicroClave and then twist mating luer away from MicroClave until loose. Do not hold catheter hub during disconnect as this may cause accidental removal of MicroClave from hub.



### Flush After Each Use

- › Flush the MicroClave with normal saline or in accordance with facility protocol. After blood use, the MicroClave can be flushed clean and does not require change-out.
- › Use routine flushing in accordance with facility protocol in order to maintain catheter patency.
- › Change MicroClave in accordance with facility protocol and CDC Guidelines.

## Functional Attributes

- › Lipid and Blood Compatible
- › Radiographic Imaging Compatible (CT Compatible)
- › Contains No Latex, Phthalate's (DEHP) or Metal Components, MRI Compatible



MicroClave's saline flush option can eliminate the risk of Heparin Induced Thrombocytopenia (HIT).

**icumedical**  
human connections