

1000 mL

List No. 7927-21

0—		—0
1—		—1
2—		—2
3—		—3
—		—
4—		—4
—		—
5—		—5
—		—
6—		—6
—		—
7—		—7
—		—
8—		—8
—		—
9—		—9

**EMPTY
CONTAINER**
with Attached
Y-Transfer Set,
Convertible Pins

*SINGLE-USE
I.V. CONTAINER
FOR ADMIXING
COMPATIBLE MEDICATIONS.*

FLUID PATH AND AREAS BENEATH UNDISTURBED PROTECTIVE SET COVERS ARE STERILE AND NONPYROGENIC IN INTACT UNIT PACKAGE. USE ASEPTIC TECHNIQUE. MIX CONTAINER CONTENTS THOROUGHLY. ALL GRADUATIONS ARE APPROXIMATE. LABEL APPROPRIATELY. BEFORE ADMINISTRATION USE, CHECK FOR MINUTE LEAKS BY SQUEEZING CONTAINER FIRMLY. DO NOT USE IF CONTAINER IS DAMAGED. ADHERE TO STORAGE REQUIREMENTS OF ADDED MEDICATIONS. MUST NOT BE USED IN SERIES CONNECTIONS.

MADE IN COSTA RICA.

PRODUCT INQUIRIES SHOULD BE DIRECTED TO
ICU MEDICAL, INC., LAKE FOREST, ILLINOIS, 60045, USA



Not made with
natural rubber latex



CONTAINS DEHP
Recyclable as Type-3 Plastic

Rx ONLY

IM-4241 (06/2017) **icumedical**