

1000 mL

RX ONLY

List No. 7951-19

EMPTY CONTAINER

*SINGLE-USE I.V. CONTAINER
FOR ADMIXING
COMPATIBLE MEDICATIONS.*

FLUID PATH AND AREAS BENEATH UNDISTURBED
PROTECTIVE COVERS ARE STERILE AND NONPYROGENIC
IN INTACT UNIT PACKAGE.

USE ASEPTIC TECHNIQUE.

MIX CONTAINER CONTENTS THOROUGHLY.
ALL GRADUATIONS ARE APPROXIMATE. LABEL
APPROPRIATELY.

BEFORE ADMINISTRATION USE, CHECK FOR MINUTE
LEAKS BY SQUEEZING CONTAINER FIRMLY.

DO NOT USE IF CONTAINER IS DAMAGED.

ADHERE TO STORAGE REQUIREMENTS OF ADDED
MEDICATIONS.

MUST NOT BE USED IN SERIES CONNECTIONS.

MADE IN COSTA RICA.

PRODUCT INQUIRIES SHOULD BE DIRECTED TO
ICU MEDICAL, INC., LAKE FOREST, ILLINOIS, 60045, USA



Not made with
natural rubber latex



CONTAINS DEHP
Recyclable as Type-3 Plastic