

HOMEPUMP*

Ambulatory Infusion System

AVANOS

HOMEPUMP* AMBULATORY INFUSION SYSTEM: A DISPOSABLE SOLUTION FOR YOUR PATIENTS



Avanos offers a full portfolio of IV Infusion Pumps -
HOMEPUMP* has been on the market as an IV
therapy solution for over 20 years.

HOMEPUMP ECLIPSE*
HOMEPUMP C-SERIES*

HOMEPUMP ECLIPSE* AND HOMEPUMP C-SERIES* AMBULATORY INFUSION SYSTEMS

The Homepump system is used for intermittent antibiotic therapies as well as continuous infusions of chemotherapy.

ECLIPSE* SYSTEM

Indicated for continuous delivery of medications through intravenous routes.

C-SERIES* SYSTEM

Indicated for continuous delivery of medications through intravenous, intra-arterial, intramuscular, subcutaneous or epidural routes.



WHY HOMEPUMP*?

Benefits for Patients and HCPs


- Disposable solution for your patients to continue intravenous infusion at home
- Independently tested for drug stability of multiple medications
- Pre-set flow rate
- +/- 15% accuracy when pump is filled to the labeled (nominal) volume
- HOMEPUMP ECLIPSE* System: Approximate delivery times between 15 minutes to 5 hours and 15 minutes
- HOMEPUMP C-SERIES* System: Approximate delivery times between 12 hours to 12 days



HOMEPUMP* DESIGN


HOMEPUMP ECLIPSE*

Ideal for Antibiotic Delivery



HOMEPUMP C-SERIES*

Ideal for Chemotherapy and other Protocols



HOME PUMP ECLIPSE*

Ideal for Antibiotic Delivery

Model #	Size/Flow Rate	Qty/Box
E050500	50 ml vol, 50 ml/hr	48
E100500	100 ml vol, 50 ml/hr	48
E101000	100 ml vol, 100 ml/hr	48
E102000	100 ml vol, 200 ml/hr	48
E251750	250 ml vol, 175 ml/hr	24
E252500	250 ml vol, 250 ml/hr	24
E401000	400 ml vol, 100 ml/hr	24
E402000	400 ml vol, 200 ml/hr	24
E502500	500 ml vol, 250 ml/hr	24

HOME PUMP C-SERIES*

Ideal for Chemotherapy and other Protocols

Model #	Size/Flow Rate	Qty/Box
C060020	60 ml, 2 ml/hr	24
C100005	100 ml, 0.5 ml/hr	24
C100020	100 ml, 2 ml/hr	24
C125050	125 ml, 5 ml/hr	24
C270010	270 ml, 1 ml/hr	24
C270020	270 ml, 2 ml/hr	24
C270050	270 ml, 5 ml/hr	24
C270100	270 ml, 10 ml/hr	24
C300060	300 ml, 6 ml/hr	24

HOME PUMP* ACCESSORIES

Model #	Description	Qty/Box
5001096	Carry Pouch, 50-200 ml pumps	12
5001097	Carry Pouch, 250-500 ml pumps	12
5001087	E-Clip	100

There are inherent risks in all medical devices. Please refer to the product labeling for **Indications, Cautions, Warnings** and **Contraindications**. Failure to follow the product labeling could directly impact patient safety. Physician is responsible for prescribing and administering medications per instructions provided by the drug manufacturer. Refer to www.avanospainmanagement.com for additional product safety Technical Bulletins.

For more information for clinicians and product code ordering information, please visit: [avanospainmanagement.com](http://www.avanospainmanagement.com)
Call 1-800-448-3569 in the United States and Canada.