



RESPIRATORY HEALTH PRODUCTS



CONTENTS

ADULT

ENDOTRACHEAL TUBES

HALYARD* MICROCUFF* ENDOTRACHEAL TUBES

Oral/Nasal Magill Endotracheal Tube page 4

HALYARD* MICROCUFF* SUBGLOTTIC SECRETION DRAINAGE ENDOTRACHEAL TUBE

Subglottic Secretion Drainage Endotracheal Tube page 4

CLOSED SUCTION SYSTEMS

HALYARD* CLOSED SUCTION SYSTEMS

Double Swivel Elbows page 5

WET PAK* Double Swivel Elbows page 5

T-Pieces page 6

WET PAK* T-Pieces page 6

Standard Dual Lumens page 7

Standard Directional Tips page 7

HALYARD* TURBO-CLEANING CLOSED SUCTION SYSTEMS

Double Swivel Elbows page 8

WET PAK* Double Swivel Elbows page 8

T-Pieces page 9

WET PAK* T-Pieces page 9

HALYARD* MULTI-ACCESS PORT CLOSED SUCTION SYSTEMS

Double Swivel Elbows page 10

Replacement Catheter page 10

WET PAK* Double Swivel Elbows page 11

Seal Cassette page 11

HALYARD* CLOSED SUCTION SYSTEMS – ACCESSORIES

HALYARD* Catheter Accessories page 15

BRONCHIAL ASPIRATE SAMPLING CATHETERS

HALYARD* Bal Cath* Bronchial Aspirate Sampling Catheters page 12

ORAL CARE SOLUTIONS

HALYARD* ORAL CARE SOLUTIONS

Oral Suction Single Patient Components and Packs page 13

DENTASWAB* Oral Care Swabs page 13

24-Hour Oral Care Kit page 14

Single Patient Suction Kit page 14

PAEDIATRIC

ENDOTRACHEAL TUBES

HALYARD* MICROCUFF* ENDOTRACHEAL TUBES

Oral/Nasal Magill Endotracheal Tube	page 16
Oral Curved Endotracheal Tube	page 16

CLOSED SUCTION SYSTEMS

HALYARD* CLOSED SUCTION SYSTEMS – NEONATES/PAEDIATRIC

Elbow Tube Adapters	page 17
Y Tube Adapters	page 17
Multi-Access Catheters	page 18

ORAL CARE SOLUTIONS

HALYARD* ORAL CARE SOLUTIONS

Oral Suction Single Patient Components	page 18
--	---------

ENDOTRACHEAL TUBES

HALYARD* MICROCUFF* ORAL/NASAL MAGILL ENDOTRACHEAL TUBE

- Advanced micro-thin polyurethane cuff material designed to reduce micro-aspiration
- Provides a superior tracheal seal and reduces leakage past the cuff
- Seals at low pressure
- Oral/Nasal Magill Curve with Murphy Eye

Code	Inner Diameter	Description	Packaging
35210	5.0mm	Oral/Nasal Magill, Murphy Eye, Endotracheal Tube	10/case
35211	5.5mm	Oral/Nasal Magill, Murphy Eye, Endotracheal Tube	10/case
35212	6.0mm	Oral/Nasal Magill, Murphy Eye, Endotracheal Tube	10/case
35213	6.5mm	Oral/Nasal Magill, Murphy Eye, Endotracheal Tube	10/case
35214	7.0mm	Oral/Nasal Magill, Murphy Eye, Endotracheal Tube	10/case
35215	7.5mm	Oral/Nasal Magill, Murphy Eye, Endotracheal Tube	10/case
35216	8.0mm	Oral/Nasal Magill, Murphy Eye, Endotracheal Tube	10/case
35217	8.5mm	Oral/Nasal Magill, Murphy Eye, Endotracheal Tube	10/case
35218	9.0mm	Oral/Nasal Magill, Murphy Eye, Endotracheal Tube	10/case
35220	10.0mm	Oral/Nasal Magill, Murphy Eye, Endotracheal Tube	10/case



HALYARD* MICROCUFF* SUBGLOTTIC SECRETION DRAINAGE ENDOTRACHEAL TUBE

Code	Inner Diameter	Outer Diameter	Description	Packaging
13220	7.0mm	10.4mm	Subglottic Secretion Drainage Endotracheal Tube	10/case
13221	7.5mm	11.2mm	Subglottic Secretion Drainage Endotracheal Tube	10/case
13222	8.0mm	11.8mm	Subglottic Secretion Drainage Endotracheal Tube	10/case
13223	8.5mm	12.3mm	Subglottic Secretion Drainage Endotracheal Tube	10/case
13224	9.0mm	12.8mm	Subglottic Secretion Drainage Endotracheal Tube	10/case



CLOSED SUCTION SYSTEMS

HALYARD* DOUBLE SWIVEL ELBOWS

- TRACH CARE* Technology
- Gamma Sterilised
- One-Way Lavage/Irrigation Port
- Dual Swivels provides flexibility without connectors
- Tactile Sleeve for Better Feel
- Numbered catheter allowing the clinician to limit the distance of catheter advance
- Tracheostomy sizes
- Multiple French sizes


Code	Size (Outer Diameter)	Length	Description	Packaging
2160-5	12 Fr / 4mm	54cm	Double Swivel Elbow	20/case
21603-5	12 Fr / 4mm	30.5cm	Double Swivel Elbow, Tracheostomy	20/case
2210-5	14 Fr / 4.6mm	54cm	Double Swivel Elbow,	20/case
22103-5	14 Fr / 4.6mm	30.5cm	Double Swivel Elbow, Tracheostomy	20/case
22108-5	14 Fr / 4.6mm	54cm	Double Swivel Elbow, MDI Built-In	20/case
2260-5	16 Fr / 5.3mm	54cm	Double Swivel Elbow	20/case
25-2210-5	14 Fr / 4.6mm	63.5cm	Double Swivel Elbow, Safety Lavage Port	20/case



HALYARD* WET PAK* DOUBLE SWIVEL ELBOWS

- Catheters packaged separate from Vials
- 12 Saline Vials per catheter
- Vials designed for respiratory use



Code	Size (Outer Diameter)	Length	Description	Packaging
21606-5	12 Fr / 4mm	54cm	WET PAK*, Double Swivel Elbow	20/case
216036-5	12 Fr / 4mm	30.5cm	WET PAK*, Double Swivel Elbow, Tracheostomy	20/case
221036-5	14 Fr / 4.6mm	30.5cm	WET PAK*, Double Swivel Elbow, Tracheostomy	20/case
22106-5	14 Fr / 4.6mm	54cm	WET PAK*, Double Swivel Elbow	20/case
22606-5	16 Fr / 5.3mm	54cm	WET PAK*, Double Swivel Elbow	20/case



HALYARD* T-PIECES

- TRACH CARE* Technology
- Gamma Sterilised
- Lavage / Irrigation Port
- Locking Thumb Port
- T-Piece
- Multiple French Sizes


Code	Size (Outer Diameter)	Length	Description	Packaging
2155-5	12 Fr / 4mm	54cm	Clear T-Piece	20/case
2205-5	14 Fr / 4.6mm	54cm	Clear T-Piece	20/case
220135-5	14 Fr / 4.6mm	30.5cm	Clear T-Piece, Swivel Adapter, Tracheostomy	20/case
2255-5	16 Fr / 5.3mm	54cm	Clear T-Piece	20/case



HALYARD* WET PAK* T-PIECES

- TRACH CARE* Technology
- Gamma Sterilised
- Lavage / Irrigation Port
- Catheters packaged separate from Vials
- 12 Saline Vials per Catheter
- Vials designed for respiratory use
- Multiple French sizes


Code	Size (Outer Diameter)	Length	Description	Packaging
22056-5	14 Fr / 4.6mm	54cm	WET PAK*, Clear T-Piece	20/case
220568-5	14 Fr / 4.6mm	54cm	WET PAK*, Clear T-Piece, MDI Adapter	20/case
22012356-5	14 Fr / 4.6mm	30.5cm	WET PAK*, Clear T-Piece, Swivel Adapter, 15mm x 22mm Adapter, Tracheostomy	20/case
220123568-5	14 Fr / 4.6mm	30.5cm	WET PAK*, Clear T-Piece, Swivel Adapter, 15mm x 22mm Adapter, MDI Adapter, Tracheostomy	20/case



HALYARD* STANDARD DUAL LUMENS

- TRACH CARE* Technology
- Gamma Sterilised
- Dual Lavage / Irrigation Ports
- Available in T-Piece and Double Swivel Elbow
- Available in Directional Tip
- Allows Lavage to take place at the patient end of the Catheter
- Multiple French Sizes



Code	Size (Outer Diameter)	Length	Description	Packaging
22107-5	14 Fr / 4.6mm	54cm	Double Swivel Elbow, Dual Lumen, Endotracheal	20/case
221037-5	14 Fr / 4.6mm	30.5cm	Double Swivel Elbow, Dual Lumen, Tracheostomy	20/case



HALYARD* STANDARD DIRECTIONAL TIPS

- TRACH CARE* Technology
- Gamma Sterilised
- Lavage / Irrigation Port
- Available in T-Piece, Elbow and Double Swivel Elbow
- Allows for Suctioning of the Right or Left Main Stem Bronchus
- Radiopaque Line provides for Observation of Catheter Path under Fluoroscopy
- Coloured Centre Line allows maintaining manipulation of the Catheter Tip to access proper location of the Right or Left Main Stem Bronchus
- Multiple French Sizes


Code	Size (Outer Diameter)	Length	Description	Packaging
21609-5	12 Fr / 4mm	56cm	Double Swivel Elbow, Directional Tip, Endotracheal	20/case
22059-5	14 Fr / 4.6mm	56cm	Clear T-Piece, Directional Tip, Endotracheal	20/case
220569-5	14 Fr / 4.6mm	56cm	WET PAK*, Clear T-Piece, Directional Tip, Endotracheal	20/case
221069-5	14 Fr / 4.6mm	56cm	WET PAK*, Double Swivel Elbow, Directional Tip	20/case
22109-5	14 Fr / 4.6mm	56cm	Double Swivel Elbow, Directional Tip, Endotracheal	20/case



CLOSED SUCTION SYSTEMS


HALYARD* TURBO-CLEANING CLOSED SUCTION SYSTEMS

- TRACH CARE* Technology
- Gamma Sterilised
- One-Way Lavage / Irrigation Port
- Dual Swivels provides flexibility without connectors
- Tactile Sleeve for better feel
- Numbered Catheter allowing the clinician to limit the distance of catheter advance
- Tracheostomy sizes
- Multiple French sizes



HALYARD* TURBO-CLEANING – DOUBLE SWIVEL ELBOWS


Code	Size (Outer Diameter)	Length	Description	Packaging
22716-5	12 Fr / 4mm	54cm	Double Swivel Elbow, Endotracheal	20/case
2271603-5	12 Fr / 4mm	30.5cm	Double Swivel Elbow, Tracheostomy	20/case
227-5	14 Fr / 4.6mm	54cm	Double Swivel Elbow, Endotracheal	20/case
22714186-5	14Fr/ 4.6mm	54cm	Double Swivel Elbow, MDI Built-In, Endotracheal	20/case
22703-5	14 Fr / 4.6mm	30.5cm	Double Swivel Elbow, Tracheostomy	20/case
22716186-5	16Fr/ 5.3mm	54cm	Double Swivel Elbow, MDI Built-In, Endotracheal	20/case



HALYARD* TURBO-CLEANING – WET PAK* DOUBLE SWIVEL ELBOWS

- Catheters packaged separate from Vials
- 12 Saline Vials per catheter
- Vials designed for respiratory use
- Multiple French sizes

Code	Size (Outer Diameter)	Length	Description	Packaging
227166-5	12 Fr / 4mm	54cm	WET PAK*, Double Swivel Elbow	20/case
22716036-5	12 Fr / 4mm	30.5cm	WET PAK*, Double Swivel Elbow, Tracheostomy	20/case
2276-5	14 Fr / 4.6mm	54cm	WET PAK*, Double Swivel Elbow	20/case
227012586	14Fr/ 4.6mm	54cm	WETPAK* Double Swivel Elbow, MDI Built-In, Endotracheal	20/case
227036-5	14 Fr / 4.6mm	30.5cm	WET PAK*, Double Swivel Elbow, Tracheostomy	20/case




NOTE: Please see package inserts for directions for use for each of these products.

HALYARD* TURBO-CLEANING – T-PIECES

- Available in T-Piece, Elbow and Double Swivel Elbow
- Lavage / Irrigation Port
- Multiple French sizes


Code	Size (Outer Diameter)	Length	Description	Packaging
22715-5	12 Fr / 4mm	54cm	T-Piece Endotracheal	20/case
227135-5	12 Fr / 4mm	30.5cm	T-Piece Tracheostomy	20/case



HALYARD* TURBO-CLEANING – WET PAK* T-PIECES


- Catheters packaged separate from Vials
- 12 saline Vials per Catheter
- Vials designed for respiratory Use
- Multiple French sizes

Code	Size (Outer Diameter)	Length	Description	Packaging
227156-5	12 Fr / 4.0mm	54cm	WET PAK* T-Piece, Standard	20/case
227056-5	14 Fr / 4.6mm	54cm	WET PAK* T-Piece, Standard	20/case
22701356-5	14 Fr / 4.6mm	30.5cm	WET PAK* T-Piece, Tracheostomy	20/case




MULTI-ACCESS PORT CLOSED SUCTION SYSTEMS

- TRACH CARE* Technology
- Gamma sterilised
- Available with Metered Dose Inhaler (MDI) Port
- Built-in Cleaning Chamber
- Rotational Manifold with multiple ports (Suctioning Port and Alternate Therapy for Bronchoscopy or BAL CATH*)
- Replaceable Catheters available separately
- DEHP Free
- Tethered Sleeve
- Tethered Cap
- Multiple French sizes



HALYARD* MULTI-ACCESS PORT – DOUBLE SWIVEL ELBOWS

Code	Size (Outer Diameter)	Length	Description	Packaging
240128-5	12 Fr / 4mm	58.5cm	Double Swivel Elbow, Endotracheal, MDI Port	20/case
24001238-5	12 Fr / 4mm	35cm	Double Swivel Elbow, Tracheostomy, MDI Port	20/case
240148-5	14 Fr / 4.6mm	58.5cm	Double Swivel Elbow, Endotracheal, MDI Port	20/case
2401438-5	14 Fr / 4.6mm	35cm	Double Swivel Elbow, Tracheostomy, MDI Port	20/case




HALYARD* MULTI-ACCESS PORT REPLACEMENT CATHETERS

Code	Size (Outer Diameter)	Length	Description	Packaging
240125-5	12 Fr / 4mm	58.5cm	Endotracheal, Replacement Catheter	20/case
2401253-5	12 Fr / 4mm	35cm	Tracheostomy, Replacement Catheter	20/case
240145-5	14 Fr / 4.6mm	58.5cm	Endotracheal, Replacement Catheter	20/case
2401453-5	14 Fr / 4.6mm	35cm	Endotracheal, Replacement Catheter	20/case



HALYARD* MULTI-ACCESS PORT – WET PAK* DOUBLE SWIVEL ELBOWS

- Catheters packaged separate from Vials
- 12 saline Vials per Catheter
- Vials designed for respiratory Use



Code	Size (Outer Diameter)	Length	Description	Packaging
2401236-5	12 Fr / 4mm	35cm	WET PAK*, Double Swivel Elbow, MDI Port, Endotracheal	20/case
24012386-5	12 Fr / 4mm	35cm	WET PAK*, Double Swivel Elbow, MDI Port, Tracheostomy	20/case
2401286-5	12 Fr / 4mm	58.5cm	WET PAK*, Double Swivel Elbow, MDI Port, Endotracheal	20/case
24014386-5	14 Fr / 4.6mm	35cm	WET PAK*, Double Swivel Elbow, MDI Port, Tracheostomy	20/case
2401486-5	14 Fr / 4.6mm	58.5cm	WET PAK, Double Swivel Elbow, MDI Port, Endotracheal	20/case

HALYARD* SEAL CASSETTE


Code	Description	Packaging
24013-5	Seal Cassette	20/case



BRONCHIAL ASPIRATE SAMPLING CATHETERS

HALYARD* BAL CATH* BRONCHIAL ASPIRATE SAMPLING CATHETERS

- Bronchial Aspirate Sampling Catheters
- HALYARD* BAL CATH* Bronchial Aspirate Sampling Catheters
- Gamma sterilised
- Atraumatic "Mushroom" Tip
- Directional Tip for R/L Insertion
- Early, accurate diagnosis of VAP
- Simple, effective non-bronchoscopic bronchoalveolar lavage at the bedside
- Includes an Access Port Elbow
- Multiple French sizes




Code	Description	Packaging
142A	BAL CATH* Bronchial Aspirate Sampling Catheter, 16 Fr Outer/12 Fr Inner Telescoping Catheter	5/case
143A	BAL CATH* Bronchial Aspirate Sampling Catheter, 13 Fr Outer/8 Fr Inner Telescoping Catheter	5/case

ORAL CARE SOLUTIONS



HALYARD* ORAL SUCTION SINGLE PATIENT COMPONENTS AND PACKS

Code	Description	Packaging
Components:		
12601	1-Oral Suction Catheter	25/case
12602	1-Suction Toothbrush	25/case
12603	1-Angled-Tip Suction Swab	25/case
12604	1-Oral Suction Tip	25/case
12606	1-Suction Valve and Adapter	25/case
99785	1-Self-Cleaning Covered Yankauer	25/case
99786	5-Suction Handle	25/case
99789	1-Oral Care Rep Pack: 1-Self-Cleaning Covered Yankauer 1-Suction Adapter 1-Y-Adapter	40/case



HALYARD* DENTASWAB* ORAL CARE SWABS

Code	Description	Packaging
DENTASWAB* – Single Wrapped		
12240	250-Single Wrapped DENTASWAB* Swabs	4/case
12245	250-Single Wrapped DENTASWAB* Swabs with Dentifrice	4/case
12246	250-Single Wrapped DENTASWAB* Swabs with Dentifrice	2/case
12254	250-Single Wrapped Polyplus DENTASWAB* Swabs	4/case
12255	250-Single Wrapped Polyplus DENTASWAB* Swabs with Dentifrice	4/case
DENTASWAB* – Bulk Packaged		
12249	20-DENTASWAB* Swabs with Dentifrice	50/case
12252	500-Polyplus DENTASWAB* Swabs	2/case




HALYARD* 24-HOUR ORAL CARE KIT

Code	Description	Packaging
25928	24-Hour Oral Care Advanced Kit	5/case
Components:		
1-Oral Care Prep Pack with Self Cleaning Covered Yankauer, Suction Handle and "Y" Adapter		
2-Toothbrush Packs with Empty Solution Cups		
4-Suction Swab Packs with Empty Solution Cups		



HALYARD* ORAL CARE SINGLE PATIENT SUCTION KIT

Code	Description	Packaging
123	Oral Care Single Patient Suction Kit	50/case
Components:		
1-Suction Valve		
2.4m Tubing		
10-Oral Suction probes		
1-Suction Line Clamp		
1-"Y" Adapter		





HALYARD* CLOSED SUCTION SYSTEMS – ACCESSORIES

HALYARD* TRACH CARE* CATHETER ACCESSORIES

- Complement the HALYARD* Closed System line of products
- Can be sold individually or with TRACH CARE*
- Designed for specific patient needs

Code	Description	Packaging
11000	Manual Percussor for Neonatal Pulmonary Use	12/case
111	Swivel Adapter 15mm Inner Diameter x 15mm Outer Diameter	50/case
112	Fixed Adapter 15mm Outer Diameter x 22mm Inner Diameter	50/case
114	Parallel "Y"	50/case
118	Unit Dose Saline .9%, Pink	144/case
119	Metered Dose Adapter	20/case
19000	Transport Kit "Y" Plug, SUC CVC	144/case
19030	HALYARD* 3mm Endotracheal Adapters	50/case
19035	HALYARD* 3.5mm Endotracheal Adapters	50/case




PAEDIATRIC ENDOTRACHEAL TUBES

HALYARD* MICROCUFF* ENDOTRACHEAL TUBES


- Ultra-low cuff pressure can reduce the risk of tracheal trauma
- Short, cylindrical cuff placed near the tracheal tube tip
- Four precision bands to facilitate and confirm optimal tube placement
- No Murphy Eye
- Oral/Nasal Magill and Pre-Curved shapes

Code	Inner Diameter	Description	Packaging
35111	3.0mm	Paediatric Oral/Nasal Magill, Endotracheal Tube	10/case
35112	3.5mm	Paediatric Oral/Nasal Magill, Endotracheal Tube	10/case
35113	4.0mm	Paediatric Oral/Nasal Magill, Endotracheal Tube	10/case
35114	4.5mm	Paediatric Oral/Nasal Magill, Endotracheal Tube	10/case
35115	5.0mm	Paediatric Oral/Nasal Magill, Endotracheal Tube	10/case
35116	5.5mm	Paediatric Oral/Nasal Magill, Endotracheal Tube	10/case
35117	6.0mm	Paediatric Oral/Nasal Magill, Endotracheal Tube	10/case
35118	6.5mm	Paediatric Oral/Nasal Magill, Endotracheal Tube	10/case
35119	7.0mm	Paediatric Oral/Nasal Magill, Endotracheal Tube	10/case



HALYARD* MICROCUFF* ORAL CURVED ENDOTRACHEAL TUBE

Code	Inner Diameter	Description	Packaging
35161	3.0mm	Paediatric Oral Curved, Endotracheal Tube	10/case
35162	3.5mm	Paediatric Oral Curved, Endotracheal Tube	10/case
35163	4.0mm	Paediatric Oral Curved, Endotracheal Tube	10/case
35164	4.5mm	Paediatric Oral Curved, Endotracheal Tube	10/case
35165	5.0mm	Paediatric Oral Curved, Endotracheal Tube	10/case
35166	5.5mm	Paediatric Oral Curved, Endotracheal Tube	10/case
35167	6.0mm	Paediatric Oral Curved, Endotracheal Tube	10/case
35168	6.5mm	Paediatric Oral Curved, Endotracheal Tube	10/case
35169	7.0mm	Paediatric Oral Curved, Endotracheal Tube	10/case



PAEDIATRIC/NEONATE CLOSED SUCTION SYSTEMS

HALYARD* CLOSED SUCTION SYSTEMS – NEONATES

- TRACH CARE* Technology
- Gamma sterilised
- Lavage/Irrigation Port
- Locking Thumb Port
- "Y" and Elbow Adapters
- Multiple French sizes



HALYARD* ELBOW TUBE ADAPTERS

Code	Size (Outer Diameter)	Length	Description	Packaging
206-5	6 Fr / 2mm	30.5cm	Elbow, Endotracheal Tube Adapters (2.5mm, 3mm, 3.5mm)	20/case
207-5	7 Fr / 2.3mm	35.5cm	Elbow, Endotracheal Tube Adapters (2.5mm, 3mm, 3.5mm, 4mm)	20/case
208-5	8 Fr / 2.6mm	35.5cm	Elbow, Endotracheal Tube Adapters (3mm, 3.5mm, 4mm)	20/case
210-5	10 Fr / 3.3mm	40.5cm	Elbow, Endotracheal Tube Adapters (4.5mm, 5mm, 5.5mm)	20/case



HALYARD* Y TUBE ADAPTERS

Code	Size (Outer Diameter)	Length	Description	Packaging
195-5	5 Fr / 1.66mm	30.5cm	"Y" Endotracheal Tube Adapters (2mm, 2.5mm)	20/case
196-5	6 Fr / 2mm	30.5cm	"Y" Endotracheal Tube Adapters (2.5mm, 3mm, 3.5mm)	20/case
197-5	7 Fr / 2.3mm	30.5cm	"Y" Endotracheal Tube Adapters (3mm, 3.5mm, 4mm)	20/case
198-5	8 Fr / 2.6mm	30.5cm	"Y" Endotracheal Tube Adapters (3.5mm, 4mm, 4.5mm)	20/case
1910-5	10 Fr / 3.3mm	40.5cm	"Y" Endotracheal Tube Adapters (4.5mm, 5mm, 5.5mm)	20/case
1912-5	12 Fr / 4mm	40.5cm	"Y" Endotracheal Tube Adapters (5.5mm, 6mm, 6.5mm)	20/case




PAEDIATRIC

HALYARD* MULTI-ACCESS CATHETERS

- TRACH CARE* Technology
- Gamma sterilised
- Multi-Access Catheter designed to deliver surfactants
- Multiple French sizes

Code	Size (Outer Diameter)	"Y" Size/ Plug Size	Description	Packaging
1900	5 Fr / 1.66mm	2.5mm / 7.5mm, 3.0mm / 7.5mm, 3.5mm / 7.5mm	Mac Catheter	5/case



HALYARD* ORAL CARE SOLUTIONS

HALYARD* ORAL SUCTION SINGLE PATIENT COMPONENTS

- Multiple French sizes

Code	Size (Outer Diameter)	Description	Packaging
1221	6 Fr / 2mm	5-Oral Care Suction Catheters for Neonates / Paediatrics	50/case
1222	8 Fr / 2.6mm	5-Oral Care Suction Catheters for Neonates / Paediatrics	50/case
1223	10 Fr / 3.3mm	5-Oral Care Suction Catheters for Neonates / Paediatrics	50/case



INDEX

Code	Page	Code	Page	Code	Page	Code	Page
ADULT				PAEDIATRIC			
HALYARD* MICROCUFF* ENDOTRACHEAL TUBES				HALYARD* MICROCUFF* ENDOTRACHEAL TUBES			
35210	4	HALYARD* TURBO-CLEANING CLOSED SUCTION SYSTEMS		HALYARD* BRONCHIAL ASPIRATE SAMPLING CATHETERS		35111	16
35211	4	227-5	8	142A	12	35112	16
35212	4	22701356-5	9	143A	12	35113	16
35213	4	22703-5	8	HALYARD* ORAL CARE SOLUTIONS			
35214	4	227036-5	8	12240	13	35114	16
35215	4	227056-5	9	12245	13	35115	16
35216	4	227135-5	9	12246	13	35116	16
35217	4	22715-5	9	12249	13	35117	16
35218	4	227156-5	9	12252	13	35118	16
35220	4	22716-5	8	12254	13	35119	16
HALYARD* CLOSED SUCTION SYSTEMS				12255	13	35161	16
2155-5	6	2271603-5	8	123	14	35162	16
2160-5	5	22716036-5	8	12601	13	35163	16
21603-5	5	227166-5	8	12602	13	35164	16
216036-5	5	2276-5	8	12603	13	35165	16
21606-5	5	HALYARD* MULTI-ACCESS PORT CLOSED SUCTION SYSTEMS		12604	13	35166	16
21609-5	7	24001238-5	10	12606	13	35167	16
22012356-5	6	2401236-5	11	25928	15	35168	16
220123568-5	6	24012386-5	11	99785	13	35169	16
220135-5	6	240125-5	10	99786	13	NEONATES/ PAEDIATRIC	
2205-5	6	2401253-5	10	99789	13	HALYARD* CLOSED SUCTION SYSTEMS	
22056-5	6	240128-5	10	1900	18	1910-5	17
220568-5	6	2401286-5	11	1912-5	17	195-5	17
220569-5	7	2401438-5	10	196-5	17	197-5	17
22059-5	7	24014386-5	11	198-5	17	198-5	17
2210-5	5	240145-5	10	206-5	17	207-5	17
22103-5	5	2401453-5	10	208-5	17	210-5	17
221036-5	5	240148-5	10	HALYARD* ORAL CARE SOLUTIONS			
221037-5	7	2401486-5	11	1221	18	1222	18
22106-5	5	HALYARD* CLOSED SUCTION SYSTEMS ACCESSORIES		1223	18		
221069-5	5	11000	15				
22107-5	7	111	15				
22108-5	5	112	15				
22109-5	7	114	15				
2255-5	6	118	15				
2260-5	5	119	15				
22606-5	5	19000	15				
25-2210-5	6	19030	15				
		19035	15				
		24013-5	11				



Infection Prevention Product Solutions
Halyard* Education Foundation
Tools & Best Practice
Commitment to Excellence and Sustainability
Experienced Customer Support
Knowledgeable Sales Force
Innovation & Technology



Halyard Australia Pty. Limited
Phone: 1800 101 021

Halyard New Zealand
Phone: 0800 482 211

www.halyardhealth.com.au

www.halyardhealth.com.au



Halyard Health makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. © 2015 HYH. All rights reserved.
RH:1506:WFHY631