HOMEPUMP*

Ambulatory Infusion System

Ordering information

Homepump Eclipse*


Ideal for Antibiotic Delivery


E050500 50 ml vol, 50 ml/hr 48 E100500 100 ml vol, 50 ml/hr 48 E101000 100 ml vol, 100 ml/hr 48 E102000 100 ml vol, 200 ml/hr 48
E101000 100 ml vol, 100 ml/hr 48
E102000 100 ml vol 200 ml/hr 48
E251750 250 ml vol, 175 ml/hr 24
E252500 250 ml vol, 250 ml/hr 24
E401000 400 ml vol, 100 ml/hr 24
E402000 400 ml vol, 200 ml/hr 24
E502500 500 ml vol, 250 ml/hr 24

Homepump C-Series*

Ideal for Chemotherapy and Other Protocols

Model #	Size/Flow Rate	Duration	Qty/Box
C060020	60 ml vol, 2 ml/hr	30 hrs	24
C100005	100 ml vol, 0.5 ml/hr	8 days	24
C100020	100 ml vol, 2 ml/hr	2 days	24
C125050	125 ml vol, 5 ml/hr	1 day	24
C270010	270 ml vol, 1 ml/hr	11 days	24
C270020	270 ml vol, 2 ml/hr	5.5 days	5 24
C270050	270 ml vol, 5 ml/hr	2 days	24
C270100	270 ml vol, 10 ml/hr	1 day	24
C300060	300 ml vol, 6 ml/hr	2 days	24
Model #	Accessories Carry Pouch		Qty/Box
5001096	50-200 ml pumps		12
5001097	250-500 ml pumps		12
5001087	E-Clip		100





For customer service call: 800.448.3569 • 949.206.2700 Fax 949.206.2600

Rx only.

*Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc. or its affiliates. © 2012 KCWW. All Rights Reserved.

US Patent Nos 5,284,481C1; 5083,741. Foreign Pat. Pending.

There are inherent risks in all medical devices. Please refer to the product labeling for Indications, Cautions, Warnings and Contraindications. Failure to follow the product labeling could directly impact patient safety. Physician is responsible for prescribing and administering medications per instructions provided by the drug manufacturer. Refer to www.iflo.com for product safety Technical Bulletins. MK-00284

Sig Kimberly-Clark

Trusted Clinical Solutions*



www.iflo.com

HOMEPUMP* Ambulatory Infusion System



Homepump Eclipse^{*} Homepump C-Series^{*}

The infusion pump that's at home anywhere.



At home or anywhere, the Homepump^{*} system is the trusted brand you know.



For more than two decades the Homepump system has delivered cost-effective, patient minded benefits, with a simplicity of use that may save time for both clinicians and patients.

Homepump Eclipse* and Homepump C-Series* Ambulatory Infusion Systems

The Homepump system offers a complete family of products for your infusion needs.

- Eclipse system the pump designed for antibiotic therapy providing delivery times from 30 minutes to 5 hours.
- C-Series system the portable pump for continuous infusion with delivery times from 1-11 days for chemotherapy and other protocols.

The Homepump system combines high quality standards with the volumes and flow rates you need for successful drug delivery. The pump's proven ease of use allows patients to receive their therapy in the comfort of home.




*Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc. or its affiliates © 2012 KCWW. All Rights Reserved

Performance you can trust.

- Comprehensive drug compatibility and stability data.
- Patented, elastomeric membrane delivers performance and dependability.
- Pre-set flow rate no complex programming or drop counting simplifies life for clinicians as well as patients.
- Overfill/partial fill capacity allows dosing flexibility.

Compact features. Huge benefits.

- Discreet, portable and patient-friendly.
- Allows patient mobility for patients to maintain lifestyle at home.
- Non-electronic: no complex programming, maintenance, power cords, batteries or IV poles.
- Disposable design lets patients discard unit after use.
- Cost-efficient without sacrificing choice of volumes and flow rates.
- Small size requires less space to store both in pharmacy and in the patient's home.



World-class performance.

At home or on the road, the Homepump system lives up to its long history of performance, dependability, convenience and safety. Learn more about the infusion pump that's at home anywhere. For customer service call:

800-448-3569

www.iflo.com

HOMEPUMP* Ambulatory Infusion System