

# BD | BARD PRODUCT OFFERING



Integrated Medical Systems, Inc. carries a broad line of medical delivery solutions from BD-Bard in order to meet all your parenteral medication preparation and administration needs. By bringing together two leading portfolios, BD-Bard offers an expanded variety of medication delivery solutions and drive standardization across your health system. IMS is proud to provide an unprecedented selection of innovative, high-quality products to meet all your vascular access management needs. Click on product names to view the full product list in the IMS Portal.

## VASCULAR MANAGEMENT SYSTEMS

### CARE & MAINTENANCE

#### Pre-Filled IV Flush Syringes & Caps

##### [BD PosiFlush™ Pre-Filled Saline Syringe](#)

Pre-Filled Saline Syringe complies with major PICC manufacturer recommendations for flushing with a 10-mL diameter syringe. 0.5mL graduations. Sterilized product. Steam autoclaved. Needleless. DEHP-free, Latex-free, Preservative free. For single use only.

##### [BD PosiFlush™ SF Saline Flush Syringe](#)

Sterile Field (SF) Saline Flush Syringe is terminally sterilized in its peel pouch, enabling it to be aseptically presented to a sterile field. Typically used in oncology, interventional radiology and critical care, the syringe supports applications such as PICC or CVC insertions. It is intended to be used only for the flushing of indwelling vascular access devices.

##### [BD PosiFlush™ Pre-Filled Heparin Lock Flush Syringe](#)

Pre-Filled Heparin Lock Flush Syringe is available in two concentrations—10 USP units per mL and 100 USP units per mL—and in a variety of sizes and fill volumes to help support your catheter maintenance practice. Specifications Include: 0.5mL graduations. Sterilized product. Steam autoclaved. Needleless. DEHP-free, Latex-free, Preservative free. For single use only.



#### Saline & SF Saline

<a href="#">306544</a>	3mL BD PosiFlush™ Normal Saline Syringe w/ 10mL diameter	480/Case
<a href="#">306545</a>	5mL BD PosiFlush™ Normal Saline Syringe w/ 10mL diameter	480/Case
<a href="#">306546</a>	10mL BD PosiFlush™ Normal Saline Syringe w/ 10mL diameter	480/Case
<a href="#">306553</a>	10 mL BD PosiFlush™ Standard Plunger Rod Sterile Field Syringe	240/Case

#### Heparin

<a href="#">306413</a>	3 mL BD PosiFlush™ Heparin Lock Flush Syringe, 10 USP units/mL	480/Case
<a href="#">306414</a>	5 mL BD PosiFlush™ Heparin Lock Flush Syringe, 10 USP units/mL	480/Case
<a href="#">306423</a>	3 mL BD PosiFlush™ Heparin Lock Flush Syringe, 100 USP units/mL	480/Case
<a href="#">306424</a>	5 mL BD PosiFlush™ Heparin Lock Flush Syringe, 100 USP units/mL	480/Case

## BD PureHub™ Disinfecting Caps

This disinfecting cap is used for the prevention of CRBSI, which can result in fewer infections, more cost savings, more beds available for new admissions and potential time savings. Designed for compatibility with market-leading needle-free connectors prior to access and maintains a physical barrier to contamination for up to seven days. Maintains security and ease of use.



<a href="#">306596</a>	BD PureHub™ Disinfecting Cap, Singles	300/Box
<a href="#">306597</a>	BD PureHub™ Disinfecting Cap, IV Pole Strips [10 Count]	30 Strips/Box

## SiteScrub® IPA Device

Site-Scrub® IPA Device is intended for use on injection ports and female Luer hubs as a disinfecting cleaner. The Site-Scrub® designed to disinfect through a patent-pending friction scrub design. Patented foam “fingers” design with 70% isopropyl alcohol for 8x back-and-forth scrub over 10-seconds, 5-second dry.



<a href="#">3882100</a>	Site-Scrub IPA Device	5 Boxes of 100/ Cases
-------------------------	-----------------------	--------------------------

## Stabilization Devices & Dressings

### GuardIVa® Antimicrobial Hemostatic IV Dressing


GuardIVa® IV Dressing is an IV site care solution that combines the antimicrobial protection of chlorhexidine gluconate (CHG) with the proprietary hemostatic agent which controls surface bleeding from percutaneous catheters and vascular access sites.



<a href="#">FP-23-AD008</a>	GuardIVa Antimicrobial IV Dressing, 4mm Hole x 1" Disk	10/Box
-----------------------------	--	--------

### StatLock® PICC/CVC Stabilization Device

The StatLock® Stabilization Devices are proprietary, mechanical, IV catheter stabilization devices designed to minimize complications resulting from movement of a catheter and may help to reduce vessel trauma. The StatLock® eliminates securement with suture, accidental needlesticks, reduces unscheduled restarts, reduces catheter complications, and increases catheter dwell times.



<a href="#">VPPCSP</a>	StatLock PICC Plus Stabilization Device, Foam Anchor Pad, Sliding Posts, Crescent	50/Case
<a href="#">VPPPSP</a>	StatLock PICC Plus Stabilization Device, Pediatric, Foam Anchor Pad, Sliding Posts	50/Case
<a href="#">PIC0220</a>	StatLock PICC Plus Stabilization Device, Tricot Anchor Pad, Sliding Posts	50/Case
<a href="#">PIC0222</a>	StatLock PICC Plus Stabilization Device, Tricot Anchor Pad, Fixed Post	50/Case
<a href="#">IV0520</a>	StatLock IV Ultra Stabilization Device, Adult	50/Case

### **BD ChloraShield™ IV Dressings with CHG Antimicrobial**

The innovative BD ChloraShield™ IV Dressing with CHG Antimicrobial is designed to address some of the limitations of standard IV catheter dressings, including antimicrobial efficacy and moisture management. BeneHold™ adhesive technology formulated with CHG provides sustained antimicrobial activity against skin flora for up to 7 days. Improved moisture management increases the dressing's total fluid handling capacity by up to 40% compared to a standard film dressing. Thin, flexible material is easy to apply and remove, and is minimally irritating to the skin



<a href="#">41011</a>	ChloraShield™ IV Film Dressing w/ CHG Antimicrobial Within the Adhesive, 3" x 2.4" [for peripheral lines]	100/Box
<a href="#">410122</a>	Bordered & Slit ChloraShield™ IV Surround Dressing with CHG Antimicrobial Within the Adhesive, 3" x 2.4" [for peripheral lines]	100/Box
<a href="#">410123</a>	Bordered & Slit ChloraShield™ Dressing Designed for BD Nexiva Closed IV Catheter System, with CHG Antimicrobial & 2 sterile, pre-cut tape strips, 3.3" x 2.4" [for peripheral lines]	100/Box
<a href="#">410131</a>	ChloraShield™ IV Film Dressing w/ CHG Antimicrobial Within the Adhesive, 4" x 4.3" [for central lines]	100/Box
<a href="#">410132</a>	Bordered & Slit ChloraShield™ IV Surround Dressing with CHG Antimicrobial Within the Adhesive, 4" x 4.3" [for central lines]	100/Box

### **ChloraPrep® Patient Preoperative Skin Prep Applicator**

All ChloraPrep™ applicators contain 2% w/v chlorhexidine gluconate (CHG) in 70% v/v isopropyl alcohol (IPA). Available in different sizes to facilitate skin preparation for various preoperative and pre-injection procedures. An antiseptic approved for one-step preoperative skin preparation.




<a href="#">930480</a>	ChloraPrep® 1mL Clear Applicator with Sterile Solution	60/Box
<a href="#">930400</a>	ChloraPrep® 3mL Clear Applicator with Sterile Solution	25/Box
<a href="#">930415</a>	ChloraPrep® 3mL Hi-Lite Orange Applicator with Sterile Solution	25/Box
<a href="#">930700</a>	ChloraPrep® 10.5mL Clear Applicator with Sterile Solution	25/Box
<a href="#">930715</a>	ChloraPrep® 10.5mL Hi-Lite Orange Applicator with Sterile Solution	25/Box
<a href="#">930725</a>	ChloraPrep® 10.5mL Scrub Teal Applicator with Sterile Solution	100/Case
<a href="#">930815</a>	ChloraPrep® 26mL Hi-Lite Orange Applicator with Sterile Solution	25/Box
<a href="#">930825</a>	ChloraPrep® 26mL Scrub Teal Applicator with Sterile Solution	25/Box

## Disposable Syringes

### BD Syringes Only

Disposable syringes fitted with conventional or safety needles. Luer-Lok, Slip or Eccentric Tips available.

<a href="#">309628</a>	BD 1mL Syringe, Luer-Lok Tip	200/Box
<a href="#">309657</a>	BD 3mL Syringe, Luer-Lok Tip	200/Box
<a href="#">309656</a>	BD 3mL Syringe, Slip Tip	200/Box
<a href="#">309646</a>	BD 5mL Syringe, Luer-Lok Tip	125/Box
<a href="#">309647</a>	BD 5mL Syringe, Slip Tip	125/Box
<a href="#">302995</a>	BD 10mL Syringe, Luer-Lok Tip	200/Box
<a href="#">303134</a>	BD 10mL Syringe, Slip Tip	200/Box
<a href="#">302830</a>	BD 20mL Syringe, Luer-Lok Tip	40/Box
<a href="#">302831</a>	BD 20mL Syringe, Slip-Tip	48/Box
<a href="#">302832</a>	BD 30mL Syringe, Luer-Lock Tip	56/Box
<a href="#">302833</a>	BD 30mL Syringe, Slip-Tip	56/Box
<a href="#">309653</a>	BD 50ml Syringe, Luer-Lok Tip	40/Box
<a href="#">309654</a>	BD 50mL Syringe, Slip Tip	40/Box
<a href="#">309620</a>	BD 50mL Catheter Tip Syringe w/ Tip Shield	40/Box



## VASCULAR ACCESS DEVICES

### PERIPHERAL ACCESS

#### Peripheral IV Catheters

##### BD Nexiva™ Closed IV Catheter System

The BD Nexiva™ system is power-injection compatible, minimizing the need to exchange catheter or tubing. Single and dual ports and 18- to 22-G catheters are indicated for use in computed tomography (CT) power injection. Integrated system prevents leaking at the catheter-extension set connection.

##### Dual Port

<a href="#">383530</a>	24G x 0.56", Yellow	80/Case
<a href="#">383531</a>	24G x 0.75", Yellow	80/Case
<a href="#">383532</a>	22G x 1.00", Blue	80/Case
<a href="#">383536</a>	20G x 1.00", Pink	80/Case
<a href="#">383537</a>	20G x 1.25", Pink	80/Case
<a href="#">383538</a>	20G x 1.75", Pink	80/Case
<a href="#">383539</a>	18G x 1.25", Green	80/Case
<a href="#">383540</a>	18G x 1.75", Green	80/Case

##### Single Port

<a href="#">383510</a>	24G x 0.56", Yellow	80/Case
<a href="#">383511</a>	24G x 0.75", Yellow	80/Case
<a href="#">383512</a>	22G x 1.00", Blue	80/Case
<a href="#">383516</a>	20G x 1.00", Pink	80/Case



## BD Insyte™ Autoguard™ Shielded IV Catheter

The Insyte™ Autoguard™ incorporates advanced safety features including a push button that instantly retracts the needle, helping reduce the risk of accidental needlestick injuries. Notched needle technology lets you easily visualize blood return along the catheter shaft, which can enhance care for fragile vessels or vessels without venous tone. Available in winged or non-winged.



### Non-Winged

<a href="#">381412</a>	24G x 0.75", Yellow	50/Box
<a href="#">381423</a>	22G x 1.00", Blue	50/Box
<a href="#">381433</a>	20G x 1.00", Pink	50/Box
<a href="#">381434</a>	20G x 1.16", Pink	50/Box
<a href="#">381437</a>	20G x 1.88", Pink	50/Box
<a href="#">381444</a>	18G x 1.16", Green	50/Box
<a href="#">381447</a>	18G x 1.88", Green	50/Box
<a href="#">381454</a>	16G x 1.16", Gray	50/Box
<a href="#">381457</a>	16G x 1.77", Gray	50/Box
<a href="#">381467</a>	14G x 1.75", Orange	50/Box

### Winged

<a href="#">381512</a>	24G x 0.75", Yellow	50/Box
<a href="#">381523</a>	22G x 1.00", Blue	50/Box
<a href="#">381533</a>	20G x 1.00", Pink	50/Box
<a href="#">381534</a>	20G x 1.16", Pink	50/Box
<a href="#">381537</a>	20G x 1.88", Pink	50/Box
<a href="#">381544</a>	18G x 1.16", Green	50/Box
<a href="#">381547</a>	18G x 1.88", Green	50/Box
<a href="#">381554</a>	16G x 1.16", Gray	50/Box
<a href="#">381557</a>	16G x 1.77", Gray	50/Box

## BD Insyte™ Autoguard™ Shielded IV Catheter w/ Blood Control Technology

Insyte™ Autoguard™ BC shielded IV catheter with BC technology is clinically demonstrated to reduce the risk of blood exposure by 95%. BD blood control technology minimizes blood exposure, cleanup time and use of supplies associated with blood spills.



### Non-Winged

<a href="#">382512</a>	24G x 0.75", Yellow	50/Box
<a href="#">382523</a>	22G x 1.00", Blue	50/Box
<a href="#">382533</a>	20G x 1.00", Pink	50/Box
<a href="#">382534</a>	20G x 1.16", Pink	50/Box
<a href="#">382537</a>	20G x 1.88", Pink	50/Box
<a href="#">382544</a>	18G x 1.16", Green	50/Box
<a href="#">382547</a>	18G x 1.88", Green	50/Box
<a href="#">382554</a>	16G x 1.16", Gray	50/Box
<a href="#">382557</a>	16G x 1.77", Gray	50/Box

### Winged

<a href="#">382612</a>	24G x 0.75", Yellow	50/Box
<a href="#">382623</a>	22G x 1.00", Blue	50/Box
<a href="#">382633</a>	20G x 1.00", Pink	50/Box
<a href="#">382634</a>	20G x 1.16", Pink	50/Box
<a href="#">382637</a>	20G x 1.88", Pink	50/Box
<a href="#">382644</a>	18G x 1.16", Green	50/Box
<a href="#">382647</a>	18G x 1.88", Green	50/Box
<a href="#">382654</a>	16G x 1.16", Gray	50/Box
<a href="#">382657</a>	16G x 1.77", Gray	50/Box

## CENTRAL ACCESS

### Port Access Needles - Power Injectable

#### PowerLoc® Safety Infusion Set

PowerLoc® Safety Infusion Set's low profile, small footprint and design make it an ideal needle for dressing the access site. The small size of the needle and the downward angle of the tubing coming off of the back of the device help to maintain dressing integrity.



#### With Y-Injection Site

<a href="#">0672034</a>	20G x .75", 5mL/s	20/Case
<a href="#">0672010</a>	20G x 1", 5mL/s	20/Case
<a href="#">0672015</a>	20G x 1", 5mL/s	20/Case

#### Without Y-Injection Site

<a href="#">0652034</a>	20G x 0.75", 5mL/s	20/Case
<a href="#">0652010</a>	20G x 0.75", 5mL/s	20/Case
<a href="#">0652015</a>	20G x 1", 5mL/s	20/Case

#### PowerLoc® EZ Safety Infusion Set

The PowerLoc® EZ Safety Infusion Set's small needle size helps to maintain dressing integrity. Integrated Bio-Bag designed to minimize contact with blood and fluids following safety activation. Padded footprint designed for patient comfort.



#### With Y-Injection Site

<a href="#">SHW20-75Y</a>	PowerLoc EZ 20G x 0.75"	25/Case
<a href="#">SHW20-100Y</a>	PowerLoc EZ 20G x 1"	25/Case
<a href="#">SHW20-150Y</a>	PowerLoc EZ 20G x 1.5"	25/Case
<a href="#">SHW22-75Y</a>	PowerLoc EZ 22G x 0.75"	25/Case
<a href="#">SHW22-100Y</a>	PowerLoc EZ 22G x 1"	25/Case
<a href="#">SHW22-150Y</a>	PowerLoc EZ 22G x 1.5"	25/Case

#### Without Y-Injection Site

<a href="#">SHW20-75</a>	PowerLoc EZ 20G x 0.75"	25/Case
<a href="#">SHW20-100</a>	PowerLoc EZ 20G x 1"	25/Case
<a href="#">SHW20-150</a>	PowerLoc EZ 20G x 1.5"	25/Case
<a href="#">SHW22-75</a>	PowerLoc EZ 22G x 0.75"	25/Case
<a href="#">SHW22-100</a>	PowerLoc EZ 22G x 1"	25/Case
<a href="#">SHW22-150</a>	PowerLoc EZ 22G x 1"	25/Case

#### PowerLoc® Max Safety Infusion Set

The textured handle is easy to grasp and rotates to enable maneuvering and positioning of tubing to enable comfortable and proper site dressing. Small footprint helps minimize dressing materials while the clear plastic base allows clinicians to view and inspect the insertion site throughout the course of therapy.



#### With Y-Injection Site

<a href="#">0132075</a>	20G x 0.75", 5mL/s	25/Case
<a href="#">0132275</a>	22G x 0.75", 2mL/s	25/Case
<a href="#">0132010</a>	20G x 1", 5mL/s	25/Case
<a href="#">0132210</a>	22G x 1", 2mL/s	25/Case
<a href="#">0132015</a>	20G x 1", 5mL/s	25/Case
<a href="#">0132215</a>	22G x 1", 2mL/s	25/Case

#### Without Y-Injection Site

<a href="#">0142075</a>	20G x 0.75", 5mL/s	25/Case
<a href="#">0142275</a>	22G x 0.75", 2mL/s	25/Case
<a href="#">0142010</a>	20G x 1", 5mL/s	25/Case
<a href="#">0142210</a>	22G x 1", 2mL/s	25/Case
<a href="#">0142015</a>	20G x 1.5", 5mL/s	25/Case
<a href="#">0142215</a>	22G x 1.5", 2mL/s	25/Case

## Port Access Needles - Safety

### HuberPlus® Safety Infusion Set

HuberPlus® reduces the incidence of needlestick injuries by giving you control in hand placement and grip and single handed needle access & safety activation.



#### With Injection Site & Needleless Injection Cap

<a href="#">011934NY</a>	Huber Plus 19G x .75"	25/Case
<a href="#">012034NY</a>	Huber Plus 20G x .75"	25/Case
<a href="#">011901NY</a>	Huber Plus 19G x 1"	25/Case
<a href="#">012001NY</a>	Huber Plus 20G x 1"	25/Case
<a href="#">012201NY</a>	Huber Plus 22G x 1"	25/Case
<a href="#">011912NY</a>	Huber Plus 19G x 1.25"	25/Case
<a href="#">012012NY</a>	Huber Plus 20G x 1.25"	25/Case
<a href="#">011915NY</a>	Huber Plus 19G x 1.5"	25/Case
<a href="#">012015NY</a>	Huber Plus 20G x 1.5"	25/Case
<a href="#">012215NY</a>	Huber Plus 22G x 1.5"	25/Case

#### Without Y-Injection Site

<a href="#">012005</a>	Huber Plus 20G x .5"	25/Case
<a href="#">012205</a>	Huber Plus 22G x .5"	25/Case
<a href="#">011934</a>	Huber Plus 19G x .75"	25/Case
<a href="#">012034</a>	Huber Plus 20G x .75"	25/Case
<a href="#">012234</a>	Huber Plus 22G x .75"	25/Case
<a href="#">011901</a>	Huber Plus 19G x 1"	25/Case
<a href="#">012001</a>	Huber Plus 20G x 1"	25/Case
<a href="#">012201</a>	Huber Plus 22G x 1"	25/Case
<a href="#">011912</a>	Huber Plus 19G x 1.25"	25/Case
<a href="#">012012</a>	Huber Plus 20G x 1.25"	25/Case
<a href="#">012212</a>	Huber Plus 22G x 1.25"	25/Case
<a href="#">011915</a>	Huber Plus 19G x 1.5"	25/Case
<a href="#">012015</a>	Huber Plus 20G x 1.5"	25/Case
<a href="#">012215</a>	Huber Plus 22G x 1.5"	25/Case

### EZ Huber® Safety Infusion Set

The low profile, small footprint and design make it an ideal needle for dressing the access site. Has a dual utility safety mechanism designed to protect clinicians from needle base stick injuries and minimize clinician contact with blood and fluids.



#### With Wings

<a href="#">SHW19-75YS</a>	EZ Huber Safety, 19G x 0.75"	25/Case
<a href="#">SHW19-100YS</a>	EZ Huber Safety, 19G x 1"	25/Case
<a href="#">SHW19-150YS</a>	EZ Huber Safety, 19G x 1.5"	25/Case
<a href="#">SHW20-75YS</a>	EZ Huber Safety, 20G x 0.75"	25/Case
<a href="#">SHW20-100YS</a>	EZ Huber Safety, 20G x 1"	25/Case
<a href="#">SHW20-150YS</a>	EZ Huber Safety, 20G x 1.5"	25/Case
<a href="#">SHW22-50YS</a>	EZ Huber Safety, 22G x .50"	25/Case
<a href="#">SHW22-75YS</a>	EZ Huber Safety, 22G x 0.75"	25/Case
<a href="#">SHW22-100YS</a>	EZ Huber Safety, 22G x 1"	25/Case
<a href="#">SHW22-150YS</a>	EZ Huber Safety, 22G x 1.5"	25/Case

#### Without Wings

<a href="#">SH19-75YS</a>	EZ Huber Safety, 19G x 0.75"	25/Case
<a href="#">SH19-100YS</a>	EZ Huber Safety, 19G x 1"	25/Case
<a href="#">SH19-150YS</a>	EZ Huber Safety, 19G x 1.5"	25/Case
<a href="#">SH20-75YS</a>	EZ Huber Safety, 20G x 0.75"	25/Case
<a href="#">SH20-100YS</a>	EZ Huber Safety, 20G x 1"	25/Case
<a href="#">SH20-150YS</a>	EZ Huber Safety, 20G x 1.5"	25/Case
<a href="#">SH22-50YS</a>	EZ Huber Safety, 22G x .50"	25/Case
<a href="#">SH22-75YS</a>	EZ Huber Safety, 22G x 0.75"	25/Case
<a href="#">SH22-100YS</a>	EZ Huber Safety, 22G x 1"	25/Case
<a href="#">SH22-150YS</a>	EZ Huber Safety, 22G x 1.5"	25/Case

## SafeStep® Huber Needle Set

The textured handle is easy to grasp and rotates to enable maneuvering and positioning of tubing to enable comfortable and proper site dressing. Also the small footprint helps minimize dressing materials while the clear plastic base allows clinicians to view and inspect the insertion site throughout the course of therapy.



### With Y-Injection & Needleless Injection Cap


<a href="#">LH-0039YN</a>	SafeStep 19G x 0.5"	25/Case
<a href="#">LH-0033YN</a>	SafeStep 19G x 0.75"	25/Case
<a href="#">LH-0034YN</a>	SafeStep 19G x 1"	25/Case
<a href="#">LH-0040YN</a>	SafeStep 19G x 1.5"	25/Case
<a href="#">LH-0037YN</a>	SafeStep 20G x 0.5"	25/Case
<a href="#">LH-0031YN</a>	SafeStep 20G x 0.75"	25/Case
<a href="#">LH-0032YN</a>	SafeStep 20G x 1"	25/Case
<a href="#">LH-0038YN</a>	SafeStep 20G x 1.5"	25/Case
<a href="#">LH-0035YN</a>	SafeStep 22G x 0.5"	25/Case
<a href="#">LH-0029YN</a>	SafeStep 22G x 0.75"	25/Case
<a href="#">LH-0030YN</a>	SafeStep 22G x 1"	25/Case
<a href="#">LH-0036YN</a>	SafeStep 22G x 1.5"	25/Case

### Without Y-Injection & Needleless Injection Cap

<a href="#">LH-0039</a>	SafeStep 19G x 0.5"	25/Case
<a href="#">LH-0033</a>	SafeStep 19G x 0.75"	25/Case
<a href="#">LH-0034</a>	SafeStep 19G x 1"	25/Case
<a href="#">LH-0040</a>	SafeStep 19G x 1.5"	25/Case
<a href="#">LH-0037</a>	SafeStep 20G x 0.5"	25/Case
<a href="#">LH-0031</a>	SafeStep 20G x 0.75"	25/Case
<a href="#">LH-0032</a>	SafeStep 20G x 1"	25/Case
<a href="#">LH-0038</a>	SafeStep 20G x 1.5"	25/Case
<a href="#">LH-0035</a>	SafeStep 22G x 0.5"	25/Case
<a href="#">LH-0029</a>	SafeStep 22G x 0.75"	25/Case
<a href="#">LH-0030</a>	SafeStep 22G x 1"	25/Case
<a href="#">LH-0036</a>	SafeStep 22G x 1.5"	25/Case

## MiniLoc® Safety Infusion Set

The MiniLoc® Safety Infusion Set's low profile and small footprint design make it the ideal for dressing access sites. The small size of the needle and the downward angle of the tubing coming off of the back of the device help to maintain dressing integrity.



### With Y-Injection Site

<a href="#">0681950</a>	MiniLoc 19G x .50"	20/Case
<a href="#">0682050</a>	MiniLoc 20G x .50"	20/Case
<a href="#">0682250</a>	MiniLoc 22G x .50"	20/Case
<a href="#">0681934</a>	MiniLoc 19G x .75"	20/Case
<a href="#">0682034</a>	MiniLoc 20G x .75"	20/Case
<a href="#">0682234</a>	MiniLoc 22G x .75"	20/Case
<a href="#">0681910</a>	MiniLoc 19G x 1"	20/Case
<a href="#">0682010</a>	MiniLoc 20G x 1"	20/Case
<a href="#">0682210</a>	MiniLoc 22G x 1"	20/Case
<a href="#">0681915</a>	MiniLoc 19G x 1"	20/Case
<a href="#">0682015</a>	MiniLoc 20G x 1"	20/Case
<a href="#">0682215</a>	MiniLoc 22G x 1"	20/Case

### Without Y-Injection Site

<a href="#">0631950</a>	MiniLoc 19G x .50"	20/Case
<a href="#">0632050</a>	MiniLoc 20G x .50"	20/Case
<a href="#">0632250</a>	MiniLoc 22G x .50"	20/Case
<a href="#">0631934</a>	MiniLoc 19G x .75"	20/Case
<a href="#">0632034</a>	MiniLoc 20G x .75"	20/Case
<a href="#">0632234</a>	MiniLoc 22G x .75"	20/Case
<a href="#">0631910</a>	MiniLoc 19G x 1"	20/Case
<a href="#">0632010</a>	MiniLoc 20G x 1"	20/Case
<a href="#">0632210</a>	MiniLoc 22G x 1"	20/Case
<a href="#">0631915</a>	MiniLoc 19G x 1"	20/Case
<a href="#">0632015</a>	MiniLoc 20G x 1"	20/Case
<a href="#">0632215</a>	MiniLoc 22G x 1"	20/Case



## INJECTION

### BD Eclipse™ Needle

The BD Eclipse™ needle is a pivoting-shield needle for efficient safety during activation. The needle is bevel oriented to the safety cover to facilitate low-angle injections.



#### BD Eclipse Needle


<a href="#">305757</a>	30G x .50" needle	100/Box
<a href="#">305758</a>	27G x .50" needle	100/Box
<a href="#">305759</a>	25G x .625" needle	100/Box
<a href="#">305761</a>	25G x 1" needle	100/Box
<a href="#">305767</a>	25G x 1.5" needle	100/Box
<a href="#">305762</a>	23G x 1" needle	100/Box
<a href="#">305768</a>	22G x 1" needle	100/Box
<a href="#">305763</a>	22G x 1.5" needle	100/Box
<a href="#">305764</a>	21G x 1" needle	100/Box
<a href="#">305765</a>	21G x 1.5" needle	100/Box
<a href="#">305766</a>	18G x 1.5" needle	100/Box

#### BD Eclipse Needle with Syringe

<a href="#">305778</a>	Syringe 1mL Luer Lock, Eclipse Needle 30G x .50"	50/Box
<a href="#">305789</a>	Syringe 1mL Luer Lock, Eclipse Needle 27G x .50"	50/Box
<a href="#">305780</a>	Syringe 1mL Luer Lock, Eclipse Needle 25G x .625"	50/Box
<a href="#">305779</a>	Syringe 3mL Luer Lock, Eclipse Needle 21G x 1"	50/Box
<a href="#">305781</a>	Syringe 3mL Luer Lock, Eclipse Needle 25G x .625"	50/Box
<a href="#">305783</a>	Syringe 3mL Luer Lock, Eclipse Needle 22G x 1.5"	50/Box
<a href="#">305784</a>	Syringe 3mL Luer Lock, Eclipse Needle 21G x 1.5"	50/Box
<a href="#">305787</a>	Syringe 3mL Luer Lock, Eclipse Needle 25G x 1"	50/Box
<a href="#">305788</a>	Syringe 3mL Luer Lock, Eclipse Needle 22G x 1"	50/Box
<a href="#">305785</a>	Syringe 5mL Luer Lock, Eclipse Needle 22G x 1.5"	50/Box
<a href="#">305786</a>	Syringe 10mL Luer Lock, Eclipse Needle 22G x 1.5"	50/Box

### BD SafetyGlide™ Shielding Hypodermic Needle

The BD SafetyGlide™ shielding hypodermic needle features BD Activation-Assist™ technology, making a safe injection easier and more efficient. Needles are sterile.



#### Needle

<a href="#">305921</a>	27G x .625" needle	50/Box
<a href="#">305901</a>	25G x .625" needle	50/Box
<a href="#">305916</a>	25G x 1" needle	50/Box
<a href="#">305902</a>	23G x 1" needle	50/Box
<a href="#">305900</a>	22G x 1.5" needle	50/Box
<a href="#">305915</a>	21G x 1" needle	50/Box
<a href="#">305917</a>	21G x 1.5" needle	50/Box
<a href="#">305918</a>	18G x 1.5" needle	50/Box

#### Needle with Syringe

<a href="#">305927</a>	1mL Syringe, 27G x .625" needle	50/Box
<a href="#">305903</a>	1mL Syringe, 25G x .625" needle	50/Box
<a href="#">305924</a>	3mL Syringe, 25G x 1" needle	50/Box
<a href="#">305904</a>	3mL Syringe, 25G x .625" needle	50/Box
<a href="#">305905</a>	3mL Syringe, 23G x 1" needle	50/Box
<a href="#">305906</a>	3mL Syringe, 22G x 1.5" needle	50/Box
<a href="#">305909</a>	3mL Syringe, 21G x 1.5" needle	50/Box
<a href="#">305907</a>	5mL Syringe, 22G x 1.5" needle	50/Box
<a href="#">305908</a>	10mL Syringe, 22G x 1.5" needle	50/Box

## BD Integra™ Retracting Syringe

The BD Integra™ Syringe combines innovative technology with the BD quality and comfort you trust. The BD Integra™ Syringe with Tru-Lok™ technology is designed to have significantly lower waste space volume as compared to a traditional BD luer connection.



### 3mL Syringe with Detachable Needle

<a href="#">305269</a>	25G x .625"	100/Box	<a href="#">305272</a>	22G x 1.5"	100/Box
<a href="#">305270</a>	25G x 1"	100/Box	<a href="#">305273</a>	21G x 1"	100/Box
<a href="#">305271</a>	23G x 1"	100/Box	<a href="#">305274</a>	21G x 1.5"	100/Box

### 3mL Syringe Only

<a href="#">305283</a>	Integra 3mL Syringe	100/Box
------------------------	---------------------	---------

## Sharps Containers

### Patient/Exam Room Sharps Collectors

All BD Patient Room Collectors have visible fill-level indicators and handles for safe transportation. Compatible with locking cabinets to provide security and prevent unwanted access to collector contents; Counterbalance door facilitates hands-free disposal and is designed to improve tamper resistance and to reduce the potential of sharps injuries.

<a href="#">305517</a>	5.4 Quart Counterbalance Door Red	20/Case
<a href="#">305426</a>	5.4 Quart Counterbalance Door Red	12/Case
<a href="#">305436</a>	3 Gallon Counterbalance Door Red	10/Case
<a href="#">305444</a>	5.4 Quart Counterbalance Side Entry Pearl	20/Case



### Multi-Use Nestable Sharps Collectors

Clear tops with multiple opening styles to accommodate different sizes of waste and environments; Handles designed for safe transport; Funnel-style lids limit accessibility to the collector contents and minimize spillage if the collector tips over; Open-top lids feature large diameter openings to accommodate disposal needs in non-patient/public areas.

<a href="#">305344</a>	8 Quart Hinge Cap with Petals	24/Case
<a href="#">305457</a>	6 Gallon Hinge Cap	12/Case
<a href="#">305465</a>	6 Gallon Hinge Cap with Petals	12/Case



### Multi-Use One Piece Sharps Collectors

Key product features include: Handles designed for safe movement and transportation; Hinge caps designed to provide temporary and final closure features.

<a href="#">305477*</a>	5 Gallon Hinge Cap with Petals * Cannot be autoclave and is not vented	8/Case
<a href="#">305488</a>	3.3 Quart Hinge Cap with Petals	24/Case
<a href="#">305577</a>	5 Gallon with Screw Cap	8/Case



## Phlebotomy Sharps Collectors

Built-in standard, safety and pen needle removal ports designed to prevent contact with sharps during disposal; Clear tops featuring a clear view of the fill line; Temporary and final closures for safety and security. Can be used with brackets for added stability.

<a href="#">305635</a>	1 Quart with Needle Ports	60/Case
<a href="#">305487</a>	1.5 Quart Needle Port	36/Case



## Chemotherapy Sharps Collectors

Gasket or liquid-absorbing materials to prevent leakage; Visual, audible, tactile perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering; Clear lids with multiple opening styles accommodate different types of waste and environments.

<a href="#">305076</a>	3 Gallon Hinge Cap One Piece	12/Case
<a href="#">305493</a>	5 Gallon Hinge Cap	8/Case
<a href="#">305603</a>	9 Gallon Hinge Top with Gasket	8/Case
<a href="#">305614</a>	17 Gallon Hinge Cap with Gasket	5/Case



### Welcome to the IMS Customer Portal!

Please select the section of the portal you would like to go to.



**Shop for Disposables!**  
Shop from our huge catalog of supplies.



**Manage My Rental Fleet!**  
View your items and available equipment.



**My Biomed Toolbox!**  
Manage Facility Owned Equipment



**Medtrakk**  
Track all of your assets in one place.

Click on the item numbers above to view in MyIMS.Cloud. **MyIMS.cloud** now makes it easier than ever to shop from our extensive catalog of medical disposables and devices, all in one place. If you are not a current customer of IMS and would like more information on BD/Bard Vascular Access Products, please call customer service at (800)755-3800 to be connected with an IMS Sales Representative in your area. With patients and providers in mind, IMS has you covered.

## Small Enough to Care, Big Enough to Deliver

If you are not a current customer of IMS and would like more information on BD/Bard Vascular Access Products, please contact us at (800) 755-3800. We will connect you with an IMS Sales Representative in your area.

**Corporate Office**  
594 Territorial Drive, Unit A  
Bolingbrook, IL 60440  
1.800.755.3800  
www.integratedmedsys.com