


Perifix®

0.2µm Flat Epidural Filter + PinPad



High performance filtration
and ergonomic securement

B|BRAUN

Perifix® 0.2µm Flat Epidural Filter + PinPad



0.2µm filter

Anti-bacterial filtration
Particulate filtration

Bright Yellow Filter Body

Assists healthcare workers in distinguishing the epidural filter from IV filters

Flat Filter Style

Promotes patient comfort


Optional PinPad

Enables optimal positioning of filter and attached epidural catheter connector
Minimizes direct manipulation of connector
Provides strain relief for catheter



Technical Specifications

Residual volume 0.45 ml
Rated burst pressure 7 bar
Filter surface area 4cm²
Porosity 0.2µm



Product Code	Item Number	Description	Case Quantity
FEF100	415000	Perifix 0.2µm flat epidural filter	25
FPP100	415004	Perifix PinPad adhesive pad for skin	25

*For more product information, contact your
B. Braun representative, or call 1-800-BBRAUN2
(1-800-227-2862)*

B | BRAUN

B. Braun Medical Inc.
824 Twelfth Avenue
Bethlehem, PA 18018

Rx only

Perifix is a registered trademark of B. Braun Melsungen AG
© 2001 B. Braun Medical Inc.
PC3-02 12/01

1-800-BBRAUN2
www.bbraunusa.com/perifix