

# MaxZero™ needleless connector

Prolonged catheter patency and effective disinfection unite to provide a proven seven-day change-out protocol, helping keep your patient's line closed longer.

## Promotes catheter patency

### Neutral clamping sequence

No clamping is required. Anti-reflux action helps prevent blood from entering the catheter upon disconnection.



### Blood reflux prevention when you need it

The connector helps maintain line patency without positive pressure flushing and clamping.



## Enables complete disinfection

### Proven three-second disinfection

A solid, sealed surface helps prevent contamination from entering the line.



### Bacteria has nowhere to hide

The visible fluid path and high flow straight fluid channels enhance flushing.



**Prepare. Protect. Preserve.**

Cat. no.	Description	Tubing	Length	Qty. (cs.)
MZ1000-07	MaxZero needleless connector			100
MZ5301	Pressure-rated small bore extension set, bonded MaxZero connector, slide clamp, spin male luer lock, .40 mL	Small bore	7" (18 cm)	50
MZ5302	Pressure-rated small bore extension set, removable MaxZero connector, slide clamp, spin male luer lock, .40 mL	Small bore	7" (18 cm)	50
MZ5303	Pressure-rated standard bore extension set, bonded MaxZero connector, slide clamp, spin male luer lock, .50 mL	Standard bore	7" (18 cm)	50
MZ5304	Pressure-rated standard bore extension set, removable MaxZero connector, slide clamp, spin male luer lock, .50 mL	Standard bore	7" (18 cm)	50
MZ5305	Pressure-rated standard bore extension set, bonded MaxZero connector, spin male luer lock (no clamp), .50 mL	Standard bore	7" (18 cm)	50
MZ5306	Pressure-rated small bore extension set, bonded MaxZero connector, spin male luer lock (no clamp), .40 mL	Small bore	7" (18 cm)	50
MZ5307	Bifuse pressure rated extension set, 2 bonded MaxZero connectors, 2 slide clamps, spin male luer lock, .80 mL	Small bore	7" (18 cm)	50
MZ5308	Bifuse pressure rated extension set, 2 bonded MaxZero connectors, 2 slide clamps, spin male luer lock, .80 mL	Standard bore	7" (18 cm)	50

A variety of extension sets are available



To request more information or to place an order, contact your sales representative at **800.482.4822** or visit [carefusion.com](http://carefusion.com).

CareFusion  
Vernon Hills, IL

[carefusion.com](http://carefusion.com)

