

Baxter

ExactaMix Pro
AUTOMATED COMPOUNDER

A Proven Leader in Compounding



Baxter has been innovating parenteral nutrition compounding since 1982

With diverse needs of different pharmacies, care settings and patients, Baxter offers a comprehensive portfolio of clinical nutrition products from ExactaMix, ExactaMix Pro and Q-NRG+ devices, to source ingredients and licensed, ready-to-use solutions – simplifying the complex so you can provide reliable, quality PN therapy efficiently.

ExactaMix Pro has been designed to meet the evolving needs of the modern pharmacy – combining the trusted features you rely on with improved processing power, ease-of-use and an intuitive interface for the efficient, reliable delivery of automated compounding, all to help meet patient needs.

BAXTER DISPOSABLES

Disposables validated for use on ExactaMix and ExactaMix Pro automated compounders




EVA Single Chamber Patient Bags

- Sterile, non-DEHP and latex free
- Single-layer clear or light protected; Multi-layer clear
- Available in 8 sizes: 125 mL, 250 mL, 500 mL, 1000 mL, 2000 mL, 3000 mL, 4000 mL and 5000 mL




EVA Dual Chamber Patient Bags

- Sterile, non-DEHP and latex free
- Single-layer clear
- Available in 4 sizes: 650+150 mL, 1300+300 mL, 2000+500 mL and 3000+500 mL



Valve Sets

- Custom design for easy connection
- Features numbered ports and large-bore outlet connector
- Sterile and non-DEHP
- Available in 24 ports for EM2400 compounder



Inlets

- PVC tubing used to connect valve set to source containers
- Latex-free and non-DEHP
- Available in:
 - Vented for bottles and vials – macro and micro volume
 - Non-vented for bags – macro and micro volume
 - Luer connection for syringes



Accessories

- Various caps and connectors for use in conjunction with ExactaMix and ExactaMix Pro disposables

Products may vary by country. Please consult your local Baxter representative for local product availability.

Rx Only: For safe and proper use of the devices mentioned herein, refer to the appropriate Instructions for Use or Operator Manual.

Q-NRG+:

CE
0476

EVA Bags:

CE CE CE CE
0123 0051 0483

Baxter

Baxter, ExactaMix and ExactaMix Pro are trademarks of Baxter International Inc. or its subsidiaries. Q-NRG+ is a trademark of COSMED.

US-CN8-220045 V1.0 12/22