



# Nutritional Delivery Product Catalog

**Kangaroo™: The complete enteral feeding solution**

Quality products from the leader in innovative enteral feeding pump technology







# Nutritional Delivery Products

The Kangaroo™ portfolio offers a safe, innovative and complete line of enteral feeding products for use across the continuum of care. The Kangaroo™ brand has a long history of providing reliable and trusted products for clinicians and caregivers to optimize clinical efficiency. Our high-performing products are preferred by both clinicians and procurement professionals.†

Transitioning to ENFit®? Be sure to note the product description throughout the catalog for ENFit® design. Where applicable, the ENFit® version of the code is listed next to the non-ENFit® version of the code for ease of reference during a conversion.

## American Society of Parenteral and Enteral Nutrition (ASPEN) 2016 Clinical Guidelines<sup>§</sup>

ASPEN developed standards of practice for nutrition support. The clinical guidelines include:

- Nutrition support therapy in the form of early enteral nutrition (EN) should be initiated within 24 to 48 hours in the critically ill patient who is unable to maintain volitional intake.
- Enteral feeding protocols should be designed and implemented to increase the overall percentage of goal calories provided.
- Feeding tubes should be flushed with water after each feeding and medication administration, in between medications and every four hours during a continuous infusion.

† Customer Driven Growth Survey; March 2016:96-98.

§ 2016 Guidelines for the Provision and Assessment of Nutrition Support Therapy in the Adult Critically Ill Patient JPEN 2016 Vol 40, Issue 2, pp. 159 - 211.

<b>Enteral Feeding Pumps and Pump Sets</b>	<b>4</b>
Kangaroo™ ePump, Joey and Connect Enteral Feeding Pumps	4
<b>Extension Sets, Gravity Sets and Irrigation Kits</b>	<b>9</b>
Kangaroo™ Extension Sets	9
Kangaroo™ Gravity Sets	9
Kangaroo™ I.V. Pouch and Irrigation Kits	9
<b>Enteral Access Devices</b>	<b>10</b>
Kangaroo™ Feeding Tube with IRIS Technology	10
Kangaroo™ Nasogastric Feeding Tubes	11
Kangaroo™ Naso-Jejunal Feeding Tubes	12
Kangaroo™ Neonatal and Pediatric Feeding Tubes	13
Kangaroo™ Neonatal and Pediatric Feeding System	14
Kangaroo™ Gastrostomy Tubes	16
Kangaroo™ Safety P.E.G. Kits	17
Kangaroo™ Skin Level Tubes	20
Accessories	21
<b>Pharmacy Solutions</b>	<b>22</b>
Monoject™ Oral/Enteral Syringes	22
Monoject™ Syringe Tip Caps	23
Monoject™ Medicine Straws	23
Monoject™ UV/Light Inhibiting Transport Bags	23
Monoject™ ENFit® Universal Bottle Adapters	23
IVA™ Syringe Seals	23
<b>GI Products</b>	<b>24</b>
Salem Sump™ Tubes	24
Stomach Tubes	26
CO <sub>2</sub> Detectors	29
<b>Index</b>	<b>30</b>

## Kangaroo™ ePump, Joey and Connect

### Kangaroo™ ePump Enteral Feeding Pump

Features feed and flush technology that automates flushing for tube patency and delivery of hydration.

- Fully programmable variable feed and flush pump
- Flow rate: 1 to 400 mL/hr in 1 mL increments
- Flush rate: 1960 mL/hr
- Flush volume: 10 to 500 mL/hr in 1 mL increments
- Flush interval: 1 to 24 hr in 1 hr increments
- Multilingual: 15 languages
- Battery life: 15 hours
- Weight without pole clamp: 2.4 lbs
- Weight with pole clamp: 2.7 lbs
- Feeding history up to 72 hours
- Attitude independent
- Intuitive user interface
- Auto-prime



382400



773656



773662



716154

Item code	Description	Ship case
<b>Pump</b>		
382400	Kangaroo™ ePump with Pole Clamp	1
<b>Pump Accessories</b>		
382491	Kangaroo™ ePump Power Adapter	1
382492	Kangaroo™ ePump Pole Clamp	1
1054715	Kangaroo™ ePump Battery	1

### Kangaroo™ ePump Enteral Feeding Sets

- Anti-free flow design
- DEHP-free
- ENFit® connection

Item code	Description	Ship case
<b>Pump sets</b>		
772055	Kangaroo™ ePump 500 mL Pump Set	30
773656	Kangaroo™ ePump 1000 mL Pump Set	30
674668	Kangaroo™ ePump 500 mL Feed Set with 500 mL Flush Bag	30
773662	Kangaroo™ ePump 1000 mL Feed Set with 1000 mL Flush Bag	30
716154	Kangaroo™ ePump Sterile 100 mL Burette Set	30
775659	Kangaroo™ ePump ENPlus Spike Set	30
775100	Kangaroo™ ePump ENPlus Spike with 1000 mL Flush Bag	30
776150	Kangaroo™ ePump 100 mL Burette Recertification Set (not ENFit®)	5

## Kangaroo™ Joey Enteral Feeding Pump

Features feed and flush technology that automates flushing for tube patency and delivery of hydration with a compact, portable design.

- Fully programmable variable feed and flush pump
- Feed rate: 1 to 400 mL/hr in 1 mL increments
- Flush volume: 10 to 500 mL/hr in 1 mL increments
- Flush interval: 1 to 24 hr in 1 hr increments
- Feeding history up to 72 hours
- Battery life: 18 hours
- Weight without pole clamp: 1.7 lbs
- Weight with pole clamp: 2.3 lbs
- Multilingual: 15 languages
- Ambulatory-specific feed set design
- Tubing length 9 in (from the bottom of the bag to the distal end of the set)
- Attitude independent
- Intuitive user interface
- Auto-prime



383400

Item code	Description	Ship case
<b>Pump</b>		
383400	Kangaroo™ Joey Enteral Feeding Pump with Pole Clamp	1
<b>Pump Accessories</b>		
383493	Kangaroo™ Joey Pole Clamp	1
383491	Kangaroo™ Joey Power Cord	1
F010506	Kangaroo™ Joey Battery	1



762055


## Kangaroo™ Joey Enteral Feeding Pump Sets

- Anti-free flow design
- DEHP-free
- ENFit® connection

Item code	Description	Ship case
<b>Pump sets</b>		
762055	Kangaroo™ Joey 500 mL Pump Set	30
763656	Kangaroo™ Joey 1000 mL Pump Set	30
763662	Kangaroo™ Joey 1000 mL Pump Set with Flush Bag	30
765100	Kangaroo™ Joey Safety Screw Spike with 1000 mL Flush Bag	30
765559	Kangaroo™ Joey ENPlus Spike Set	30



763662



765559

## Kangaroo™ Joey Backpacks

Kangaroo™ Joey Backpacks are specially designed and tested to hold Kangaroo™ Joey Enteral Feeding Pumps.

Item code	Description	Ship case
<b>Backpacks</b>		
770026	Kangaroo™ Joey Blue Mini Backpack	1
770028	Kangaroo™ Joey Blue Super Mini Backpack	1
770034	Kangaroo™ Joey Pink Mini Backpack	1
770029	Kangaroo™ Joey Pink Super Mini Backpack	1
770027	Kangaroo™ Joey Green Mini Backpack	1
770032	Kangaroo™ Joey Green Super Mini Backpack	1
770025	Kangaroo™ Joey Black Mini Backpack	1
770031	Kangaroo™ Joey Black Super Mini Backpack	1

## Kangaroo™ Connect Enteral Feeding Pump and Accessories

Designed for the ambulatory patient, the Kangaroo™ Connect Enteral Feeding Pump is the first wirelessly connected enteral feeding pump on the market and is the smallest, lightest and quietest Kangaroo™ pump available.

- Most accurate Kangaroo™ enteral feeding pump at +/- 5% accuracy
- Expanded feed rate range
- Intuitive color screen with clear visual notifications
- Up to 24-hour battery life
- 72 hours of feeding history
- Travel accessories including car charger
- Water resistance rating of IPX6 to enable easier cleaning
- Simplified recertification
- Auto-prime



384400

Item code	Description	Ship case
384400	Kangaroo™ Connect Enteral Feeding Pump with Pole Clamp and Power Adapter	1
584400	Kangaroo™ Connect Enteral Feeding Pump with Pole Clamp and International Power Adapter	1
384491	Kangaroo™ Connect Power Cord with Adapter	1
584491	Kangaroo™ Connect Power Cord with International Adapter	1
384492	Kangaroo™ Connect Pole Clamp	1
384494	Kangaroo™ Connect Car Charger	1

## Kangaroo™ Connect Enteral Feeding Sets for Regular Pump — Not for Use with Wireless Connectivity

- Anti-free flow design
- DEHP-free
- ENFit® connection
- Cassette pump set design with auto prime feature

Item code	Description	Ship case
77000FD	Kangaroo™ Connect ENPlus Spike Set, Non-sterile	30
77100FD	Kangaroo™ Connect 1000 mL Bag Set, Non-sterile	30
77500FD	Kangaroo™ Connect 500 mL Bag Set, Non-sterile	30



77100FD

## Kangaroo™ Connect Enteral Feeding Sets — For Use with Wireless Connectivity

Item code	Description	Ship case
1KCP	Kangaroo™ Connect ENPlus Spike ENTelliSet with ENFit® Connector, Non-sterile	30
1000KCP	Kangaroo™ Connect ENTelliset with ENFit® Connector, 1000 mL, Non-sterile	30
500KCP	Kangaroo™ Connect ENTelliset with ENFit® Connector, 500 mL, Non-sterile	30



1000KCP

## Kangaroo™ Connect Hubs for Wireless Connectivity

Item code	Description	Ship case
384500CBW	Kangaroo™ Connect WiFi/Ethernet Communications Hub	1
384500CBC	Kangaroo™ Connect Cellular Communications Hub	1



384500CBW



384500CBC





770035M shown with 384400

## Kangaroo™ Connect Backpacks

Order code	Description	Ship case
770035L	Kangaroo™ Connect Backpack, Purple, Large	1
770037L	Kangaroo™ Connect Backpack, Black, Large	1
770035M	Kangaroo™ Connect Backpack, Purple, Medium	1
770037M	Kangaroo™ Connect Backpack, Black, Medium	1
770035S	Kangaroo™ Connect Backpack, Purple, Small	1
770037S	Kangaroo™ Connect Backpack, Black, Small	1





# Extension Sets, Gravity Sets and Irrigation Kits

## Kangaroo™ Extension Sets

- Non-sterile

Non-ENFit® code	ENFit® code	Description	Ship case
8884705008	8884705008E	Y-site Extension Set	30
8884700207	8884700207E	Pump Set 4 ft Tubing Extension	30

## Kangaroo™ Gravity Sets

- Non-sterile
- ENFit® connection

Item code	Description	Ship case
702505	1000 mL Large Bore Gravity Set	30
8884702500	1000 mL Gravity Bag	30



702505

## Kangaroo™ ENFit® Transition Connector

- Non-sterile

Item code	Description	Ship case
765000	Kangaroo™ ENFit® Transition Connector	200
F00106	ENFit® Funnel Connector	30

## Kangaroo™ I.V. Pouch and Irrigation Kits

- Convenient tear-to-open packaging, resealable flap
- Transparent, polyethylene bag hangs easily, clearly labeled
- Non-sterile

Non-ENFit® code	ENFit® code	Description	Ship case
84065	84065E	I.V. Pouch with 60 cc Piston Syringe	30
8884700108	8884700108E	Irrigation Kit Syringe with 500 cc Container and Basin	20
8884700116	8884700116E	Irrigation Kit Syringe with 500 cc Container	30



8884700108E



386100

## IRIS Technology

### Kangaroo™ Feeding Tube with IRIS Technology

**The first disposable feeding tube with an integrated real-time imaging system (IRIS) to visually aid in the placement of small-bore feeding tubes.**

- Features a camera integrated with a small-bore feeding tube to allow clinicians to identify anatomical markers during the placement of the tube

Item code	Description	Ship case
<b>Console</b>		
386100	Kangaroo™ IRIS Console with Accessories (Global Power Adapter, Battery, Interface Cable and Pole Clamp)	1
<b>Feeding Tubes</b>		
460836E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 36 in (91 cm)	5
460843E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 43 in (109 cm)	5
460855E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 55 in (140 cm)	5
461036E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 36 in (91 cm)	5
461043E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 43 in (109 cm)	5
461055E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 55 in (140 cm)	5
461236E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 36 in (91 cm)	5
461243E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 43 in (109 cm)	5
461255E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 55 in (140 cm)	5
<b>Accessories</b>		
388101	Kangaroo™ IRIS Power Adapter	1
388102	Kangaroo™ IRIS Replacement Battery	1
388103	Kangaroo™ IRIS Interface Cable	1
388104	Kangaroo™ IRIS Mounting Clamp	1
388105E	Kangaroo™ IRIS Insufflation Device	5
388106	Kangaroo™ IRIS Console Mounting Cart	1
388107	Kangaroo™ IRIS Battery Charging Station	1
388108	Kangaroo™ IRIS Console Carrying Case	1

## Nasogastric Feeding Tubes

### Kangaroo™ Feeding Tubes with Dobbhoff™ Tip

- Weighted feeding tubes with and without flow-through stylets
- Radiopaque polyurethane material
- Centimeter markings
- MRI conditional once stylet is removed
- DEHP-free
- ENFit® connection

Non-ENFit® code	ENFit® code	Includes stylet	Weight	French size	Length	Ship case
8884710842	8884710842E	No	7 g	8	43 in (109 cm)	10
8884710859	8884710859E	Yes	7 g	8	43 in (109 cm)	10
8884710826	8884710826E	Yes	7 g	8	55 in (140 cm)	10
8884710867	8884710867E	Yes	5 g	8	43 in (109 cm)	10
8884711006	8884711006E	Yes	7 g	10	43 in (109 cm)	10
8884711246	8884711246E	No	7 g	12	43 in (109 cm)	10
8884711253	8884711253E	Yes	7 g	12	43 in (109 cm)	10




8884711246E

### Kangaroo™ Feeding Tubes with Standard Tip

- Weighted feeding tubes with and without flow-through stylets
- Radiopaque polyurethane material
- Centimeter markings
- MRI conditional once stylet is removed
- DEHP-free
- ENFit® connection

Non-ENFit® code	ENFit® code	Includes stylet	Weight	French size	Length	Ship case
8884720817	8884720817E	No	3 g	8	36 in (91 cm)	10
8884720841	8884720841E	No	3 g	8	43 in (109 cm)	10
8884720825	8884720825E	Yes	3 g	8	36 in (91 cm)	10
8884720858	8884720858E	Yes	3 g	8	43 in (109 cm)	10
8884720855	8884720855E	Yes	3 g	8	55 in (140 cm)	10
8884721088	8884721088E	Yes	5 g	10	43 in (109 cm)	10
8884721096	8884721096E	No	5 g	10	43 in (109 cm)	10
8884721055	8884721055E	Yes	5 g	10	55 in (140 cm)	10
8884721211	8884721211E	No	5 g	12	36 in (91 cm)	10
8884721237	8884721237E	No	5 g	12	43 in (109 cm)	10
8884721252	8884721252E	Yes	5 g	12	43 in (109 cm)	10
8884721255	8884721255E	Yes	5 g	12	55 in (140 cm)	10



8884721211



## Kangaroo™ Non-Weighted Feeding Tubes


- Clear polyurethane material
- Proximal Y-site connector with catheter tip adapter (not available in ENFit®)
- Large exit holes and open distal end
- Radiopaque stripe
- DEHP-free

Item code	French size	Length	Ship case
8884711568	8	36 in (91 cm)	10
8884711550	10	36 in (91 cm)	10
8884711501	12	36 in (91 cm)	10
8884711519	14	36 in (91 cm)	10

## Specialty Tubes

### Kangaroo™ Naso-Jejunal Feeding Tubes

- Designed for precision placement in the small bowel using either endoscopy and/or fluoroscopy (includes guidewire and transfer cannula)
- Unique hollow-bolus weighted tip
- Radiopaque polyurethane material
- MRI conditional once stylet is removed
- ENFit® connection




8884752505

Non-ENFit® code	ENFit® code	Weight	French size	Length	Ship case
8884752505	8884752505E	7 g	12	60 in (152 cm)	2

### Kangaroo™ Naso-Jejunal Feeding and Gastric Decompression Tubes

- For “push” insertion over guidewire (kit included)
- Radiopaque polyurethane material
- Unique adjustable-length feeding tube design
- Spiral feeding outlet holes
- Jejunal tube can be retracted within gastric portion
- MRI conditional once stylet is removed
- ENFit® connection



8884752513

Non-ENFit® code	ENFit® code	Overall diameter	Jejunal tube diameter	Overall length	Ship case
8884752513	8884752513E	16 Fr	9 Fr	67 in (170 cm)	2


## Neonatal and Pediatric Feeding Tubes

**Kangaroo™ and Argyle™ neonatal and pediatric feeding tubes with ENFit® connection help address the 2014 JCAHO Tubing Misconnections Sentinel Event Alert. The ENFit® feeding tubes:**

- Are compatible with ENFit® syringes and extension sets
- Are DEHP- and latex-free
- Feature markings from 5-25 cm, which provide a visual indicator to assist in placement and migration checks
- Are transparent with Sentinel Line™ and Sentinel Eye™ radiopaque markings to assist in proper placement

### Kangaroo™ Polyurethane Feeding Tubes

Non-ENFit® code	ENFit® code	French size	Length	Ship case
461355	461355E	3.5 (1.2 mm)	12 in (31 cm)	10
461412	461412E	5 (1.7 mm)	20 in (51 cm)	10
461503	461503E	5 (1.7 mm)	36 in (91 cm)	10
461420	461420E	6.5 (2.2 mm)	20 in (51 cm)	10
461438	461438E	6.5 (2.2 mm)	36 in (91 cm)	10
461701	461701E	8 (2.7 mm)	20 in (51 cm)	10
461800	461800E	8 (2.7 mm)	42 in (107 cm)	10
461042	461042E	10 (3.3 mm)	42 in (107 cm)	10



461412E

## Argyle™ PVC Feeding Tubes

Non-ENFit® code	ENFit® code	French size	Length	Ship case
461206	461206E	3.5 (1.2 mm)	12 in (31 cm)	10
460802	460802E	5 (1.7 mm)	16 in (41 cm)	50
461008	461008E	5 (1.7 mm)	36 in (91 cm)	50
460711	460711E	6.5 (2.2 mm)	16 in (41 cm)	50
460604	460604E	8 (2.7 mm)	16 in (41 cm)	50
460406	460406E	8 (2.7 mm)	42 in (107 cm)	50
460208	460208E	10 (3.3 mm)	42 in (107 cm)	50

## Neonatal and Pediatric Feeding System

### Kangaroo™ Feeding Tube with ENFit® Connection Extension Sets

- PVC tubing
- DEHP-free
- With or without bifurcation (second port for irrigation and medicine delivery)
- Latex-free



60VENY

Non-ENFit® code	ENFit® code	Description	Length	Ship case
35ES	35ENS	Extension Set	35 in (89 cm)	50
60ES	60ENS	Extension Set	60 in (152 cm)	50
35EY	35ENY	Bifurcated Extension Set	35 in (89 cm)	50
60EY	60ENY	Bifurcated Extension Set	60 in (152 cm)	50

### Kangaroo™ Gravity Feeding System with ENFit® Connection

Item code	Description	Length	Ship case
4VENS	Gravity Extension Set	4 in (10.2 cm)	50
60VENS	Gravity Extension Set	60 in (152 cm)	50
60VENY	Gravity Extension Set, Bifurcated	60 in (152 cm)	50



## Kangaroo™ Milk Straw with ENFit® Connection

Item code	Description	Length	Ship case
461000E	Kangaroo™ Milk Straw	5 in (12.7 cm)	100



## Gastrostomy Tubes

### Kangaroo™ Gastrostomy Feeding Tubes with Y-Ports

- Medical-grade silicone tubing
- Rounded skin disc for external retention with feet
- Graduated shaft
- Gastric balloon for internal retention
- Rounded flush tip with open distal end
- DEHP-free



8884715148

Non-ENFit® code	ENFit® code	French size	Balloon size	Ship case
8884715122	8884715122E	12	5 cc	5
8884715148	8884715148E	14	5 cc	5
8884720163	8884720163E	16	20 cc	5
8884720189	8884720189E	18	20 cc	5
8884720205	8884720205E	20	20 cc	5
8884720221	8884720221E	22	20 cc	5
8884720247	8884720247E	24	20 cc	5
8884720265	8884720265E	26	20 cc	5
8884720285	8884720285E	28	20 cc	5

### Kangaroo™ Y-Port-G-Tube/P.E.G. Adapter

- For use with the Kangaroo™ P.E.G. kits

Non-ENFit® code	ENFit® code	Description	Ship case
8884751622	8884751622E	16 Fr Y-port/P.E.G. Adapter, Non-sterile	20
8884752042	8884752042E	20 Fr Y-port/P.E.G. Adapter, Non-sterile	20

## Safety P.E.G. Kits

### Kangaroo™ Safety P.E.G. Kits

- Featuring the Grant 90° bend, which straightens easily for P.E.G. placement and insertion of accessory feeding tubes
- Endoscopic removal
- Radiopaque polyurethane material
- Accessory jejunal feeding/gastric decompression tube available for 20 Fr P.E.G.
- Sterile
- DEHP-free

Non-ENFit® code	ENFit® code	French size	Ship case	Description
<b>Safety Percutaneous Endoscopic Gastrostomy “Pull” Kits</b>				
761610	761610E	16	2	<b>Kit Contains:</b> 21G x 1 1/2” Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6 cc Syringe, Pinch Clamp, #11 Safety Scalpel, Catheter Assembly w/90° Grant Bend, 18G Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare
762010	762010E	20	2	<b>Kit Contains:</b> 21G x 1 1/2” Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6 cc Syringe, Pinch Clamp, #11 Safety Scalpel, Catheter Assembly w/90° Grant Bend, 18G Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare



761610



## Kangaroo™ Safety Single Pass P.E.G. Kits

- Radiopaque internal retention bolster
- Non-endoscopic removal
- Easy-to-read placement marks
- Hydromer™ hydrophilic grafting on the “Pull” dilator tip
- Accessory jejunal feeding/gastric decompression tube available for 20 Fr P.E.G.
- Available in both “Pull” and “Push” over the guidewire methods
- Sterile
- DEHP-free

Non-ENFit® code	ENFit® code	French size	Ship case	Description
<b>Safety Percutaneous Endoscopic Gastrostomy “Pull” Kits</b>				
761600	761600E	16	2	<b>Kit Contains:</b> 21G x 1 1/2” Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6 cc Syringe, Pinch Clamp, #11 Safety Scalpel, Catheter Assembly w/90° Grant Bend, 18G Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Obturator Removal Kit, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare
762000	762000E	20	2	<b>Kit Contains:</b> 21G x 1 1/2” Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6 cc Syringe, Pinch Clamp, #11 Safety Scalpel, Catheter Assembly w/90° Grant Bend, 18G Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Obturator Removal Kit, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare





## Skin Level Tubes

### Kangaroo™ Skin Level Balloon Gastrostomy Kits

- 100% silicone material
- Raised feet design
- ClipStar™ enteral security clip
- DEHP-free

#### Contents of Kit:

Skin Level Gastrostomy Device, Step-Up Shim, Bolus Feeding Set, Continuous Feeding Set with ClipStar™ Enteral Security Clip, 60 cc Irrigation Syringe, 6 cc Syringe, Gauze Pads

Item code	Description	Ship case	Item code	Description	Ship case
712080	Balloon SLG 12 Fr x 0.8 cm	1	714350	Balloon SLG 14 Fr x 3.5 cm	1
712100	Balloon SLG 12 Fr x 1.0 cm	1	714400	Balloon SLG 14 Fr x 4.0 cm	1
712120	Balloon SLG 12 Fr x 1.2 cm	1	714450	Balloon SLG 14 Fr x 4.5 cm	1
712150	Balloon SLG 12 Fr x 1.5 cm	1	714500	Balloon SLG 14 Fr x 5.0 cm	1
712170	Balloon SLG 12 Fr x 1.7 cm	1	716080	Balloon SLG 16 Fr x 0.8 cm	1
712200	Balloon SLG 12 Fr x 2.0 cm	1	716100	Balloon SLG 16 Fr x 1.0 cm	1
712230	Balloon SLG 12 Fr x 2.3 cm	1	716120	Balloon SLG 16 Fr x 1.2 cm	1
712250	Balloon SLG 12 Fr x 2.5 cm	1	716150	Balloon SLG 16 Fr x 1.5 cm	1
712270	Balloon SLG 12 Fr x 2.7 cm	1	716170	Balloon SLG 16 Fr x 1.7 cm	1
712300	Balloon SLG 12 Fr x 3.0 cm	1	716200	Balloon SLG 16 Fr x 2.0 cm	1
712350	Balloon SLG 12 Fr x 3.5 cm	1	716230	Balloon SLG 16 Fr x 2.3 cm	1
712400	Balloon SLG 12 Fr x 4.0 cm	1	716250	Balloon SLG 16 Fr x 2.5 cm	1
712450	Balloon SLG 12 Fr x 4.5 cm	1	716270	Balloon SLG 16 Fr x 2.7 cm	1
712500	Balloon SLG 12 Fr x 5.0 cm	1	716300	Balloon SLG 16 Fr x 3.0 cm	1
714080	Balloon SLG 14 Fr x 0.8 cm	1	716350	Balloon SLG 16 Fr x 3.5 cm	1
714100	Balloon SLG 14 Fr x 1.0 cm	1	716400	Balloon SLG 16 Fr x 4.0 cm	1
714120	Balloon SLG 14 Fr x 1.2 cm	1	716450	Balloon SLG 16 Fr x 5.0 cm	1
714150	Balloon SLG 14 Fr x 1.5 cm	1	716500	Balloon SLG 16 Fr x 5.0 cm	1
714170	Balloon SLG 14 Fr x 1.7 cm	1	718080	Balloon SLG 18 Fr x 0.8 cm	1
714200	Balloon SLG 14 Fr x 2.0 cm	1	718100	Balloon SLG 18 Fr x 1.0 cm	1
714230	Balloon SLG 14 Fr x 2.3 cm	1	718120	Balloon SLG 18 Fr x 1.2 cm	1
714250	Balloon SLG 14 Fr x 2.5 cm	1	718150	Balloon SLG 18 Fr x 1.5 cm	1
714270	Balloon SLG 14 Fr x 2.7 cm	1	718170	Balloon SLG 18 Fr x 1.7 cm	1
714300	Balloon SLG 14 Fr x 3.0 cm	1	718200	Balloon SLG 18 Fr x 2.0 cm	1

## Kangaroo™ Skin Level Balloon Gastrostomy Kits (continued)

Item code	Description	Ship case	Item code	Description	Ship case
718230	Balloon SLG 18 Fr x 2.3 cm	1	720350	Balloon SLG 20 Fr x 3.5 cm	1
718250	Balloon SLG 18 Fr x 2.5 cm	1	720400	Balloon SLG 20 Fr x 4.0 cm	1
718270	Balloon SLG 18 Fr x 2.7 cm	1	720450	Balloon SLG 20 Fr x 4.5 cm	1
718300	Balloon SLG 18 Fr x 3.0 cm	1	720500	Balloon SLG 20 Fr x 5.0 cm	1
718350	Balloon SLG 18 Fr x 3.5 cm	1	724080	Balloon SLG 24 Fr x 0.8 cm	1
718400	Balloon SLG 18 Fr x 4.0 cm	1	724100	Balloon SLG 24 Fr x 1.0 cm	1
718450	Balloon SLG 18 Fr x 4.5 cm	1	724120	Balloon SLG 24 Fr x 1.2 cm	1
718500	Balloon SLG 18 Fr x 5.0 cm	1	724150	Balloon SLG 24 Fr x 1.5 cm	1
720080	Balloon SLG 20 Fr x 0.8 cm	1	724170	Balloon SLG 24 Fr x 1.7 cm	1
720100	Balloon SLG 20 Fr x 1.0 cm	1	724200	Balloon SLG 24 Fr x 2.0 cm	1
720120	Balloon SLG 20 Fr x 1.2 cm	1	724230	Balloon SLG 24 Fr x 2.3 cm	1
720150	Balloon SLG 20 Fr x 1.5 cm	1	724250	Balloon SLG 24 Fr x 2.5 cm	1
720170	Balloon SLG 20 Fr x 1.7 cm	1	724270	Balloon SLG 24 Fr x 2.7 cm	1
720200	Balloon SLG 20 Fr x 2.0 cm	1	724300	Balloon SLG 24 Fr x 3.0 cm	1
720230	Balloon SLG 20 Fr x 2.3 cm	1	724350	Balloon SLG 24 Fr x 3.5 cm	1
720250	Balloon SLG 20 Fr x 2.5 cm	1	724400	Balloon SLG 24 Fr x 4.0 cm	1
720270	Balloon SLG 20 Fr x 2.7 cm	1	724450	Balloon SLG 24 Fr x 4.5 cm	1
720300	Balloon SLG 20 Fr x 3.0 cm	1	724500	Balloon SLG 24 Fr x 5.0 cm	1

## Accessories

Item code	Description	Ship case
8884741805	Stoma Measuring Device, Sterile	5
8884741813	Bolus Feeding Set, Non-sterile	10
8884741821	Continuous Feeding Set, Non-sterile	10



Standard Syringes

## Monoject™ Oral Syringes — Standard Tip

Item code	Description	Ship case
401S	1 mL Non-sterile Purple Oral Syringe	500
403S	3 mL Non-sterile Purple Oral Syringe	500
406S	6 mL Non-sterile Purple Oral Syringe	500
410S	10 mL Non-sterile Purple Oral Syringe	500
401SG	1 mL Sterile Purple Oral Syringe	240
403SG	3 mL Sterile Purple Oral Syringe	480
406SG	6 mL Sterile Purple Oral Syringe	400
412SG	12 mL Sterile Purple Oral Syringe	192
435SG	35 mL Sterile Purple Oral Syringe	160
460SG	60 mL Sterile Purple Oral Syringe	120



ENFit® Syringes

## Monoject™ Enteral Syringes — ENFit® Connection

Item code	Description	Ship case
8881101015	1 mL Non-sterile Enteral Syringe	3200
8881103015	3 mL Non-sterile Enteral Syringe	2000
8881106015	6 mL Non-sterile Enteral Syringe	1200
8881112015	12 mL Non-sterile Enteral Syringe	1000
8881135015	35 mL Non-sterile Enteral Syringe	300
8881160015	60 mL Non-sterile Enteral Syringe	175
401SE	1 mL Sterile Enteral Syringe	240
403SE	3 mL Sterile Enteral Syringe	400
406SE	6 mL Sterile Enteral Syringe	400
412SE	12 mL Sterile Enteral Syringe	480
435SE	35 mL Sterile Enteral Syringe	160
460SE	60 mL Sterile Enteral Syringe	120



## Monoject™ Enteral Syringe Tip Caps with ENFit® Connection

Item code	Description	Ship case
400SE	ENFit® Syringe Tip Cap, Sterile	1000
400BE	ENFit® Syringe Tip Cap, Non-sterile	1000



400BE

## Monoject™ Medicine Straw with ENFit® Connection

Item code	Description	Length	Ship case
462000E	Monoject™ Medicine Straw 6 cm	2 in (6 cm)	100
464000E	Monoject™ Medicine Straw 10 cm	4 in (10 cm)	100

## Monoject™ UV/Light Inhibiting Transfer Bag

Item code	Description	Ship case
46312	UV/Light Inhibiting Transfer Bag	250



46312

## Monoject™ Universal Bottle Adapter with ENFit® Connection

Item code	Description	Ship case
4699E	Universal Bottle Adapter	50



4699E

## IVA™ Syringe Seals

Item code	Description	Ship case
CP3020R	Syringe Seal, Large, Red	1000
CP3020Y	Syringe Seal, Large, Yellow	1000
CP3030R	Syringe Seal, Small, Red	1000
CP3720A	Syringe Seal, Large, "Don't Use if Seal's Broken"	1000
CP3730A	Syringe Seal, Small, "Don't Use if Seal's Broken"	1000



CP3720A

## Salem Sump™ Tubes

### Salem Sump™ PVC Tubes

- Double lumen PVC tube—one for suction drainage, one for sump vent
- Centimeter markings assist the clinician with tube placement by providing reference points on the tube and allow for tracking of tube migration after tube securement
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Funnel end of vent lumen can be used to cap 5-in-1 adapter
- Sterile




8888268060

Item code	French size	Length	Ship case
8888268060	6 (2.0 mm)	24 in (61 cm)	10
8888268086	8 (2.7 mm)	24 in (61 cm)	10
8888264911	10 (3.3 mm)	36 in (91 cm)	50
8888264929	12 (4.0 mm)	48 in (122 cm)	50
8888264945	14 (4.7 mm)	48 in (122 cm)	50
8888264960	16 (5.3 mm)	48 in (122 cm)	50
8888264986	18 (6.0 mm)	48 in (122 cm)	50

### Salem Sump™ PVC Tubes with Anti-Reflux Valve

- Double lumen PVC tube—one for suction drainage, one for sump vent
- Gastric reflux prevention valve
- Centimeter markings assist the clinician with tube placement by providing reference points on the tube and allow for tracking of tube migration after tube securement
- Sterile
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Funnel end of vent lumen can be used to cap 5-in-1 adapter



8888266106

Item code	French size	Length	Ship case
8888266106	10 (3.3 mm)	36 in (91 cm)	50
8888266114	12 (4.0 mm)	48 in (122 cm)	50
8888266122	14 (4.7 mm)	48 in (122 cm)	50
8888266130	16 (5.3 mm)	48 in (122 cm)	50
8888266148	18 (6.0 mm)	48 in (122 cm)	50

## Salem Sump™ Anti-Reflux Valves

- Gastric reflux prevention valve
- Allows for maximum air flow to assist with stomach ulcer prevention
- Sterile

Item code	Description	Ship case
8888266197	Valve only	10



8888266197

## Salem Sump™ Silicone Tubes

- Double lumen silicone tube—one for suction drainage, one for venting action
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Removable 5-in-1 connector
- Capping system with sump vent lumen cap, suction lumen adapter and suction lumen adapter cap
- Sterile

Item code	French size	Length	Ship case
8888265108	10 (3.3 mm)	36 in (91 cm)	10
8888265116	12 (4.0 mm)	48 in (122 cm)	10
8888265124	14 (4.7 mm)	48 in (122 cm)	10
8888265132	16 (5.3 mm)	48 in (122 cm)	10
8888265140	18 (6.0 mm)	48 in (122 cm)	10



8888265108

## Salem Sump™ Silicone Tubes with Anti-Reflux Valve

- Double lumen silicone tube—one for suction drainage, one for venting action
- Gastric reflux prevention valve
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Removable 5-in-1 connector
- Capping system with sump vent lumen cap, suction lumen adapter and suction lumen adapter cap
- Sterile

Item code	French size	Length	Ship case
8888266213	10 (3.3 mm)	36 in (91 cm)	10
8888266221	12 (4.0 mm)	48 in (122 cm)	10
8888266247	14 (4.7 mm)	48 in (122 cm)	10
8888266262	16 (5.3 mm)	48 in (122 cm)	10
8888266288	18 (6.0 mm)	48 in (122 cm)	10



8888266213

## Stomach Tubes

### Argyle™ Stomach Tubes (Levin Type)

- Clear, thermosensitive PVC tube with four eyes
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Open distal end
- Integral funnel connector
- Markings 18 in (46 cm), 22 in (56 cm), 26 in (66 cm) and 30 in (76 cm) from distal end
- Non-sterile



8888265314

Item code	French size	Length	Ship case
8888265314	12 (4.0 mm)	48 in (122 cm)	50
8888265322	16 (5.3 mm)	48 in (122 cm)	50
8888265330	18 (6.0 mm)	48 in (122 cm)	50

### Stomach Tubes (Levin Type)

- Tube has four staggered eyes and center opening to reduce trauma
- Depth marks begin 20 in from distal tip with three additional marks at 4 in intervals
- Smooth, rounded tip reduces insertion trauma
- Non-sterile



155712

Item code	French size	Length	Ship case
155710	12 (4.0 mm)	48 in (122 cm)	50
155711	14 (4.7 mm)	48 in (122 cm)	50
155712	16 (5.3 mm)	48 in (122 cm)	50
155713	18 (6.0 mm)	48 in (122 cm)	50

## Argyle™ Edlich Gastric Lavage Trays

- Large 140 cc (5 oz) graduated Monoject™ syringes
- Depth marks at 18 in (46 cm), 22 in (56 cm), 26 in (66 cm) and 30 in (76 cm)
- Clamp
- Transparent PVC gastric tube with closed end and four eyes



8888750018

Item code	French size	Length	Ship case
8888750018	Large, 34 (11.3 mm)	36 in (91 cm)	20
8888750513	Medium/Pediatric, 24 (8.0 mm)	36 in (91 cm)	20






## Kangaroo™ Multifunctional Ports with Safe Enteral Connections

- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- All-in-one closed port system
- Allows for feeding, suction, irrigation and medication
- Internal ARV design
- DEHP-free
- Radiopaque tubing

Non-ENFit® code	ENFit® code	Description	Ship case
771111	771111E	Kangaroo™ Multifunctional Port with Single Lumen Adapter	10

## Salem Sump™ Tubes

Non-ENFit® code	ENFit® code	Description	Ship case
7771010	7771010E	10 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	10
7771210	7771210E	12 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	10
7771410	7771410E	14 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	10
7771610	7771610E	16 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	10
7771810	7771810E	18 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	10



7771810E

## Silicone Salem Sump™ Tubes with Anti-Reflux Valve Irrigation


Item code	Description	Ship case
7771011	10 Fr Silicone Salem Sump™ Tube with Anti-Reflux Valve Irrigation and Kangaroo™ Multifunctional Port	10
7771211	12 Fr Silicone Salem Sump™ Tube with Anti-Reflux Valve Irrigation and Kangaroo™ Multifunctional Port	10
7771411	14 Fr Silicone Salem Sump™ Tube with Anti-Reflux Valve Irrigation and Kangaroo™ Multifunctional Port	10
7771611	16 Fr Silicone Salem Sump™ Tube with Anti-Reflux Valve Irrigation and Kangaroo™ Multifunctional Port	10
7771811	18 Fr Silicone Salem Sump™ Tube with Anti-Reflux Valve Irrigation and Kangaroo™ Multifunctional Port	10

## CO<sub>2</sub> Detectors

### Kangaroo™ CO<sub>2</sub> Detector

- Proprietary colorimetric paper technology
- Detects presence of CO<sub>2</sub> within seconds
- Easily attaches to either a bellows or syringe
- Requires no change in insertion technique
- Convenient bellows device:
  - Clears obstructions from distal tip
  - Prevents CO<sub>2</sub> from entering tube prematurely
  - Facilitates airflow for quick results

Non-ENFit® code	ENFit® code	Description	Ship case
777702	777702E	Kangaroo™ CO <sub>2</sub> Detector	10
710859CN	710859CNE	8 Fr Kangaroo™ Dobbhoff™ Tip; 43 in with Stylet and CO <sub>2</sub> Detector	10
711006CN	711006CNE	10 Fr Kangaroo™ Dobbhoff™ Tip; 43 in with Stylet and CO <sub>2</sub> Detector	10
711253CN	711253CNE	12 Fr Kangaroo™ Dobbhoff™ Tip; 43 in with Stylet and CO <sub>2</sub> Detector	10
720858CN	720858CNE	8 Fr Kangaroo™ Standard Tip; 43 in with Stylet and CO <sub>2</sub> Detector	10
721088CN	721088CNE	10 Fr Kangaroo™ Standard Tip; 43 in with Stylet and CO <sub>2</sub> Detector	10
721252CN	721252CNE	12 Fr Kangaroo™ Standard Tip; 43 in with Stylet and CO <sub>2</sub> Detector	10
Item code		Description	Ship case
266130CN		16 Fr Salem Sump™ Tube with Anti-Reflux Valve and CO <sub>2</sub> Detector	50
266148CN		18 Fr Salem Sump™ Tube with Anti-Reflux Valve and CO <sub>2</sub> Detector	50
7771610CN		16 Fr Salem Sump™ Tube with Multifunctional Port and CO <sub>2</sub> Detector	10
7771810CN		18 Fr Salem Sump™ Tube with Multifunctional Port and CO <sub>2</sub> Detector	10



# Index

Code	Product	Page	Code	Product	Page
1KCP	Kangaroo™ Connect ENPlus Spike ENTelliSet with ENFit® Connector, Non-sterile	7	388102	Kangaroo™ IRIS Replacement Battery	10
1000KCP	Kangaroo™ Connect ENTelliset with ENFit® Connector, 1000 mL, Non-sterile	7	388103	Kangaroo™ IRIS Interface Cable	10
1041411	Kangaroo™ ePump Battery	4	388104	Kangaroo™ IRIS Mounting Clamp	10
155710	12 Fr, 48 in (122 cm) Stomach Tubes, Levin Type	26	388105E	Kangaroo™ IRIS Insufflation Device	10
155711	14 Fr, 48 in (122 cm) Stomach Tubes, Levin Type	26	388106	Kangaroo™ IRIS Console Mounting Cart	10
155712	16 Fr, 48 in (122 cm) Stomach Tubes, Levin Type	26	388107	Kangaroo™ IRIS Battery Charging Station	10
155713	18 Fr, 48 in (122 cm) Stomach Tubes, Levin Type	26	388108	Kangaroo™ IRIS Console Carrying Case	10
266130CN	16 Fr Salem Sump™ Tube with Anti-Reflux Valve and CO <sub>2</sub> Detector	29	4VENS	4" Gravity Extension Set	14
266148CN	18 Fr Salem Sump™ Tube with Anti-Reflux Valve and CO <sub>2</sub> Detector	29	400SE	ENFit® Syringe Tip Cap, Sterile	23
35ES/35ENS	35" Extension Set	14	400BE	ENFit® Syringe Tip Cap, Non-sterile	23
35EY/35ENY	35" Bifurcated Extension Set	14	401S	1 mL Non-sterile Purple Oral Syringe	22
382400	Kangaroo™ ePump with Pole Clamp	4	401SE	1 mL Sterile Enteral Syringe	22
382491	Kangaroo™ ePump Power Adapter	4	401SG	1 mL Sterile Purple Oral Syringe	22
382492	Kangaroo™ ePump Pole Clamp	4	403S	3 mL Non-sterile Purple Oral Syringe	22
383400	Kangaroo™ Joey Enteral Feeding Pump with Pole Clamp	5	403SE	3 mL Sterile Enteral Syringe	22
383491	Kangaroo™ Joey Power Cord	5	403SG	3 mL Sterile Purple Oral Syringe	22
383493	Kangaroo™ Joey Pole Clamp	5	406S	6 mL Non-sterile Purple Oral Syringe	22
384400	Kangaroo™ Connect Enteral Feeding Pump with Pole Clamp and Power Adapter	6	406SE	6 mL Sterile Enteral Syringe	22
384491	Kangaroo™ Connect Power Cord with Adapter	6	406SG	6 mL Non-sterile Purple Oral Syringe	22
384492	Kangaroo™ Connect Pole Clamp	6	410S	10 mL Non-sterile Purple Oral Syringe	22
384494	Kangaroo™ Connect Car Charger	6	412SG	12 mL Sterile Purple Oral Syringe	22
384500CBW	Kangaroo™ Connect WiFi/Ethernet Communications Hub	7	412SE	12 mL Sterile Enteral Syringe	22
384500CBC	Kangaroo™ Connect Cellular Communications Hub	7	435SG	35 mL Sterile Purple Oral Syringe	22
386100	Kangaroo™ IRIS Console with Accessories (Global Power Adapter, Battery, Interface Cable and Pole Clamp)	10	435SE	35 mL Sterile Enteral Syringe	22
388101	Kangaroo™ IRIS Power Adapter	10	460SE	60 mL Sterile Enteral Syringe	22
			460SG	60 mL Sterile Purple Oral Syringe	22
			460208/ 460208E	10 Fr, 42 in (107 cm) Argyle™ PVC Feeding Tube	14
			460406/ 460406E	8 Fr, 42 in (107 cm) Argyle™ PVC Feeding Tube	14
			460604/ 460604E	8 Fr, 16 in (41 cm) Argyle™ PVC Feeding Tube	14

Code	Product	Page	Code	Product	Page
460711/ 460711E	6.5 Fr, 16 in (41 cm) Argyle™ PVC Feeding Tube	14	461701/ 461701E	8 Fr, 20 in (51 cm) Kangaroo™ Polyurethane Feeding Tube	13
460802/ 460802E	5 Fr, 16 in (41 cm) Argyle™ PVC Feeding Tube	14	461800/ 461800E	8 Fr, 42 in (107 cm) Kangaroo™ Polyurethane Feeding Tube	13
460836E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 36 in (91 cm)	10	462000E	Monoject™ Medicine Straw 6 cm	23
460843E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 43 in (109 cm)	10	46312	UV/Light Inhibiting Transfer Bag	23
460855E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 55 in (140 cm)	10	464000E	Monoject™ Medicine Straw 10 cm	23
461000E	Kangaroo™ Milk Straw	15	4699E	Universal Bottle Adapter	23
461008/ 461008E	5 Fr, 36 in (91 cm) Argyle™ PVC Feeding Tube	14	500KCP	Kangaroo™ Connect ENtellisets with ENFit® Connector, 500 mL, Non-sterile	7
461036E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 36 in (91 cm)	10	584400	Kangaroo™ Connect Enteral Feeding Pump with Pole Clamp and International Power Adapter	6
461043E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 43 in (109 cm)	10	584491	Kangaroo™ Connect Power Cord with International Adapter	6
461055E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 55 in (140 cm)	10	60ES/60ENS	60" Extension Set	14
461236E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 36 in (91 cm)	10	60EY/60ENY	60" Bifurcated Extension Set	14
461243E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 43 in (109 cm)	10	60VENS	60" Gravity Extension Set	14
461255E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 55 in (140 cm)	10	60VENY	60" Gravity Extension Set, Bifurcated	14
461042/ 461042E	10 Fr, 42 in (107 cm) Kangaroo™ Polyurethane Feeding Tube	13	674668	Kangaroo™ ePump 500 mL Feed Set with 500 mL Flush Bag	14
461206/ 461206E	3.5 Fr, 12 in (31 cm) Argyle™ PVC Feeding Tube	14	702505	1000 mL Large Bore Gravity Set	9
461355/ 461355E	3.5 Fr, 12 in (31 cm) Kangaroo™ Polyurethane Feeding Tube	13	710859CN/ 710859CNE	Kangaroo™ Dobbhoff™ Tip 8 Fr; 43 in with Stylet and CO <sub>2</sub> Detector	29
461412/ 461412E	5 Fr, 20 in (51 cm) Kangaroo™ Polyurethane Feeding Tube	13	711006CN/ 711006CNE	Kangaroo™ Dobbhoff™ Tip 10 Fr; 43 in with Stylet and CO <sub>2</sub> Detector	29
461420/ 461420E	6.5 Fr, 20 in (51 cm) Kangaroo™ Polyurethane Feeding Tube	13	711253CN/ 711253CNE	Kangaroo™ Dobbhoff™ Tip 12 Fr; 43 in with Stylet and CO <sub>2</sub> Detector	29
461438/ 461438E	6.5 Fr, 36 in (91 cm) Kangaroo™ Polyurethane Feeding Tube	13	712080	Balloon SLG 12 Fr x 0.8 cm	20
461503/ 461503E	5 Fr, 36 in (91 cm) Kangaroo™ Polyurethane Feeding Tube	13	712100	Balloon SLG 12 Fr x 1.0 cm	20
			712120	Balloon SLG 12 Fr x 1.2 cm	20
			712150	Balloon SLG 12 Fr x 1.5 cm	20
			712170	Balloon SLG 12 Fr x 1.7 cm	20
			712200	Balloon SLG 12 Fr x 2.0 cm	20
			712230	Balloon SLG 12 Fr x 2.3 cm	20
			712250	Balloon SLG 12 Fr x 2.5 cm	20
			712270	Balloon SLG 12 Fr x 2.7 cm	20

<b>Code</b>	<b>Product</b>	<b>Page</b>	<b>Code</b>	<b>Product</b>	<b>Page</b>
712300	Balloon SLG 12 Fr x 3.0 cm	20	716500	Balloon SLG 16 Fr x 5.0 cm	20
712350	Balloon SLG 12 Fr x 3.5 cm	20	718080	Balloon SLG 18 Fr x 0.8 cm	20
712400	Balloon SLG 12 Fr x 4.0 cm	20	718100	Balloon SLG 18 Fr x 1.0 cm	20
712450	Balloon SLG 12 Fr x 4.5 cm	20	718120	Balloon SLG 18 Fr x 1.2 cm	20
712500	Balloon SLG 12 Fr x 5.0 cm	20	718150	Balloon SLG 18 Fr x 1.5 cm	20
714080	Balloon SLG 14 Fr x 0.8 cm	20	718170	Balloon SLG 18 Fr x 1.7 cm	20
714100	Balloon SLG 14 Fr x 1.0 cm	20	718200	Balloon SLG 18 Fr x 2.0 cm	20
714120	Balloon SLG 14 Fr x 1.2 cm	20	718230	Balloon SLG 18 Fr x 2.3 cm	21
714150	Balloon SLG 14 Fr x 1.5 cm	20	718250	Balloon SLG 18 Fr x 2.5 cm	21
714170	Balloon SLG 14 Fr x 1.7 cm	20	718270	Balloon SLG 18 Fr x 2.7 cm	21
714200	Balloon SLG 14 Fr x 2.0 cm	20	718300	Balloon SLG 18 Fr x 3.0 cm	21
714230	Balloon SLG 14 Fr x 2.3 cm	20	718350	Balloon SLG 18 Fr x 3.5 cm	21
714250	Balloon SLG 14 Fr x 2.5 cm	20	718400	Balloon SLG 18 Fr x 4.0 cm	21
714270	Balloon SLG 14 Fr x 2.7 cm	20	718450	Balloon SLG 18 Fr x 4.5 cm	21
714300	Balloon SLG 14 Fr x 3.0 cm	20	718500	Balloon SLG 18 Fr x 5.0 cm	21
714350	Balloon SLG 14 Fr x 3.5 cm	20	720080	Balloon SLG 20 Fr x 0.8 cm	21
714400	Balloon SLG 14 Fr x 4.0 cm	20	720100	Balloon SLG 20 Fr x 1.0 cm	21
714450	Balloon SLG 14 Fr x 4.5 cm	20	720120	Balloon SLG 20 Fr x 1.2 cm	21
714500	Balloon SLG 14 Fr x 5.0 cm	20	720150	Balloon SLG 20 Fr x 1.5 cm	21
716080	Balloon SLG 16 Fr x 0.8 cm	20	720170	Balloon SLG 20 Fr x 1.7 cm	21
716100	Balloon SLG 16 Fr x 1.0 cm	20	720200	Balloon SLG 20 Fr x 2.0 cm	21
716120	Balloon SLG 16 Fr x 1.2 cm	20	720230	Balloon SLG 20 Fr x 2.3 cm	21
716150	Balloon SLG 16 Fr x 1.5 cm	20	720250	Balloon SLG 20 Fr x 2.5 cm	21
716154	Kangaroo™ ePump Sterile 100 mL Burette Set	4	720270	Balloon SLG 20 Fr x 2.7 cm	21
716170	Balloon SLG 16 Fr x 1.7 cm	20	720300	Balloon SLG 20 Fr x 3.0 cm	21
716200	Balloon SLG 16 Fr x 2.0 cm	20	720350	Balloon SLG 20 Fr x 3.5 cm	21
716230	Balloon SLG 16 Fr x 2.3 cm	20	720400	Balloon SLG 20 Fr x 4.0 cm	21
716250	Balloon SLG 16 Fr x 2.5 cm	20	720450	Balloon SLG 20 Fr x 4.5 cm	21
716270	Balloon SLG 16 Fr x 2.7 cm	20	720500	Balloon SLG 20 Fr x 5.0 cm	21
716300	Balloon SLG 16 Fr x 3.0 cm	20	720858CN/ 720858CNE	Kangaroo™ Standard Tip 8 Fr; 43 in with Stylet and CO <sub>2</sub> Detector	29
716350	Balloon SLG 16 Fr x 3.5 cm	20	721088CN/ 721088CNE	Kangaroo™ Standard Tip 10 Fr; 43 in with Stylet and CO <sub>2</sub> Detector	29
716400	Balloon SLG 16 Fr x 4.0 cm	20	721252CN/ 721252CNE	Kangaroo™ Standard Tip 12 Fr; 43 in with Stylet and CO <sub>2</sub> Detector	29
716450	Balloon SLG 16 Fr x 4.5 cm	20			



Code	Product	Page	Code	Product	Page
724080	Balloon SLG 24 Fr x 0.8 cm	21	770032	Kangaroo™ Joey Green Super Mini Backpack	6
724100	Balloon SLG 24 Fr x 1.0 cm	21	770034	Kangaroo™ Joey Pink Mini Backpack	6
724120	Balloon SLG 24 Fr x 1.2 cm	21	770035S	Kangaroo™ Connect Backpack Purple, Small	8
724150	Balloon SLG 24 Fr x 1.5 cm	21	770035M	Kangaroo™ Connect Backpack Purple, Medium	8
724170	Balloon SLG 24 Fr x 1.7 cm	21	770035L	Kangaroo™ Connect Backpack Purple, Large	8
724200	Balloon SLG 24 Fr x 2.0 cm	21	770037S	Kangaroo™ Connect Backpack Black, Small	8
724230	Balloon SLG 24 Fr x 2.3 cm	21	770037M	Kangaroo™ Connect Backpack Black, Medium	8
724250	Balloon SLG 24 Fr x 2.5 cm	21	770037L	Kangaroo™ Connect Backpack Black, Large	8
724270	Balloon SLG 24 Fr x 2.7 cm	21	77100FD	Kangaroo™ Connect 1000 mL Bag Set, Non-sterile	7
724300	Balloon SLG 24 Fr x 3.0 cm	21	771111/ 771111E	Kangaroo™ Multifunctional Port with Single Lumen Adapter	28
724350	Balloon SLG 24 Fr x 3.5 cm	21	772055	Kangaroo™ ePump 500 mL Pump Set	4
724400	Balloon SLG 24 Fr x 4.0 cm	21	773656	Kangaroo™ ePump 1000 mL Pump Set	4
724450	Balloon SLG 24 Fr x 4.5 cm	21	773662	Kangaroo™ ePump 1000 mL Feed Set with 1000 mL Flush Bag	4
724500	Balloon SLG 24 Fr x 5.0 cm	21	77500FD	Kangaroo™ Connect 500 mL Bag Set, Non-sterile	7
761600/ 761600E	16 Fr, Single Pass Safety P.E.G. "Pull" Kit	18	775100	Kangaroo™ ePump ENPlus Spike with 1000 mL Flush Bag	4
762000/ 762000E	20 Fr, Single Pass Safety P.E.G. "Pull" Kit	18	775659	Kangaroo™ ePump ENPlus Spike Set	4
761610/ 761610E	16 Fr, Safety P.E.G. "Pull" Kit	17	776150	Kangaroo™ ePump 100 mL Burette Recertification Set	4
762010/ 762010E	20 Fr, Safety P.E.G. "Pull" Kit	17	7771010/ 7771010E	10 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	28
762055	Kangaroo™ Joey 500 mL Pump Set	5	7771011	10 Fr Silicone Salem Sump™ Tube with Anti-Reflux Valve Irrigation and Kangaroo™ Multifunctional Port	28
763656	Kangaroo™ Joey 1000 mL Pump Set	5	7771210/ 7771210E	12 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	28
763662	Kangaroo™ Joey 1000 mL Pump Set with Flush Bag	5	7771211	12 Fr Silicone Salem Sump™ Tube with Anti-Reflux Valve Irrigation and Kangaroo™ Multifunctional Port	28
765100	Kangaroo™ Joey Safety Screw Spike 1000 mL Flush Bag	5	7771410/ 7771410E	14 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	28
765000	Kangaroo™ ENFit® Transition Connector	9	7771411	14 Fr Silicone Salem Sump™ Tube with Anti-Reflux Valve Irrigation and Kangaroo™ Multifunctional Port	28
765559	Kangaroo™ Joey ENPlus Spike Set	5	7771610/ 7771610E	16 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	28
77000FD	Kangaroo™ Connect ENPlus Spike Set, Non-sterile	7	7771610CN	16 Fr Salem Sump™ Tube Multifunctional Port and CO <sub>2</sub> Detector	29
770025	Kangaroo™ Joey Black Mini Backpack	6	7771611	16 Fr Silicone Salem Sump™ Tube with Anti-Reflux Valve Irrigation and Kangaroo™ Multifunctional Port	28
770026	Kangaroo™ Joey Blue Mini Backpack	6			
770027	Kangaroo™ Joey Green Mini Backpack	6			
770028	Kangaroo™ Joey Blue Super Mini Backpack	6			
770029	Kangaroo™ Joey Pink Super Mini Backpack	6			
770031	Kangaroo™ Joey Black Super Mini Backpack	6			

Code	Product	Page	Code	Product	Page
7771810/ 7771810E	18 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	28	8884711519	14 Fr, 36 in (91 cm) Kangaroo™ Non-weighted Feeding Tube	12
7771810CN	18 Fr Salem Sump™ Tube with Multifunctional Port and CO <sub>2</sub> Detector	29	8884711550	10 Fr, 36 in (91 cm) Kangaroo™ Non-weighted Feeding Tube	12
7771811	18 Fr Silicone Salem Sump™ with Anti-Reflux Valve Irrigation and Kangaroo™ Multifunctional Port	28	8884711568	8 Fr, 36 in (91 cm) Kangaroo™ Non-weighted Feeding Tube	12
777702/ 777702E	Kangaroo™ CO <sub>2</sub> Detector	29	8884715122/ 8884715122E	12 Fr, 5 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
84065/84065E	I.V. Pouch with 60 cc Piston Syringe	9	8884715148/ 8884715148E	14 Fr, 5 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
8881101015	1 mL Non-sterile Enteral Syringe	22	8884720163/ 8884720163E	16 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
8881103015	3 mL Non-sterile Enteral Syringe	22	8884720189/ 8884720189E	18 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
8881106015	6 mL Non-sterile Enteral Syringe	22	8884720205/ 8884720205E	20 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
8881112015	12 mL Non-sterile Enteral Syringe	22	8884720221/ 8884720221E	22 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
8881135015	35 mL Non-sterile Enteral Syringe	22	8884720247/ 8884720247E	24 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
8881160015	60 mL Non-sterile Enteral Syringe	22	8884720265/ 8884720265E	26 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
8884700108/ 8884700108E	Irrigation Kit Syringe with 500 cc Container and Basin	9	8884720285/ 8884720285E	28 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
8884700116/ 8884700116E	Irrigation Kit Syringe with 500 cc Container	9	8884720817/ 8884720817E	8 Fr, 36 in (91 cm) Kangaroo Feeding Tube with Standard Tip, No Stylet	11
8884700207/ 8884700207E	Pump Set 4 ft Tubing Extension	9	8884720825/ 8884720825E	8 Fr, 36 in (91 cm) Kangaroo™ Feeding Tube with Standard Tip	11
8884702500	1000 mL Gravity Bag	9	8884720841/ 8884720841E	8 Fr, 43 in (109 cm) Kangaroo™ Feeding Tube with Standard Tip, No Stylet	11
8884705008/ 8884705008E	Y-site Extension Set	9	8884720855/ 8884720855E	8 Fr, 55 in (140 cm) Kangaroo™ Feeding Tube with Standard Tip	11
8884710826/ 8884710826E	8 Fr, 55 in (140 cm) Kangaroo™ Feeding Tube with Dobbhoff™ Tip	11	8884720858/ 8884720858E	8 Fr, 43 in (109 cm) Kangaroo™ Feeding Tube with Standard Tip	11
8884710842/ 8884710842E	8 Fr, 43 in (109 cm) Kangaroo™ Feeding Tube with Dobbhoff™ Tip, No Stylet	11	8884721055/ 8884721055E	10 Fr, 55 in (140 cm) Kangaroo™ Feeding Tube with Standard Tip	11
8884710859/ 8884710859E	8 Fr, 43 in (109 cm), Kangaroo™ Feeding Tube with Dobbhoff™ Tip	11	8884721088/ 8884721088E	10 Fr, 43 in (109 cm) Kangaroo™ Feeding Tube with Standard Tip	11
8884710867/ 8884710867E	8 Fr, 43 in (109 cm) Kangaroo™ Feeding Tube with Standard Tip	11	8884721096/ 8884721096E	10 Fr, 43 in (109 cm) Kangaroo™ Feeding Tube with Standard Tip, No Stylet	11
8884711006/ 8884711006E	10 Fr, 43 in (109 cm) Kangaroo™ Feeding Tube with Dobbhoff™ Tip	11	8884721211/ 8884721211E	12 Fr, 36 in (91 cm) Kangaroo™ Feeding Tube with Standard Tip, No Stylet	11
8884711246/ 8884711246E	12 Fr, 43 in (109 cm), Kangaroo™ Feeding Tube with Dobbhoff™ Tip, No Stylet	11			
8884711253/ 8884711253E	12 Fr, 43 in (109 cm) Kangaroo™ Feeding Tube with Dobbhoff™ Tip	11			
8884711501	12 Fr, 36 in (91 cm) Kangaroo™ Non-weighted Feeding Tube	12			

Code	Product	Page	Code	Product	Page
8884721237/ 8884721237E	12 Fr, 43 in (109 cm) Kangaroo™ Feeding Tube with Standard Tip, No Stylet	11	8888266122	14 Fr, 48 in (122 cm) Salem Sump™ PVC Tube with Anti-Reflux Valve	24
8884721252/ 8884721252E	12 Fr, 43 in (109 cm) Kangaroo™ Feeding Tube with Standard Tip	11	8888266130	16 Fr, 48 in (122 cm) Salem Sump™ PVC Tube with Anti-Reflux Valve	24
8884721255/ 8884721255E	12 Fr, 55 in (140 cm) Kangaroo™ Feeding Tube with Standard Tip	11	8888266148	18 Fr, 48 in (122 cm) Salem Sump™ PVC Tube with Anti-Reflux Valve	24
8884741805	Stoma Measuring Device, Sterile	21	8888266197	Salem Sump™ Anti-Reflux Valve	25
8884741813	Bolus Feeding Set, Non-sterile	21	8888266213	10 Fr, 36 in (91 cm) Salem Sump™ Silicone Tube with Anti-Reflux Valve	25
8884741821	Continuous Feeding Set, Non-sterile	21	8888266221	12 Fr, 48 in (122 cm) Salem Sump™ Silicone Tube with Anti-Reflux Valve	25
8884751622/ 8884751622E	16 Fr Y-port/P.E.G. Adapter, Non-sterile	16	8888266247	14 Fr, 48 in (122 cm) Salem Sump™ Silicone Tube with Anti-Reflux Valve	25
8884752042/ 8884752042E	20 Fr Y-port/P.E.G. Adapter, Non-sterile	16	8888266262	16 Fr, 48 in (122 cm) Salem Sump™ Silicone Tube with Anti-Reflux Valve	25
8884752505/ 8884752505E	12 Fr, 60 in (152 cm) Kangaroo™ Naso-Jejunal Feeding Tube	12	8888266288	18 Fr, 48 in (122 cm) Salem Sump™ Silicone Tube with Anti-Reflux Valve	25
8884752513/ 8884752513E	16 Fr, 67 in (170 cm) Kangaroo™ Naso-Jejunal Feeding and Gastric Decompression Tube	12	8888268060	6 Fr, 24 in (61 cm) Salem Sump™ PVC Tube	24
8888264911	10 Fr, 36 in (91 cm) Salem Sump™ PVC Tube	24	8888268086	8 Fr, 24 in (61 cm) Salem Sump™ PVC Tube	24
8888264929	12 Fr, 48 in (122 cm) Salem Sump™ PVC Tube	24	8888750018	Argyle™ Edlich Gastric Lavage Trays, Large	27
8888264945	14 Fr, 48 in (122 cm) Salem Sump™ PVC Tube	24	8888750513	Argyle™ Edlich Gastric Lavage Trays, Medium	27
8888264960	16 Fr, 48 in (122 cm) Salem Sump™ PVC Tube	24	CP3020R	Syringe Seal, Large, Red	23
8888264986	18 Fr, 48 in (122 cm) Salem Sump™ PVC Tube	24	CP3020Y	Syringe Seal, Large, Yellow	23
8888265108	10 Fr, 36 in (91 cm) Salem Sump™ Silicone Tube	25	CP3030R	Syringe Seal, Small, Red	23
8888265116	12 Fr, 48 in (122 cm) Salem Sump™ Silicone Tube	25	CP3720A	Syringe Seal, Large, "Don't Use if Seal's Broken"	23
8888265124	14 Fr, 48 in (122 cm) Salem Sump™ Silicone Tube	25	CP3730A	Syringe Seal, Small, "Don't Use if Seal's Broken"	23
8888265132	16 Fr, 48 in (122 cm) Salem Sump™ Silicone Tube	25	F010506	ENFit® Funnel Connector	9
8888265140	18 Fr, 48 in (122 cm) Salem Sump™ Silicone Tube	25			
8888265314	12 Fr, 48 in (122 cm) Argyle™ Stomach Tube, Levin Type	26			
8888265322	16 Fr, 48 in (122 cm) Argyle™ Stomach Tube, Levin Type	26			
8888265330	18 Fr, 48 in (122 cm) Argyle™ Stomach Tube, Levin Type	26			
8888266106	10 Fr, 36 in (91 cm) Salem Sump™ PVC Tube with Anti-Reflux Valve	24			
8888266114	12 Fr, 48 in (122 cm) Salem Sump™ PVC Tube with Anti-Reflux Valve	24			

**Cardinal Health Customer Service**  
**800.964.5227**

**[cardinalhealth.com/kangaroo](https://www.cardinalhealth.com/kangaroo)**

The International Organization for Standardization (ISO) and the Global Enteral Device Supplier Association (GEDSA) have created a new standard design for enteral connectors known as ENFit® connectors. ENFit® connectors are intended to improve patient safety and decrease the risk of placing IV medicines in an enteral feeding device.

© 2019 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, KANGAROO, Salem Sump, Monoject and Argyle are trademarks or registered trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Patent [www.cardinalhealth.com/patents](https://www.cardinalhealth.com/patents). Lit. No. 2MS19-1048032 (11/2019)



Cardinal Health  
777 West Street  
Mansfield, MA 02048