

**Cardinal Health™ FLEXAL™**  
**Powder-free nitrile examination gloves**



# Best-selling<sup>1</sup> brand for protection and dependability

Cardinal Health™ FLEXAL™ offers dependable protection against chemicals and chemotherapy drugs, while providing a desired fit and feel.

## Features and benefits:

- FLEXAL™ level protection
  - Standard level protection within our High Use Nitrile portfolio
- Tested for use with 9 chemotherapy drugs<sup>2</sup>
- Tested against 56 chemicals<sup>3</sup>
- USP 800 compliant\*



\* This Cardinal Health product helps healthcare facilities comply with USP <800> standard for the safe handling of hazardous drugs.

# Cardinal Health™ FLEXAL™ Powder-free nitrile examination gloves



Cat. No.	Size	Pkg.
<b>88TN01XS</b>	X-Small	200/bx, 10 bx/cs
<b>88TN02S</b>	Small	200/bx, 10 bx/cs
<b>88TN03M</b>	Medium	200/bx, 10 bx/cs
<b>88TN04L</b>	Large	200/bx, 10 bx/cs
<b>88TN05XL</b>	X-Large	200/bx, 10 bx/cs
<b>88TN06XXL</b>	XX-Large	200/bx, 10 bx/cs

- Powder-free
- Not made with natural rubber latex
- Ambidextrous
- Color: Cornflower blue
- Nominal fingertip thickness: 3.7 mil (.10 mm)

For more information or to order, contact your Cardinal Health sales representative. Call 877.CARDINAL or visit [cardinalhealth.com](http://cardinalhealth.com).

1 "Best-selling brand" based on FLEXAL™ sales in Acute and Non-Acute Healthcare channels, GHX Market Intelligence, CY 2018 Report

2 This glove has been tested for resistance to permeation of various chemotherapy drugs per ASTM D6978, "Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs." Caution: Gloves used for protection against chemotherapy drug and laboratory chemical exposure should be selected specifically for the type of chemicals being used.

Note: Review Safety Data Sheets (SDS) for the chemicals/drugs being used to determine required level of protection.

3 ASTM F739 - Standard Test Method for Permeation of Liquids and Gases through Protective Clothing Materials under Conditions of Continuous Contact