

New gowns for enhanced clinician safety

Reduce exposure risk when working with chemotherapy drugs

As the nature of chemotherapy treatment evolves, so do the exposure risks to healthcare personnel. Our newly developed, chemo-tested AAMI Level 3 gown helps protect against exposure to 16 types of chemotherapy drugs during preparation, handling, administration and more — wherever therapy is delivered.

Chemotherapy Drug Permeation Resistance Testing	FDA 510(k) Clearance Number: K193327 16 chemotherapy drugs tested per ASTM F739-20
AAMI 3 Barrier Protection	AATCC TM42, Test Method for Water Resistance: Impact Penetration Test requirement: spray impact ≤ 1.0 g AATCC TM127, Test Method for Water Resistance: Hydrostatic Pressure Test Test requirement: hydrostatic pressure ≥ 50 cm See AAMI Levels Specification Sheet for more information
Materials	Gown: polyethylene coated blue spunbond polypropylene Cuffs: polyester knit
Setting usage	This gown is intended to protect healthcare personnel from exposure to chemotherapy drugs during preparation, handling, and administration
Properties	Found to be non-cytotoxic, non-irritating and non-sensitizing after being tested for biological clearance per ISO 10993 Not made with natural rubber latex
Packaging	Non-sterile: Polyolefin bags contain the gowns. Sterile: Individual pouches made from Tyvek™ protective circle vent material All gowns are packed out in a shipping carton that is recyclable with corrugated products
Manufacturing	Made in the Dominican Republic of materials primarily from the United States



Cardinal Health™ ChemoPlus™ Chemotherapy Products

Ordering information

Sterile				
Product code	Size	Description	Qty.	Case dimensions (L x W x H)
CT5502TS	M	ChemoPlus™ Full Coverage Gown, closed back	1 gown/bag; 14 bags/cs	14.875 x 11.125 x 9.5 in.
CT5503TS	L	ChemoPlus™ Full Coverage Gown, closed back	1 gown/bag; 14 bags/cs	14.875 x 11.125 x 9.5 in.
CT5504TS	XL	ChemoPlus™ Full Coverage Gown, closed back	1 gown/bag; 14 bags/cs	15.625 x 11.625 x 12.375 in.
CT5505TS	XXL	ChemoPlus™ Full Coverage Gown, closed back	1 gown/bag; 14 bags/cs	15.625 x 11.625 x 12.375 in.
Non-sterile				
Product code	Size	Description	Qty.	Case dimensions (L x W x H)
CT5502T	M	ChemoPlus™ Full Coverage Gown, closed back	11 gowns/polybag; 2 polybag/cs	15.625 x 11.625 x 12.375 in.
CT5503T	L	ChemoPlus™ Full Coverage Gown, closed back	11 gowns/polybag; 2 polybag/cs	15.625 x 11.625 x 12.375 in.
CT5504T	XL	ChemoPlus™ Full Coverage Gown, closed back	11 gowns/polybag; 2 polybag/cs	15.625 x 11.625 x 12.375 in.
CT5505T	XXL	ChemoPlus™ Full Coverage Gown, closed back	11 gowns/polybag; 2 polybag/cs	15.625 x 11.625 x 12.375 in.
DP5001GT	L	ChemoPlus™ Full Coverage Gown, open back	14 gowns/polybag; 2 polybag/cs	15.625 x 11.625 x 12.375 in.
DP5002GT	XL	ChemoPlus™ Full Coverage Gown, open back	14 gowns/polybag; 2 polybag/cs	15.625 x 11.625 x 12.375 in.
DP5003GT	M	ChemoPlus™ Full Coverage Gown, open back	14 gowns/polybag; 2 polybag/cs	15.625 x 11.625 x 12.375 in.
DP5004GT	XXL	ChemoPlus™ Full Coverage Gown, open back	14 gowns/polybag; 2 polybag/cs	15.625 x 11.625 x 12.375 in.
Sleeves				
Product code	Size	Description	Qty.	Case dimensions (L x W x H)
CT5500TS	NA	ChemoPlus™ Full Coverage Sleeve, sterile	2 sleeves/bag; 60 bags/cs	14.875 x 11.062 x 14.5 in.
CT5500T	NA	ChemoPlus™ Full Coverage Sleeve, non-sterile	120 sleeves/cs	15.625 x 11.625 x 12.375 in.



Open back gown

	DP5003GT	DP5001GT	DP5002GT	DP5004GT
Size	Medium	Large	X-Large	XX-Large
Chest	24.75 in.	27.5 in.	29.125 in.	32.75 in.
Length	46.5 in.	46.5 in.	46.5 in.	46.5 in.



Closed back gown

	CT5502T/TS	CT5503T/TS	CT5504T/TS	CT5505T/TS
Size	Medium	Large	X-Large	XX-Large
Chest	25.6 in.	29.3 in.	31.5 in.	33.5 in.
Length	43.3 in.	47.3 in.	49.3 in.	51.1 in.

Invest in clinician safety. To develop your comprehensive infection control plan, contact your Cardinal Health sales representative.

