



## Cardinal Health™ Monoject™ Sharps Safety Products



Medication preparation



Medication delivery



Sharps disposal



# Engineered for safety





## Medication preparation



## Medication delivery



## Sharps and waste disposal



# CardinalHealth™ Monoject™

## Sharps Safety Products

# Engineered for safety

You're prepared, focused and committed to creating a culture of safety throughout your facility. For this, you need tools that are intuitively designed, easy to activate and require little to no training to use. When you standardize with Cardinal Health™ Monoject™ Products, you get a line of safety-engineered products designed to decrease the risk of needlestick injuries, minimize exposure to blood and bodily fluid and comply with OSHA regulations.

In fact, the CDC estimates that up to **75%** of sharps injuries can be prevented simply by choosing safer medical devices.<sup>1</sup>

With fewer types of safety mechanisms to master, you'll reduce training hours and enjoy the confidence that comes from simplifying your sharps products. From clinicians and patients to your downstream team, Monoject™ Sharps Safety Products are **designed to keep everyone safe.**


Count on **Monoject™ Sharps Safety resources** to support you with tools for product selection and utilization to keep you safe and ready your supply.

#### Reference:

1. OSHA. Bloodborne pathogens and needlestick prevention: Evaluating and controlling exposure. <https://www.osha.gov/SLTC/bloodbornepathogens/evaluation.html>. Accessed August 12, 2020.

Click product title  
to jump to page.

# Contents



## Medication preparation

[SmarTip™ Vial Access Cannulas](#)  
[Blunt fill needles](#)  
[Blunt fill cannulas](#)  
[Blunt cannulas with aluminum hubs](#)  
[Blunt cannulas and syringes](#)  
[Medication transfer needles](#)  
[Blunt Tip Plastic Cannulas](#)



## Medication delivery

**Safety needles and syringes:**  
[Magellan™ Hypodermic Safety Needles](#)  
[Magellan™ Hypodermic Safety Needles with Syringes](#)  
[Magellan™ Insulin and Tuberculin Safety Syringes](#)  
[Monoject™ Hypodermic Safety Needles](#)  
[Monoject™ Hypodermic Safety Needles with Syringes](#)  
[Monoject™ Insulin and Tuberculin Safety Syringes](#)

**Standard needles and syringes:**  
[Hypodermic standard needles](#)  
[Hypodermic standard needles with syringes](#)  
[Insulin standard syringes](#)  
[Tuberculin standard syringes](#)  
[Allergy trays](#)  
[Syringes](#)  
[Pharmacy trays](#)  
[Control syringes](#)  
[Non-sterile syringes](#)  
[Syringe tip caps](#)

## Dental needles and syringes:

Dental needles:  
• [Monoject™ Dental Needles – Plastic Hubs with Metal Inserts](#)  
• [Monoject™ Dental Needles – Metal Hubs](#)  
• [Monoject™ Endodontic Needles](#)  
[Monoject™ Endodontic Needles with Irrigation Syringes](#)  
[Monoject™ Curved Tip Syringes](#)

## IV sets and accessories:

[ChemoPlus™ IVA Security Seals](#)  
[Intermittent infusion plugs](#)  
[IV extension tubing](#)



## Sharps and waste disposal

**Biohazardous containers:**  
In-room sharps containers:  
• [SharpStar™ In-Room Sharps Containers](#)  
• [Safety in-room sharps containers](#)  
• [Mailbox in-room sharps containers](#)  
• [Always-open in-room sharps containers](#)  
[Multi-purpose sharps containers](#)  
[Phlebotomy sharps containers](#)  
[Large volume sharps containers](#)  
[Transportable sharps containers](#)

**Trace chemotherapy containers**  
**Hazardous waste pharmaceutical (RCRA) and universal containers**

## Non-hazardous pharmaceutical waste containers

## Non-infectious waste containers

## PGII Rated containers

## Container hardware:


[Brackets](#)  
[Wall enclosures](#)  
[Glove box dispensers](#)  
[Table top holders](#)  
[Large volume carts and foot pedal carts](#)

## Hardware and container cross reference guide




**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.


■ Biohazardous ■ Trace chemotherapy ■ Hazardous (RCRA) and universal ■ Non-hazardous ■ Non-infectious



## Medication preparation

Cardinal Health offers a line of safety-engineered products to prepare medication while avoiding needlesticks — like Monoject™ SmarTip™ Vial Access Cannula. Developed to replace a variety of medication preparation needles in your sharps line with little training, our SmarTip™ Cannula is the needle-free solution that's superior by design. Other products in our portfolio include the Monoject™ Blunt Fill Needle, a safe and reliable alternative to standard needles, as well as the Monoject™ Blunt Fill Needle with Filter and Monoject™ Medication Transfer Needle.





## Common medication preparation products

Standard

Safest



**Standard needles**  
18 G – 20 G



**Large bore safety hypodermic**  
18 G – 20 G



**Blunt fill**  
18 G




**Safety vial access cannula**

“Avoid the use of needles where safe and effective alternatives are available.”\*

- National Institute of Occupational Safety and Health (NIOSH)

**Reference:**

\*OSHA [1999]. How to prevent needlestick injuries: answers to some important questions. Washington, DC: U.S. Department of Labor, Occupational Safety and Health Administration. Accessed 3/4/2019 at: <https://www.cdc.gov/niosh/docs/2000-108/pdfs/2000-108.pdf>





## Monoject™ SmarTip™ Vial Access Cannulas

# No needle means no needlestick

Patient and clinician safety starts with smart design — like Monoject™ SmarTip™ needleless med prep technology from Cardinal Health. Our plastic SmarTip™ Vial Access Cannula is designed to avoid needlesticks.

Developed to replace a variety of medication preparation needles in your sharps line with little training, the SmarTip™ Cannula is the needle-free solution that's **superior by design**.



**Needleless**  
for clinician safety  
Plastic cannula engineered to help prevent accidental needlesticks

**Anti-coring**  
for patient safety  
Eliminates coring and particulate by piercing, not puncturing, the vial<sup>2</sup>

**Efficient**  
medication draws  
16-gauge equivalent flow rate and dual eyelets draw up even the most viscous medication from a vial easily and completely

- Features:**
- Soft pack packaging
  - Not manufactured with DEHP\* or natural rubber latex
  - Sterile, single-use and individually packaged

Compared with blunt and standard needles, the Monoject™ SmarTip™ Cannula:

135%

Requires 135% more force to penetrate than the BD™ Blunt Fill Needle<sup>1</sup>



Is less likely to penetrate the skin than the BD™ Blunt Fill Needle<sup>2</sup>

489%

Requires 489% more force to penetrate than a standard needle: 18 G x 1.5 in.<sup>3</sup>

2x

Requires twice as much force to penetrate than a blunt fill needle<sup>1</sup>


| Order code | Description                             | Box qty. | Case qty. |
|------------|-----------------------------------------|----------|-----------|
| 8881540111 | Needleless vial access cannula          | 100      | 1,000     |
| 8881540133 | 3 mL needleless vial access combination | 100      | 800       |
| 8881540166 | 6 mL needleless vial access combination | 100      | 400       |

**References:**

1. BD™ Blunt Fill needle reorder#305180. DOF. 2016. 2. Based on internal test results using simulated skin. Data on file. Cardinal Health, 2016. 3. Monoject 18 Gage needle reorder # 11811022. DOF. 2016.

\*Di(2-ethylhexyl) phthalate.






## Monoject™ Blunt Fill Needles

The Monoject™ Blunt Fill Needle provides clinicians with a safe and reliable alternative to standard needles during medication preparation.


Blunt fill needles are attached to a syringe for aspiration from multi-dose medication vials or injection into IV systems.


### Features:

- Standard luer hub – Fits all standard luer lock syringes
- 45° bevel – Helps reduce needlestick injuries
- 5 micron filter – Effectively filters contaminant particles for safe ampule medication preparation
- Red identification cap – For quick reference as blunt fill needle
- Available with and without a 5 micron filter
- Soft pack packaging
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



| Order code | Description                              | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------------------|--------------|---------------|----------|-----------|
| 11811022   | 18 G x 1½ in. (1.27 x 38 mm)             | 18 G         | 1½ in.        | 100      | 1,000     |
| 11811022F  | 18 G x 1½ in. (1.27 x 38 mm) with filter | 18 G         | 1½ in.        | 100      | 1,000     |





## Monoject™ Blunt Fill Cannulas

### Features:

- Polypropylene luer-lock hub
- Rigid pack packaging with heat stake, cap-to-sheath seal
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



| Order code | Description                  | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------|--------------|---------------|----------|-----------|
| 8881202018 | 18 G x 1½ in. (1.27 x 38 mm) | 18 G         | 1½ in.        | 100      | 1,000     |

## Monoject™ Blunt Cannulas with Aluminum Hubs


### Features:

- Stainless steel, non-lubricated cannula
- Rounded blunt end
- Aluminum luer-lock hub
- Rigid pack packaging with heat stake, cap-to-sheath seal
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



| Order code | Description                  | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------|--------------|---------------|----------|-----------|
| 8881202314 | 15 G x 1½ in. (1.83 x 38 mm) | 15 G         | 1½ in.        | 25       | 100       |
| 8881202322 | 16 G x 1½ in. (1.65 x 38 mm) | 16 G         | 1½ in.        | 25       | 100       |
| 8881202330 | 17 G x 1½ in. (1.47 x 38 mm) | 17 G         | 1½ in.        | 25       | 100       |
| 8881202348 | 18 G x 1 in. (1.27 x 25 mm)  | 18 G         | 1 in.         | 25       | 100       |
| 8881202355 | 19 G x 1½ in. (1.07 x 38 mm) | 19 G         | 1½ in.        | 25       | 100       |
| 8881202363 | 20 G x 1½ in. (0.90 x 38 mm) | 20 G         | 1½ in.        | 25       | 100       |
| 8881202371 | 21 G x 1 in. (0.81 x 25 mm)  | 21 G         | 1 in.         | 25       | 100       |
| 8881202389 | 22 G x 1½ in. (0.71 x 38 mm) | 22 G         | 1½ in.        | 25       | 100       |
| 8881202397 | 23 G x 1 in. (0.64 x 25 mm)  | 23 G         | 1 in.         | 25       | 100       |





## Monoject™ Blunt Cannulas and Syringes

### Features:

- Stainless steel cannula
- Rounded blunt end
- For use with LifeShield™\* System
- Polypropylene luer-lock hub
- Rigid pack packaging
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



| Order code | Description                          | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|--------------------------------------|----------------|--------------|---------------|----------|-----------|
| 8881202017 | 18 G x 1 in.<br>(1.27 x 25 mm)       | –              | 18 G         | 1 in.         | 100      | 1,000     |
| 8881202033 | 3 mL, 18 G x 1 in.<br>(1.27 x 25 mm) | 3 mL           | 18 G         | 1 in.         | 100      | 1,000     |

## Monoject™ Medication Transfer Needles


### Features:

- Double-ended stainless steel needle with anti-coring bevels
- Smooth, rounded aluminum hub
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



| Order code | Description                 | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|-----------------------------|--------------|---------------|----------|-----------|
| 8881204005 | 20 G x 1 in. (0.90 x 25 mm) | 20 G         | 1 in.         | 500      | 1,000     |

\*LifeShield is a trademark of ICU Medical




## Monoject™ Blunt Tip Plastic Cannulas

### Features:

- Universal cannula fits a wide variety of pre-pierced ports
- 15 gauge x 1/2 in. plastic cannula
- Rigid pack packaging
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged

| Order code | Description                                 | Box qty. | Case qty. |
|------------|---------------------------------------------|----------|-----------|
| 8881540101 | Blunt tip plastic cannula                   | 100      | 1,000     |
| 8881541034 | 3 mL syringe with blunt tip plastic cannula | 100      | 1,000     |
| 8881541067 | 6 mL syringe with blunt tip plastic cannula | 50       | 500       |





## Medication delivery

Cardinal Health offers a line of needles and syringes, both standard and safety, to support medication delivery. Our Monoject™ and Magellan™ Needles and Syringes are engineered to keep patients and clinicians safe at every step, from use to disposal. We also offer a variety of ChemoPlus™ IVA Security Seals, a solution to support your facility's USP <797> and <800> compliance practices designed with patient and clinician safety in mind.






## Tips for handling needles and syringes\*

When you handle needles and syringes, you must maintain aseptic technique at all times to prevent contamination and potential for infection for your patient.


- The shaft, bevel and hub of the needle and tip of syringe (Fig. A) should be kept sterile.
- Immediately after use, place used needles into a puncture-proof sharps container.
- Do not recap used needles to avoid risk of contaminated needlesticks.



**Figure A.**  
**Magellan™ Hypodermic Safety technology**

**Reference:**  
\*Burton, M., Smith, D. W., & Ludwig, L. J. (2019). *Fundamentals of nursing care: Concepts, connections & skills* (3rd ed.). Philadelphia: F.A. Davis Company.

## Types of needles and injections



| Injection type                         | Intramuscular                                                | Subcutaneous                                       | Intradermal                                     |
|----------------------------------------|--------------------------------------------------------------|----------------------------------------------------|-------------------------------------------------|
| <b>Common medical examples</b>         | Vaccines & toxoids, analgesics, antibiotics, hormonal agents | Insulin, opioids, heparin, epinephrine             | Allergens PPD (Mantoux)                         |
| <b>Usual maximum injection volume</b>  | <= 2 mL                                                      | <= 1 mL                                            | <= 0.1 mL                                       |
| <b>Injection site examples</b>         | Dorsogluteal, ventrogluteal, deltoid, vastus, lateralis      | Upper outer arm, anterior thigh, abdomen, buttocks | Inner (volar) forearm, upper back (subscapular) |
| <b>Typical needle lengths (adults)</b> | 1 in.                                                        | 5/8 in.                                            | 1/2 in. or 5/8 in.                              |
| <b>Typical needle gauges</b>           | 21 G - 25 G                                                  | 25 G - 27 G                                        | 27 G - 28 G                                     |



## Safety needles and syringes:

### Monoject™ Magellan™ Safety Needles and Syringes

# Clinically preferred safety needles designed to protect

Because more than 40% of needlestick injuries happen after use<sup>1</sup>, we designed Cardinal Health™ Monoject™ Magellan™ Safety Needles to keep patients and clinicians safe at every step. Magellan™ Safety Needles are engineered to comply with OSHA regulations and help you standardize solutions across every job with one safety technology. Intuitive, easy-to-use design offers protection with just the push of a finger, improving clinician comfort and confidence with less follow-up training.


#### Features:

- Easy to use, practical and intuitive one-handed activation via multiple methods (finger, thumb or surface) for clinical flexibility
- Even if bent inadvertently, safety shield provides 360-degree protection after activation to improve safety during use and downstream workers after disposal
- Audible and tactile click confirms safety activation and ready to be discarded in an approved sharps container
- Safety shield is an integral part of the needle hub, reducing the possibility of removal in accordance with OSHA standards
- Bevel orientation facilitates injection technique; if the safety shield faces up, the needle bevel also faces up
- Inspected multiple times throughout the manufacturing process to ensure optimum quality
- Designed to simplify safety technology standardization throughout your facility
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged

## One-handed activation



Thumb



Finger



Flat surface

Although combinations come preattached, ensure that the needle is seated.




#### Reference:

1. NIOSH. Alert: preventing needlestick injuries in health care settings. Cincinnati, OH: US Dept of Health and Human Services; 1999. <https://www.cdc.gov/niosh/docs/2000-108/pdfs/2000-108.pdf?id=10.26616/NIOSH PUB2000108>.




### Monoject™ Magellan™ Hypodermic Safety Needles

| Order code | Description                 | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|-----------------------------|--------------|---------------|----------|-----------|
| 8881850810 | 18 G x 1 in. (1.2 x 25 mm)  | 18 G         | 1 in.         | 50       | 500       |
| 8881850815 | 18 G x 1½ in. (1.2 x 38 mm) | 18 G         | 1½ in.        | 50       | 500       |
| 8881850910 | 19 G x 1 in. (1.1 x 25 mm)  | 19 G         | 1 in.         | 50       | 500       |
| 8881850915 | 19 G x 1½ in. (1.1 x 38 mm) | 19 G         | 1½ in.        | 50       | 500       |
| 8881850010 | 20 G x 1 in. (0.9 x 25 mm)  | 20 G         | 1 in.         | 50       | 500       |
| 8881850015 | 20 G x 1½ in. (0.9 x 38 mm) | 20 G         | 1½ in.        | 50       | 500       |
| 8881850110 | 21 G x 1 in. (0.8 x 25 mm)  | 21 G         | 1 in.         | 50       | 500       |
| 8881850115 | 21 G x 1½ in. (0.8 x 38 mm) | 21 G         | 1½ in.        | 50       | 500       |
| 8881850210 | 22 G x 1 in. (0.7 x 25 mm)  | 22 G         | 1 in.         | 50       | 500       |
| 8881850215 | 22 G x 1½ in. (0.7 x 38 mm) | 22 G         | 1½ in.        | 50       | 500       |
| 8881850358 | 23 G x ⅝ in. (0.6 x 16 mm)  | 23 G         | ⅝ in.         | 50       | 500       |
| 8881850310 | 23 G x 1 in. (0.6 x 25 mm)  | 23 G         | 1 in.         | 50       | 500       |
| 8881850558 | 25 G x ⅝ in. (0.5 x 16 mm)  | 25 G         | ⅝ in.         | 50       | 500       |
| 8881850510 | 25 G x 1 in. (0.5 x 25 mm)  | 25 G         | 1 in.         | 50       | 500       |
| 8881850515 | 25 G x 1½ in. (0.5 x 38 mm) | 25 G         | 1½ in.        | 50       | 500       |



### Monoject™ Magellan™ Hypodermic Safety Needles with Syringes


| Order code | Description                      | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|----------------------------------|----------------|--------------|---------------|----------|-----------|
| 8881811558 | 1 mL 25 G x ⅝ in. (0.5 x 16 mm)  | 1 mL           | 25 G         | ⅝ in.         | 50       | 200       |
| 8881811510 | 1 mL 25 G x 1 in. (0.5 x 25 mm)  | 1 mL           | 25 G         | 1 in.         | 50       | 200       |
| 8881811310 | 1 mL 23 G x 1 in. (0.6 x 25 mm)  | 1 mL           | 23 G         | 1 in.         | 50       | 200       |
| 8881833810 | 3 mL 18 G x 1 in. (1.2 x 25 mm)  | 3 mL           | 18 G         | 1 in.         | 50       | 400       |
| 8881833010 | 3 mL 20 G x 1 in. (0.9 x 25 mm)  | 3 mL           | 20 G         | 1 in.         | 50       | 400       |
| 8881833015 | 3 mL 20 G x 1½ in. (0.9 x 38 mm) | 3 mL           | 20 G         | 1½ in.        | 50       | 400       |
| 8881833110 | 3 mL 21 G x 1 in. (0.8 x 25 mm)  | 3 mL           | 21 G         | 1 in.         | 50       | 400       |
| 8881833115 | 3 mL 21 G x 1½ in. (0.8 x 38 mm) | 3 mL           | 21 G         | 1½ in.        | 50       | 400       |
| 8881833210 | 3 mL 22 G x 1 in. (0.7 x 25 mm)  | 3 mL           | 22 G         | 1 in.         | 50       | 400       |
| 8881833215 | 3 mL 22 G x 1½ in. (0.7 x 38 mm) | 3 mL           | 22 G         | 1½ in.        | 50       | 400       |
| 8881833310 | 3 mL 23 G x 1 in. (0.6 x 25 mm)  | 3 mL           | 23 G         | 1 in.         | 50       | 400       |
| 8881833510 | 3 mL 25 G x 1 in. (0.5 x 25 mm)  | 3 mL           | 25 G         | 1 in.         | 50       | 400       |
| 8881833558 | 3 mL 25 G x ⅝ in. (0.5 x 16 mm)  | 3 mL           | 25 G         | ⅝ in.         | 50       | 400       |



All hypodermic safety needles and syringes are ISO color coded.



## Safety needles and syringes: Routine diabetic care



## Rotating injection sites

An important part of routine diabetic care includes rotating injection sites. Injecting in the same tissue repeatedly has a negative impact on proper absorption and action of insulin. There are numerous sites on the body where insulin is routinely injected, with the idea of moving from site-to-site to maintain the viability of all appropriate injection sites.



## Safety needles and syringes:

### Monoject™ Magellan™ Insulin and Tuberculin Safety Syringes

Monoject™ Magellan™ Insulin and Tuberculin Safety Syringes are integral needle-based safety solutions. Intuitive, easy-to-use designs offer protection with just the push of a finger, improving clinician comfort and confidence with less follow-up training.

#### Features:

- Supporting diabetes care
- Easy to use, practical and intuitive one-handed activation via finger or thumb
- Safety shield provides 360-degree protection after activation to improve safety during use and downstream workers after disposal
- Permanently attached needle design provides low dead space and reduces medication loss
- Audible and tactile click confirms safety activation and ready to be discarded in an approved sharps container
- Safety shield is an integral part of the needle hub, reducing the possibility of removal in accordance with OSHA standards
- Bold, unobstructed graduations ensure visibility for precise measurements
- Designed to simplify safety technology standardization throughout your facility
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



### Monoject™ Magellan™ Insulin Safety Syringes with Permanent Needles

| Order code | Description                       | Syringe volume | Needle gauge | Needle length   | Box qty. | Case qty. |
|------------|-----------------------------------|----------------|--------------|-----------------|----------|-----------|
| 8881892910 | 1 mL 29 G x ½ in. (0.33 x 13 mm)  | 1 mL           | 29 G         | ½ in. (13 mm)   | 50       | 500       |
| 8881892950 | ½ mL 29 G x ½ in. (0.33 x 13 mm)  | ½ mL           | 29 G         | ½ in. (13 mm)   | 50       | 500       |
| 8881892930 | ¾ mL 29 G x ½ in. (0.33 x 13 mm)  | ¾ mL           | 29 G         | ½ in. (13 mm)   | 50       | 500       |
| 8881893010 | 1 mL 30 G x 5/16 in. (0.3 x 8 mm) | 1 mL           | 30 G         | 5/16 in. (8 mm) | 50       | 500       |
| 8881893050 | ½ mL 30 G x 5/16 in. (0.3 x 8 mm) | ½ mL           | 30 G         | 5/16 in. (8 mm) | 50       | 500       |
| 8881893030 | ¾ mL 30 G x 5/16 in. (0.3 x 8 mm) | ¾ mL           | 30 G         | 5/16 in. (8 mm) | 50       | 500       |

**Note:** 3/10 mL = 30 units, 1/2 mL = 50 units, 1 mL = 100 units

### Monoject™ Magellan™ Tuberculin Safety Syringes with Permanent Needles

| Order code | Description                      | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|----------------------------------|----------------|--------------|---------------|----------|-----------|
| 8881882712 | 1 mL 27 G x ½ in. (0.4 x 13 mm)  | 1 mL           | 27 G         | ½ in. (13 mm) | 50       | 500       |
| 8881882812 | 1 mL 28 G x ½ in. (0.36 x 13 mm) | 1 mL           | 28 G         | ½ in. (13 mm) | 50       | 500       |

### Monoject™ Magellan™ Tuberculin Safety Syringe Trays

| Order code  | Description                      | Syringe volume | Needle gauge | Needle length | Box qty.     | Case qty.     |
|-------------|----------------------------------|----------------|--------------|---------------|--------------|---------------|
| 8881882712T | 1 mL 27 G x ½ in. (0.4 x 13 mm)  | 1 mL           | 27 G         | ½ in. (13 mm) | 20 each/tray | 40 trays/case |
| 8881882812T | 1 mL 28 G x ½ in. (0.36 x 13 mm) | 1 mL           | 28 G         | ½ in. (13 mm) | 20 each/tray | 40 trays/case |





## Safety needles and syringes:


### **Monoject™ Hypodermic Safety Needles and Syringes**

# Safety needles designed to protect for strong performance and value

The Monoject™ Hypodermic Safety Needle is an intuitive, easy-to-use hinged safety device offering protection with just the push of a finger, improving clinician comfort and confidence with less follow-up training.

#### **Features:**

- Easy to use, practical and intuitive one-handed activation via multiple methods (finger, thumb or surface) for clinical flexibility
- Double-lock, hinged arm technology designed to encompass the needle after activation to improve safety during use and downstream workers after disposal
- Audible and tactile click confirms safety activation and ready to be discarded in an approved sharps container
- Safety shield sizes are consistent with needle length, ensuring full needle coverage
- Designed to simplify safety technology standardization throughout your facility
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



## **One-handed activation**



### **Flat surface activation**

To activate using a flat surface, gently press the safety shield against a flat surface until an audible and tactile click verifies the safety device has been activated.



### **Thumb activation**

To activate using your thumb, gently push down on the safety shield. An audible and tactile click will verify the safety device has been activated. Your thumb will remain above the needle at all times.



## Safety needles and syringes:

# Monoject™ Hypodermic Safety Needles and Syringes

## Monoject™ Hypodermic Safety Needles

| Order code | Description                 | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|-----------------------------|--------------|---------------|----------|-----------|
| 1181810    | 18 G x 1 in. (1.2 x 25 mm)  | 18 G         | 1 in.         | 100      | 800       |
| 1181815    | 18 G x 1½ in. (1.2 x 38 mm) | 18 G         | 1½ in.        | 100      | 800       |
| 1181910    | 19 G x 1 in. (1.1 x 25 mm)  | 19 G         | 1 in.         | 100      | 800       |
| 1181915    | 19 G x 1½ in. (1.1 x 38 mm) | 19 G         | 1½ in.        | 100      | 800       |
| 1182010    | 20 G x 1 in. (0.9 x 25 mm)  | 20 G         | 1 in.         | 100      | 800       |
| 1182015    | 20 G x 1½ in. (0.9 x 38 mm) | 20 G         | 1½ in.        | 100      | 800       |
| 1182110    | 21 G x 1 in. (0.8 x 25 mm)  | 21 G         | 1 in.         | 100      | 800       |
| 1182115    | 21 G x 1½ in. (0.8 x 38 mm) | 21 G         | 1½ in.        | 100      | 800       |
| 1182210    | 22 G x 1 in. (0.7 x 25 mm)  | 22 G         | 1 in.         | 100      | 800       |
| 1182215    | 22 G x 1½ in. (0.7 x 38 mm) | 22 G         | 1½ in.        | 100      | 800       |
| 1182310    | 23 G x 1 in. (0.6 x 25 mm)  | 23 G         | 1 in.         | 100      | 800       |
| 1182315    | 23 G x 1½ in. (0.6 x 38 mm) | 23 G         | 1½ in.        | 100      | 800       |
| 1182510    | 25 G x 1 in. (0.5 x 25 mm)  | 25 G         | 1 in.         | 100      | 800       |
| 1182515    | 25 G x 1½ in. (0.5 x 38 mm) | 25 G         | 1½ in.        | 100      | 800       |
| 1182558    | 25 G x ⅝ in. (0.5 x 16 mm)  | 25 G         | ⅝ in.         | 100      | 800       |
| 1182705    | 27 G x ½ in. (0.4 x 13 mm)  | 27 G         | ½ in.         | 100      | 800       |
| 1182734    | 27 G x ¾ in. (0.4 x 19 mm)  | 27 G         | ¾ in.         | 100      | 800       |
| 1183005    | 30 G x ½ in. (0.3 x 13 mm)  | 30 G         | ½ in.         | 100      | 800       |



## Monoject™ Hypodermic Safety Needles with Syringes

| Order code | Description                      | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|----------------------------------|----------------|--------------|---------------|----------|-----------|
| 11812558   | 1 mL 25 G x ⅝ in. (0.5 x 16 mm)  | 1 mL           | 25 G         | ⅝ in.         | 100      | 1,200     |
| 11812705   | 1 mL 27 G x ½ in. (0.4 x 13 mm)  | 1 mL           | 27 G         | ½ in.         | 100      | 1,200     |
| 11813005   | 1 mL 30 G x ½ in. (0.3 x 13 mm)  | 1 mL           | 30 G         | ½ in.         | 100      | 1,200     |
| 11832010   | 3 mL 20 G x 1 in. (0.9 x 25 mm)  | 3 mL           | 20 G         | 1 in.         | 100      | 1,200     |
| 11832015   | 3 mL 20 G x 1½ in. (0.9 x 38 mm) | 3 mL           | 20 G         | 1½ in.        | 100      | 1,200     |
| 11832110   | 3 mL 21 G x 1 in. (0.8 x 25 mm)  | 3 mL           | 21 G         | 1 in.         | 100      | 1,200     |
| 11832115   | 3 mL 21 G x 1½ in. (0.8 x 38 mm) | 3 mL           | 21 G         | 1½ in.        | 100      | 1,200     |
| 11832210   | 3 mL 22 G x 1 in. (0.7 x 25 mm)  | 3 mL           | 22 G         | 1 in.         | 100      | 1,200     |
| 11832215   | 3 mL 22 G x 1½ in. (0.7 x 38 mm) | 3 mL           | 22 G         | 1½ in.        | 100      | 1,200     |
| 11832310   | 3 mL 23 G x 1 in. (0.6 x 25 mm)  | 3 mL           | 23 G         | 1 in.         | 100      | 1,200     |
| 11832315   | 3 mL 23 G x 1½ in. (0.6 x 38 mm) | 3 mL           | 23 G         | 1½ in.        | 100      | 1,200     |
| 11832510   | 3 mL 25 G x 1 in. (0.5 x 25 mm)  | 3 mL           | 25 G         | 1 in.         | 100      | 1,200     |
| 11832558   | 3 mL 25 G x ⅝ in. (0.5 x 16 mm)  | 3 mL           | 25 G         | ⅝ in.         | 100      | 1,200     |



All hypodermic safety needles and syringes are ISO color coded.



## Safety needles and syringes:

### Monoject™ Insulin and Tuberculin Safety Syringes

Monoject™ Insulin and Tuberculin Safety Syringes utilize a barrel-based safety technology, designed to keep patients and clinicians safe with increased flexibility.

#### Features:

- Supporting diabetes care
- Safety shield provides 360-degree protection after activation to improve safety during use and downstream workers after disposal
- Transport position designed for safe transportation of medication without recapping the needle
- Permanently attached needle design provides low dead space and reduces medication loss
- Bold graduations ensure visibility for precise measurements
- Designed to simplify safety technology standardization throughout your facility
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



### Monoject™ Insulin Safety Syringes with Permanent Needles

| Order code | Description                           | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|---------------------------------------|----------------|--------------|---------------|----------|-----------|
| 8881511144 | 3/10 mL 29 G x 1/2 in. (0.33 x 13 mm) | 3/10 mL        | 29 G         | 1/2 in.       | 100      | 500       |
| 8881511136 | 1/2 mL 29 G x 1/2 in. (0.33 x 13 mm)  | 1/2 mL         | 29 G         | 1/2 in.       | 100      | 500       |
| 8881511110 | 1 mL 29 G x 1/2 in. (0.33 x 13 mm)    | 1 mL           | 29 G         | 1/2 in.       | 100      | 500       |
| 8881511344 | 3/10 mL 30 G x 5/16 in. (0.3 x 8 mm)  | 3/10 mL        | 30 G         | 5/16 in.      | 100      | 500       |
| 8881511336 | 1/2 mL 30 G x 5/16 in. (0.3 x 8 mm)   | 1/2 mL         | 30 G         | 5/16 in.      | 100      | 500       |
| 8881511310 | 1 mL 30 G x 5/16 in. (0.3 x 8 mm)     | 1 mL           | 30 G         | 5/16 in.      | 100      | 500       |

### Monoject™ Tuberculin Safety Syringes with Permanent Needles

| Order code | Description                        | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------------|----------------|--------------|---------------|----------|-----------|
| 8881511235 | 1 mL 25 G x 5/8 in. (0.5 x 16 mm)  | 1 mL           | 25 G         | 5/8 in.       | 100      | 500       |
| 8881511201 | 1 mL 28 G x 1/2 in. (0.36 x 13 mm) | 1 mL           | 28 G         | 1/2 in.       | 100      | 500       |



## Standard needles and syringes:

### Monject™ Soft Pack Hypodermic Standard Needles – Polypropylene Hubs

#### Features:


- Stainless steel lubricated cannula
- Color-coded needle hub
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



### Monject™ Soft Pack Hypodermic Standard Needles – Polypropylene Hubs

| Order code | Description                        | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------------|--------------|---------------|----------|-----------|
| 1188818100 | 18 G x 1 in. (1.27 x 25 mm) A      | 18 G         | 1 in.         | 100      | 1,000     |
| 1188818112 | 18 G x 1½ in. (1.27 x 38 mm) A     | 18 G         | 1½ in.        | 100      | 1,000     |
| 1188818113 | 18 G x 1½ in. (1.27 x 38 mm) B     | 18 G         | 1½ in.        | 100      | 1,000     |
| 1188819100 | 19 G TW* x 1 in. (1.07 x 25 mm) A  | 19 G         | 1 in.         | 100      | 1,000     |
| 1188819112 | 19 G TW* x 1½ in. (1.07 x 38 mm) A | 19 G         | 1½ in.        | 100      | 1,000     |
| 1188820100 | 20 G x 1 in. (0.9 x 25 mm) A       | 20 G         | 1 in.         | 100      | 1,000     |
| 1188820112 | 20 G x 1½ in. (0.9 x 38 mm) A      | 20 G         | 1½ in.        | 100      | 1,000     |
| 1188822100 | 22 G x 1 in. (0.711 x 25 mm) A     | 22 G         | 1 in.         | 100      | 1,000     |
| 1188822112 | 22 G x 1½ in. (0.71 x 38 mm) A     | 22 G         | 1½ in.        | 100      | 1,000     |
| 1188823100 | 23 G x 1 in. (0.63 x 25 mm) A      | 23 G         | 1 in.         | 100      | 1,000     |
| 1188825058 | 25 G x ⅝ in. (0.51 x 16 mm) A      | 25 G         | ⅝ in.         | 100      | 1,000     |
| 1188825100 | 25 G x 1 in. (0.51 x 25 mm) A      | 25 G         | 1 in.         | 100      | 1,000     |
| 1188825112 | 25 G x 1½ in. (0.51 x 38 mm) A     | 25 G         | 1½ in.        | 100      | 1,000     |
| 1188826112 | 26 G x 1½ in. (0.46 x 38 mm) A     | 26 G         | 1½ in.        | 100      | 1,000     |
| 1188827012 | 27 G x ½ in. (0.41 x 13 mm) A      | 27 G         | ½ in.         | 100      | 1,000     |
| 1188827112 | 27 G x 1½ in. (0.41 x 38 mm) A     | 27 G         | 1½ in.        | 100      | 1,000     |
| 1188830340 | 30 G x ¾ in. (0.31 x 19 mm) A      | 30 G         | ¾ in.         | 100      | 1,000     |

A Bevel - Regular bevel for intramuscular and subcutaneous injections  
 B Bevel - Short bevel for intravenous or intra-arterial applications and nerve blocks  
 IN Bevel - Intermediate bevel for intradermal injections/skin testing  
 \*Thin wall needle






## Standard needles and syringes:

### Monoject™ Rigid Pack Hypodermic Standard Needles – Polypropylene Hubs

#### Features:

- Stainless steel lubricated cannula
- Color-coded needle hub
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



### Monoject™ Rigid Pack Hypodermic Standard Needles – Polypropylene Hubs

| Order code | Description                        | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------------|--------------|---------------|----------|-----------|
| 8881250040 | 18 G x 1 in. (1.27 x 25 mm) A      | 18 G         | 1 in.         | 100      | 1,000     |
| 8881250016 | 18 G x 1½ in. (1.27 x 38 mm) A     | 18 G         | 1½ in.        | 100      | 1,000     |
| 8881250024 | 18 G x 1½ in. (1.27 x 38 mm) B     | 18 G         | 1½ in.        | 100      | 1,000     |
| 8881250081 | 19 G TW* x 1 in. (1.07 x 25 mm) A  | 19 G         | 1 in.         | 100      | 1,000     |
| 8881250099 | 19 G TW* x 1 in. (1.07 x 25 mm) B  | 19 G         | 1 in.         | 100      | 1,000     |
| 8881250065 | 19 G TW* x 1½ in. (1.07 x 38 mm) A | 19 G         | 1½ in.        | 100      | 1,000     |
| 8881250073 | 19 G TW* x 1½ in. (1.07 x 38 mm) B | 19 G         | 1½ in.        | 100      | 1,000     |
| 8881250057 | 18 G x 1 in. (1.27 x 25 mm) B      | 18 G         | 1 in.         | 100      | 1,000     |
| 8881250123 | 20 G x 1 in. (0.9 x 25 mm) A       | 20 G         | 1 in.         | 100      | 1,000     |
| 8881250107 | 20 G x 1½ in. (0.9 x 38 mm) A      | 20 G         | 1½ in.        | 100      | 1,000     |
| 8881250115 | 20 G x 1½ in. (0.9 x 38 mm) B      | 20 G         | 1½ in.        | 100      | 1,000     |
| 8881250339 | 20 G x ½ in. (0.9 x 13 mm) A       | 20 G         | ½ in.         | 100      | 1,000     |
| 8881250172 | 21 G x 1 in. (0.81 x 25 mm) A      | 21 G         | 1 in.         | 100      | 1,000     |
| 8881250149 | 21 G x 1½ in. (0.81 x 38 mm) A     | 21 G         | 1½ in.        | 100      | 1,000     |
| 8881250222 | 22 G x 1 in. (0.71 x 25 mm) A      | 22 G         | 1 in.         | 100      | 1,000     |
| 8881250230 | 22 G x 1 in. (0.71 x 25 mm) B      | 22 G         | 1 in.         | 100      | 1,000     |
| 8881250206 | 22 G x 1½ in. (0.71 x 38 mm) A     | 22 G         | 1½ in.        | 100      | 1,000     |
| 8881250214 | 22 G x 1½ in. (0.71 x 38 mm) B     | 22 G         | 1½ in.        | 100      | 1,000     |
| 8881250248 | 22 G x ¾ in. (0.71 x 19 mm) A      | 22 G         | ¾ in.         | 100      | 1,000     |
| 8881250289 | 23 G x ¾ in. (0.64 x 19 mm) A      | 23 G         | ¾ in.         | 100      | 1,000     |
| 8881250255 | 23 G x 1 in. (0.64 x 25 mm) A      | 23 G         | 1 in.         | 100      | 1,000     |
| 8881250313 | 25 G x ⅝ in. (0.51 x 16 mm) A      | 25 G         | ⅝ in.         | 100      | 1,000     |
| 8881250305 | 25 G x 1 in. (0.51 x 25 mm) A      | 25 G         | 1 in.         | 100      | 1,000     |
| 8881250545 | 25 G x 1½ in. (0.51 x 38 mm) A     | 25 G         | 1½ in.        | 100      | 1,000     |
| 8881250321 | 26 G x ½ in. (0.46 x 13 mm) IN     | 26 G         | ½ in.         | 100      | 1,000     |
| 8881250362 | 27 G x ½ in. (0.41 x 13 mm) IN     | 27 G         | ½ in.         | 100      | 1,000     |
| 8881250032 | 30 G x ¾ in. (0.31 x 19 mm) A      | 30 G         | ¾ in.         | 100      | 1,000     |

A Bevel - Regular bevel for intramuscular and subcutaneous injections

B Bevel - Short bevel for intravenous or intra-arterial applications and nerve blocks

IN Bevel - Intermediate bevel for intradermal injections/skin testing

\*Thin wall needle



## Standard needles and syringes:

### Monoject™ Rigid Pack Hypodermic Standard Needles – Aluminum Hubs

#### Features:


- Stainless steel lubricated cannula
- Color-coded sheath
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



### Monoject™ Rigid Pack Hypodermic Standard Needles – Aluminum Hubs

| Order code | Description                        | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------------|--------------|---------------|----------|-----------|
| 8881200805 | 14 G x 1 in. (2.11 x 25 mm) IN     | 14 G         | 1 in.         | 100      | 1,000     |
| 8881200011 | 14 G x 1½ in. (2.11 x 38 mm) IN    | 14 G         | 1½ in.        | 100      | 1,000     |
| 8881200573 | 14 G x 2 in. (2.11 x 51 mm) IN     | 14 G         | 2 in.         | 100      | 1,000     |
| 8881200029 | 15 G x 1½ in. (1.83 x 38 mm) B     | 15 G         | 1½ in.        | 100      | 1,000     |
| 8881200755 | 16 G x 5⁄8 in. (1.65 x 16 mm) B    | 16 G         | 5⁄8 in.       | 100      | 1,000     |
| 8881200060 | 16 G x ¾ in. (1.65 x 19 mm) IN     | 16 G         | ¾ in.         | 100      | 1,000     |
| 8881200052 | 16 G x 1 in. (1.65 x 25 mm) A      | 16 G         | 1 in.         | 100      | 1,000     |
| 8881200037 | 16 G x 1½ in. (1.65 x 38 mm) A     | 16 G         | 1½ in.        | 100      | 1,000     |
| 8881200045 | 16 G x 1½ in. (1.65 x 38 mm) B     | 16 G         | 1½ in.        | 100      | 1,000     |
| 8881200714 | 18 G x 1 in. (1.27 x 25 mm) A      | 18 G         | 1 in.         | 100      | 1,000     |
| 8881200078 | 18 G x 1½ in. (1.27 x 38 mm) A     | 18 G         | 1½ in.        | 100      | 1,000     |
| 8881200169 | 19 G TW* x 1 in. (1.07 x 25 mm) A  | 19 G         | 1 in.         | 100      | 1,000     |
| 8881200136 | 19 G TW* x 1½ in. (1.07 x 38 mm) A | 19 G         | 1½ in.        | 100      | 1,000     |
| 8881200219 | 20 G x 1 in. (0.9 x 25 mm) A       | 20 G         | 1 in.         | 100      | 1,000     |
| 8881200185 | 20 G x 1½ in. (0.9 x 38 mm) A      | 20 G         | 1½ in.        | 100      | 1,000     |
| 8881200318 | 21 G x 2 in. (0.81 x 51 mm) A      | 21 G         | 2 in.         | 100      | 1,000     |
| 8881200342 | 22 G x 1 in. (0.71 x 25 mm) A      | 22 G         | 1 in.         | 100      | 1,000     |
| 8881200326 | 22 G x 1½ in. (0.71 x 38 mm) A     | 22 G         | 1½ in.        | 100      | 1,000     |
| 8881200466 | 25 G x 5⁄8 in. (0.51 x 16 mm) A    | 25 G         | 5⁄8 in.       | 100      | 1,000     |
| 8881200433 | 25 G x 1¼ in. (0.51 x 32 mm) A     | 25 G         | 1¼ in.        | 100      | 1,000     |
| 8881200441 | 25 G x 2 in. (0.51 x 51 mm) A      | 25 G         | 2 in.         | 100      | 1,000     |
| 8881200516 | 27 G x ½ in. (0.41 x 13 mm) IN     | 27 G         | ½ in.         | 100      | 1,000     |
| 8881200508 | 27 G x 1¼ in. (0.41 x 32 mm) A     | 27 G         | 1¼ in.        | 100      | 1,000     |

A Bevel - Regular bevel for intramuscular and subcutaneous injections  
 B Bevel - Short bevel for intravenous or intra-arterial applications and nerve blocks  
 IN Bevel - Intermediate bevel for intradermal injections/skin testing  
 \*Thin wall needle





## Standard needles and syringes: Monoject™ Soft Pack Syringes with Hypodermic Standard Needles – Polypropylene Hubs

Monoject™ Hypodermic Needles are available as needle only or needle and syringe pre-attached syringe combinations.

### Features:

- Stainless steel lubricated cannula
- Color-coded needle hub
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



## Monoject™ Soft Pack Syringes with Hypodermic Standard Needles – Polypropylene Hubs

| Order code | Description                        | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------------|----------------|--------------|---------------|----------|-----------|
| 1180320100 | 3 mL 20 G x 1 in. (0.9 x 25 mm)    | 3 mL           | 20 G         | 1 in.         | 100      | 800       |
| 1180320112 | 3 mL 20 G x 1½ in. (0.9 x 38 mm)   | 3 mL           | 20 G         | 1½ in.        | 100      | 800       |
| 1180321100 | 3 mL 21 G x 1 in. (0.81 x 25 mm)   | 3 mL           | 21 G         | 1 in.         | 100      | 800       |
| 1180321112 | 3 mL 21 G x 1½ in. (0.81 x 38 mm)  | 3 mL           | 21 G         | 1½ in.        | 100      | 800       |
| 1180322112 | 3 mL 22 G x 1½ in. (0.71 x 38 mm)  | 3 mL           | 22 G         | 1½ in.        | 100      | 800       |
| 1180323100 | 3 mL 23 G x 1 in. (0.64 x 25 mm)   | 3 mL           | 23 G         | 1 in.         | 100      | 800       |
| 1180322100 | 3 mL 22 G x 1 in. (0.71 x 25 mm)   | 3 mL           | 22 G         | 1 in.         | 100      | 800       |
| 1180325058 | 3 mL 25 G x 5/8 in. (0.51 x 16 mm) | 3 mL           | 25 G         | 5/8 in.       | 100      | 800       |
| 1180325100 | 3 mL 25 G x 1 in. (0.51 x 25 mm)   | 3 mL           | 25 G         | 1 in.         | 100      | 800       |
| 1180325114 | 3 mL 25 G x 1¼ in. (0.51 x 32 mm)  | 3 mL           | 25 G         | 1¼ in.        | 100      | 800       |
| 1180327114 | 3 mL 27 G x 1¼ in. (0.41 x 32 mm)  | 3 mL           | 27 G         | 1¼ in.        | 100      | 800       |
| 1181620112 | 6 mL 20 G x 1½ in. (0.9 x 38 mm)   | 6 mL           | 20 G         | 1½ in.        | 100      | 400       |
| 1181621100 | 6 mL 21 G x 1 in. (0.81 x 25 mm)   | 6 mL           | 21 G         | 1 in.         | 100      | 400       |
| 1181621112 | 6 mL 21 G x 1½ in. (0.81 x 38 mm)  | 6 mL           | 21 G         | 1½ in.        | 100      | 400       |
| 1181622112 | 6 mL 22 G x 1½ in. (0.71 x 38 mm)  | 6 mL           | 22 G         | 1½ in.        | 100      | 400       |





## Standard needles and syringes:

### Monoject™ Rigid Pack Syringes with Hypodermic Standard Needles – Polypropylene Hubs

Monoject™ Hypodermic Needles are available as needle only or needle and syringe combinations.


#### Features:

- Stainless steel lubricated cannula
- Color-coded needle hub
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



### Monoject™ Rigid Pack Syringes with Hypodermic Standard Needles – Polypropylene Hubs

| Order code | Description                       | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|-----------------------------------|----------------|--------------|---------------|----------|-----------|
| 8881513025 | 3 mL 20 G x ¾ in. (0.9 x 19 mm)   | 3 mL           | 20 G         | ¾ in.         | 100      | 1,000     |
| 8881513033 | 3 mL 20 G x 1 in. (0.9 x 25 mm)   | 3 mL           | 20 G         | 1 in.         | 100      | 1,000     |
| 8881513058 | 3 mL 20 G x 1½ in. (0.9 x 38 mm)  | 3 mL           | 20 G         | 1½ in.        | 100      | 1,000     |
| 8881513132 | 3 mL 21 G x 1 in. (0.81 x 25 mm)  | 3 mL           | 21 G         | 1 in.         | 100      | 1,000     |
| 8881513157 | 3 mL 21 G x 1½ in. (0.81 x 38 mm) | 3 mL           | 21 G         | 1½ in.        | 100      | 1,000     |
| 8881513231 | 3 mL 22 G x 1 in. (0.71 x 25 mm)  | 3 mL           | 22 G         | 1 in.         | 100      | 1,000     |
| 8881513256 | 3 mL 22 G x 1½ in. (0.71 x 38 mm) | 3 mL           | 22 G         | 1½ in.        | 100      | 1,000     |
| 8881513330 | 3 mL 23 G x 1 in. (0.64 x 25 mm)  | 3 mL           | 23 G         | 1 in.         | 100      | 1,000     |
| 8881513512 | 3 mL 25 G x ⅝ in. (0.51 x 16 mm)  | 3 mL           | 25 G         | ⅝ in.         | 100      | 1,000     |
| 8881513538 | 3 mL 25 G x 1 in. (0.51 x 25 mm)  | 3 mL           | 25 G         | 1 in.         | 100      | 1,000     |
| 8881513744 | 3 mL 27 G x 1¼ in. (0.41 x 32 mm) | 3 mL           | 27 G         | 1¼ in.        | 100      | 1,000     |
| 8881516051 | 6 mL 20 G x 1½ in. (0.90 x 38 mm) | 6 mL           | 20 G         | 1½ in.        | 50       | 500       |
| 8881516135 | 6 mL 21 G x 1 in. (0.81 x 25 mm)  | 6 mL           | 21 G         | 1 in.         | 50       | 500       |







## Standard needles and syringes: Monoject™ Insulin Standard Syringes

### Features:

- Supporting diabetes care
- For U-100 insulin
- Stainless steel cannula
- Polypropylene hub
- ISO color coded
- Single unit graduations
- Available in soft and rigid pack packaging
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged

### Monoject™ Soft Pack Insulin Standard Syringes

| Order code | Description                        | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------------|----------------|--------------|---------------|----------|-----------|
| 8881600004 | ½ mL 28 G x ½ in. (0.36 x 13 mm)   | ½ mL           | 28 G         | ½ in.         | 100      | 300       |
| 1188528012 | ½ mL 28 G x ½ in. (0.36 x 13 mm)   | ½ mL           | 28 G         | ½ in.         | 100      | 500       |
| 8881601101 | 1 mL 28 G x ½ in. (0.36 x 13 mm)   | 1 mL           | 28 G         | ½ in.         | 100      | 300       |
| 1188128012 | 1 mL 28 G x ½ in. (0.36 x 13 mm)   | 1 mL           | 28 G         | ½ in.         | 100      | 500       |
| 8881600145 | ¾ mL 29 G x ½ in. (0.33 x 13 mm)   | ¾ mL           | 29 G         | ½ in.         | 100      | 300       |
| 8881600350 | ½ mL 29 G x ½ in. (0.33 x 13 mm)   | ½ mL           | 29 G         | ½ in.         | 100      | 300       |
| 8881601358 | 1 mL 29 G x ½ in. (0.33 x 13 mm)   | 1 mL           | 29 G         | ½ in.         | 100      | 300       |
| 8881600800 | ¾ mL 30 G x 5/16 in. (0.31 x 8 mm) | ¾ mL           | 30 G         | 5/16 in.      | 100      | 300       |
| 8881600700 | ½ mL 30 G x 5/16 in. (0.31 x 8 mm) | ½ mL           | 30 G         | 5/16 in.      | 100      | 300       |
| 8881601600 | 1 mL 30 G x 5/16 in. (0.31 x 8 mm) | 1 mL           | 30 G         | 5/16 in.      | 100      | 300       |
| 1188100555 | 1 mL insulin syringe luer-slip tip | 1 mL           | -            | -             | 100      | 500       |
| 1188100777 | 1 mL insulin syringe luer-lock tip | 1 mL           | -            | -             | 60       | 240       |



### Monoject™ Rigid Pack Insulin Standard Syringes

| Order code | Description                        | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------------|----------------|--------------|---------------|----------|-----------|
| 8881501210 | 1 mL 28 G x ½ in. (0.36 x 13 mm)   | 1 mL           | 28 G         | ½ in.         | 100      | 500       |
| 8881500014 | ½ mL 28 G x ½ in. (0.36 x 13 mm)   | ½ mL           | 28 G         | ½ in.         | 100      | 500       |
| 8881501384 | 1 mL insulin syringe luer-slip tip | 1 mL           | -            | -             | 100      | 500       |





## Standard needles and syringes: Monoject™ Tuberculin Standard Syringes

### Features:

- Stainless steel cannula
- 0.01 mL graduations
- Available in soft and rigid pack packaging
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged

### Monoject™ Soft Pack Tuberculin Standard Syringes

| Order code        | Description                          | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|-------------------|--------------------------------------|----------------|--------------|---------------|----------|-----------|
| <b>1180125158</b> | 1 mL 25 G x 5/8 in. (0.51 x 16 mm)   | 1 mL           | 25 G         | 5/8 in.       | 100      | 500       |
| <b>1180125058</b> | 1 mL 25 G x 5/8 in. (0.51 x 16 mm) * | 1 mL           | 25 G         | 5/8 in.       | 100      | 500       |
| <b>1180126038</b> | 1 mL 26 G x 3/8 in. (0.46 x 10 mm) * | 1 mL           | 26 G         | 3/8 in.       | 100      | 500       |
| <b>1180127012</b> | 1 mL 27 G x 1/2 in. (0.41 x 13 mm) * | 1 mL           | 27 G         | 1/2 in.       | 100      | 500       |
| <b>1180128012</b> | 1 mL 28 G x 1/2 in. (0.36 x 13 mm)   | 1 mL           | 28 G         | 1/2 in.       | 100      | 500       |

\*Detachable needle



### Monoject™ Rigid Pack Tuberculin Standard Syringes

| Order code        | Description                          | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|-------------------|--------------------------------------|----------------|--------------|---------------|----------|-----------|
| <b>8881501640</b> | 1 mL 25 G x 5/8 in. (0.51 x 16 mm)   | 1 mL           | 25 G         | 5/8 in.       | 100      | 500       |
| <b>8881500105</b> | 1/2 mL 28 G x 1/2 in. (0.36 x 13 mm) | 1/2 mL         | 28 G         | 1/2 in.       | 100      | 500       |
| <b>8881501608</b> | 1 mL 28 G x 1/2 in. (0.36 x 13 mm)   | 1 mL           | 28 G         | 1/2 in.       | 100      | 500       |
| <b>8881501160</b> | 1 mL 25 G x 5/8 in. (0.51 x 16 mm)*  | 1 mL           | 25 G         | 5/8 in.       | 100      | 500       |
| <b>8881501178</b> | 1 mL 26 G x 3/8 in. (0.46 x 10 mm)*  | 1 mL           | 26 G         | 3/8 in.       | 100      | 500       |
| <b>8881501368</b> | 1 mL 27 G x 1/2 in. (0.41 x 13 mm)*  | 1 mL           | 27 G         | 1/2 in.       | 100      | 500       |

\*Detachable needle





## Standard needles and syringes: Monoject™ Standard Allergy Trays

### Features:

- Bulk packaged needle and syringe combinations for pre-filling and batch filling procedures
- Graduations 0.01 mL
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use



### Monoject™ Standard Allergy Trays

| Order code        | Description                                        | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|-------------------|----------------------------------------------------|----------------|--------------|---------------|----------|-----------|
| <b>8881500501</b> | ½ mL permanent needle, 28 G x ½ in. (0.36 x 13 mm) | ½ mL           | 28 G         | ½ in.         | 25       | 1,000     |
| <b>8881501236</b> | 1 mL permanent needle, 28 G x ½ in. (0.36 x 13 mm) | 1 mL           | 28 G         | ½ in.         | 25       | 1,000     |
| <b>8881501962</b> | 1 mL 27 G x ½ in. (0.41 x 13 mm)*                  | 1 mL           | 27 G         | ½ in.         | 25       | 1,000     |

\*Detachable needle



## Standard needles and syringes:

### Monoject™ Syringes

Monoject™ Syringes are available sterile or non-sterile, in 1, 3, 6, 10, 12, 20, 35, 60 or 140 mL sizes, with several different tip designs for convenience.

#### Features:

- Polypropylene barrel and plunger rod
- Graduated barrel markings
- Available in soft and rigid pack packaging
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use, and individually packaged or in a bulk pack
- Available in luer-slip, Toomey, luer-lock, catheter and eccentric tip configurations

**Luer-slip:** For preparation and injection of medications and fluids. Friction fit attachment to needles and other female luer devices.

**Toomey:** For irrigation. Mates with cystoscopes or resectoscopes.

**Luer-lock:** For preparation and injection of medications and fluids. Secure twist-on attachment to needles and other female luer devices.

**Catheter:** For irrigation and bolus enteral feeding. Connects to compatible urinary and gastric systems.

**Eccentric:** For intravenous aspiration or injection. Friction fit attachment to needles and other luer devices. Facilitates low angle use.





## Standard needles and syringes:

### Monoject™ Syringes

#### Monoject™ Syringes – Soft Pack

| Order code  | Syringe volume | Tip type  | Box qty. | Case qty. |
|-------------|----------------|-----------|----------|-----------|
| 1180100555  | 1 mL           | Luer-slip | 100      | 500       |
| 1180100777  | 1 mL           | Luer-lock | 60       | 240       |
| 1180300555  | 3 mL           | Luer-slip | 100      | 800       |
| 1180300777  | 3 mL           | Luer-lock | 100      | 800       |
| 1180600555  | 6 mL           | Luer-slip | 100      | 400       |
| 1180600777  | 6 mL           | Luer-lock | 100      | 400       |
| 1181200555K | 10 mL          | Luer-slip | 100      | 1,000     |
| 1181200777K | 10 mL          | Luer-lock | 100      | 1,000     |
| 1181200777T | 12 mL          | Luer-lock | 100      | 1,000     |
| 1182000555  | 20 mL          | Luer-slip | 40       | 160       |
| 1182000777  | 20 mL          | Luer-lock | 40       | 160       |
| 1183500555  | 35 mL          | Luer-slip | 40       | 160       |
| 1183500777  | 35 mL          | Luer-lock | 40       | 160       |
| 1183500888  | 35 mL          | Catheter  | 40       | 160       |
| 1186000555  | 60 mL          | Luer-slip | 30       | 120       |
| 1186000555T | 60 mL          | Luer-slip | 30       | 360       |
| 1186000777  | 60 mL          | Luer-lock | 30       | 120       |
| 1186000777T | 60 mL          | Luer-lock | 30       | 360       |
| 1186000444  | 60 mL          | Catheter  | 30       | 120       |
| 1186060265  | 60 mL          | Toomey    | 30       | 120       |



#### Monoject™ Syringes – Rigid Pack

| Order code | Syringe volume | Tip type  | Box qty. | Case qty. |
|------------|----------------|-----------|----------|-----------|
| 8881501400 | 1 mL           | Luer-slip | 100      | 500       |
| 8881513918 | 3 mL           | Luer-slip | 100      | 1,000     |
| 8881513934 | 3 mL           | Luer-lock | 100      | 1,000     |
| 8881516911 | 6 mL           | Luer-slip | 50       | 500       |
| 8881516937 | 6 mL           | Luer-lock | 50       | 500       |
| 8881512852 | 12 mL          | Luer-slip | 80       | 480       |
| 8881512878 | 12 mL          | Luer-lock | 80       | 480       |
| 8881520673 | 20 mL          | Luer-slip | 50       | 300       |
| 8881520657 | 20 mL          | Luer-lock | 50       | 300       |
| 8881520665 | 20 mL          | Eccentric | 50       | 300       |
| 8881535796 | 35 mL          | Luer-slip | 30       | 180       |
| 8881535762 | 35 mL          | Luer-lock | 30       | 180       |
| 8881535788 | 35 mL          | Eccentric | 30       | 180       |
| 8881535770 | 35 mL          | Catheter  | 30       | 180       |
| 8881560224 | 60 mL          | Luer-slip | 20       | 100       |
| 8881560125 | 60 mL          | Luer-lock | 20       | 100       |
| 8881560182 | 60 mL          | Eccentric | 20       | 100       |
| 8881560141 | 60 mL          | Catheter  | 20       | 100       |
| 8881560265 | 60 mL          | Toomey    | 20       | 100       |






## Standard needles and syringes: Monoject™ Pharmacy Trays

### Features:

- Sterile, single-use and bulk packaged syringes for pre-filling and batch filling procedures
- Not manufactured with DEHP or natural rubber latex

| Order code | Syringe volume | Tip type  | Tray qty. | Case qty. |
|------------|----------------|-----------|-----------|-----------|
| 8881501459 | 1 mL           | Luer-slip | 25        | 1,000     |
| 8881513207 | 3 mL           | Luer-lock | 25        | 200       |
| 8881516200 | 6 mL           | Luer-lock | 25        | 200       |
| 8881512258 | 12 mL          | Luer-lock | 25        | 200       |
| 8881520251 | 20 mL          | Luer-lock | 25        | 200       |
| 8881535101 | 35 mL          | Luer-lock | 10        | 120       |
| 8881560232 | 60 mL          | Luer-lock | 10        | 120       |



## Standard needles and syringes: Monoject™ Control Syringes

### Features:

- Polypropylene barrel and plunger rod
- Graduated barrel markings
- Control rings permanently bonded to syringe
- Soft pack packaging
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged

| Order code | Syringe volume | Tip type  | Box qty. | Case qty. |
|------------|----------------|-----------|----------|-----------|
| 8881512977 | 12 mL          | Luer-lock | 40       | 160       |
| 8881520178 | 20 mL          | Luer-lock | 40       | 160       |





## Standard needles and syringes: Monoject™ Non-Sterile Syringes

### Features:

- Polypropylene barrel and plunger rod
- Not manufactured with DEHP or natural rubber latex
- Non-sterile, bulk packaged


| Order code | Syringe volume | Tip type  | Case qty. |
|------------|----------------|-----------|-----------|
| 8881103025 | 3 mL           | Luer-slip | 1,000     |
| 8881103066 | 3 mL           | Luer-lock | 2,000     |
| 8881106010 | 6 mL           | Luer-lock | 500       |
| 8881106028 | 6 mL           | Luer-slip | 500       |
| 8881112059 | 12 mL          | Luer-slip | 500       |
| 8881112083 | 12 mL          | Luer-lock | 800       |
| 8881120037 | 20 mL          | Luer-slip | 250       |
| 8881120029 | 20 mL          | Eccentric | 250       |
| 8881135084 | 35 mL          | Luer-slip | 100       |
| 8881160405 | 60 mL          | Luer-slip | 100       |
| 8881160629 | 60 mL          | Luer-lock | 155       |
| 8881114055 | 140 mL         | Catheter  | 20        |



## Standard needles and syringes: Monoject™ Syringe Tip Caps

### Features:

- Caps positioned upright in tray wells
- Designed to fit luer-lock or luer-slip tip syringes
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use, multiple-packaged



## Monoject™ Syringe Tip Caps

| Order code  | Description                 | Carton qty. | Case qty. |
|-------------|-----------------------------|-------------|-----------|
| 8881682101- | Syringe tip cap, tray of 10 | 1,000       | 12,000    |
| 8881682119  | Syringe tip cap, tray of 25 | 1,000       | 12,000    |
| 8881682085  | Syringe tip cap, bag of 100 | 1,000       | 12,000    |

## Monoject™ Safety Syringe Tip Caps

| Order code | Description            | Box qty. | Case qty. |
|------------|------------------------|----------|-----------|
| 8881682010 | Safety syringe tip cap | 50       | 250       |





# Monoject™ Dental Needles and Syringes

Designed to meet your dental needs



## Dental needles and syringes: Monoject™ Dental Needles – Plastic Hubs with Metal Inserts

**Features:**

- Stainless steel cannula
- Polypropylene hub with metal insert
- Anti-coring butt end point design
- Rigid polypropylene hard pack packaging
- Tamper-evident heat stake
- Sterile, single-use, DEHP free
- Not made with natural rubber latex



| Order code | Description                  | Color            | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------|------------------|--------------|---------------|----------|-----------|
| 8881400033 | 25 G long x 1¼ in. (32 mm)   | Red              | 25 G         | 1¼ in.        | 100      | 1,000     |
| 8881400041 | 25 G short x ¾ in. (19 mm)   | Red              | 25 G         | ¾ in.         | 100      | 1,000     |
| 8881400058 | 27 G long x 1¼ in. (32 mm)   | Yellow           | 27 G         | 1¼ in.        | 100      | 1,000     |
| 8881400066 | 27 G short x ¾ in. (19 mm)   | Yellow           | 27 G         | ¾ in.         | 100      | 1,000     |
| 8881400074 | 30 G short x ¾ in. (19 mm)   | Blue             | 30 G         | ¾ in.         | 100      | 1,000     |
| 8881400173 | 30 G X-short x ½ in. (13 mm) | Blue (white cap) | 30 G         | ½ in.         | 100      | 1,000     |





## Dental needles and syringes: Monoject™ Dental Needles – Metal Hubs

### Features:

- Stainless steel cannula
- Metal hub
- Anti-coring butt end point design
- Rigid polypropylene hard pack packaging
- Tamper-evident heat stake
- Sterile, single-use, DEHP free
- Not made with natural rubber latex



| Order code | Description                  | Color  | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|------------------------------|--------|--------------|---------------|----------|-----------|
| 8881401031 | 25 G long x 1¼ in. (32 mm)   | Red    | 25 G         | 1¼ in.        | 100      | 1,000     |
| 8881401049 | 25 G short x 1 in. (25 mm)   | Red    | 25 G         | 1 in.         | 100      | 1,000     |
| 8881401056 | 27 G long x 1⅜ in. (35 mm)   | Yellow | 27 G         | 1⅜ in.        | 100      | 1,000     |
| 8881401064 | 27 G short x 1 in. (25 mm)   | Yellow | 27 G         | 1 in.         | 100      | 1,000     |
| 8881401072 | 30 G short x 1 in. (25 mm)   | Blue   | 30 G         | 1 in.         | 100      | 1,000     |
| 8881401171 | 30 G X-short x ½ in. (13 mm) | Blue   | 30 G         | ½ in.         | 100      | 1,000     |

## Dental needles and syringes: Monoject™ Endodontic Needles

### Features:

- Stainless steel cannula
- Notched tip relieves pressure when the cannula tip is pressed against tissue
- Polypropylene hub
- Rigid polypropylene hard pack packaging
- Tamper-evident heat stake
- Sterile, single-use, DEHP free
- Not made with natural rubber latex



| Order code | Description           | Color  | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|-----------------------|--------|--------------|---------------|----------|-----------|
| 8881471232 | 23 G x 1¼ in. (32 mm) | Orange | 23 G         | 1¼ in.        | 25       | 100       |
| 8881471273 | 27 G x 1¼ in. (32 mm) | Yellow | 27 G         | 1¼ in.        | 25       | 100       |




## Dental needles and syringes:

### Monoject™ Endodontic Needles with Irrigation Syringes

**Features:**

- Stainless steel cannula
- Notched tip relieves pressure when the cannula tip is pressed against tissue
- Polypropylene hub
- Polypropylene syringe barrel and plunger rod
- Heat-etched, permanent graduations
- Available in soft and rigid pack packaging
- Tamper-evident heat stake
- Sterile, single-use, DEHP free
- Not made with natural rubber latex



### Monoject™ Endodontic Needles with Irrigation Syringes – Soft Pack

| Order code  | Description           | Color  | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|-------------|-----------------------|--------|----------------|--------------|---------------|----------|-----------|
| 1180323114E | 23 G x 1¼ in. (32 mm) | Orange | 3 mL           | 23 G         | 1¼ in.        | 100      | 800       |
| 1180327114E | 27 G x 1¼ in. (32 mm) | Yellow | 3 mL           | 27 G         | 1¼ in.        | 100      | 800       |

### Monoject™ Endodontic Needles with Irrigation Syringes – Rigid Pack

| Order code | Description           | Color  | Syringe volume | Needle gauge | Needle length | Box qty. | Case qty. |
|------------|-----------------------|--------|----------------|--------------|---------------|----------|-----------|
| 8881513843 | 23 G x 1¼ in. (32 mm) | Orange | 3 mL           | 23 G         | 1¼ in.        | 100      | 1,000     |
| 8881513850 | 27 G x 1¼ in. (32 mm) | Yellow | 3 mL           | 27 G         | 1¼ in.        | 100      | 1,000     |

## Dental needles and syringes:

### Monoject™ Curved Tip Syringes

**Features:**

- Polypropylene barrel and plunger rod
- Ungraduated, 12 mL capacity
- Non-sterile, bulk packaged
- Not manufactured with DEHP or natural rubber latex
- Curved tapered tip provides better access for irrigating gums



| Order code | Syringe volume | Tip type | Box qty. | Case qty. |
|------------|----------------|----------|----------|-----------|
| 8881412012 | 12 mL          | Curved   | 50       | 500       |



## IV sets and accessories:

### ChemoPlus™ IVA Security Seals

#### Promote safety and security of IV admixtures

These tamper-evident seals are a solution to support your facility's USP <797> and <800> compliance practices. Available in a variety of sizes and colors for use with vials, IV bags and medication containers.

#### Sterile

- Protects point of entry


#### Dual warning system

- Ensures reconstituted medication is not further diluted or enhanced prior to administering to a patient
- Tamper-evident seal self destructs upon attempted removal
- Needle punctures are highly visible

#### Visual indicator

- Aids in the identification of contents, routes of administration, beyond-use dates and other alerts

#### Made in the U.S.



#### For vial tops

| Order code | Description | Color  | Carton qty. |
|------------|-------------|--------|-------------|
| CP3005A-   | 9 mm        | Silver | 1,100       |
| CP3005R-   | 9 mm        | Red    | 1,100       |
| CP3003A-   | 13 mm       | Silver | 1,100       |
| CP3003G-   | 13 mm       | Green  | 1,100       |
| CP3003B-   | 13 mm       | Blue   | 1,100       |
| CP3003R-   | 13 mm       | Red    | 1,100       |
| CP3003Y    | 13 mm       | Yellow | 1,100       |
| CP3002Y-   | 20 mm       | Yellow | 1,000       |
| CP3002A-   | 20 mm       | Silver | 1,000       |
| CP3002G-   | 20 mm       | Green  | 1,000       |
| CP3002B-   | 20 mm       | Blue   | 1,000       |
| CP3002O-   | 20 mm       | Orange | 1,000       |
| CP3002R-   | 20 mm       | Red    | 1,000       |




#### For bottle tops

| Order code | Description | Color  | Carton qty. |
|------------|-------------|--------|-------------|
| CP3001A-   | 28 mm       | Silver | 1,000       |
| CP3001G-   | 28 mm       | Green  | 1,000       |
| CP3001B-   | 28 mm       | Blue   | 1,000       |
| CP3001Y-   | 28 mm       | Yellow | 1,000       |
| CP3001O-   | 28 mm       | Orange | 1,000       |
| CP3001R-   | 28 mm       | Red    | 1,000       |
| CP3001T-   | 28 mm       | Tan    | 1,000       |
| CP3009A-   | 36 mm       | Silver | 1,000       |
| CP3009G-   | 36 mm       | Green  | 1,000       |
| CP3009B-   | 36 mm       | Blue   | 1,000       |
| CP3009R-   | 36 mm       | Red    | 1,000       |



## IV sets and accessories: ChemoPlus™ IVA Security Seals



For Pfizer Large Bags and ADD-Vantage™ System

| Order code | Color  | Carton qty. |
|------------|--------|-------------|
| CP3004A    | Silver | 1,000       |
| CP3004G-   | Green  | 1,000       |
| CP3004B-   | Blue   | 1,000       |
| CP3004O    | Orange | 1,000       |
| CP3004R-   | Red    | 1,000       |



For B. Braun Excel™ Bags and Small Bags

| Order code | Color  | Carton qty. |
|------------|--------|-------------|
| CP3013A-   | Silver | 1,000       |
| CP3013G-   | Green  | 1,000       |
| CP3013B-   | Blue   | 1,000       |
| CP3013R-   | Red    | 1,000       |



For Baxter Viaflex™ Bags and Mini-Bag Plus Containers

| Order code | Color                           | Carton qty. |
|------------|---------------------------------|-------------|
| CP3011A-   | Silver                          | 1,000       |
| CP3011G-   | Green                           | 1,000       |
| CP3011B-   | Blue                            | 1,000       |
| CP3011Y-   | Yellow                          | 1,000       |
| CP3011R-   | Red                             | 1,000       |
| CP3111Y-   | Yellow/black ink "Chemotherapy" | 1,000       |



For smaller syringes, 3.625 in. length

| Order code | Color                                       | Carton qty. |
|------------|---------------------------------------------|-------------|
| CP3030R-   | Red                                         | 1,000       |
| CP3730A-   | Silver/red ink "Don't use if seal's broken" | 1,000       |



For larger syringes, 4.25 in. length

| Order code | Color                                       | Carton qty. |
|------------|---------------------------------------------|-------------|
| CP3020B-   | Blue                                        | 1,000       |
| CP3020Y-   | Yellow                                      | 1,000       |
| CP3020R-   | Red                                         | 1,000       |
| CP3720A-   | Silver/red ink "Don't use if seal's broken" | 1,000       |



## IV sets and accessories:

### Monoject™ Intermittent Infusion Plug

**Features:**

- Male luer-lock tip design also accommodates luer-slip tip fitting
- Special synthetic diaphragm retained by secure plastic shrink sleeve
- Protective sleeve over connector end
- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged



| Order code | Description                | Box qty. | Case qty. |
|------------|----------------------------|----------|-----------|
| 8888115006 | Heparin lock for IV device | 50       | 250       |

### Monoject™ IV Extension Tubing

- Not manufactured with DEHP or natural rubber latex
- Sterile, single-use and individually packaged

| Order code | Description   | Case qty. |
|------------|---------------|-----------|
| 155922     | 20 in. length | 50        |
| 155923     | 30 in. length | 50        |



**Reference:**

1. Op. Cit. Poster presented at the 2010 annual meeting of AVA.



# Sharps and waste disposal

Count on Cardinal Health to protect through the final steps of disposal with superior safety and value. Developed with safety-engineered controls to meet OSHA standards for ease of use and needlestick prevention, our complete line of disposal containers is designed to be safe for clinicians, patients and the environment.

Choose from a variety of single-use containers designed for the disposal of waste, including biohazardous, trace chemotherapy, hazardous pharmaceutical and universal, non-hazardous pharmaceutical and non-infectious. We also offer a complete line of hardware for use with our sharps containers.

## Select your waste type:



Biohazardous



Trace  
Chemotherapy



Hazardous  
(RCRA) and  
Universal



Non-Hazardous  
Pharmaceutical  
Waste



Non-Infectious

### Disposal and processing recommendations:

- All non-gasket sharps disposal containers manufactured by Cardinal Health are autoclavable<sup>1</sup>.
- All trace chemotherapy, hazardous and non-hazardous pharmaceutical waste must be incinerated.
- Material Safety Data Sheets (MSDS) are not needed for Cardinal Health™ Sharps Disposal Containers. They contain no heavy metals or DEHP and they are latex free.



**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

■ Biohazardous ■ Trace chemotherapy ■ Hazardous (RCRA) and universal ■ Non-hazardous ■ Non-infectious

### References:

1. It is recommended that customers perform their own autoclave testing to determine safety and efficacy. Our testing consists of autoclaving at 275 degrees Fahrenheit for one hour to assure the container does not lose integrity that the contents will spill. Puncture resistance has not been tested after autoclaving.



## Monoject™ Renewables Sharps Containers


### Reduce. Reuse. Recycle.

Cardinal Health™ Monoject™ Renewables Sharps Containers are an environmentally responsible option made from post-industrial recycled resin. These containers all maintain the same high level of safety and clinical acceptability. While the containers may differ visually, they retain the same overall robust design and puncture resistance for safe usage and handling that has made Cardinal Health an industry leader in the disposable sharps containers category.

- Manufactured using recycled plastics
- Reduces waste into landfills and incinerators
- Saves fossil fuels, water and other precious resources

### Look for the icon

Renewable sharps containers can be found under the corresponding use type categories.





## Monoject™ In-Room Sharps Containers

### Lid design matters

Select the lids that align with the safety standards and disposal needs of your facility

Level of safety:

**Standard**

**Safest**



**Always-open lid**


**Mailbox lid**


**Safety In-Room and SharpStar™ lid**

**Safety engineering features:**

|                          |   |   |   |
|--------------------------|---|---|---|
| Forced horizontal drop   | • | • | • |
| "Hands-out" locking tabs |   | • | • |
| One-handed disposal      | • |   | • |
| Counterbalanced lid      |   |   | • |
| "Hands-out" technology   |   |   | • |
| FULL indication          |   |   | • |
| Automatic shut-off       |   |   | • |
| Anti-kickback technology |   |   | • |



 **Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

 Biohazardous






## Monoject™ SharpStar™ Safety In-Room Sharps Containers

### Disposal designed to protect

Your goal is for every sharp to pass through your facility safely without injury or exposure. Count on Monoject™ SharpStar™ In-Room Containers from Cardinal Health to protect through the final steps of disposal with superior safety and value.

We engineered our SharpStar™ In-Room Containers to meet OSHA standards for ease of use and needlestick prevention, including one-handed disposal, secure closures, locking tabs positioned to keep hands out and an automatic shut-off mechanism to avoid overfilling. Safety-engineered controls work behind the scenes to improve clinician confidence, while a forced horizontal drop stores waste flat to help make the most of your container space and use fewer containers.



**SharpStar™ In-Room Container**



**Sharps container without a counterbalanced lid**

The counterbalanced lid forces needles and syringes to drop horizontally, maximizing space and ultimately allowing a user to fit **10%\* more needles and syringes** than a non-counterbalanced lid.

**Reference:**

\* Data on file. Internal Study. Cardinal Health, 2020



## Monoject™ SharpStar™ In-Room Sharps Containers

### Engineering features include:

- Horizontal drop maximizes container space
- Bold red FULL indicator alerts user to replace container
- Restricts hand accessibility. Directs sharps away from user
- Safety mechanism prevents over-filling
- Prevents sharps from exiting the lid
- Permanently locks to keep hands away from contents



8541SA

8541SA

| Order code | Description                     | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W x H in. | Case qty. |
|------------|---------------------------------|---------------------------------------|-----------------------------------------|-----------|
| 8506SA     | 5 quart, clear                  | 12.5 x 5.5 x 10.75                    | 8.75 x 1.25 x 1.25                      | 20        |
| 8507SA     | 5 quart, transparent red        | 12.5 x 5.5 x 10.75                    | 8.75 x 1.25 x 1.25                      | 20        |
| 8507MW     | Renewable, 5 quart, opaque red  | 12.5 x 5.5 x 10.75                    | 8.75 x 1.25 x 1.25                      | 20        |
| 8533SA     | 2 gallon, clear                 | 11.75 x 6 x 13.75                     | 8.75 x 3.75 x 1.75                      | 10        |
| 8534SA     | 2 gallon, transparent red       | 11.75 x 6 x 13.75                     | 8.75 x 3.75 x 1.75                      | 10        |
| 8536SA     | 3 gallon, clear                 | 16.5 x 6 x 13.75                      | 8.75 x 3.75 x 1.75                      | 10        |
| 8537SA     | 3 gallon, transparent red       | 16.5 x 6 x 13.75                      | 8.75 x 3.75 x 1.75                      | 10        |
| 8537MW     | Renewable, 3 gallon, opaque red | 16.5 x 6 x 13.75                      | 8.75 x 3.75 x 1.75                      | 10        |
| 8540SA     | 4 gallon, clear                 | 18.5 x 6 x 16.5                       |                                         | 10        |
| 8541SA     | 4 gallon, transparent red       | 18.5 x 6 x 16.5                       |                                         | 10        |



**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

Biohazardous



## Monoject™ Safety In-Room Sharps Containers

### Engineering features include:

- Wider opening\* to accommodate larger waste or devices
- Easy visibility of recommended fill line
- Signals when container's capacity is reached with a full indicator
- Horizontal drop system maximizes full use of container
- Thumb tabs lock the container closed for final disposal when full
- Complete line of brackets and accessories available



| Order code | Description               | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W x H in. | Case qty. |
|------------|---------------------------|---------------------------------------|-----------------------------------------|-----------|
| 31353603   | 5 quart, transparent red  | 14 x 6 x 13                           | 11.5 x 5 x 1.75                         | 14        |
| 31323333   | 2 gallon, transparent red | 17.5 x 6 x 14                         | 11.5 x 5 x 1.75                         | 12        |
| 31314886   | 3 gallon, transparent red | 20.5 x 6 x 14                         | 11.5 x 5 x 1.75                         | 12        |

### Reference:

\*2.75 in. L x 3.75 in. W x .5 in. H wider lid opening dimensions compared to 5 quart SharpStar™ In-room container



**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

■ Biohazardous



## Monoject™ Mailbox In-Room Sharps Containers

### Engineering features include:

- Mailbox-style lid can be kept closed when container is not in use
- Economical horizontal drop maximizes container volume
- Containers lock for final disposal
- Wall enclosures, brackets and adhesive-backed holders are available

| Order code | Description              | Container dimensions<br>H x D x W in. | Case<br>qty. |
|------------|--------------------------|---------------------------------------|--------------|
| 85121      | 5 quart, clear           | 11 x 4.75 x 10.75                     | 20           |
| 85131      | 5 quart, transparent red | 11 x 4.75 x 10.75                     | 20           |



85121

## Monoject™ Always-Open In-Room Sharps Containers

### Engineering features include:

- Always-open lids provide convenient, quick and easy point-of-use sharps disposal
- Accommodates larger syringes
- Economical horizontal drop maximizes container volume
- Containers easily lock for final disposal
- 8.75 in. wide opening will accept up to 60mL syringe
- Wall enclosures, brackets and adhesive-backed table top holders are available

| Order code | Description               | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W x H in. | Case<br>qty. |
|------------|---------------------------|---------------------------------------|-----------------------------------------|--------------|
| 851201     | 5 quart, clear            | 11 x 4.75 x 10.75                     | 8.75 x 1.75 x 1.75                      | 20           |
| 851301     | 5 quart, transparent red  | 11 x 4.75 x 10.75                     | 8.75 x 1.75 x 1.75                      | 20           |
| 85221      | 3 gallon, clear           | 16.25 x 6 x 13.75                     | 8.75 x 2.25 x 2.5                       | 10           |
| 85221R     | 3 gallon, transparent red | 16.25 x 6 x 13.75                     | 8.75 x 2.25 x 2.5                       | 10           |



851201



851301



**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

■ Biohazardous



## Monoject™ Multi-Purpose Sharps Containers

### Engineering features include:

- Adjustable rotor opening or hinged lids accommodate a variety of sharps sizes and provide temporary and permanent closure
- Easy, simple assembly
- Containers lock for final disposal
- Nestable containers save valuable storage space
- Available holders and brackets secure containers wherever needed
- Plastics and pigments are safe to autoclave or incinerate
- For use on countertops or mounted on the wall
- Economical horizontal drop maximizes container volume
- Large horizontal openings accommodate larger syringes

### Monoject™ Multi-Purpose Sharps Containers with Sliding Lids

| Order code | Description   | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in. | Case<br>qty. |
|------------|---------------|---------------------------------------|-------------------------------------|--------------|
| 31143699   | 1 gallon, red | 6.25 x 7 x 10.5                       | 8.5 x 1.7                           | 32           |
| 31142222   | 2 gallon, red | 10.25 x 7 x 10.5                      | 8.5 x 1.7                           | 20           |



### Monoject™ Multi-Purpose Sharps Containers with Horizontal-Drop Opening Lids

| Order code | Description               | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W x H in. | Case<br>qty. |
|------------|---------------------------|---------------------------------------|-----------------------------------------|--------------|
| 89671      | 2 gallon, transparent red | 12.75 x 7.25 x 10.5                   | 9 x 2.5 x 2.125                         | 20           |
| 89651      | 2 gallon, red             | 12.75 x 7.25 x 10.5                   | 9 x 2.5 x 2.125                         | 20           |



**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

■ Biohazardous




## Monoject™ Multi-Purpose Sharps Containers

### Monoject™ Multi-Purpose Sharps Containers with Hinged Lids

| Order code | Description   | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in.  | Case<br>qty. |
|------------|---------------|---------------------------------------|--------------------------------------|--------------|
| 8990SA     | 2 gallon, red | 10 x 7.25 x 10.5                      | Hinged: 4 x 6<br>Snap cap: 2.25 dia. | 20           |



### Monoject™ Multi-Purpose Sharps Containers with Rotor Opening Lids

| Order code | Description                                                                                                       | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in. | Case<br>qty. |
|------------|-------------------------------------------------------------------------------------------------------------------|---------------------------------------|-------------------------------------|--------------|
| 8950SA     | Round, 5 quart, red                                                                                               | 6.75 x 8.75                           | 3.5 x 1.5                           | 40           |
| 8970       | 2 gallon, red                                                                                                     | 10 x 7.25 x 10.5                      | 3.5 x 1.5                           | 20           |
| 8979MW     |  Renewable, 2 gallon, opaque red | 10 x 7.25 x 10.5                      | 3.5 x 1.5                           | 20           |



8970

### Monoject™ Multi-Purpose Sharps Containers with Rotor Opening and Hinged Lids

| Order code | Description               | Lid opening dimensions<br>L x W in.  | Case<br>qty. |
|------------|---------------------------|--------------------------------------|--------------|
| 8527R      | 3 gallon, transparent red | Rotor: 3.5 x 1.5<br>Hinged: 6 x 5.25 | 10           |




### Monoject™ Multi-Purpose Sharps Containers with Vertical Drop Lids

| Order code | Description   | Container dimensions<br>L x W x H in. | Case<br>qty. |
|------------|---------------|---------------------------------------|--------------|
| 8881676434 | 14 quart, red | 6.75 x 10.56 x 17.77                  | 10           |
| 8881676285 | 8 quart, red  | 6.75 x 10.56 x 10.89                  | 20           |
| 8881676236 | 4 quart, red  | 6.75 x 10.56 x 7.08                   | 40           |



8881676434

 **Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

 Biohazardous



## Monoject™ Phlebotomy Sharps Containers

### Engineering features include:

- Fit in most common caddies and trays
- Available holders and brackets
- A variety of container sizes and openings allow for disposal of other sharps such as lancets, butterfly tubing and small syringes
- Containers lock for final disposal

| Order code | Description                    | Container dimensions<br>H x D x W in. | Case<br>qty. |
|------------|--------------------------------|---------------------------------------|--------------|
| 8900SA     | 1 quart, red                   | 6.25 x 4.25 x 4.25                    | 100          |
| 1522SA     | 2.2 quart, red                 | 7.25 x 4.75 x 6.5                     | 60           |
| 8900MW     | Renewable, 1 quart, opaque red | 6.25 x 4.25 x 4.5                     | 100          |



## Monoject™ Large Volume Sharps Containers

### Engineering features include:

- Sliding, hinged and sealing gasket lid
- Adjustable sliding and hinged lids accommodate large sharps
- Lids can provide temporary closure when container is not in use to reduce exposure to contents
- Carrying handles and carts with wheels are available for easy transport and lifting of heavy, full containers
- Optional gasket lids provide leak-resistant containment of contents
- Easy, simple assembly
- Containers lock for final disposal
- Nestable containers save valuable storage space
- Sliding gasket lids have a round 6 in. opening to accommodate up to 1,000 mL bottles



## Monoject™ Safety Large Volume Sharps Containers

| Order code | Description               | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in. | Case<br>qty. |
|------------|---------------------------|---------------------------------------|-------------------------------------|--------------|
| 31143665   | 10 gallon, red, split lid | 21.5 x 12 x 15.5                      | Rectangular, 8.5 x 1.75             | 6            |



**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

Biohazardous



## Monoject™ Large Volume Sharps Containers

### Monoject™ Large Volume Containers with Hinged Lids

| Order code     | Description                                        | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in.                    | Case qty. |
|----------------|----------------------------------------------------|---------------------------------------|--------------------------------------------------------|-----------|
| <b>8980-</b>   | 8 gallon, red                                      | 17.75 x 11 x 15.5                     | Hinged opening: 5 x 9<br>Rotor opening: 3.5 x 1.5      | 10        |
| <b>8997</b>    | 8 gallon, red, with sealing gasket lid             | 17.75 x 11 x 15.5                     | 5.25 x 5.75                                            | 10        |
| <b>8997PG2</b> | 8 gallon, red, PGII rated with sealing gasket lid  | 17.75 x 11 x 15.5                     | 5.25 x 5.75                                            | 10        |
| <b>8933</b>    | 12 gallon, red                                     | 18.75 x 12.75 x 18.25                 | Hinged opening: 5.8 x 12.7<br>Rotor opening: 3.5 x 1.5 | 10        |
| <b>8932</b>    | 12 gallon, red, with sealing gasket lid            | 18.75 x 12.75 x 18.25                 | Hinged opening: 5.8 x 12.7<br>Rotor opening: 3.5 x 1.5 | 10        |
| <b>8932PG2</b> | 12 gallon, red, PGII rated with sealing gasket lid | 18.75 x 12.75 x 18.25                 | 5.875 x 12.7                                           | 10        |
| <b>8991</b>    | 18 gallon, red                                     | 26 x 12.75 x 18.25                    | Rotor opening: 3.5 x 1.5                               | 5         |
| <b>8998</b>    | 18 gallon, red, with sealing gasket lid            | 26 x 12.75 x 18.25                    | 5.875 x 12.7                                           | 5         |
| <b>8998PG2</b> | 18 gallon, red, PGII rated with sealing gasket lid | 26 x 12.75 x 18.25                    | 5.875 x 12.7                                           | 5         |
| <b>8930SA</b>  | 30 gallon, red, PGII rated with sealing gasket lid | 27.5 x 15.25 x 21.25                  |                                                        | 3         |



### Monoject™ Large Volume Containers with Sliding Lids

| Order code      | Description                                        | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in. | Case qty. |
|-----------------|----------------------------------------------------|---------------------------------------|-------------------------------------|-----------|
| <b>8980S</b>    | 8 gallon, red                                      | 17.75 x 11 x 15.5                     | 5.5 x 6.5                           | 10        |
| <b>8997S</b>    | 8 gallon, red, with sealing gasket lid             | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10        |
| <b>8997SPG2</b> | 8 gallon, red, PGII rated with sealing gasket lid  | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10        |
| <b>8935</b>     | 12 gallon, red                                     | 18.75 x 12.75 x 18.25                 | 7.25 x 6                            | 10        |
| <b>8936SA</b>   | 12 gallon, red with sealing gasket lid             | 18.75 x 12.75 x 18.25                 | 5.25 x 5.75                         | 10        |
| <b>8936PG2</b>  | 12 gallon, red, PGII rated with sealing gasket lid | 18.75 x 12.75 x 18.25                 | 5.25 x 5.75                         | 10        |
| <b>8938</b>     | 18 gallon, red                                     | 26 x 12.75 x 18.25                    | 7.25 x 6                            | 5         |
| <b>8998S</b>    | 18 gallon, red with sealing gasket lid             | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5         |
| <b>8998SPG2</b> | 18 gallon, red, PGII rated with sealing gasket lid | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5         |



**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

■ Biohazardous





## Monoject™ Transportable Sharps Containers

### Engineering features include:

- Simple, one-handed operation opens and closes the one-quart disposal opening
- Convenient size fits into carry bags for nursing, EMS and fire/rescue
- Large disposal opening easily accommodates small syringes and other sharps waste
- Sharps shuttle snap closure for one-handed opening and closing

| Order code | Description  | Container dimensions<br>H x D x W in. | Case<br>qty. |
|------------|--------------|---------------------------------------|--------------|
| 8303SA     | 1 quart, red | 8.75 x 2.5 x 4.5                      | 20           |



8303SA



**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

■ Biohazardous



## Monoject™ Trace Chemotherapy Containers

- For the disposal of trace chemotherapy waste
- Sealing gasket lids provide leak-resistant containment of fluid
- Absorbent pad in bottom of container absorbs excess liquid
- Adjustable sliding lid reduces the possibility of aerosolization and accommodates chemotherapy waste of various types and sizes
- Available carts with wheels and built-in carrying handles for easy transport and lifting of heavy, full containers
- Containers provide temporary and final locking closure for disposal
- Nestable containers save valuable storage space
- Sliding gasket lids have a round 6 in. opening to accommodate up to 1,000 mL bottles

### Monoject™ Trace Chemotherapy Containers with Sliding Lids


| Order code      | Description                                           | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in. | Case qty. |
|-----------------|-------------------------------------------------------|---------------------------------------|-------------------------------------|-----------|
| <b>8985S</b>    | 8 gallon, yellow                                      | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10        |
| <b>8985SPG2</b> | 8 gallon, yellow, PGII rated with sealing gasket lid  | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10        |
| <b>8934</b>     | 12 gallon, yellow                                     | 18.75 x 12.75 x 18.25                 | 5.25 x 5.75                         | 10        |
| <b>8934PG2</b>  | 12 gallon, yellow, PGII rated with sealing gasket lid | 18.75 x 12.75 x 18.25                 | 5.25 x 5.75                         | 10        |
| <b>8939</b>     | 18 gallon, yellow                                     | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5         |
| <b>8939PG2</b>  | 18 gallon, yellow, PGII rated with sealing gasket lid | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5         |



### Monoject™ Trace Chemotherapy Containers with Hinged Lids

| Order code     | Description                                           | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in. | Case qty. |
|----------------|-------------------------------------------------------|---------------------------------------|-------------------------------------|-----------|
| <b>8982</b>    | 2 gallon, yellow                                      | 10 x 7.25 x 10.5                      | 4 x 6                               | 20        |
| <b>8985</b>    | 8 gallon, yellow                                      | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10        |
| <b>8985PG2</b> | 8 gallon, yellow, PGII rated with sealing gasket lid  | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10        |
| <b>8931</b>    | 12 gallon, yellow                                     | 18.75 x 12.75 x 18.25                 | 5.875 x 12.7                        | 10        |
| <b>8931PG2</b> | 12 gallon, yellow, PGII rated with sealing gasket lid | 18.75 x 12.75 x 18.25                 | 5.875 x 12.7                        | 10        |
| <b>8989</b>    | 18 gallon, yellow                                     | 26 x 12.75 x 18.25                    | 5.875 x 12.7                        | 5         |
| <b>8989PG2</b> | 18 gallon, yellow, PGII rated with sealing gasket lid | 26 x 12.75 x 18.25                    | 5.875 x 12.7                        | 5         |



 **Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

 Trace chemotherapy



## Hazardous Waste Pharmaceutical (RCRA) and Universal Containers

- For the disposal of hazardous pharmaceutical waste
- RCRA hazardous waste label meets EPA satellite accumulation labeling requirements
- White and black color easily distinguished from sharps, trace chemo and other disposal containers
- Tamper-evident label helps deter tampering during final disposal
- Hazardous waste pharmaceutical label meets EPA 2019 guidelines
- Most are PGII rated, meet D.O.T. requirements for transport without secondary packaging
- Leak-resistant gasket and absorbent pad help contain liquid contents should a primary container break
- Locking floor carts

| Order code    | Description                                              | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in. | Case<br>qty. |
|---------------|----------------------------------------------------------|---------------------------------------|-------------------------------------|--------------|
| <b>8605RC</b> | RCRA, safety in-room counterbalanced lid, 5 quart, black | 11 x 6 x 13                           | 11.5 x 5 x 1.75                     | 14           |
| <b>8607U</b>  | Universal, hinged lid, 8 gallon, black                   | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8611U</b>  | Universal, hinged lid, 12 gallon, black                  | 18.75 x 12.75 x 18.25                 | 5.875 x 12.7                        | 10           |
| <b>8617U</b>  | Universal, hinged lid, 18 gallon, black                  | 26 x 12.75 x 18.25                    | 5.875 x 12.7                        | 5            |
| <b>8602RC</b> | RCRA, hinged lid, 2 gallon, black                        | 10 x 7.25 x 10.5                      | 4 x 6                               | 20           |
| <b>8607RC</b> | RCRA, hinged lid, 8 gallon, black                        | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8611RC</b> | RCRA, hinged lid, 12 gallon, black                       | 18.75 x 12.75 x 18.25                 | 5.875 x 12.7                        | 10           |
| <b>8617RC</b> | RCRA, hinged lid, 18 gallon, black                       | 26 x 12.75 x 18.25                    | 5.875 x 12.7                        | 5            |
| <b>8608U</b>  | Universal, sliding lid, 8 gallon, black                  | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8612U</b>  | Universal, sliding lid, 12 gallon, black                 | 18.75 x 12.75 x 18.25                 | 5.25 x 5.75                         | 10           |
| <b>8618U</b>  | Universal, sliding lid, 18 gallon, black                 | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5            |
| <b>8608RC</b> | RCRA, sliding lid, 8 gallon, black                       | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8612RC</b> | RCRA, sliding lid, 12 gallon, black                      | 18.75 x 12.75 x 18.25                 | 5.25 x 5.75                         | 10           |
| <b>8618RC</b> | RCRA, sliding lid, 18 gallon, black                      | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5            |



8618RC



8602RC



8605RC



**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

■ Hazardous (RCRA) and universal





## Non-Hazardous Pharmaceutical Waste Containers

- For the disposal of non-hazardous pharmaceutical waste
- White and blue color easily distinguished from sharps, chemotherapy and other disposal containers
- Leak-resistant gasket and absorbent pad help contain liquid contents
- Tamper-evident label helps deter tampering during final disposal
- Sturdy, rigid polypropylene construction provides excellent impact and leak resistance (unlike cardboard or fiber) to contain heavy vials, IV bags and more
- Two-gallon, three-gallon and eight-gallon capacity sizes are suitable for multiple locations — nursing units, med carts, anesthesia, pharmacies, etc.
- Secure, hinged lid with temporary and final closure accepts larger, bulkier waste; closes temporarily between uses; locks for final disposal
- Nestable containers save valuable storage space
- Locking wall-mounting brackets and floor carts available allow for point-of-use placement and one-handed disposal; increases mobility, stability and security



| Order code | Description                                            | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in. | Case<br>qty. |
|------------|--------------------------------------------------------|---------------------------------------|-------------------------------------|--------------|
| 8820       | 2 gallon, white, gasketed hinged lid                   | 10 x 7.25 x 10.5                      | 4 x 6                               | 20           |
| 8836SA     | 3 gallon, white, counterbalanced lid                   | 16.5 x 6 x 13.75                      |                                     | 10           |
| 8850       | 8 gallon, white, gasketed hinged lid                   | 17.75 x 11 x 15.5                     |                                     | 10           |
| 8851       | 8 gallon, white, gasketed sliding lid                  | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| 8860-      | 12 gallon, white gasketed hinged lid                   | 18.75 x 12.75 x 18.25                 |                                     | 10           |
| 8870-      | 18 gallon, white, gasketed hinged lid                  | 26 x 12.75 x 18.25                    |                                     | 5            |
| 8871       | 18 gallon, white, PGII rated with gasketed sliding lid | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5            |

 **Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

 Non-hazardous



## Non-Infectious Waste Containers

EnviroStar Sharps Disposal Containers:

- Environmental green in color differentiate from other containers
- Designed to help segregate waste to reduce waste costs



| Order code | Description                 | Container dimensions<br>H x D x W in. | Case<br>qty. |
|------------|-----------------------------|---------------------------------------|--------------|
| 8781       | Hinged lid, 8 gallon, green | 17.75 x 15.5 x 11                     | 10           |



**Reminder:** Containers may vary in color, but are most generally indicated by colors in construction.

■ Non-infectious



## PG II Rated Containers

The PG II Rated containers allow for transport of sharps, cultures and stocks. This rating eliminates the need for costly secondary packaging and provides a solution for safe disposal of sharps, chemotherapy waste, cultures and stocks.



| Order code      | Description                                                            | Container dimensions<br>H x D x W in. | Lid opening dimensions<br>L x W in. | Case<br>qty. |
|-----------------|------------------------------------------------------------------------|---------------------------------------|-------------------------------------|--------------|
| <b>8985PG2</b>  | Chemotherapy, 8 gallon, yellow, PGII rated with hinged gasketed lid    | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8985SPG2</b> | Chemotherapy, 8 gallon, yellow, PGII rated with sliding gasketed lid   | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8931PG2</b>  | Chemotherapy, 12 gallon, yellow, PGII rated with hinged gasketed lid   | 18.75 x 12.75 x 18.25                 | 5.875 x 12.7                        | 10           |
| <b>8934PG2</b>  | Chemotherapy, 12 gallon, yellow, PGII rated with sliding gasketed lid  | 18.75 x 12.75 x 18.25                 | 5.25 x 5.75                         | 10           |
| <b>8989PG2</b>  | Chemotherapy, 18 gallon, yellow, PGII rated with hinged gasketed lid   | 26 x 12.75 x 18.25                    | 5.875 x 12.7                        | 5            |
| <b>8939PG2</b>  | Chemotherapy, 18 gallon, yellow, PGII rated with sliding gasketed lid  | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5            |
| <b>8997PG2</b>  | Biohazardous, 8 gallon, red, PGII rated with hinged gasketed lid       | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8997SPG2</b> | Biohazardous, 8 gallon, red, PGII rated with sliding gasketed lid      | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8932PG2</b>  | Biohazardous, 12 gallon, red, PGII rated with hinged gasketed lid      | 18.75 x 12.75 x 18.25                 | 5.875 x 12.7                        | 10           |
| <b>8936PG2</b>  | Biohazardous, 12 gallon, red, PGII rated with sliding gasketed lid     | 18.75 x 12.75 x 18.25                 | 5.25 x 5.75                         | 10           |
| <b>8998PG2</b>  | Biohazardous, 18 gallon, red, PGII rated with hinged gasketed lid      | 26 x 12.75 x 18.25                    | 5.875 x 12.7                        | 5            |
| <b>8998SPG2</b> | Biohazardous, 18 gallon, red, PGII rated with sliding gasketed lid     | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5            |
| <b>8930SA</b>   | Biohazardous, 30 gallon, red, PGII rated with hinged gasketed lid      | 27.5 x 15.25 x 21.25                  |                                     | 3            |
| <b>8871</b>     | Pharmaceutical, 18 gallon, white, PGII rated with sliding gasketed lid | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5            |
| <b>8602RC</b>   | RCRA, 2 gallon, black, PGII rated with hinged gasketed lid             | 10 x 7.25 x 10.5                      | 4 x 6                               | 20           |
| <b>8607RC</b>   | RCRA, 8 gallon, black, PGII rated with hinged gasketed lid             | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8608RC</b>   | RCRA, 8 gallon, black, PGII rated with sliding gasketed lid            | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8611RC</b>   | RCRA, 12 gallon, black, PGII rated with hinged gasketed lid            | 18.75 x 12.75 x 18.25                 | 5.875 x 12.7                        | 10           |
| <b>8612RC</b>   | RCRA, 12 gallon, black, PGII rated with sliding gasketed lid           | 18.75 x 12.75 x 18.25                 | 5.25 x 5.75                         | 10           |
| <b>8617RC</b>   | RCRA, 18 gallon, black, PGII rated with hinged gasketed lid            | 26 x 12.75 x 18.25                    | 5.875 x 12.7                        | 5            |
| <b>8618RC</b>   | RCRA, 18 gallon, black, PGII rated with sliding gasketed lid           | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5            |
| <b>8607U</b>    | Universal, 8 gallon, black, PGII rated with hinged gasketed lid        | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8608U</b>    | Universal, 8 gallon, black, PGII rated with sliding gasketed lid       | 17.75 x 11 x 15.5                     | 5.25 x 5.75                         | 10           |
| <b>8611U</b>    | Universal, 12 gallon, black, PGII rated with hinged gasketed lid       | 18.75 x 12.75 x 18.25                 | 5.875 x 12.7                        | 10           |
| <b>8612U</b>    | Universal, 12 gallon, black, PGII rated with sliding gasketed lid      | 18.75 x 12.75 x 18.25                 | 5.25 x 5.75                         | 10           |
| <b>8617U</b>    | Universal, 18 gallon, black, PGII rated with hinged gasketed lid       | 26 x 12.75 x 18.25                    | 5.875 x 12.7                        | 5            |
| <b>8618U</b>    | Universal, 18 gallon, black, PGII rated with sliding gasketed lid      | 26 x 12.75 x 18.25                    | 5.25 x 5.75                         | 5            |




## Monoject™ Container Hardware Brackets

Brackets secure sharps containers to walls and med carts by maintaining the container upright throughout use. Locking brackets are available to help prevent access to contents in high traffic areas.

- Locking wall brackets are unobtrusive and securely hold container in place
- According to NIOSH's guidance, containers should only be used with hardware that is manufactured specifically for them<sup>1</sup>
- OSHA Regulation 29 CFR 1910.1030(d)(4)(iii), "During use, containers for contaminated sharps shall be maintained upright throughout use..."
- Cardinal Health offers a complete line of brackets and holders to ensure your safety and compliance



8528C  
Container not included



8518X



31307062



8963

| Order code | Description                                                                                                   | Dimensions<br>H x D x W in. | Case<br>qty. |
|------------|---------------------------------------------------------------------------------------------------------------|-----------------------------|--------------|
| 31307054   | Locking wire bracket, 5 quart, 2 and 3 gallon                                                                 | -                           | 10           |
| 31307062   | Non-locking 5 quart, 2 and 3 gallon safety in-room                                                            | -                           | 10           |
| 31139747   | Extra keys for 31158960 (extra keys for hidden bracket set)                                                   | -                           | 10           |
| 31158960   | Hidden bracket and key for safety in-room container                                                           | -                           | 20           |
| 8518X      | Locking bracket 2 and 5 quart in-room                                                                         | 8.25 x 6.5 x 11.25          | 5            |
| 8524C      | Non-locking bracket for 2 and 3 gallon in-room                                                                | 11.75 x 6.5 x 13.75         | 5            |
| 8528C      | Locking bracket for 2 and 3 gallon in-room                                                                    | 11.75 x 6.5 x 13.75         | 5            |
| 89125A     | Non-locking bracket for 1 quart phlebotomy container                                                          | 5 x 4.25 x 5                | 5            |
| 8926       | Locking wire bracket for 2 quart multi-purpose with rotor lid                                                 | -                           | 5            |
| 8963       | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical and chemotherapy containers | 8.25 x 7.5 x 11.5           | 5            |
| 8519C      | Non-locking bracket for in-room sharps container, 2 quart and 5 quart                                         | 8.25 x 6.5 x 11.25          | 5            |
| 8881676012 | Aluminum bracket for 4, quart, 8 quart and 14 quart vertical drop containers                                  | -                           | 10           |
| 8975       | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical, chemotherapy containers                       | 7.5 x 6.5 x 10              | 5            |
| 8984-      | Locking wire wall bracket for 8 gallon containers                                                             | 18 x 9.75 x 15.5            | 5            |

**Reference:**

1. <https://www.cdc.gov/niosh/docs/97-111/pdfs/97-111.pdf?id=10.26616/NIOSH PUB97111>



## Monoject™ Container Hardware Wall Enclosures

Mounted wall enclosures secure sharps containers and help prevent access to contents in high traffic areas. Contents can be viewed by hospital personnel, but unsightly items are not fully visible to patients and visitors. Optional glove box holders keep gloves readily available.

- Sturdy, locked sharps container mounting allows for secure point-of-use disposal of sharps
- For convenience, all cabinets are universally keyed
- Durable, long lasting ABS plastic provides a secondary barrier of protection
- Easily monitored fill level through viewing windows on front and sides
- Keep gloves readily available with optional glove dispensers
- Beige, unobtrusive style fits any room decor



| Order code        | Description                                                                                                                      | Dimensions<br>H x D x W in. | Case<br>qty. |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------|-----------------------------|--------------|
| <b>85161H</b>     | Wall enclosure for in-room sharps container, 5 quart                                                                             | 11.75 x 5.5 x 13            | 1            |
| <b>85165H</b>     | Wall enclosure for in-room sharps container, 5 quart                                                                             | 11.75 x 5.5 x 13            | 5            |
| <b>8556H</b>      | Wall enclosure with glove box and container for 5 quart in-room container                                                        | 24.25 x 5.5 x 13            | 1            |
| <b>85301H</b>     | Wall enclosure for in-room sharps container, 2 and 3 gallon                                                                      | 17.5 x 6.5 x 14.5           | 1            |
| <b>85401H-</b>    | Wall enclosure for 3 and 4 gallon in-room container                                                                              | 26.5 x 7.25 x 17.75         | 1            |
| <b>8554B</b>      | Wall mounting bracket for glove dispensers, 8550LG and 8555SA                                                                    | 8 x 4.5 x 12.5              | 5            |
| <b>8516HDL</b>    | Replacement wall enclosure door with window and lock for 2 and 5 quart in-room sharps container                                  | -                           | 1            |
| <b>8550B</b>      | Glove box with attached wall mounting brackets (8554b)                                                                           | 7.25 x 4.5 x 12.25          | 5            |
| <b>31307005</b>   | Wall enclosure for 2 and 3 gallon safety in-room sharps container                                                                | 16 x 6.5 x 14.25            | 1            |
| <b>31353553</b>   | Wall enclosure for 5 quart safety in-room sharps container                                                                       | 9.75 x 6.75 x 13.5          | 1            |
| <b>8556H-</b>     | Wall enclosure for 2 and 5 quart safety in-room sharps container with glove dispenser (8550LG) and SharpStar™ container (85075A) | 24.25 x 5.5 x 13            | 1            |
| <b>88361H</b>     | Wall enclosure for 3 gallon pharmaceutical container                                                                             | 17.5 x 6.5 x 14.5           | 1            |
| <b>8881676624</b> | Locking wall bracket for 4 and 8 quart vertical drop containers                                                                  | 12 x 13 x 7                 | 2            |





## Monoject™ Container Hardware Glove Box Dispensers

- Glove boxes attach to wall enclosures as an added convenience
- Sturdy, durable, long lasting ABS plastic
- Keep gloves readily available
- Beige, unobtrusive style fits any room decor



| Order code | Description                                            | Dimensions<br>H x D x W in. | Case<br>qty. |
|------------|--------------------------------------------------------|-----------------------------|--------------|
| 31353553   | Wall cabinet for safety in-room container, 5 quart     | 9.75 x 6.75 x 13.5          | 1            |
| 31307005   | Wall cabinet, safety in-room container, 2 and 3 gallon | 16 x 6.5 x 14.25            | 1            |
| 31320032   | Glove box holder for safety in-room cabinet            | 6.5 x 4.5 x 12              | 10           |
| 31139762   | Glove box holder for safety in-room container          | 5 x 4 x 3.5                 | 20           |
| 8550B      | Glove box and mounting brackets                        | 7.25 x 4.5 x 12.25          | 5            |
| 8550LG     | Universal glove dispenser, top loading                 | 6.75 x 4.5 x 12             | 10           |
| 8554B      | Double-up wall mounting brackets                       | 8 x 4.5 x 12.5              | 5            |
| 8555SA     | Multi-glove dispenser, stackable, side loading         | 8 x 4.5 x 12.5              | 10           |



## Monoject™ Container Hardware

### Table Top Holders

Table top holders provide stability by maintaining the sharps container upright throughout use. Ideal for sharps containers used on counter tops or med carts. Adhesive strip is provided on the bottom of each tray holder to prevent movement of container.



| Order code | Description                                                                                                   | Dimensions<br>H x D x W in. | Case<br>qty. |
|------------|---------------------------------------------------------------------------------------------------------------|-----------------------------|--------------|
| 8910-      | Table top holder for 1 quart phlebotomy container                                                             | 2.75 x 3.75 x 3.75          | 10           |
| 1523SA     | Table top holder for 2.2 quart phlebotomy container                                                           | 2.25 x 4.25 x 5.75          | 5            |
| 8539SA     | Table top holder for 2 gallon safety in-room container                                                        | 3.25 x 7.75 x 11            | 5            |
| 8974       | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety container, almond, 2 gallon | 3.25 x 7.75 x 11            | 5            |
| 8529       | Table top holder, multi-purpose SharpStar™, and in-room container, almond, 3 gallon                           | 2.5 x 6.5 x 13              | 5            |
| 8515-      | Sharps container holder 5 quart                                                                               | 3.25 x 7.75 x 14.25         | 5            |



## Monoject™ Container Hardware

### Large Volume Carts and Foot Pedal Carts

- One-handed disposal
- Lids lift or slide open easily with foot-operated pedal
- Lids may remain closed when not in use to reduce exposure to contents
- Ergonomic handle is telescopic when transporting and retractable when stationary
- Heavy containers can be removed from the side with minimal lifting
- OSHA Regulation 29 CFR 1910.1030(d)(4)(iii), "During use, containers for contaminated sharps shall be maintained upright throughout use..."
- Cardinal Health offers a complete line of brackets and carts for large volume containers to ensure your safety and compliance



8991FP



8992H



8938FP - container not included

| Order code      | Description                                                | Dimensions<br>H x D x W in. | Case<br>qty. |
|-----------------|------------------------------------------------------------|-----------------------------|--------------|
| <b>8992H</b>    | Floor cart for large volume container, 8, 12 and 18 gallon | 22 x 12.75 x 22.5           | 1            |
| <b>8980FP</b>   | Foot pedal cart for hinged lid container, 8 gallon         | 22 x 21.5 x 16              | 1            |
| <b>8981FP</b>   | Foot pedal cart for slide lid container, 8 gallon          | 22 x 21.5 x 16              | 1            |
| <b>8991FP</b>   | Foot pedal cart for hinged lid container, 12 and 18 gallon | 31 x 22.75 x 15.25          | 1            |
| <b>8938FP</b>   | Foot pedal cart for slide lid container, 12 and 18 gallon  | 31 x 22.75 x 15.25          | 1            |
| <b>8993H</b>    | Rolling cart for large volume container, 30 gallon         | 35.5 x 22.25 x 15.25        | 1            |
| <b>31140109</b> | 4-castered cart for 7 and 10 gallon containers             |                             | 1            |

# Hardware and container cross reference guide

Search by container order code to identify corresponding hardware options.

| Container order code | Container size | Hardware order code | Hardware description                                                                                                                   | Locking Bracket? | Hardware case qty. |
|----------------------|----------------|---------------------|----------------------------------------------------------------------------------------------------------------------------------------|------------------|--------------------|
| <b>Biohazardous</b>  |                |                     |                                                                                                                                        |                  |                    |
| 1522SA               | 2.2 quart      | 1523SA              | Table top holder for 2.2 quart phlebotomy container                                                                                    |                  | 5 each             |
| 1525SA               | 2 gallon       | 8974                | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety container, almond, 2 gallon                          |                  | 5 each             |
| 31142222             | 2 gallon       | 8963                | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical and chemotherapy containers                          | ✓                | 5 each             |
| 31142222             | 2 gallon       | 8974                | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety container, almond, 2 gallon                          |                  | 5 each             |
| 31142222             | 2 gallon       | 8975                | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical and chemotherapy containers                                             |                  | 5 each             |
| 31143665             | 10 gallon      | 31140109            | 4-castered cart for 7 and 10 gallon containers                                                                                         |                  | 1 each             |
| 31143699             | 1 gallon       | 8963                | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical and chemotherapy containers                          | ✓                | 5 each             |
| 31143897             | 5 quart        | 31158960            | Hidden bracket and key for safety in-room container                                                                                    | ✓                | 20 each            |
| 31144010             | 5 quart        |                     | Wall cabinet, almond, for 5 quart tortuous path sharps container                                                                       | ✓                | 10 each            |
| 31144010             | 5 quart        | 31158960            | Hidden bracket and key for safety in-room container                                                                                    | ✓                | 20 each            |
| 31158960             |                | 31139747            | Extra keys for hidden bracket                                                                                                          |                  | 10 each            |
| 31314886             | 3 gallon       | 31307005            | Wall enclosure for 2 gallon and 3 gallon safety in-room sharps container                                                               | ✓                | 1 each             |
| 31314886             | 3 gallon       | 31307054            | Locking wire bracket, 5 quart, 2 and 3 gallon                                                                                          | ✓                | 10 each            |
| 31323333             | 2 gallon       | 31307005            | Wall enclosure for 2 gallon and 3 gallon safety in-room sharps container                                                               | ✓                | 1 each             |
| 31323333             | 2 gallon       | 31307054            | Locking wire bracket, 5 quart, 2 and 3 gallon                                                                                          | ✓                | 10 each            |
| 31323333             | 2 gallon       | 31307062            | Non-locking 5 quart, 2 and 3 gallon safety in-room                                                                                     |                  | 10 each            |
| 31353603             | 5 quart        | 31307054            | Locking wire bracket, 5 quart, 2 and 3 gallon                                                                                          | ✓                | 10 each            |
| 31353603             | 5 quart        | 31307062            | Non-locking 5 quart, 2 and 3 gallon safety in-room                                                                                     |                  | 10 each            |
| 31353603             | 5 quart        | 31353553            | Wall enclosure for 5 quart safety in-room sharps container                                                                             | ✓                | 1 each             |
| 8506SA               | 5 quart        | 8515-               | Sharps container holder 5 quart                                                                                                        |                  | 5 each             |
| 8506SA               | 5 quart        | 85161H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 1 each             |
| 8506SA               | 5 quart        | 85165H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 5 each             |
| 8506SA               | 5 quart        | 8518X               | Locking bracket 2 quart and 5 quart in-room                                                                                            | ✓                | 5 each             |
| 8506SA               | 5 quart        | 8519C               | Non-locking bracket for in-room sharps container, 2 quart and 5 quart                                                                  |                  | 5 each             |
| 8506SA               | 5 quart        | 8556H-              | Wall enclosure for 2 quart and 5 quart safety in-room sharps container with glove dispenser (8550lg) and SharpStar™ container (8507sa) | ✓                | 1 each             |
| 8507MW               | 5 quart        | 8515-               | Sharps container holder 5 quart                                                                                                        |                  | 5 each             |
| 8507MW               | 5 quart        | 85161H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 1 each             |
| 8507MW               | 5 quart        | 85165H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 5 each             |
| 8507MW               | 5 quart        | 8518X               | Locking bracket 2 quart and 5 quart in-room                                                                                            | ✓                | 5 each             |
| 8507MW               | 5 quart        | 8519C               | Non-locking bracket for in-room sharps container, 2 quart and 5 quart                                                                  |                  | 5 each             |
| 8507MW               | 5 quart        | 8556H-              | Wall enclosure for 2 quart and 5 quart safety in-room sharps container with glove dispenser (8550lg) and SharpStar™ container (8507sa) | ✓                | 1 each             |
| 8507SA               | 5 quart        | 8515-               | Sharps container holder 5 quart                                                                                                        |                  | 5 each             |
| 8507SA               | 5 quart        | 85161H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 1 each             |
| 8507SA               | 5 quart        | 85165H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 5 each             |
| 8507SA               | 5 quart        | 8518X               | Locking bracket 2 quart and 5 quart in-room                                                                                            | ✓                | 5 each             |
| 8507SA               | 5 quart        | 8519C               | Non-locking bracket for in-room sharps container, 2 quart and 5 quart                                                                  |                  | 5 each             |
| 8507SA               | 5 quart        | 8556H-              | Wall enclosure for 2 quart and 5 quart safety in-room sharps container with glove dispenser (8550lg) and SharpStar™ container (8507sa) | ✓                | 1 each             |
| 851201               | 5 quart        | 8515-               | Sharps container holder 5 quart                                                                                                        |                  | 5 each             |
| 851201               | 5 quart        | 85161H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 1 each             |
| 851201               | 5 quart        | 85165H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 5 each             |
| 851201               | 5 quart        | 8518X               | Locking bracket 2 quart and 5 quart in-room                                                                                            | ✓                | 5 each             |
| 851201               | 5 quart        | 8519C               | Non-locking bracket for in-room sharps container, 2 quart and 5 quart                                                                  |                  | 5 each             |
| 851201               | 5 quart        | 8556H-              | Wall enclosure for 2 quart and 5 quart safety in-room sharps container with glove dispenser (8550lg) and SharpStar™ container (8507sa) | ✓                | 1 each             |
| 85121                | 5 quart        | 8515-               | Sharps container holder 5 quart                                                                                                        |                  | 5 each             |
| 85121                | 5 quart        | 85161H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 1 each             |
| 85121                | 5 quart        | 85165H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 5 each             |
| 85121                | 5 quart        | 8518X               | Locking bracket 2 quart and 5 quart in-room                                                                                            | ✓                | 5 each             |
| 85121                | 5 quart        | 8519C               | Non-locking bracket for in-room sharps container, 2 quart and 5 quart                                                                  |                  | 5 each             |

# Hardware and container cross reference guide (continued)

Search by container order code to identify corresponding hardware options.

| Container order code | Container size | Hardware order code | Hardware description                                                                                                                   | Locking Bracket? | Hardware case qty. |
|----------------------|----------------|---------------------|----------------------------------------------------------------------------------------------------------------------------------------|------------------|--------------------|
| <b>Biohazardous</b>  |                |                     |                                                                                                                                        |                  |                    |
| 85121                | 5 quart        | 8556H-              | Wall enclosure for 2 quart and 5 quart safety in-room sharps container with glove dispenser (8550lg) and SharpStar™ container (8507sa) | ✓                | 1 each             |
| 851301               | 5 quart        | 8515-               | Sharps container holder 5 quart                                                                                                        |                  | 5 each             |
| 851301               | 5 quart        | 85161H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 1 each             |
| 851301               | 5 quart        | 85165H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 5 each             |
| 851301               | 5 quart        | 8518X               | Locking bracket 2 quart and 5 quart in-room                                                                                            | ✓                | 5 each             |
| 851301               | 5 quart        | 8519C               | Non-locking bracket for in-room sharps container, 2 quart and 5 quart                                                                  |                  | 5 each             |
| 851301               | 5 quart        | 8556H-              | Wall enclosure for 2 quart and 5 quart safety in-room sharps container with glove dispenser (8550lg) and SharpStar™ container (8507sa) | ✓                | 1 each             |
| 85131                | 5 quart        | 8515-               | Sharps container holder 5 quart                                                                                                        |                  | 5 each             |
| 85131                | 5 quart        | 85161H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 1 each             |
| 85131                | 5 quart        | 85165H              | Wall enclosure for in-room sharps container, 2 quart and 5 quart                                                                       | ✓                | 5 each             |
| 85131                | 5 quart        | 8518X               | Locking bracket 2 quart and 5 quart in-room                                                                                            | ✓                | 5 each             |
| 85131                | 5 quart        | 8519C               | Non-locking bracket for in-room sharps container, 2 quart and 5 quart                                                                  |                  | 5 each             |
| 85131                | 5 quart        | 8556H-              | Wall enclosure for 2 quart and 5 quart safety in-room sharps container with glove dispenser (8550lg) and SharpStar™ container (8507sa) | ✓                | 1 each             |
| 85161H               |                | 8516HDL             | Replacement wall enclosure door with window and lock for 2 quart and 5 quart in room sharps container                                  |                  | 1 each             |
| 85221                | 3 gallon       | 8524C               | Non-locking bracket for 2 gallon and 3 gallon in-room                                                                                  |                  | 5 each             |
| 85221                | 3 gallon       | 8528C               | Locking bracket for 2 and 3 gallon in-room                                                                                             | ✓                | 5 each             |
| 85221                | 3 gallon       | 8529                | Table top holder, Multi-Purpose SharpStar, and In-Room Container, almond, 3 gallon                                                     |                  | 5 each             |
| 85221                | 3 gallon       | 85301H              | Wall enclosure for in-room sharps container, 2 and 3 gallon                                                                            | ✓                | 1 each             |
| 85221R               | 3 gallon       | 8524C               | Non-locking bracket for 2 gallon and 3 gallon in-room                                                                                  |                  | 5 each             |
| 85221R               | 3 gallon       | 8528C               | Locking bracket for 2 and 3 gallon in-room                                                                                             | ✓                | 5 each             |
| 85221R               | 3 gallon       | 8529                | Table top holder, Multi-Purpose SharpStar, and In-Room container, almond, 3 gallon                                                     |                  | 5 each             |
| 85221R               | 3 gallon       | 85301H              | Wall enclosure for in-room sharps container, 2 and 3 gallon                                                                            | ✓                | 1 each             |
| 8527R                | 3 gallon       | 8524C               | Non-locking bracket for 2 gallon and 3 gallon in-room                                                                                  |                  | 5 each             |
| 8527R                | 3 gallon       | 8528C               | Locking bracket for 2 and 3 gallon in-room                                                                                             | ✓                | 5 each             |
| 8527R                | 3 gallon       | 8529                | Table top holder, Multi-Purpose SharpStar, and In-Room container, almond, 3 gallon                                                     |                  | 5 each             |
| 8533SA               | 2 gallon       | 8524C               | Non-locking bracket for 2 gallon and 3 gallon in-room                                                                                  |                  | 5 each             |
| 8533SA               | 2 gallon       | 8528C               | Locking bracket for 2 and 3 gallon in-room                                                                                             | ✓                | 5 each             |
| 8533SA               | 2 gallon       | 85301H              | Wall enclosure for in-room sharps container, 2 and 3 gallon                                                                            | ✓                | 1 each             |
| 8533SA               | 2 gallon       | 8539SA              | Table top holder for 2 gallon safety in-room container                                                                                 |                  | 5 each             |
| 8534SA               | 2 gallon       | 8524C               | Non-locking bracket for 2 gallon and 3 gallon in-room                                                                                  |                  | 5 each             |
| 8534SA               | 2 gallon       | 8528C               | Locking bracket for 2 and 3 gallon in-room                                                                                             | ✓                | 5 each             |
| 8534SA               | 2 gallon       | 85301H              | Wall enclosure for in-room sharps container, 2 and 3 gallon                                                                            | ✓                | 1 each             |
| 8534SA               | 2 gallon       | 8539SA              | Table top holder for 2 gallon safety in-room container                                                                                 |                  | 5 each             |
| 8536SA               | 3 gallon       | 8524C               | Non-locking bracket for 2 gallon and 3 gallon in-room                                                                                  |                  | 5 each             |
| 8536SA               | 3 gallon       | 8528C               | Locking bracket for 2 and 3 gallon in-room                                                                                             | ✓                | 5 each             |
| 8536SA               | 3 gallon       | 8529                | Table top holder, Multi-Purpose SharpStar™, and In-Room container, almond, 3 gallon                                                    |                  | 5 each             |
| 8536SA               | 3 gallon       | 85301H              | Wall enclosure for in-room sharps container, 2 and 3 gallon                                                                            | ✓                | 1 each             |
| 8537MW               | 3 gallon       | 8524C               | Non-locking bracket for 2 gallon and 3 gallon in-room                                                                                  |                  | 5 each             |
| 8537MW               | 3 gallon       | 8528C               | Locking bracket for 2 and 3 gallon in-room                                                                                             | ✓                | 5 each             |
| 8537MW               | 3 gallon       | 8529                | Table top holder, Multi-Purpose SharpStar™, and In-Room container, almond, 3 gallon                                                    |                  | 5 each             |
| 8537MW               | 3 gallon       | 85301H              | Wall enclosure for in-room sharps container, 2 and 3 gallon                                                                            | ✓                | 1 each             |
| 8537SA               | 3 gallon       | 8524C               | Non-locking bracket for 2 gallon and 3 gallon in-room                                                                                  |                  | 5 each             |
| 8537SA               | 3 gallon       | 8528C               | Locking bracket for 2 and 3 gallon in-room                                                                                             | ✓                | 5 each             |
| 8537SA               | 3 gallon       | 8529                | Table top holder, Multi-Purpose SharpStar™, and In-Room container, almond, 3 gallon                                                    |                  | 5 each             |
| 8537SA               | 3 gallon       | 85301H              | Wall enclosure for in-room sharps container, 2 and 3 gallon                                                                            | ✓                | 1 each             |
| 8540SA               | 4 gallon       | 85401H-             | Wall enclosure for 3 gallon and 4 gallon in-room container                                                                             | ✓                | 1 each             |
| 8541SA               | 4 gallon       | 85401H-             | Wall enclosure for 3 gallon and 4 gallon in-room container                                                                             | ✓                | 1 each             |
| 8881676236           | 4 quart        | 8881676012          | Aluminum bracket for 4, quart, 8 quart and 14 quart vertical drop containers                                                           | ✓                | 1 each             |
| 8881676236           | 4 quart        | 8881676624          | Locking wall bracket for 4 quart and 8 quart vertical drop containers                                                                  | ✓                | 1 each             |
| 8881676285           | 8 quart        | 8881676012          | Aluminum bracket for 4, quart, 8 quart and 14 quart vertical drop containers                                                           | ✓                | 10 each            |

# Hardware and container cross reference guide (continued)

Search by container order code to identify corresponding hardware options.

| Container order code | Container size | Hardware order code | Hardware description                                                                                           | Locking Bracket? | Hardware case qty. |
|----------------------|----------------|---------------------|----------------------------------------------------------------------------------------------------------------|------------------|--------------------|
| <b>Biohazardous</b>  |                |                     |                                                                                                                |                  |                    |
| <b>8881676285</b>    | 8 quart        | 8881676624          | Locking wall bracket for 4 quart and 8 quart vertical drop containers                                          | ✓                | 5 each             |
| <b>8881676434</b>    | 14 quart       | 8881676012          | Aluminum bracket for 4, quart, 8 quart and 14 quart vertical drop containers                                   | ✓                | 10 each            |
| <b>8900MW</b>        | 1 quart        | 8910-               | Table top holder for 1 quart phlebotomy container                                                              |                  | 10 each            |
| <b>8900MW</b>        | 1 quart        | 8912SA              | Non-locking bracket for 1 quart phlebotomy container                                                           |                  | 5 each             |
| <b>8900SA</b>        | 1 quart        | 8910-               | Table top holder for 1 quart phlebotomy container                                                              |                  | 10 each            |
| <b>8900SA</b>        | 1 quart        | 8912SA              | Non-locking bracket for 1 quart phlebotomy container                                                           |                  | 5 each             |
| <b>8930SA</b>        | 30 gallon      | 8993H               | Rolling cart for large volume container, 30 gallon                                                             |                  | 1 each             |
| <b>8932</b>          | 12 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8932</b>          | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8932PG2</b>       | 12 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8932PG2</b>       | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8933</b>          | 12 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8933</b>          | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8935</b>          | 12 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                      |                  | 1 each             |
| <b>8935</b>          | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8936PG2</b>       | 12 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                      |                  | 1 each             |
| <b>8936PG2</b>       | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8936SA</b>        | 12 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                      |                  | 1 each             |
| <b>8936SA</b>        | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8938</b>          | 18 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                      |                  | 1 each             |
| <b>8938</b>          | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>89651</b>         | 2 gallon       | 8963                | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical and chemotherapy containers  | ✓                | 5 each             |
| <b>89651</b>         | 2 gallon       | 8974                | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety container, almond, 2 gallon  |                  | 5 each             |
| <b>89651</b>         | 2 gallon       | 8975                | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical and chemotherapy containers                     |                  | 5 each             |
| <b>89671</b>         | 2 gallon       | 8963                | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical and chemotherapy containers  | ✓                | 5 each             |
| <b>89671</b>         | 2 gallon       | 8974                | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety container, almond, 2 gallon  |                  | 5 each             |
| <b>89671</b>         | 2 gallon       | 8975                | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical and chemotherapy containers                     |                  | 5 each             |
| <b>8970</b>          | 2 gallon       | 8963                | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical and chemotherapy containers  | ✓                | 5 each             |
| <b>8970</b>          | 2 gallon       | 8974                | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety container, almond, 2 gallon  |                  | 5 each             |
| <b>8970</b>          | 2 gallon       | 8975                | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical and chemotherapy containers                     |                  | 5 each             |
| <b>8979MW</b>        | 2 gallon       | 8963                | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical and chemotherapy containers  | ✓                | 5 each             |
| <b>8979MW</b>        | 2 gallon       | 8974                | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety container, almond, 2 gallon  |                  | 5 each             |
| <b>8979MW</b>        | 2 gallon       | 8975                | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical and chemotherapy containers                     |                  | 5 each             |
| <b>8980-</b>         | 8 gallon       | 8980FP              | Foot pedal cart for hinged lid container, 8 gallon                                                             |                  | 1 each             |
| <b>8980-</b>         | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                              | ✓                | 5 each             |
| <b>8980-</b>         | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8980S</b>         | 8 gallon       | 8981FP              | Foot pedal cart for slide lid container, 8 gallon                                                              |                  | 1 each             |
| <b>8980S</b>         | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                              | ✓                | 5 each             |
| <b>8980S</b>         | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8990SA</b>        | 2 gallon       | 8963                | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical, and chemotherapy containers | ✓                | 5 each             |
| <b>8990SA</b>        | 2 gallon       | 8974                | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety container, almond, 2 gallon  |                  | 5 each             |
| <b>8990SA</b>        | 2 gallon       | 8975                | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical and chemotherapy containers                     |                  | 5 each             |

# Hardware and container cross reference guide (continued)

Search by container order code to identify corresponding hardware options.

| Container order code | Container size | Hardware order code | Hardware description                                       | Locking Bracket? | Hardware case qty. |
|----------------------|----------------|---------------------|------------------------------------------------------------|------------------|--------------------|
| <b>Biohazardous</b>  |                |                     |                                                            |                  |                    |
| 8991                 | 18 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon |                  | 1 each             |
| 8991                 | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon |                  | 1 each             |
| 8997                 | 8 gallon       | 8980FP              | Foot pedal cart for hinged lid container, 8 gallon         |                  | 1 each             |
| 8997                 | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers          | ✓                | 5 each             |
| 8997                 | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon |                  | 1 each             |
| 8997PG2              | 8 gallon       | 8980FP              | Foot pedal cart for hinged lid container, 8 gallon         |                  | 1 each             |
| 8997PG2              | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers          | ✓                | 5 each             |
| 8997PG2              | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon |                  | 1 each             |
| 8997S                | 8 gallon       | 8981FP              | Foot pedal cart for slide lid container, 8 gallon          |                  | 1 each             |
| 8997S                | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers          | ✓                | 5 each             |
| 8997S                | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon |                  | 1 each             |
| 8997SPG2             | 8 gallon       | 8981FP              | Foot pedal cart for slide lid container, 8 gallon          |                  | 1 each             |
| 8997SPG2             | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers          | ✓                | 5 each             |
| 8997SPG2             | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon |                  | 1 each             |
| 8998                 | 18 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon |                  | 1 each             |
| 8998                 | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon |                  | 1 each             |
| 8998PG2              | 18 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon |                  | 1 each             |
| 8998PG2              | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon |                  | 1 each             |
| 8998S                | 18 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon  |                  | 1 each             |
| 8998S                | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon |                  | 1 each             |
| 8998SPG2             | 18 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon  |                  | 1 each             |
| 8998SPG2             | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon |                  | 1 each             |

# Hardware and container cross reference guide (continued)

Search by container order code to identify corresponding hardware options.

| Container order code      | Container size | Hardware order code | Hardware description                                                                                          | Locking Bracket? | Hardware case qty. |
|---------------------------|----------------|---------------------|---------------------------------------------------------------------------------------------------------------|------------------|--------------------|
| <b>Trace chemotherapy</b> |                |                     |                                                                                                               |                  |                    |
| <b>8931</b>               | 12 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8931</b>               | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8931PG2</b>            | 12 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8931PG2</b>            | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8934</b>               | 12 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8934</b>               | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8934PG2</b>            | 12 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8934PG2</b>            | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8939</b>               | 18 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8939</b>               | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8939PG2</b>            | 18 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8939PG2</b>            | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8982</b>               | 2 gallon       | 8963                | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical and chemotherapy containers | ✓                | 5 each             |
| <b>8982</b>               | 2 gallon       | 8974                | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety container, almond, 2 gallon |                  | 5 each             |
| <b>8982</b>               | 2 gallon       | 8975                | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical and chemotherapy containers                    |                  | 5 each             |
| <b>8985</b>               | 8 gallon       | 8980FP              | Foot pedal cart for hinged lid container, 8 gallon                                                            |                  | 1 each             |
| <b>8985</b>               | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 5 each             |
| <b>8985</b>               | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8985PG2</b>            | 8 gallon       | 8980FP              | Foot pedal cart for hinged lid container, 8 gallon                                                            |                  | 1 each             |
| <b>8985PG2</b>            | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 5 each             |
| <b>8985PG2</b>            | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8985S</b>              | 8 gallon       | 8981FP              | Foot pedal cart for slide lid container, 8 gallon                                                             |                  | 1 each             |
| <b>8985S</b>              | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 5 each             |
| <b>8985S</b>              | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8985SPG2</b>           | 8 gallon       | 8981FP              | Foot pedal cart for slide lid container, 8 gallon                                                             |                  | 1 each             |
| <b>8985SPG2</b>           | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 5 each             |
| <b>8985SPG2</b>           | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8989</b>               | 18 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8989</b>               | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8989PG2</b>            | 18 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8989PG2</b>            | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |



# Hardware and container cross reference guide (continued)

Search by container order code to identify corresponding hardware options.

| Container order code    | Container size | Hardware order code | Hardware description                                                                                          | Locking Bracket? | Hardware case qty. |
|-------------------------|----------------|---------------------|---------------------------------------------------------------------------------------------------------------|------------------|--------------------|
| <b>Hazardous (RCRA)</b> |                |                     |                                                                                                               |                  |                    |
| <b>31158960</b>         |                | 31139747            | Extra keys for hidden bracket                                                                                 |                  | 10 each            |
| <b>8602RC</b>           | 2 gallon       | 8963                | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical and chemotherapy containers | ✓                | 5 each             |
| <b>8602RC</b>           | 2 gallon       | 8974                | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety container, almond, 2 gallon |                  | 5 each             |
| <b>8602RC</b>           | 2 gallon       | 8975                | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical and chemotherapy containers                    |                  | 5 each             |
| <b>8605RC</b>           | 5 quart        | 31307054            | Locking wire bracket, 5 quart, 2 and 3 gallon                                                                 | ✓                | 10 each            |
| <b>8605RC</b>           | 5 quart        | 31307062            | Non-locking 5 quart, 2 and 3 gallon safety in-room                                                            |                  | 10 each            |
| <b>8605RC</b>           | 5 quart        | 31353553            | Wall enclosure for 5 quart safety in-room sharps container                                                    | ✓                | 1 each             |
| <b>8607RC</b>           | 8 gallon       | 8980FP              | Foot pedal cart for hinged lid container, 8 gallon                                                            |                  | 1 each             |
| <b>8607RC</b>           | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 5 each             |
| <b>8607RC</b>           | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8608RC</b>           | 8 gallon       | 8981FP              | Foot pedal cart for slide lid container, 8 gallon                                                             |                  | 1 each             |
| <b>8608RC</b>           | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 5 each             |
| <b>8608RC</b>           | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8611RC</b>           | 12 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8611RC</b>           | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8612RC</b>           | 12 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8612RC</b>           | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8617RC</b>           | 18 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8617RC</b>           | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8618RC</b>           | 18 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8618RC</b>           | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>Universal</b>        |                |                     |                                                                                                               |                  |                    |
| <b>8602U</b>            | 2 gallon       | 8975                | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical and chemotherapy containers                    |                  | 5 each             |
| <b>8605U</b>            | 5 quart        | 31307054            | Locking wire bracket, 5 quart, 2 and 3 gallon                                                                 | ✓                | 10 each            |
| <b>8605U</b>            | 5 quart        | 31307062            | Non-locking 5 quart, 2 and 3 gallon safety in-room                                                            |                  | 10 each            |
| <b>8605U</b>            | 5 quart        | 31353553            | Wall enclosure for 5 quart safety in-room sharps container                                                    | ✓                | 1 each             |
| <b>8607U</b>            | 8 gallon       | 8980FP              | Foot pedal cart for hinged lid container, 8 gallon                                                            |                  | 1 each             |
| <b>8607U</b>            | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 5 each             |
| <b>8607U</b>            | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8608U</b>            | 8 gallon       | 8981FP              | Foot pedal cart for slide lid container, 8 gallon                                                             |                  | 1 each             |
| <b>8608U</b>            | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 5 each             |
| <b>8608U</b>            | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8611U</b>            | 12 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8611U</b>            | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8612U</b>            | 12 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8612U</b>            | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8617U</b>            | 18 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8617U</b>            | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8618U</b>            | 18 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8618U</b>            | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |

# Hardware and container cross reference guide (continued)

Search by container order code to identify corresponding hardware options.

| Container order code  | Container size | Hardware order code | Hardware description                                                                                          | Locking Bracket? | Hardware case qty. |
|-----------------------|----------------|---------------------|---------------------------------------------------------------------------------------------------------------|------------------|--------------------|
| <b>Non-hazardous</b>  |                |                     |                                                                                                               |                  |                    |
| <b>8820</b>           | 2 gallon       | 8963                | Locking wire bracket for 2 gallon, multi-purpose, RCRA, universal, pharmaceutical and chemotherapy containers | ✓                | 1 each             |
| <b>8820</b>           | 2 gallon       | 8974                | Table top holder, phlebotomy, multi-purpose, RCRA, pharmaceutical and ChemoSafety Container, almond, 2 gallon |                  | 1 each             |
| <b>8820</b>           | 2 gallon       | 8975                | Wire wall bracket 2 gallon multi-purpose, RCRA, pharmaceutical and chemotherapy containers                    |                  | 1 each             |
| <b>8836SA</b>         | 3 gallon       | 8524C               | Non-locking bracket for 2 gallon and 3 gallon in-room                                                         |                  | 1 each             |
| <b>8836SA</b>         | 3 gallon       | 8528C               | Locking bracket for 2 and 3 gallon in-room                                                                    | ✓                | 1 each             |
| <b>8836SA</b>         | 3 gallon       | 8529                | Table top holder, Multi-Purpose SharpStar™, and In-Room container, almond, 3 gallon                           |                  | 1 each             |
| <b>8836SA</b>         | 3 gallon       | 88361H              | Wall enclosure for 3 gallon pharmaceutical container                                                          | ✓                | 1 each             |
| <b>8850</b>           | 8 gallon       | 8980FP              | Foot pedal cart for hinged lid container, 8 gallon                                                            |                  | 1 each             |
| <b>8850</b>           | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 5 each             |
| <b>8850</b>           | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8851</b>           | 8 gallon       | 8981FP              | Foot pedal cart for slide lid container, 8 gallon                                                             |                  | 1 each             |
| <b>8851</b>           | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 5 each             |
| <b>8851</b>           | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8860-</b>          | 12 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8860-</b>          | 12 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8870-</b>          | 18 gallon      | 8991FP              | Foot pedal cart for hinged lid container, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8870-</b>          | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>8871</b>           | 18 gallon      | 8938FP              | Foot pedal cart for slide lid container, 12 and 18 gallon                                                     |                  | 1 each             |
| <b>8871</b>           | 18 gallon      | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |
| <b>Non-infectious</b> |                |                     |                                                                                                               |                  |                    |
| <b>8781</b>           | 8 gallon       | 8980FP              | Foot pedal cart for hinged lid container, 8 gallon                                                            |                  | 1 each             |
| <b>8781</b>           | 8 gallon       | 8984-               | Locking wire wall bracket for 8 gallon containers                                                             | ✓                | 1 each             |
| <b>8781</b>           | 8 gallon       | 8992H               | Floor cart for large volume container, 8, 12 and 18 gallon                                                    |                  | 1 each             |

# Index

Find your item by order code and page number

|             |    |             |        |          |        |
|-------------|----|-------------|--------|----------|--------|
| 1180100555  | 27 | 11832010    | 16     | 31323333 | 41     |
| 1180100777  | 27 | 11832015    | 16     | 31353553 | 54, 55 |
| 1180125058  | 24 | 11832110    | 16     | 31353603 | 41     |
| 1180125158  | 24 | 11832115    | 16     | 8303SA   | 47     |
| 1180126038  | 24 | 11832210    | 16     | 8506SA   | 40     |
| 1180127012  | 24 | 11832215    | 16     | 8507MW   | 40     |
| 1180128012  | 24 | 11832310    | 16     | 8507SA   | 40     |
| 1180300555  | 27 | 11832315    | 16     | 851201   | 42     |
| 1180300777  | 27 | 11832510    | 16     | 85121    | 42     |
| 1180320100  | 21 | 11832558    | 16     | 851301   | 42     |
| 1180320112  | 21 | 1183500555  | 27     | 85131    | 42     |
| 1180321100  | 21 | 1183500777  | 27     | 8515-    | 56     |
| 1180321112  | 21 | 1183500888  | 27     | 85161H   | 54     |
| 1180322100  | 21 | 1186000444  | 27     | 85165H   | 54     |
| 1180322112  | 21 | 1186000555  | 27     | 8516HDL  | 54     |
| 1180323100  | 21 | 1186000555T | 27     | 8518X    | 53     |
| 1180323114E | 32 | 1186000777  | 27     | 8519C    | 53     |
| 1180325058  | 21 | 1186000777T | 27     | 85221    | 42     |
| 1180325100  | 21 | 1186060265  | 27     | 85221R   | 42     |
| 1180325114  | 21 | 1188100555  | 23     | 8524C    | 53     |
| 1180327114  | 21 | 1188100777  | 23     | 8527R    | 44     |
| 1180327114E | 32 | 1188128012  | 23     | 8528C    | 53     |
| 1180600555  | 27 | 1188528012  | 23     | 8529     | 56     |
| 1180600777  | 27 | 1188818100  | 18     | 85301H   | 54     |
| 11811022    | 5  | 1188818112  | 18     | 8533SA   | 40     |
| 11811022F   | 5  | 1188818113  | 18     | 8534SA   | 40     |
| 1181200555K | 27 | 1188819100  | 18     | 8536SA   | 40     |
| 1181200777K | 27 | 1188819112  | 18     | 8537MW   | 40     |
| 1181200777T | 27 | 1188820100  | 18     | 8537SA   | 40     |
| 11812558    | 16 | 1188820112  | 18     | 8539SA   | 56     |
| 11812705    | 16 | 1188822100  | 18     | 85401H-  | 54     |
| 11813005    | 16 | 1188822112  | 18     | 8540SA   | 40     |
| 1181620112  | 21 | 1188823100  | 18     | 8541SA   | 40     |
| 1181621100  | 21 | 1188825058  | 18     | 8550B    | 54, 55 |
| 1181621112  | 21 | 1188825100  | 18     | 8550LG   | 55     |
| 1181622112  | 21 | 1188825112  | 18     | 8554B    | 54, 55 |
| 1181810     | 16 | 1188826112  | 18     | 8555SA   | 55     |
| 1181815     | 16 | 1188827012  | 18     | 8556H-   | 54     |
| 1181910     | 16 | 1188827112  | 18     | 8556H    | 54     |
| 1181915     | 16 | 1188830340  | 18     | 8602RC   | 49, 52 |
| 1182000555  | 27 | 1522SA      | 45     | 8605RC   | 49     |
| 1182000777  | 27 | 1523SA      | 56     | 8607RC   | 49, 52 |
| 1182010     | 16 | 155922      | 35     | 8607U    | 49, 52 |
| 1182015     | 16 | 155923      | 35     | 8608RC   | 49, 52 |
| 1182110     | 16 | 31139747    | 53     | 8608U    | 49, 52 |
| 1182115     | 16 | 31139762    | 55     | 8611RC   | 49, 52 |
| 1182210     | 16 | 31140109    | 57     | 8611U    | 49, 52 |
| 1182215     | 16 | 31142222    | 43     | 8612RC   | 49, 52 |
| 1182310     | 16 | 31143665    | 45     | 8612U    | 49, 52 |
| 1182315     | 16 | 31143699    | 43     | 8617RC   | 49, 52 |
| 1182510     | 16 | 31158960    | 53     | 8617U    | 49, 52 |
| 1182515     | 16 | 31307005    | 54, 55 | 8618RC   | 49, 52 |
| 1182558     | 16 | 31307054    | 53     | 8618U    | 49, 52 |
| 1182705     | 16 | 31307062    | 53     | 8781     | 51     |
| 1182734     | 16 | 31314886    | 41     | 8820     | 50     |
| 1183005     | 16 | 31320032    | 55     | 88361H   | 54     |

# Index

Find your item by order code and page number

|            |        |            |    |            |    |
|------------|--------|------------|----|------------|----|
| 8836SA     | 50     | 8881250024 | 19 | 8881511110 | 17 |
| 8850       | 50     | 8881250032 | 19 | 8881511136 | 17 |
| 8851       | 50     | 8881250040 | 19 | 8881511144 | 17 |
| 8860-      | 50     | 8881250057 | 19 | 8881511201 | 17 |
| 8870-      | 50     | 8881250065 | 19 | 8881511235 | 17 |
| 8871       | 50, 52 | 8881250073 | 19 | 8881511310 | 17 |
| 8881103025 | 29     | 8881250081 | 19 | 8881511336 | 17 |
| 8881103066 | 29     | 8881250099 | 19 | 8881511344 | 17 |
| 8881106010 | 29     | 8881250107 | 19 | 8881512258 | 28 |
| 8881106028 | 29     | 8881250115 | 19 | 8881512852 | 27 |
| 8881112059 | 29     | 8881250123 | 19 | 8881512878 | 27 |
| 8881112083 | 29     | 8881250149 | 19 | 8881512977 | 28 |
| 8881114055 | 29     | 8881250172 | 19 | 8881513025 | 22 |
| 8881120029 | 29     | 8881250206 | 19 | 8881513033 | 22 |
| 8881120037 | 29     | 8881250214 | 19 | 8881513058 | 22 |
| 8881135084 | 29     | 8881250222 | 19 | 8881513132 | 22 |
| 8881160405 | 29     | 8881250230 | 19 | 8881513157 | 22 |
| 8881160629 | 29     | 8881250248 | 19 | 8881513207 | 28 |
| 8881200011 | 20     | 8881250255 | 19 | 8881513231 | 22 |
| 8881200029 | 20     | 8881250289 | 19 | 8881513256 | 22 |
| 8881200037 | 20     | 8881250305 | 19 | 8881513330 | 22 |
| 8881200045 | 20     | 8881250313 | 19 | 8881513512 | 22 |
| 8881200052 | 20     | 8881250321 | 19 | 8881513538 | 22 |
| 8881200060 | 20     | 8881250339 | 19 | 8881513744 | 22 |
| 8881200078 | 20     | 8881250362 | 19 | 8881513843 | 32 |
| 8881200136 | 20     | 8881250545 | 19 | 8881513850 | 32 |
| 8881200169 | 20     | 8881400033 | 30 | 8881513918 | 27 |
| 8881200185 | 20     | 8881400041 | 30 | 8881513934 | 27 |
| 8881200219 | 20     | 8881400058 | 30 | 8881516051 | 22 |
| 8881200318 | 20     | 8881400066 | 30 | 8881516135 | 22 |
| 8881200326 | 20     | 8881400074 | 30 | 8881516200 | 28 |
| 8881200342 | 20     | 8881400173 | 30 | 8881516911 | 27 |
| 8881200433 | 20     | 8881401031 | 31 | 8881516937 | 27 |
| 8881200441 | 20     | 8881401049 | 31 | 8881520178 | 28 |
| 8881200466 | 20     | 8881401056 | 31 | 8881520251 | 28 |
| 8881200508 | 20     | 8881401064 | 31 | 8881520657 | 27 |
| 8881200516 | 20     | 8881401072 | 31 | 8881520665 | 27 |
| 8881200573 | 20     | 8881401171 | 31 | 8881520673 | 27 |
| 8881200714 | 20     | 8881412012 | 32 | 8881535101 | 28 |
| 8881200755 | 20     | 8881471232 | 31 | 8881535762 | 27 |
| 8881200805 | 20     | 8881471273 | 31 | 8881535770 | 27 |
| 8881202017 | 7      | 8881500014 | 23 | 8881535788 | 27 |
| 8881202018 | 6      | 8881500105 | 24 | 8881535796 | 27 |
| 8881202033 | 7      | 8881500501 | 25 | 8881540101 | 8  |
| 8881202314 | 6      | 8881501160 | 24 | 8881540111 | 4  |
| 8881202322 | 6      | 8881501178 | 24 | 8881540133 | 4  |
| 8881202330 | 6      | 8881501210 | 23 | 8881540166 | 4  |
| 8881202348 | 6      | 8881501236 | 25 | 8881541034 | 8  |
| 8881202355 | 6      | 8881501368 | 24 | 8881541067 | 8  |
| 8881202363 | 6      | 8881501384 | 23 | 8881560125 | 27 |
| 8881202371 | 6      | 8881501400 | 27 | 8881560141 | 27 |
| 8881202389 | 6      | 8881501459 | 28 | 8881560182 | 27 |
| 8881202397 | 6      | 8881501608 | 24 | 8881560224 | 27 |
| 8881204005 | 7      | 8881501640 | 24 | 8881560232 | 28 |
| 8881250016 | 19     | 8881501962 | 25 | 8881560265 | 27 |

# Index

Find your item by order code and page number

|             |    |          |        |          |    |
|-------------|----|----------|--------|----------|----|
| 8881600004  | 23 | 8900MW   | 45     | CP3001O- | 33 |
| 8881600145  | 23 | 8900SA   | 45     | CP3001R- | 33 |
| 8881600350  | 23 | 8910-    | 56     | CP3001T- | 33 |
| 8881600700  | 23 | 8912SA   | 53     | CP3001Y- | 33 |
| 8881600800  | 23 | 8926     | 53     | CP3002A- | 33 |
| 8881601101  | 23 | 8930SA   | 46, 52 | CP3002B- | 33 |
| 8881601358  | 23 | 8931     | 48     | CP3002G- | 33 |
| 8881601600  | 23 | 8931PG2  | 48, 52 | CP3002O- | 33 |
| 8881676012  | 53 | 8932     | 46     | CP3002R- | 33 |
| 8881676236  | 44 | 8932PG2  | 46, 52 | CP3002Y- | 33 |
| 8881676285  | 44 | 8933     | 46     | CP3003A- | 33 |
| 8881676434  | 44 | 8934     | 48     | CP3003B- | 33 |
| 8881676624  | 54 | 8934PG2  | 48, 52 | CP3003G- | 33 |
| 8881682010  | 29 | 8935     | 46     | CP3003R- | 33 |
| 8881682085  | 29 | 8936PG2  | 46, 52 | CP3003Y  | 33 |
| 8881682101  | 29 | 8936SA   | 46     | CP3004A  | 34 |
| 8881682119  | 29 | 8938     | 46     | CP3004B- | 34 |
| 8881811310  | 12 | 8938FP   | 57     | CP3004G- | 34 |
| 8881811510  | 12 | 8939     | 48     | CP3004O  | 34 |
| 8881811558  | 12 | 8939PG2  | 48, 52 | CP3004R- | 34 |
| 8881833010  | 12 | 8950SA   | 44     | CP3005A- | 33 |
| 8881833015  | 12 | 8963     | 53     | CP3005R- | 33 |
| 8881833110  | 12 | 89651    | 43     | CP3009A- | 33 |
| 8881833115  | 12 | 89671    | 43     | CP3009B- | 33 |
| 8881833210  | 12 | 8970     | 44     | CP3009G- | 33 |
| 8881833215  | 12 | 8974     | 56     | CP3009R- | 33 |
| 8881833310  | 12 | 8975     | 53     | CP3011A- | 34 |
| 8881833510  | 12 | 8979MW   | 44     | CP3011B- | 34 |
| 8881833558  | 12 | 8980-    | 46     | CP3011G- | 34 |
| 8881833810  | 12 | 8980FP   | 57     | CP3011R- | 34 |
| 8881850010  | 12 | 8980S    | 46     | CP3011Y- | 34 |
| 8881850015  | 12 | 8981FP   | 57     | CP3013A- | 34 |
| 8881850110  | 12 | 8982     | 48     | CP3013B- | 34 |
| 8881850115  | 12 | 8984-    | 53     | CP3013G- | 34 |
| 8881850210  | 12 | 8985     | 48     | CP3013R- | 34 |
| 8881850215  | 12 | 8985PG2  | 48, 52 | CP3020B- | 34 |
| 8881850310  | 12 | 8985S    | 48     | CP3020R- | 34 |
| 8881850358  | 12 | 8985SPG2 | 48, 52 | CP3020Y- | 34 |
| 8881850510  | 12 | 8989     | 48     | CP3030R- | 34 |
| 8881850515  | 12 | 8989PG2  | 48, 52 | CP3111Y- | 34 |
| 8881850558  | 12 | 8990SA   | 44     | CP3720A- | 34 |
| 8881850810  | 12 | 8991     | 46     | CP3730A- | 34 |
| 8881850815  | 12 | 8991FP   | 57     |          |    |
| 8881850910  | 12 | 8992H    | 57     |          |    |
| 8881850915  | 12 | 8993H    | 57     |          |    |
| 8881882712  | 14 | 8997     | 46     |          |    |
| 8881882712T | 14 | 8997PG2  | 46, 52 |          |    |
| 8881882812  | 14 | 8997S    | 46     |          |    |
| 8881882812T | 14 | 8997SPG2 | 46, 52 |          |    |
| 8881892910  | 14 | 8998     | 46     |          |    |
| 8881892930  | 14 | 8998PG2  | 46, 52 |          |    |
| 8881892950  | 14 | 8998S    | 46     |          |    |
| 8881893010  | 14 | 8998SPG2 | 46, 52 |          |    |
| 8881893030  | 14 | CP3001A- | 33     |          |    |
| 8881893050  | 14 | CP3001B- | 33     |          |    |
| 8888115006  | 35 | CP3001G- | 33     |          |    |

For more information or to order Monoject™ Sharps Safety Products, contact your sales representative, call **800.964.5227** or visit **[cardinalhealth.com/monoject](https://www.cardinalhealth.com/monoject)**



**[cardinalhealth.com/monoject](https://www.cardinalhealth.com/monoject)**

© 2021 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, MONOJECT, MAGELLAN, KENDALL, SMARTIP, SHARPSTAR and CHEMOPLUS are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other trademarks are the property of their respective owners. Lit. No. 2GM21-1319108 (09/2021)

Cardinal Health  
7000 Cardinal Place  
Dublin, Ohio 43017