



# AquaPur<sup>TM</sup>ST

## Technical Data Sheet

### - BRIEF DESCRIPTION -

*AquaPur<sup>TM</sup>ST is Purified Water that has been sterilized and packaged in a double bag for use in the cleanroom areas within pharmaceutical, biotechnology, and medical device facilities or anywhere the use of sterile, USP Purified water is required.*

*AquaPur<sup>TM</sup>ST is prepared and assayed to specifications for sterile USP Purified Water at time of manufacture. AquaPur<sup>TM</sup>ST contains no preservatives. It is gamma-irradiated to a sterility assurance level (SAL) of 10<sup>-6</sup>. Each lot is tested for Endotoxin via USP LAL testing and sterility via USP 14-day sterility test. A copy of lot specific documentation, which details conductivity, TOC, irradiation, endotoxin, sterility, and water quality limits for USP Purified Water, is attached to every case of AquaPur<sup>TM</sup>ST.*

### - PRECAUTIONARY STATEMENT -

Not to be used for Parenteral Administration. Not for use in sterile irrigation, injection, infusion, or as a potable product. AquaPur<sup>TM</sup>ST is intended for manufacturing and general laboratory use.

### - STORAGE -

Store in cool area. Sterility can not be assured if the plastic bags are not sealed. AquaPur<sup>TM</sup>ST has an expiration date of two years and the date is listed on the bottle.

### - PRODUCT SPECIFICATIONS -

AquaPur<sup>TM</sup> Sterile, USP grade Purified Water, is filtered through a 0.2 µm filter. AquaPur<sup>TM</sup>ST is tested to meet USP specifications for Purified Water.

USP Water Conductivity USP <645>	≤ 1.3 umhos/cm
USP Total Organic Carbon USP <643>	< 0.5 PPM

### - PACKAGING -

1 Gallon (3.8 Liters) / 4 per case

### - CATALOG NUMBER -

<u>Manufactured for:</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	1 gallon	8904

---

**For Additional Technical Information call  
Decon Customer Service  
(800) 332-6647, or  
See our website at [www.deconlabs.com](http://www.deconlabs.com)**

---

©Decon Labs, Inc.

Revised 3/26/2018