



# SUCTION CATHETERS MINI-KITS

Dynarex Resp-O<sub>2</sub>™ Suction Catheter Mini-Kits are designed to safely and effectively remove secretions from the respiratory tract of tracheostomy patients. The soft, flexible catheter with beveled tip is gentle on delicate mucosal tissue yet firm enough to allow for catheter rotation. The catheter includes a vented cap control valve for enhanced suction control and features two eyelet openings for increased suctioning effectiveness with minimal patient trauma. Each Kit also features a mini-tray to hold irrigation solutions and a pair of powder-free vinyl exam gloves.

- Convenient Mini-Kit includes a pair of vinyl exam gloves and mini-tray to hold irrigation fluids
- Soft, pliable catheter with beveled tip provides optimal patient comfort
- Vented cap with control valve allows greater procedural control
- Catheter tip is designed with 2 staggered eyelet openings for maximum suctioning efficiency
- Catheter features easy-to-read depth-number markings
- Sterile, single-patient use
- Not made with natural rubber latex or DEHP



## CONTENTS:

Graduated Catheter with Control Valve  
Pair of Powder-Free Vinyl Exam Gloves  
Mini-Tray



8 Fr

Pediatric

REORDER NO. 32080

PACKAGE QUANTITY  
100/Case

10 Fr

Adult/Pediatric

REORDER NO. 32081

PACKAGE QUANTITY  
100/Case

12 Fr

Adult

REORDER NO. 32082

PACKAGE QUANTITY  
100/Case

14 Fr

Adult

REORDER NO. 32083

PACKAGE QUANTITY  
100/Case

16 Fr

Adult

REORDER NO. 32084

PACKAGE QUANTITY  
100/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.

©2021 Dynarex Corporation



R210823