



DISPOSABLE SUCTION CANISTERS

Dynarex Resp-O₂[™] Disposable Suction Canisters are designed to safely and effectively collect and contain aspirated secretions from patients in medical settings. The Canisters feature an ultrasonically welded shut-off valve to help prevent overflow and cross-contamination. These versatile Canisters are adaptable for use with floor, cabinet, or wall-mounted units. Each disposable, single-use Canister comes with a convenient lid equipped with a pour spout, vacuum, and patient port.

- Constructed from high-quality, durable polymer materials
- Port caps are engineered to create a snug seal to prevent vacuum loss
- Transparent canister features easy-to-read gradations for accurate monitoring
- Not made with natural rubber latex or DEHP



Disposable Suction Canister
800 cc

FLOAT VALVE SHUTOFF

REORDER NO. 32200

PACKAGE QUANTITY
24/Case

Disposable Suction Canister
800 cc

FLOAT VALVE SHUTOFF
3/16" x 6' AND 3/16" x 18" TUBING
WITH MALE CONNECTOR

REORDER NO. 32201

PACKAGE QUANTITY
12/Case

Disposable Suction Canister
1200 cc

FLOAT VALVE SHUTOFF

REORDER NO. 32220

PACKAGE QUANTITY
40/Case

Disposable Suction Canister
1200 cc

FLOAT VALVE SHUTOFF
3/16" x 6' AND 3/16" x 18" TUBING
WITH MALE CONNECTOR

REORDER NO. 32221

PACKAGE QUANTITY
30/Case

Disposable EMS Suction Canister
1200 cc

FLOAT VALVE SHUTOFF

REORDER NO. 32224

PACKAGE QUANTITY
40/Case



WE CARE LIKE FAMILY[™]

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



©2021 Dynarex Corporation

R210823