



YANKAUER SUCTION HANDLE

Dynarex Resp-O₂™ Yankauer Suction Handles are designed to safely and effectively remove particulates and fluids from patients' airways. The curved, one-piece Handle is constructed from transparent materials for clear visualization, and features a slip-resistant, ribbed design to ensure a secure connection to suction tubing.

- Designed to clear airways to help improve ventilation and oxygenation
- Constructed from shatter-resistant plastic materials
- Available in vented, non-vented, bulb-tip, and straight-tip designs
- Sterile, single-use, and disposable for optimal hygiene
- Not made with natural rubber latex or DEHP



**Bulb Tip
Vented**

REORDER NO. 32120

ALL AGES
.....
PACKAGE QUANTITY
50/Case

**Bulb Tip
Non-Vented**

REORDER NO. 32121

ALL AGES
.....
PACKAGE QUANTITY
50/Case

**Straight Tip
Vented**

REORDER NO. 32122

ALL AGES
.....
PACKAGE QUANTITY
50/Case

**Straight Tip
Non-Vented**

REORDER NO. 32123

ALL AGES
.....
PACKAGE QUANTITY
50/Case

**Bulb Tip • Non-Vented
with Pre-Attached
3/8" × 6' Suction Tubing**

REORDER NO. 32125

ALL AGES
.....
PACKAGE QUANTITY
50/Case

**Bulb Tip • Non-Vented
with Pre-Attached
1/4" × 6' Suction Tubing**

REORDER NO. 32126

ALL AGES
.....
PACKAGE QUANTITY
50/Case

**Bulb Tip • Non-Vented
with Pre-Attached
3/8" × 10' Suction Tubing**

REORDER NO. 32127

ALL AGES
.....
PACKAGE QUANTITY
20/Case

**Bulb Tip • Non-Vented
with Pre-Attached
1/4" × 10' Suction Tubing**

REORDER NO. 32128

ALL AGES
.....
PACKAGE QUANTITY
20/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



©2021 Dynarex Corporation

R210823