



SINGLE SUCTION CATHETERS

Graduated with Control Valve • Coiled Packed



Dynarex Resp-O₂™ Single Suction Catheters are designed to safely and effectively remove secretions from the respiratory tract of tracheostomy patients. The soft, flexible Catheter with beveled tip is gentle on delicate mucosal tissue yet firm enough to allow for catheter rotation. The Catheter includes a vented cap port for enhanced suction control and two staggered eyelet openings to increase suctioning effectiveness and reduce patient trauma.

- Soft, pliable catheter with beveled tip provides optimal patient comfort
- Vented cap with control valve allows for greater procedural control
- Catheter tip is designed with 2 staggered eyelet openings for maximum suctioning efficiency
- Catheter features easy-to-read depth-number markings
- Coil packed for compact storage
- Sterile, single-patient use
- Not made with natural rubber latex or DEHP



6 Fr

Pediatric

REORDER NO. 32000

PACKAGE QUANTITY
50/Case

8 Fr

Pediatric

REORDER NO. 32001

PACKAGE QUANTITY
50/Case

10 Fr

Adult/Pediatric

REORDER NO. 32002

PACKAGE QUANTITY
50/Case

12 Fr

Adult

REORDER NO. 32003

PACKAGE QUANTITY
50/Case

14 Fr

Adult

REORDER NO. 32004

PACKAGE QUANTITY
50/Case

16 Fr

Adult

REORDER NO. 32005

PACKAGE QUANTITY
50/Case

18 Fr

Adult

REORDER NO. 32006

PACKAGE QUANTITY
50/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



©2021 Dynarex Corporation

R210823