




# SUCTION CATHETERS KITS

Dynarex Resp-O<sub>2</sub>™ Single Suction Catheter Kits are designed to safely and effectively remove secretions from the respiratory tract of tracheostomy patients. The soft, flexible catheter with beveled tip is gentle on delicate mucosal tissue yet firm enough to allow for catheter rotation. The catheter features a vented cap control valve for enhanced suction control and two staggered eyelet openings to increase suctioning effectiveness and reduce patient trauma. Each Kit includes a pop-up basin for hold irrigation fluids and a pair of powder-free vinyl exam gloves.

- Kit includes pop-up basin and exam gloves for added convenience
- Soft, pliable catheter with beveled tip provides optimal patient comfort
- Vented cap with control valve allows for greater procedural control
- Catheter tip is designed with 2 staggered eyelet openings for maximum suctioning efficiency
- Catheter features easy-to-read depth-number markings
- Sterile, single-patient use
- Not made with natural rubber latex or DEHP



**CONTENTS:**  
 Graduated Catheter with Control Valve  
 Pair of Powder-Free Vinyl Exam Gloves  
 Pop-Up Basin



**6 Fr**

Pediatric

**REORDER NO. 32050**

PACKAGE QUANTITY  
50/Case

**8 Fr**

Pediatric

**REORDER NO. 32051**

PACKAGE QUANTITY  
50/Case

**10 Fr**

Adult/Pediatric

**REORDER NO. 32052**

PACKAGE QUANTITY  
50/Case

**12 Fr**

Adult

**REORDER NO. 32053**

PACKAGE QUANTITY  
50/Case

**14 Fr**

Adult

**REORDER NO. 32054**

PACKAGE QUANTITY  
50/Case

**16 Fr**

Adult

**REORDER NO. 32055**

PACKAGE QUANTITY  
50/Case

**18 Fr**

Adult

**REORDER NO. 32056**

PACKAGE QUANTITY  
50/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



©2021 Dynarex Corporation

R210823