

# SUCTION UNITS

## Vacuum Levels up to 560 mmHg

Dynarex Resp-O<sub>2</sub>™ Suction Units are premium-grade medical devices designed to safely and effectively remove excess secretions and fluids from the airways of respiratory patients. These powerful suction units are easily portable and feature high-quality vacuum regulators and anti-vibration vacuum gauge settings to ensure quiet and efficient operation.

- Designed to provide vacuum levels up to 560 mmHg
- Fiber filter guards against bacteria and contaminants
- Clearly marked gauge allows for easy pressure monitoring
- Oil-less lubrication provides consistent functionality
- Includes an 800 cc disposable suction canister, tubing and bacterial filter



**32300**



**32301**



Suction Cup Legs

**REORDER NO. 32300**

PACKAGE QUANTITY  
1/Case

Plastic Base

**REORDER NO. 32301**

PACKAGE QUANTITY  
1/Case

