


OXYGEN MASKS Medium-Concentration

Dynarex Resp-O₂[™] Medium-Concentration Masks are designed to effectively deliver oxygen at a concentration rate of 40% to 50%. These masks come with light green oxygen supply tubing, designed to help improve tube visibility.

- Transparent vinyl mask features smooth edges
- Elastic strap allows for easy adjustment
- Adult and pediatric masks feature an adjustable nose clip
- Includes 7' of crush-resistant, star lumen oxygen supply tubing
- Not made with natural rubber latex or DEHP







WE CARE LIKE FAMILY