

OXYGEN MASKS Non-Rebreather & 3-in-1


Dynarex Resp- O_2^{∞} Non-Rebreather Masks are designed to effectively deliver oxygen at high concentration rates in emergency situations. The masks feature a safety valve and a reservoir bag to ensure sufficient oxygen supply. Dynarex also offers a 3-in-1 mask that can be used as a mediumconcentration, high-concentration, or non-rebreathing mask.

- Transparent vinyl mask features smooth edges
- Mask features an adjustable nose clip and elastic strap
- Includes a vinyl reservoir bag and 7' of crush-resistant, star lumen oxygen supply tubing
- Not made with natural rubber latex or DEHP



L 1.888 DYNAREX

WE CARE LIKE FAMILY



3-in-1

ADULT

.....

50/Case

Elongated Mask

1000 ML RESERVOIR BAG

STANDARD CONNECTOR

REORDER NO. 33055

R210423

PACKAGE QUANTITY