

OXYGEN SUPPLY TUBING

Dynarex Resp-O₂™ Oxygen Supply Tubing is designed to extend a patient's reach when accessing a concentrated oxygen source for respiratory therapy. The tubing features a light green color, designed to help improve tube visibility. It comes with a standard connector on one end and either a standard or universal connector on the other.

- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP





STANDARD CONNECTOR

REORDER NO. 33200

PACKAGE QUANTITY 50/Case

UNIVERSAL CONNECTOR

REORDER NO. 33201

PACKAGE QUANTITY 50/Case

14' (4.3 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33202

PACKAGE QUANTITY 50/Case

14' (4.3 m) Crush **Resistant Tubing**

UNIVERSAL CONNECTOR

REORDER NO. 33203

PACKAGE QUANTITY 50/Case

21' (6.4 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33204

PACKAGE QUANTITY 25/Case

25' (7.6 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33205

PACKAGE QUANTITY 20/Case

30' (9.1 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33206

PACKAGE QUANTITY 20/Case

40' (12.1 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33207

PACKAGE QUANTITY 15/Case

50' (15.2 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33208

PACKAGE QUANTITY 15/Case

WE CARE LIKE FAMILY

1.888 DYNAREX



DYNAREX.COM



sales representative or call to find your nearest distributor ®2021 Dynarex Corporation

Contact your local







R210426