

TRACHEOSTOMY MASK

with Swivel Tubing Connector

Dynarex Resp- $O_2^{\text{\tiny TM}}$ Tracheostomy Masks are designed to help administer aerosol medications and oxygen to tracheostomy patients during respiratory therapy. The collar-style masks feature swivel connectors that allow for access from either side of the patient.

- Collar-style mask features smooth edges
- Non-constricting, adjustable neckband features a snap-closure for easy removal

R210427

- Swivel connector accepts 22 mm aerosol corrugated tubing
- Not made with natural rubber latex or DEHP

