

# **SMALL VOLUME NEBULIZER 6 CC CUPS**

Dynarex Resp-O<sub>2</sub>™ Small Volume Nebulizer Cups are an integral part of the nebulizing process, in which liquid medication is transformed into an inhaled mist to aid in respiratory therapy treatments. The cups feature an anti-spill design that allows patients to receive treatments upright or at an angle.


- Effectively delivers medication at angles up to 90°
- Features a threaded cap to help prevent leaks
- Disposable cup capacity: 6 mL
- Suitable for hand-held and table-top nebulizers
- Not made with natural rubber latex or DEHP



T-Piece












Mouthpiece





Small Volume Nebulizer 6 cc Cup





®2021 Dynarex Corporation







# **SMALL VOLUME NEBULIZER 6 CC CUPS**

Dynarex Resp-O<sub>2</sub>™ Small Volume Nebulizer Cups are an integral part of the nebulizing process, in which liquid medication is transformed into an inhaled mist to aid in respiratory therapy treatments. The cups feature an anti-spill design that allows patients to receive treatments upright or at an angle.

- Effectively delivers medication at angles up to 90°
- Features a threaded cap to help prevent leaks
- Disposable cup capacity: 6 mL
- Suitable for hand-held and table-top nebulizers
- Not made with natural rubber latex or DEHP

# REORDER NO. 34100

......

.....

# ACCESSORIES

6 cc cup only

## ALL AGES

PACKAGE QUANTITY

# REORDER NO. 34101

•••••

## **ACCESSORIES**

7' (2.1 m) Tubing Standard Connector

## ALL AGES

PACKAGE QUANTITY

50/Case 50/Case

# REORDER NO. 34102

## ACCESSORIES

7' (2.1 m) Tubing **Standard Connector** Spring Loaded T (22 mm OD  $\times$ 22 mm OD)

.....

.....

## ALL AGES

PACKAGE QUANTITY

50/Case

# REORDER NO. 34103

## **ACCESSORIES**

7' (2.1 m) Tubing **Standard Connector** Spring Loaded T (22 mm ID  $\times$ 22 mm OD)

## ALL AGES

# PACKAGE QUANTITY

50/Case

.....

.....

# REORDER NO. 34104

.....

.....

## ACCESSORIES

7' (2.1 m) Tubing **Standard Connector** T-Piece Mouthpiece

#### ALL AGES

# PACKAGE QUANTITY

50/Case

# REORDER NO. 34105

## **ACCESSORIES**

7' (2.1 m) Tubing **Standard Connector** T-Piece Mouthpiece 6" Flex Tube

#### ALL AGES

# PACKAGE QUANTITY

50/Case

# REORDER NO. 34106

.....

## **ACCESSORIES**

7' (2.1 m) Tubing **Universal Connector** T-Piece Mouthpiece 6" Flex Tube

#### ALL AGES ..........

PACKAGE QUANTITY

50/Case

# REORDER NO. 34107

## **ACCESSORIES**

7' (2.1 m) Tubing Standard Connector 45 Degree Pacifier Elbow

•••••

# INFANT/TODDLER

••••• PACKAGE QUANTITY 12/Case

# REORDER NO. 34120

## **ACCESSORIES**

Aerosol Elongated Mask 7' (2.1 m) Tubing Standard Connector

•••••

#### **ADULT** .........

PACKAGE QUANTITY 50/Case

# REORDER NO. 34121

# **ACCESSORIES**

Aerosol Elongated Mask 7' (2.1 m) Tubing **Universal Connector** 

.....

#### ADULT .......

PACKAGE QUANTITY 50/Case

# REORDER NO. 34122

## **ACCESSORIES**

Aerosol Elongated Mask 7' (2.1 m) Tubing Standard Connector

•••••

#### PEDIATRIC ..........

PACKAGE QUANTITY 50/Case

# REORDER NO. 34123

## ACCESSORIES

Aerosol Elongated Mask 7' (2.1 m) Tubing **Universal Connector** 

•••••

#### PEDIATRIC ..........

PACKAGE QUANTITY

50/Case



1.888 DYNAREX



DYNAREX.COM



Contact your local sales representative or call to find your nearest distributor.



