


## **SPRING LOADED T-ADAPTORS**

Dynarex Resp-O2<sup>™</sup> Spring Loaded T-Adaptors are designed to be used in conjunction with a small volume nebulizer (SVN) to help deliver aerosol medication to respiratory vent patients. The T-Adaptor features a spring-loaded valve that opens and closes the nebulizer circuit during insertion and removal.

- Eliminates the need to break the ventilator circuit between treatments
- Ventless design prevents medication leakage
- Accepts most standard aerosol tubing/ventilator circuit
- Not made with natural rubber latex or DEHP



Contact your local sales representative or call to find your nearest distributor.

®2021 Dynarex Corporation



WE CARE LIKE FAMILY