



PULMONARY FUNCTION MEASUREMENT

Spirometers

Dynarex Resp-O₂™ Incentive Spirometers are hand-held devices designed to help improve lung function and minimize the risk of post-operative respiratory complications. The compact, durable spirometer helps motivate patients to strengthen their breathing muscles, particularly when recovering from prolonged anesthesia, surgery, or extended bed rest. The device features clear, easy-to-read markings to provide visual feedback to help patients monitor their progress.

- Adjustable indicator provides target guidelines to indicate sufficient lung exercise
- Flexible mouthpiece allows for easy adjustment and optimal comfort
- Features dual-sided, clearly marked calibrations for effortless visualization
- Compact design provides enhanced efficiency and ease of handling
- Not made with natural rubber latex or DEHP



Volumetric Incentive Spirometer

2500 ML

REORDER NO. 36100

ALL AGES

PACKAGE QUANTITY
12/Case



Volumetric Incentive Spirometer

5000 ML

REORDER NO. 36101

ALL AGES

PACKAGE QUANTITY
12/Case



Incentive Spirometer

TRI-FLOW CHAMBER

REORDER NO. 36102

ALL AGES

PACKAGE QUANTITY
12/Case

Peak Flow Meters

Dynarex Resp-O₂™ Peak Flow Meters are hand-held monitoring devices that measure peak expiratory flow (PEF). This lightweight, durable Meter is designed to help patients track daily changes in their breathing patterns. The Dynarex Resp-O₂™ Peak Flow Meter is available in Adult size (standard range 50-800 L/min) and Pediatric size (low range 50-300 L/min).

- Color-coded indicators help to track daily progress
- Easy-grip handle for steady handling
- Not made with natural rubber latex or DEHP
- Includes peak flow rate diary to track daily progress
- Disposable mouthpieces are available for added protection against bacteria and cross-contamination



Peak Flow Meter

LOW RANGE 50-300 L/MIN

REORDER NO. 36110

PEDIATRIC

PACKAGE QUANTITY
10/Case



Peak Flow Meter

FULL RANGE 50-800 L/MIN

REORDER NO. 36111

ADULT

PACKAGE QUANTITY
10/Case



Durable Peak Flow Mouthpiece

WITH ONE-WAY VALVE

REORDER NO. 36116

ADULT

PACKAGE QUANTITY
100/Case



Cardboard Peak Flow Mouthpiece

DISPOSABLE

REORDER NO. 36117

PEDIATRIC

PACKAGE QUANTITY
100/Case



Cardboard Peak Flow Mouthpiece

DISPOSABLE

REORDER NO. 36118

ADULT

PACKAGE QUANTITY
100/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.

©2021 Dynarex Corporation



R210823