



NASAL OXYGEN CANNULAS Cushion Tip

Dynarex Resp-O₂™ Nasal Oxygen Cannulas with Cushion Tips are designed to efficiently deliver supplemental oxygen to respiratory patients. The cannulas come with either Standard or Universal connectors and feature light green oxygen supply tubing, designed to help improve tube visibility.

- Lightweight cushioned nasal prongs provide added comfort
- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP



7' (2.1 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33120

NEONATAL

PACKAGE QUANTITY 50/Case

7' (2.1 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33121

INFANT

PACKAGE QUANTITY 50/Case

7' (2.1 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33122

PEDIATRIC

PACKAGE QUANTITY 50/Case

4' (1.2 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33123

ADULT

PACKAGE QUANTITY 50/Case

4' (1.2 m) Crush **Resistant Tubing**

UNIVERSAL CONNECTOR

REORDER NO. 33124

ADULT

PACKAGE QUANTITY 50/Case

7' (2.1 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33125

PACKAGE QUANTITY 50/Case

14' (4.3 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33126

PACKAGE QUANTITY 50/Case

21' (6.4 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33127

ADULT

PACKAGE QUANTITY 25/Case

25' (7.6 m) Crush **Resistant Tubing**

STANDARD CONNECTOR

REORDER NO. 33128

ADULT

PACKAGE QUANTITY 25/Case

50' (15.2 m) Crush **Resistant Tubing**


UNIVERSAL CONNECTOR

REORDER NO. 33129

ADULT

PACKAGE QUANTITY 15/Case







DYNAREX.COM



Contact your local sales representative or call to find your nearest distributor

®2021 Dynarex Corporation







