

# OXYGEN THERAPY ACCESSORIES

Dynarex Resp-O<sub>2</sub>™ Oxygen Therapy Accessories are designed to improve the efficacy, convenience, and overall comfort of respiratory equipment. They are not made with natural rubber latex or DEHP.

### Nasal Oxygen Cannula Ear Foam Covers

- Lightweight foam construction
- Suitable for adult and pediatric patients

### Inline Water Traps

- Traps and collects excess moisture from oxygen tubing
- Keeps liquid away from patients' nasal passages
- Tapered fittings ensure a secure connection to the tubing
- Transparent body allows for clear monitoring of fluid

### Oxygen Tubing Adaptors

- Designed to extend adaptability of oxygen tubing

### Portable Oxygen Liter Meter

- Designed to measure gas-flow from oxygen source
- Convenient pen-style design
- Not to be used as a calibration device



Nasal Oxygen Cannula Deluxe Ear Foam Covers

**REORDER NO. 33300**

**PACKAGE QUANTITY**  
50 pairs/Case



Nasal Oxygen Cannula Advantage Ear Foam Covers

**REORDER NO. 33301**

**PACKAGE QUANTITY**  
50 pairs/Case



Oxygen Tubing Inline Water Trap

**REORDER NO. 33302**

**PACKAGE QUANTITY**  
10/Case



Oxygen Tubing 22 mm I.D. x 5-7 mm Tapered Nipple Adaptor

**REORDER NO. 33303**

**PACKAGE QUANTITY**  
50/Case



Oxygen Tubing 5-7 mm End Adaptor

**REORDER NO. 33304**

**PACKAGE QUANTITY**  
50/Case



Oxygen Tubing Swivel Male/Male Adaptor

**REORDER NO. 33305**

**PACKAGE QUANTITY**  
50/Case



Oxygen Tubing Y Adaptor

**REORDER NO. 33306**

**PACKAGE QUANTITY**  
50/Case



Oxygen Tubing Nut Nipple Adaptor

**REORDER NO. 33307**

**PACKAGE QUANTITY**  
50/Case



Oxygen Tubing Nut Nipple Swivel Adaptor

**REORDER NO. 33308**

**PACKAGE QUANTITY**  
50/Case



Portable Oxygen Liter Meter

**REORDER NO. 33309**

**PACKAGE QUANTITY**  
15/Case

