



# TRACHEOSTOMY CARE KITS



Dynarex Resp-O<sub>2</sub>™ Tracheostomy Care Kits are designed to provide the essential components needed for the preparation, cleaning, and care of tracheostomy tubes. These time-saving Kits are disposable, portable, and compact for easy use and storage. The components in these fully sealed, single-use Kits are sterile and come packed in sequence of use for maximum efficiency.

- Fully sealed and sterile package and contents
- Components are packed in order of procedural use
- Disposable, single-patient use
- Not made with natural rubber latex or DEHP

REORDER NUMBER	35000	35001	35002	35003	35004	35005	35006	35008	35009
TRAY TYPE	ONE COMPARTMENT	THREE COMPARTMENT	TWO COMPARTMENT	TWO COMPARTMENT	TWO COMPARTMENT	THREE COMPARTMENT SALINE	THREE COMPARTMENT SALINE	ONE COMPARTMENT	MINI
REMOVABLE PLASTIC BASIN				1	1				
PAIR OF VINYL STRETCH GLOVES, M	1	1	1	1			1		
PAIR OF NITRILE GLOVES, M					1	1			
3-PLY WATERPROOF DRAPE, 45 x 34 CM	1	1	1	1	1	1	1	1	
NON-WOVEN GAUZE, 4" x 4"	4	4	4	4	4	4	4	4	2
NON-WOVEN DRESSING, 4" x 4"	1	1	1	1	1	1	1	1	2
COTTON TIP APPLICATOR, 5"	2	2	2	2	2	2	2	2	2
TWILL TAPE, 34"	1	1	1	1	1	1	1	1	
PIPE CLEANER, 5.5"	2	2	2	2	2	2	2	2	2
TRACHEOSTOMY BRUSH	1	1	1	1	1	1	1	1	
14 FR. SUCTION CATHETER			1						
STERILE SALINE 0.9%, 100 ML						1			
APLICARE HYDROGEN PEROXIDE 1 OZ. PACK						2	2		1
PACKAGING QUANTITY	20/CASE	20/CASE	20/CASE	20/CASE	20/CASE	20/CASE	20/CASE	20/CASE	40/CASE



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



©2021 Dynarex Corporation

R210823

# TRACHEOSTOMY CARE KIT COMPONENTS



One-Compartment Tray



Mini-Tray



Pair of Vinyl Stretch Gloves Size M



Pair of Nitrile Gloves Size M



3-Ply Waterproof Drape Size 45 x 34 cm



Two-Compartment Tray



Three-Compartment Tray



Sterile Saline 0.9%, 100 mL



14 Fr. Suction Catheter



Non-Woven Gauze Size 4" x 4"



Three-Compartment Saline Tray



Removable Plastic Basin



Aplicare Hydrogen Peroxide 1 oz. Pack



Cotton Tip Applicator Size 5"



Tracheostomy Brush



Pipe Cleaner Size 5.5"



Tracheostomy Dressing, Non-Woven Size 4" x 4"



Twill Tape Size 34"

