



SURGICAL GOWNS

- Fluid-resistant, sterile gowns designed to protect patients and medical providers
- Offers ANSI/AAMI PB70, Level 3 protection
- Non-woven spunbond materials help protect against microorganisms, body fluids, and particulates
- Constructed from breathable, low-linting material for optimal performance
- Elastic cuffs offer improved protection and easy donning
- Integrated belt ties at the back for secure closure
- Hook-and-loop closures at the neck provide added security
- Contains gown, hand towels, and wrap
- Included wrap helps provide sterile field
- Not made with natural rubber latex

8192	Medium	1 / pouch	20 pouches / cs
8193	Large	1 / pouch	20 pouches / cs
8194	X-Large	1 / pouch	20 pouches / cs
8195	2X-Large	1 / pouch	20 pouches / cs



| 1.888 DYNAREX | DYNAREX.COM |

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |