


# OPTFLOW® SUBCUTANEOUS ADMINISTRATION SETS

Integrated Medical Systems, Inc. is proud to carry EMEDs Soft-Glide® technology, which has been incorporated in a new family of SCIG needle sets branded under OPTflow®. Click on the item numbers to view in MyIMS.Cloud.

OPTflow® provides a solution for the increased market demands for faster flow rates and low penetration force to provide optimal comfort and safety to patients and clinicians. With comparable flow rates to 24 gauge SAF-Q® needle sets, OPTflow® needle sets provide optimized comfort. Patented easy-to-use wing design has a stabilizer lip for wing alignment during insertion.



## Features

- Soft-Glide® coating technology to minimize needle discomfort and anxiety
- The widest selection of single and multi-needle site infusion sets
- Needle lengths: 4, 6, 9, 12, and 14mm to accommodate all age groups and skin types
- Stabilizer lip for wing alignment to facilitate 90 degree insertion
- Proprietary ultra-flexible polyethylene tubing
- Sterile EMED hypoallergenic site dressing for each needle, enhanced with dry window
- Recessed safety feature for optimal comfort
- Gripping lip to aid needle removal from skin

## Benefits

- Facilitates 90 degree insertion for reduction of pain
- Requires significantly less penetration force when inserting the needle into the skin
- Maximizes patient comfort throughout the infusion
- Great for skin-sensitive patients (thin-skinned, older patients and children)

## Advantages

- Maximum flow rate meets and exceeds market requirements
- Easier needle insertion with minimization of tissue damage
- Manufactured with superior polymers to minimize chemical interactions with pharmaceuticals
- Optimized to enhance patient experience

### Single Needle Set

Item #	Description	UOM
<a href="#">OPT12606</a>	Single Needle Set, 70cm Tubing Length, 6mm Needle, 26G, Priming Volume: 0.28	(50/Box)
<a href="#">OPT12609</a>	Single Needle Set, 70cm Tubing Length, 9mm Needle, 26G, Priming Volume: 0.28	(50/Box)
<a href="#">OPT12612</a>	Single Needle Set, 70cm Tubing Length, 12mm Needle, 26G, Priming Volume: 0.28	(50/Box)
<a href="#">OPT12614</a>	Single Needle Set, 70cm Tubing Length, 14mm Needle, 26G, Priming Volume: 0.28	(50/Box)

### Bi-Furcated Needle Set

Item #	Description	UOM
<a href="#">OPT22606</a>	Bi-Furcated Needle Set, 70cm Tubing Length, 6mm Needle, 26G, Priming Volume: 0.53	(10/Box)
<a href="#">OPT22609</a>	Bi-Furcated Needle Set, 70cm Tubing Length, 9mm Needle, 26G, Priming Volume: 0.53	(10/Box)
<a href="#">OPT22612</a>	Bi-Furcated Needle Set, 70cm Tubing Length, 12mm Needle, 26G, Priming Volume: 0.53	(10/Box)
<a href="#">OPT22614</a>	Bi-Furcated Needle Set, 70cm Tubing Length, 14mm Needle, 26G, Priming Volume: 0.53	(10/Box)

### Tri-Furcated Needle Set

Item #	Description	UOM
<a href="#">OPT32609</a>	Tri-Furcated Needle Set, 70cm Tubing Length, 9mm Needle, 26G, Priming Volume: 0.77	(10/Box)
<a href="#">OPT32614</a>	Tri-Furcated Needle Set, 70cm Tubing Length, 14mm Needle, 26G, Priming Volume: 0.77	(10/Box)

### Quad-Furcated Needle Set

Item #	Description	UOM
<a href="#">OPT42609</a>	Quad-Furcated Needle Set, 70cm Tubing Length, 9mm Needle, 26G, Priming Volume: 1.00	(10/Box)
<a href="#">OPT42612</a>	Quad-Furcated Needle Set, 70cm Tubing Length, 12mm Needle, 26G, Priming Volume: 1.00	(10/Box)
<a href="#">OPT42614</a>	Quad-Furcated Needle Set, 70cm Tubing Length, 14mm Needle, 26G, Priming Volume: 1.00	(10/Box)

### Penta-Furcated Needle Set

Item #	Description	UOM
<a href="#">OPT52609</a>	Penta-Furcated Needle Set, 70cm Tubing Length, 9mm Needle, 26G, Priming Volume: 1.24	(10/Box)

### Hexa-Furcated Needle Set

Item #	Description	UOM
<a href="#">OPT62609</a>	Hexa-Furcated Needle Set, 70cm Tubing Length, 9mm Needle, 26G, Priming Volume: 1.47	(10/Box)

### Small Enough to Care, Big Enough to Deliver

If you are not a current customer of IMS and would like more information on EMED OPTFlow® Administration Sets, please contact us at (800) 755-3800. We will connect you with an IMS Sales Representative in your area.

**Corporate Office**  
 594 Territorial Drive, Unit A  
 Bolingbrook, IL 60440  
 1.800.755.3800  
[www.integratedmedsys.com](http://www.integratedmedsys.com)