



Out of sight. Out of mind

AccuSert® is a reusable insertion device for subcutaneous needles specifically designed to be compatible with OPTflow®. AccuSert received FDA clearance in 2023.



Advantages

AccuSert® is designed to simplify insertion and conceal the needle from view, helping patients manage self-insertion confidently, thereby facilitating patient self-insertion of subcutaneous needles.





Ergonomic design.
Simple & Intuitive to use



Supports OPTFlow sets, various infusion therapies, and allows flexibility depending on patient preference



Reusable. No need to repurchase after each use



Visual and Audible Feedback



Reliable access to subcutaneous tissue for optimized access and drug absorption



Cost effective. Low cost per infusion



Protects users from accidental needle sticks and improves device safety compliance



Streamlines multi-site infusions

Features

Mechanical Needle Inserter. Spring-loaded (1–3 Newtons)

Compatibility With OPTFlow

Insertion Angle: 90°

Benefits

Ensures optimal medication absorption

Enhances comfort and control, making self-administration easier

Reusable, low cost per infusion, and validated for two years of SCIg use

Advantages

Saves time for high volume treatments

Enhances patient confidence. Easier setup encourages therapy adherence

Confirms successful loading and insertion with a "click" noise, reducing uncertainty

Ordering Information

For orders and inquiries, contact your local EMED representative or EMED's customer service department at 916-932-0071. Orders can also be placed via email at orders@emedtc.com.

EMED Reorder Number

Item Description

Qty Per Box

FP-0010033

AccuSert®

1