

EZ-FLOW Elastomeric Pumps are a convenient alternative to electronic infusion pumps. The EZ-FLOW offers a safe and effective way to deliver medications while allowing patients freedom in their day-to-day lives without the worry of programming and battery power.

FEATURES & BENEFITS

- No battery or electronics required enabling patient freedom
- Sterile, single-use design with tethered end cap and pinch clamp for controlled, reliable infusion
- Long duration and chemotherapy pumps include a priming cap, maximizing efficiency
- Color-coded flow rates enable quick identification and help reduce the risk of administration errors
- USFDA certified for a wide range of applications including infusion and pain management
- Not made from natural rubber latex
- DEHP-free



DESIGNED FOR COMFORT. BUILT FOR LIVING.

The EZ-FLOW elastomeric pump is designed with patient comfort at the forefront, offering a lightweight, quiet alternative to electronic infusion devices. Its portable design allows patients to comfortably carry the pump in a backpack or pouch, keeping therapy discreet and unobtrusive. Without batteries, alarms, or programming, elastomerics support a more natural day-to-day experience during treatment. This freedom of movement helps patients stay active, comfortable, and engaged in daily life while receiving reliable, continuous infusion therapy.



ORDER DETAILS: EZ-FLOW ELASTOMERICS

ITEM #	FILL VOLUME	COLOR CODE	FLOW RATE	UOM (CASE)
6-PACK EZ-FLOW ELASTOMERICS				
IM100100S-6	100mL		100 mL/Hr	48/Case
IM100200S-6	100mL		200 mL/Hr	48/Case
IM270010L-6	270mL		10 mL/Hr	48/Case
SHORT DURATION				
IM050050S	50mL		50 mL/Hr	24/Case
IM100050S	100mL		50 mL/Hr	24/Case
IM100100S	100mL		100 mL/Hr	24/Case
IM100200S	100mL		200 mL/Hr	24/Case
IM200100S	200mL		100 mL/Hr	24/Case
IM200200S	200mL		200 mL/Hr	24/Case
IM250050S	250mL		50 mL/Hr	24/Case
IM250100S	250mL		100 mL/Hr	24/Case
IM250125S	250mL		125 mL/Hr	24/Case
IM250175S	250mL		175 mL/Hr	24/Case
IM250250S	250mL		250 mL/Hr	24/Case
IM250500S	250mL		500 mL/Hr	24/Case
IM270175S	270mL		175 mL/Hr	24/Case
IM400100S	400mL		100 mL/Hr	24/Case
IM400200S	400mL		200 mL/Hr	24/Case
IM500250S	500mL		250 mL/Hr	24/Case

BUILT FOR WORKFLOW. BUNDLED FOR VALUE.

IMed EZ-FLOW Elastomerics are available in convenient 6-packs designed to simplify purchasing and inventory management. The multi-pack format reduces ordering frequency, minimizes stockouts, and delivers better overall value for busy infusion settings. With ready-to-use pumps bundled together, teams save time, reduce handling, and keep care moving efficiently.



ORDER DETAILS: EZ-FLOW ELASTOMERICS

ITEM #	FILL VOLUME	COLOR CODE	FLOW RATE	UOM (CASE)
LONG DURATION				
IM060005L	60mL		5 mL/Hr	24/Case
IM100002L	100mL		2 mL/Hr	24/Case
IM100005L	100mL		5 mL/Hr	24/Case
IM125005L	125mL		5 mL/Hr	24/Case
IM270001L	270mL		1 mL/Hr	24/Case
IM270002L	270mL		2 mL/Hr	24/Case
IM270004L	270mL		4 mL/Hr	24/Case
IM270005L	270mL		5 mL/Hr	24/Case
IM270010L	270mL		10 mL/Hr	24/Case
IM500020L	500mL		20 mL/Hr	24/Case
CHEMOTHERAPY				
IM100002C	100mL		2 mL/Hr	12/Case
IM125005C	125mL		5 mL/Hr	12/Case
IM270002C	270mL		2 mL/Hr	12/Case
IM270005C	270mL		5 mL/Hr	12/Case
IM270010C	270mL		10 mL/Hr	12/Case
IM300006C	300mL		6 mL/Hr	12/Case

PRIMING MADE PRECISE

The IMed EZ-FLOW Elastomeric priming cap (available only on Long Duration and Chemotherapy pumps) is designed to simplify and standardize pump preparation for long-duration and chemotherapy infusions. Its integrated design supports efficient, controlled priming while helping reduce setup time and variability across users. By streamlining the priming process, EZ-FLOW pumps enhance workflow efficiency and support reliable therapy delivery—so clinicians can focus on patient care, not pump setup.



Strengthened by Nationwide Distribution and Dedicated Customer Support

Proudly
Distributed by

IMS
INTEGRATED MEDICAL SYSTEMS



IMed[®]
Products

IMed Products delivers reliable, high-quality medical support solutions at value-driven pricing—designed to ensure consistent supply, trusted clinical performance, and operational ease for outpatient centers, clinics, hospitals, and pharmacies.

The product details provided herein, including any images or supplements, reflect information from various sources regarding the product; this information has not been independently verified by Integrated Medical Systems ("IMS") or its affiliates. IMS is not liable for any inaccuracies or omissions in the product information, which is provided solely for your convenience.

The information presented here is subject to change without prior notice. Customers should review product packaging and labeling before use or purchase to ensure accuracy, especially when referencing information from IMS's online or printed materials. IMS disclaims any responsibility for the accuracy or completeness of the information reflected herein, regardless of whether any inaccuracies or omissions are due to errors or other factors.

© 2026 IMed® Products. All rights reserved. 02/26 IMS25095 IMed Sell Sheets_Brand Refresh EZ-FLOW 6-Pack_01