

NOW AVAILABLE IN
6-PACKS
to maximize convenience

EZ FLOW™

ELASTOMERIC PUMP

THE EFFICIENT SOLUTION FOR YOUR CLEANROOM SPACE

EZ-FLOW™ Elastomeric Pumps are a proven, safe, and effective alternative to electronic infusion pumps. Now, with the NEW 6-pack format, we offer enhanced efficiency and convenience—perfect for cleanrooms and clinical settings.

This new packaging provides a streamlined solution for better inventory management and reduced handling time.

Patient Freedom

No battery or electronics required, leaving patients with freedom to live their day-to-day lives.

Sleek Design

These sterile, single-use pumps include a tethered end cap and pinch clamp to aid in effective infusion.

Proven Safe

USFDA certifies these pumps safe for a wide range of applications, including infusion and pain management.

Color-Coded Flow Rates

Easily identify flow rate, reducing the risk of administration error with our color-coded labels.

Priming Cap Efficiencies

Long duration pumps include a priming cap, which reduces the time to prime and maximizes efficiencies.

NEW! Bulk Convenience

The NEW 6-pack format streamlines high-volume workflows, reducing handling time and maximizing efficiency in your cleanroom environment.

EZ-FLOW™ ELASTOMERICS

NEW! 6-PACK EZ-FLOW™ ELASTOMERICS

Item #	Fill Volume	Flow Rate	Color Code	UOM
IM100100S-6	100mL	100 mL/Hr		48/Case
IM100200S-6	100mL	200 mL/Hr		48/Case
IM270010L-6	270mL	10 mL/Hr		48/Case

SHORT DURATION

IM050050S	50mL	50 mL/Hr		24/Case
IM100050S	100mL	50 mL/Hr		24/Case
IM100100S	100mL	100 mL/Hr		24/Case
IM100200S	100mL	200 mL/Hr		24/Case
IM200100S	200mL	100 mL/Hr		24/Case
IM200200S	200mL	200 mL/Hr		24/Case
IM250050S	250mL	50 mL/Hr		24/Case
IM250100S	250mL	100 mL/Hr		24/Case
IM250125S	250mL	125 mL/Hr		24/Case
IM250175S	250mL	175 mL/Hr		24/Case
IM250250S	250mL	250 mL/Hr		24/Case
IM250500S	250mL	500 mL/Hr		24/Case
IM270175S	270mL	175 mL/Hr		24/Case
IM400100S	400mL	100 mL/Hr		24/Case
IM400200S	400mL	200 mL/Hr		24/Case
IM500250S	500mL	250 mL/Hr		24/Case

LONG DURATION

IM060005L	60mL	5 mL/Hr		24/Case
IM100002L	100mL	2 mL/Hr		24/Case
IM100005L	100mL	5 mL/Hr		24/Case
IM125005L	125mL	5 mL/Hr		24/Case
IM270001L	270mL	1 mL/Hr		24/Case
IM270002L	270mL	2 mL/Hr		24/Case
IM270004L	270mL	4 mL/Hr		24/Case
IM270005L	270mL	5 mL/Hr		24/Case
IM270010L	270mL	10 mL/Hr		24/Case
IM500020L	500mL	20 mL/Hr		24/Case

CHEMOTHERAPY

IM100002C	100mL	2 mL/Hr		12/Case
IM125005C	125mL	5 mL/Hr		12/Case
IM270002C	270mL	2 mL/Hr		12/Case
IM270005C	270mL	5 mL/Hr		12/Case
IM270010C	270mL	10 mL/Hr		12/Case
IM300006C	300mL	6 mL/Hr		12/Case



Leading the healthcare industry towards a more efficient and effective future through our commitment to excellence in developing cutting-edge medical products and solutions.

Proudly Distributed by:



To order, contact your IMS Territory Manager or Customer Service at 800.755.3800

The product details provided herein, including any images or supplements, reflect information from various sources regarding the product; this information has not been independently verified by Integrated Medical Systems ("IMS") or its affiliates. IMS is not liable for any inaccuracies or omissions in the product information, which is provided solely for your convenience.

The information presented here is subject to change without prior notice. Customers should review product packaging and labeling before use or purchase to ensure accuracy, especially when referencing information from IMS's online or printed materials. IMS disclaims any responsibility for the accuracy or completeness of the information reflected herein, regardless of whether any inaccuracies or omissions are due to errors or other factors.