

eitan  
medical  
sapphire

# sapphire™

Infusion Platform  
US Catalog 2021



care & beyond™

## Contents

●	Infusion Pump Kits	4
●	Accessories	8
●	Administration Sets	12
●	PC Applications	15

**Sapphire™ Multi-therapy infusion pump kit | SW Rev 15.00.0**

**Sapphire™ Epidural infusion pump kit | SW Rev 15.00.0**



**Includes:**

- Sapphire™ Multi-therapy pump
- Power supply
- Mini cradle
- Mini cradle knob key
- User manual

**Power supply type:**

- Fixed head power supply Kit P/N 17000-028-0079
- Integrated power supply & AC Power Cord – Type B Kit P/N 17000-028-0077
- Interchangeable head power supply Kit P/N 17000-028-0078

**Includes:**

- Sapphire™ Epidural pump
- Power supply
- Bolus handle

- Mini cradle with splitter
- Mini cradle knob key
- User manual

Kit P/N 17000-031-0035

Sapphire™ Multi-therapy infusion pump kit | SW Rev 13.23.2



**Includes:**

- Sapphire™ Multi-therapy pump
- Power supply
- Mini cradle
- Mini cradle knob key
- User manual

Kit P/N 17000-028-0061

Sapphire™ Epidural infusion pump kit | SW Rev 13.23.2



**Includes:**

- Sapphire™ Epidural pump
- Power supply
- Bolus handle
- Mini cradle with splitter
- Mini cradle knob key
- User manual

Kit P/N 17000-031-0027

**Sapphire™**  
Mini Cradle

**Sapphire™**  
Mini Cradle

**With Splitter**

P/N 17000-026-0007

**Without Splitter**

P/N 15042-000-0010



**Sapphire™**  
Multi Pump Mounting System



P/N 17000-020-0008

**Sapphire™**  
Splitter



P/N 15043-000-0005

**Sapphire™**  
Integrated Power Supply



P/N 15072-000-0010

**Sapphire™**  
Power Supply



P/N 15072-000-0007

**Sapphire™**  
External Battery Pack



P/N 15079-111-0001

**Sapphire™**  
AC Power Cord - Type B



P/N 15070-112-0005

**Sapphire™**  
Power Cord - compatible  
with IPS



P/N 15171-112-0002

Accessories | Bags

**Sapphire™**  
Infusion Pouch 500 mL



P/N 05080-605-0001

**Sapphire™**  
Travel Case



P/N 15020-010-0002

**Sapphire™**  
Homecare Large Backpack



P/N 15088-000-0001

Accessories | Pain Management

**Sapphire™**  
PCA lockbox 250mL



P/N 15063-000-0001

**Sapphire™**  
Lockbox 250mL strap and pouch



P/N 15063-000-0004

**Sapphire™**  
PCA lockbox 100mL








P/N 15061-131-0004

**Sapphire™**  
Bolus handle







P/N 15071-000-0005






## Administration Sets | Standard Sets

NAME	P/N	DESCRIPTION	LENGTH (CM)	PRIME VOLUME	UNITS/BOX
AP 403-01	12003-000-0013	Microbore infusion set with non-vented spike	290 cm 114 in	~6 mL	30
					
*AP 409-01	12005-000-0003	Macro bore primary infusion set with vented/non-vented drip chamber, back check valve, roller clamp and 2 needleless Y-site	277 cm 109 in	~23 mL	30
					
* Alternative set available, BP 450-01, P/N - 12005-000-0010, Units /Box: 20					
AP 423-01	12003-000-0023	Microbore infusion set with non-vented spike and back check valve	292 cm 115 in	~7 mL	30
					
AP 433-01	12000-000-0044	Macro bore infusion set with vented/non-vented drip chamber, roller clamp and needleless Y-site	286 cm 113 in	~22 mL	20
					
AP 416-01	12003-000-0012	Microbore infusion set with female luer lock, slide clamp and needleless Y-site	237 cm 93.3 in	~5 mL	30
					

## Administration Sets | Filter Sets

NAME	P/N	DESCRIPTION	LENGTH (CM)	PRIME VOLUME	UNITS/BOX
AP 204-01	12003-000-0070	Microbore infusion set with non-vented spike, 0.2µ filter and back check valve	295 cm 116 in	~8 mL	30
					
AP 206-01	12000-000-0022	Macro bore infusion set with non-vented spike, 1.2µ filter	295 cm 116 in	~23 mL	20
					
AP 214-01	12004-000-0002	Microbore yellow marked infusion set with vented/non-vented spike and 0.2µ filter	295 cm 116 in	~8 mL	30
					
AP 240-01	12003-000-0069	Microbore infusion set with non-vented spike, 1.2µ filter and back check valve	295 cm 116 in	~8 mL	20
					

## Administration Sets | Specialty Sets

NAME	P/N	DESCRIPTION	LENGTH (CM)	PRIME VOLUME	UNITS/BOX
AP 133-01	12006-000-0010	Macrobore blood infusion set with 2 non-vented spikes, 200µ filter, 3 roller clamps and needleless Y-site	292 cm 115 in	30 mL	20
					
AP 402-01	12004-000-0004	Microbore yellow marked infusion set with non-vented spike	290 cm 114 in	~6 mL	30
					
AP 410-01	12005-000-0007	Macrobore secondary infusion set with vented/non-vented drip chamber. Includes an infusion bag hanger	71 cm 28 in	~9 mL	30
					
AP 424-01	12003-000-0024	Microbore infusion set with non-vented spike, side arm with female luer lock and back check valve	294 cm 116 in	~7 mL	30
					
AP 504-01	12000-000-0045	Macrobore UV-Resistant infusion set with vented/non-vented drip chamber, roller clamp and needleless Y-site	287 cm 113 in	~26 mL	20
					

## PC Applications

The Sapphire modern, software based platform expands beyond the pump to the PC applications suite that complements the pumps;



### Pump loader

Enables upgrading or reloading software to the Sapphire pump via a wire connection.



### Drug library

Enables to set configuration settings, drug specific limit and more. When complete, the library can be uploaded to multiple Sapphire pumps.



### Fastest pm

(SW + disposable kit [16024-000-0002])  
 FasTest PM is Preventative Maintenance Servicing Tool for Sapphire Infusion Pump platform. FasTest PM allows end-users to quickly and economically conduct full maintenance tests, including occlusion detection, air detection and pump accuracy checks. FasTest PM provides end-users with an immediate, downloadable authentic manufacturer. The tool is designed to significantly improve maintenance costs and runtime.



### Annual certification kit

(SW + disposable kit [15075-020-0001])  
 Enables to run performance test (air detection, pressure detection and accuracy). Results are compared to manufacturing database per pump.







### Event viewer

Enables downloading and viewing the event log of the Sapphire infusion pump on a computer.














## PC Applications | Accessories

PC APPLICATION	ITEM	P/N	PICTURE
FasTest PM	FasTest PM kit	16024-000-0002	
Annual Certification Test	Annual certification kit	15075-020-0001	
All	Communication cable	05020-110-0213	
All	RS232 to USB 4 port adapter	15077-000-0001	

## Accessories | Spare Parts

SPARE PART FOR	ITEM	P/N	PICTURE
Sapphire pump	Li-Ion battery	05020-160-0001	
Mini cradle	connection knob flange Key	05040-380-0019	
Mini cradle	connection knob key	05040-380-0018	
Mini cradle	MC- AVP Connection Bolt Knob Tip	05040-380-0027	
Mini cradle	MC- AVP Press Plane	05040-380-0026	
PCA lockbox 100	Key	05070-190-0015	
PCA lockbox 250	Key	10010-122-0046	
PCA lockbox 250	Release clip (eliminates lockbox locking onto cradle)	15063-121-0029	
PCA Lockbox 500	Key	15062-132-0001	

REF 10050-003-0009 | Rev.05 / 01.2021

**U.S. Contact information:**

[customerservice@eitanmedical.com](mailto:customerservice@eitanmedical.com)

877-541-9944 | 949-303-7845

949-288-5482

**Eitan Group North America**

65 Enterprise 3rd Floor Aliso Viejo, CA 92656, USA

**Technical Service Contact Information:**

[service@eitanmedical.com](mailto:service@eitanmedical.com)



Q Core Medical Ltd. is the manufacturer of Sapphire products.