

Meeting the Evolving Needs of Specialty Infusion



Delivering cost-effective care whilst optimizing patient safety, patient outcomes and patient experience is at the heart of specialty infusion management.

With the trend of patients moving away from acute care forecast to continue, how can Specialty Infusion Providers ensure that they can effectively meet the needs of this growing population amid current increasing pressures?

The Patient-Care Revolution Starts Here

Introducing Avoset™ and Eitan Insights™ – the connected infusion management platform designed to transform specialty infusion, improve patient safety, enhance quality of life and support proactive data-driven care.



Key Specifications

| | | | | | |
|-------------------------------|---|--------------------------------|--------------------------------------|-----------------------------|---|
| Dimensions (W x H x D) | 97 x 97 x 47 mm / 3.8 x 3.8 x 1.85 inches | Infusion delivery modes | Continuous, PCA, Intermittent, Taper | Downstream occlusion | Low (1.6 Bar/23.2 PSI) Medium (1.2 Bar/17.4 PSI) High (0.5 Bar /7.25 PSI) |
| Accuracy | ± 5% | Priming | Automatic or manual | Pumping mechanism | Volumetric |
| Weight | 450 g / 15.9 oz including batteries | Volume (VTBI) | 0.1 – 9,999 mL | Connectivity | Bluetooth®, NFC |
| Battery type | Rechargeable or alkaline 3 AA | Flow rate | 0.1 – 300 mL/h | Ingress protection | IP34 |
| Battery life | Up to 280 hrs (1 mL/h) | KVO rate | 0.1 – 10 mL/h | | |



About Eitan Medical

Patient safety and care are only the starting point, as Eitan Medical goes beyond— delivering connected, intuitive drug delivery and infusion solutions designed to improve patient quality of life across the continuum of care. For over a decade, Eitan Medical has provided solutions that meet evolving drug delivery needs.

Contact information

Eitan Medical UK
customersupport.uk@eitanmedical.com
+44.1865.504289

Eitan Medical North America
customerservice@eitanmedical.com
+1.877.541.9944

15132-041-1048-MMTL | Rev.03 | 07.2023



www.eitanmedical.com



avoset™
Connected Infusion System



Designed to meet the needs of Specialty Infusion



- Compact, lightweight device with long battery life allows patient mobility and may improve quality of life
- Versatile ambulatory pump with several modes: Continuous, Intermittent, PCA and Taper
- Simple web-based PC programming allows programming precision and supports medication safety
- Delivering from a variety of medication containers, Avoset™ allows a streamlined pump inventory and effortless clinician onboarding
- Effortless set up with easy one-click administration set designed to reduce risk of error
- Connected by design to support informed and proactive care, which can lead to improved patient outcomes
- Reliable and secure design to minimize maintenance and support cost savings



Effortless to Set Up Simple to Operate



Connected by Design to Support Proactive Care

Smart, Secure Connectivity

Eitan Insights™ is built on cutting-edge cloud architecture that utilizes the latest security features and best practices with state-of-the-art encryption, multi-factor authentication, and is HIPAA and GDPR compliant.

ISO 27001:2013 • HIPAA Compliant • GDPR Compliant



Advanced Pump Programming


Our web-based application supports PC programming, thus reducing potential for user error and allows the patient to receive a pump that is ready to go!

Avoset™ pumps are within reach, transmitting infusion data through Bluetooth®, technology and designated mobile app to Eitan Insights™.



Transforming Specialty Infusion

Eitan Insights™ can help optimize revenue through accurate billing, auto-documentation, and can support drive operational efficiency with GPS location capability. Locate and track your infusion pump inventory and help avoid costly capital equipment loss.



Data-Driven Insights Help to Improve Patient Outcomes

Eitan Insights™ provides access to remote infusion data, which can support fast patient triage, and drive proactive care planning - supporting patients **when they need it most.**

Eitan Insights™ eliminates subjectivity that can help avoid anxiety-provoking troubleshooting, thus helping improve both patient and care-giver experience.

