

Avoset™ Catalog

Redefining
Specialty Infusion



eitan
medical
avoset™



The connected Avoset™ Infusion System is designed to transform specialty infusion with compact and simplified technology that can monitor infusion treatment data remotely, aiming to enhance patient safety, and improve the user experience.

Together, the Avoset™ Infusion Pump, and Eitan Insights™, our cloud-based infusion management system, support proactive care planning and remote troubleshooting.



avoset™
Infusion Pump

P/N 17000-036-0004




Medication Reservoirs & Administration Sets



Avoset Medication Reservoir – 100mL

100 ml medication reservoir with free-flow protection, slide clamp and female luer lock. Nonvented stopper included.




Length: ~330 mm (tube) **Priming volume:** ~1 mL **Units per box:** 10

P/N 12222-000-0002

Avoset Microbore Extension Set

Microbore infusion set with spin male luer lock, slide clamp and one-way check valve with male luer lock.




Length: ~1,687 mm **Priming volume:** ~3 mL **Units per box:** 20

P/N 12222-000-0016

Avoset Microbore Extension Set with 0.2 Micron Filter

Microbore infusion set with spin male luer lock, slide clamp, vented 0.2 µm filter and one-way check valve with male luer lock.





Length: ~1,697 mm **Priming volume:** ~3 mL **Units per box:** 20

P/N 12222-000-1000

Avoset Microbore Set with Backcheck Valve

Microbore infusion set with non-vented bag spike, administration cassette with free-flow protection, slide clamp and one-way check valve with male luer lock.





Length: ~2,911 mm **Priming volume:** ~6 mL **Units per box:** 20

P/N 12222-000-0001



Avoset Microbore Set with 0.2 Micron Filter

Microbore infusion set with non-vented bag spike, administration cassette with free-flow protection, slide clamp, 0.2 µm air elimination filter and one-way check valve with male luer lock.




Length: ~2,951 mm **Priming volume:** ~8 mL **Units per box:** 20

P/N 12222-000-0013

Avoset Microbore Set with 1.2 Micron Filter

Microbore infusion set with non-vented bag spike, administration cassette with free-flow protection, slide clamp, 1.2 µm air elimination filter and one-way check valve with male luer lock.




Length: ~2,916 mm **Priming volume:** ~8 mL **Units per box:** 20

P/N 12222-000-0010

Avoset Microbore Set with Female Luer

Microbore infusion set with female luer lock, administration cassette with free-flow protection, slide clamp and one-way check valve with male luer lock.



Length: ~400 mm **Priming volume:** ~1 mL **Units per box:** 20

P/N 12222-000-0004



AvosetPad

The AvosetPad uses Bluetooth® and near-field communication to connect the pump to a PC. It enables wireless connection between the Avoset Infusion Pump and the web-based Avoset Programming Tool.




Avoset Cradle

The Avoset Cradle offers flexible positioning of the Avoset Infusion Pump at any angle or height on a pole or bed. It can be attached vertically or horizontally and can be rotated to keep the pump display conveniently oriented.

Avoset Cassette Lock

Designed for safety, the Avoset Cassette Lock is a transparent secure cover that prevents unauthorised users from disconnecting the medication reservoir from the Avoset Infusion Pump.





Our cloud-based software platform, Eitan Insights™, transforms specialty infusion with a web-based programming tool and a mobile app that provides remote visibility of treatment data.



For more information about Eitan Insights™, please contact your sales representative.

15132-041-1076-MMTL | Rev.01 | 09.2023



Eitan Medical
877.541.9944
customerservice@eitanmedical.com



www.eitanmedical.com