

Agilia[®] Connect

Volumetric Pump

Part of the **Agilia Connect** Infusion System



Your patient-side ally

Designed with usability engineering and smart technology in mind, and manufactured to be both durable and serviceable, the Agilia® Connect Volumetric Pump fits easily into every user's everyday workflow. It's your ally in safely and efficiently caring for patients in your healthcare facility.

Clinically Focused on patient safety

In addition to the safety features of the Vigilant Master Med drug library, the Agilia Connect Volumetric Pump has built-in features and functions that help to keep safety front and center for your team.

Technically Packed with powerful technology

Though small in size, the Agilia Connect Volumetric Pump enhanced microprocessor and memory are big on capacity – able to handle the complex libraries of Vigilant Master Med and regular data transfer to Vigilant Insight.

Functionally Refined for intuitive operation

Human factors design has been applied to every aspect of the Agilia Connect Volumetric Pump – from form to function – to help ease your team's training and everyday use.

Practically Engineered for reliable performance

The Agilia Connect Volumetric Pump is built to perform with accuracy and to maintain and repair with efficiency.

Lightweight.

Less than

5 lbs

Flexible.

Infusion rate from

0.1 to
1500 mL/h

Reliable.

3 year

Preventative
Maintenance cycle

"Focus-friendly", incremental keypad helps avoid key bounce and programming errors

Visually-cued elements like color-coded parts and pictogram button labels make set-up speedy

Auto-restart minimizes patient-side occlusion alarms caused by momentary kinking of tubing, IV pushes, etc.

Ramp, taper and sequential functionality safely and efficiently deliver medications

Stackable design occupies a smaller footprint making it easier to transport patient

Swiveling pole clamp attaches easily to IV pole or bed rail

Wide range of parameters

- for flow rates, infusion volumes, modes and times, as well as pressure modes
- can be programmed to launch as automatic defaults

Vigilant Master Med dose error reduction software

- can define soft & hard limits for each drug.
- establish pressure & air-in-line settings per each drug to help minimize nuisance alarms
- create a data set with up to 19 profiles that may include up to 3,800 different drug/therapies entries
- contains 10,000 drugs/therapies, 50 different configurations, drug categories, drug libraries, profiles & data sets

3-year preventative maintenance cycle helps to minimize maintenance activities

20-minute diagnostic test using Agilia Partner, helps to optimize your resources.

High-performance battery

- 3-year replacement cycle
- long-lasting between charges:
 - Wireless: at 25 mL/h > 5 hrs and at 1500 mL/h > 4 hrs;
 - Non-wireless: at 25 mL/h > 8 hrs; at 1500 mL/h > 5 hrs

Rugged ABS shell provides shock-protection

Removable universal power cord for the entire Agilia Connect family enhances patient mobility and avoids tripping hazard





Dedicated Experts

Account Manager
Project Manager
Clinical Pharmacy Services Liaison
Clinical Nurse Liaison
IT Service Liaison

SCP certification, 10 years running

Our technical team's commitment to annual certification by the global standard in field service quality – Service Capability & Performance (SCP) – whose audits cover a broad scope of practices, ensures installation and training, as well as ongoing service, are provided by individuals with deep knowledge.

Digital training tools

Learning Management System resources bring Fresenius Kabi subject matter expertise to your teams, anytime, anywhere.

Exceptionally Supported to drive best practices

From initial implementation to everyday application, Fresenius Kabi clinical and technical experts work alongside your teams to put Agilia® Connect Volumetric Pump to best use, getting everyone started quickly and ready to reach your infusion practice goals.

Globally Recognized as a proven partner

Fresenius Kabi's medical device division is a leader throughout the world, with more than 35 years of infusion technology experience and an installed fleet of more than a million infusion pumps. Fresenius Kabi is also a steadily-growing force in the North American generic injectables space, having worked shoulder-to-shoulder with acute- and alternate-care institutions over the last decade to help ensure the availability and affordability of essential IV drugs for anesthesia and analgesia, critical care, and oncology.

SCP is a registered trademark of Service Strategies.



**FRESENIUS
KABI**

caring for life

Fresenius Kabi, USA, LLC
Three Corporate Drive
Lake Zurich, Illinois 60047
www.fresenius-kabi.us/agiliconnect

©2022 Fresenius Kabi, USA, LLC. All Rights Reserved.
2574-AG-10-03/22