



New **TapTight™**
Locking Lid Design!



HARDY DIAGNOSTICS

Contact Plates



HARDY DIAGNOSTICS

Contact Plates

Hardy Diagnostics' contact plates are available in a variety of formats and formulas to meet your needs for surface monitoring.

Irradiated, triple-bagged plates are available when an added level of sterility assurance is needed.

All contact plates feature the new TapTight™ locking lid to secure the lid after inoculation.



800.266.2222
805.614.9274 (fax)
sales@hardydiagnostics.com
www.hardydiagnostics.com

Santa Maria, CA Phoenix, AZ Salt Lake City, UT Tacoma, WA Lima, OH



Contact Plates

Features & Benefits

New TapTight™ Locking Lids

Simply tap the contact plate lid to secure it into place.

Molded Grid

Each plate has a grid molded into the bottom of the plate to aid in the counting of microorganisms.

Certificates of Analysis

Available with each order!

Many formats

Hardy Diagnostics contact plates are available in many formats, including triple-bagged, and irradiated.

SabDex with Lecithin and Tween, triple-bagged and irradiated.



TapTight™ Locking Lids



1. Place lid on contact plate.
2. Tap lid into place.
3. Lid is secure.

Contact Plates

Product/Quantity	Cat. no.
TSA with Lecithin & Tween 80, 10/bag	P34
TSA with Lecithin & Tween 80, triple-bagged, irradiated, 20/box	P74
SabDex Agar, 10/bag	P36
SabDex with Lecithin & Tween 80, triple-bagged, irradiated, 20/box	P75
Blood Agar 5%, 10/bag	P33
Baird Parker Agar, 10/bag	P16
Eosin Methylene Blue (EMB), 10/bag	P09
Mannitol Salt Agar, 10/bag	P98
Rose Bengal Agar, 10/bag	P42
R2A Agar, 10/bag	P32
D/E Neutralizing Agar, 10/bag	P99
D/E Neutralizing Agar, triple-bagged, irradiated, 20/box	P13